## CERTEST ADENOVIRUS POSITIVE CONTROL

### INTENDED USE

The CerTest Adenovirus Positive Control is an external quality control for CerTest Adenovirus Card test and Adenovirus test line of CerTest Rotavirus-Adenovirus card test.  
Use this control swab to test that the extraction reagents and the test are working properly.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The reagent is stable through the expiration date printed on the sealed pouch. The reagent must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not use the swab if the aluminium pouch has been damaged during storage.
- The swab should be discarded in a proper biohazard container after testing.

### MATERIALS PROVIDED

- 20 Adenovirus Positive Control swabs
- 20 disposable testing tubes
- 20 disposable plastic pipettes
- 1 instructions for use:


### MATERIALS REQUIRED BUT NO PROVIDED

- CerTest Adenovirus Card test (cat. CER A 820001F) or CerTest Rota-Adeno Card (cat. CER X 820001F)
- Sample diluent (stool collection vials included in the kit: cat. CER A820001F or cat. CER X820001F)
- Disposable gloves
- Timer


### QUALITY CONTROL TESTING PROCEDURE

Allow the tests and swab control to reach to room temperature (15-30°C) prior to testing. Do not open the package until ready to perform the assay. Only bring to room temperature the number of tests required to assay before opening it.

1. Add 15 drops (1) of Sample diluent and immediately put the Positive Control Swab into the tube.
2. Mix the solution by rotating the Positive Control Swab forcefully against the side of the tube at least 10 times. Extract as much liquid as possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
3. Remove the CerTest Adenovirus Card or CerTest Rota-Adeno Card devices from its sealed bag just before using.
4. Place the test on a flat surface. Use a separate pipette and device for each control. Dispense exactly 4 drops or 100 µL from the testing tube, into the circular window marked with an arrow (3).
5. Read the result at 10 minutes.

### INTERPRETATION OF RESULTS (please refer to the illustration below)

**Control + Adenovirus POSITIVE:** A GREEN control band (control line) in the control region (C), a BLUE band in the result region (T) appear at the same time due to presence of antigens in the Control + Swab.

**INVALID:** A total absence of the control coloured band (GREEN) regardless the appearance of the result line (BLUE) or only the control (C) coloured band (GREEN). Incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new control. If the problem persists, discontinue using the test kit and contact your local distributor.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Batch code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		For IVD performance evaluation only

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