CERTEST INFLUENZA A+B BLISTER TEST
One step Influenza A+B Blister Test for the detection of Influenza type A and type B from nasal swabs, nasal wash or nasal aspirate specimens

INTENDED USE

The CerTest Influenza A+B Blister Test is a one step coloured chromatographic immunoassay for the qualitative detection of influenza type A and type B antigens. It can be used directly with nasal swabs or nasal wash or nasal aspirate specimens.

INTRODUCTION

Influenza is caused by a virus that attacks mainly the upper respiratory tract – the nose, throat and bronchi and rarely also the lungs. The infection usually lasts for about a week. It is characterized by sudden onset of high fever, myalgia, headache and severe malaise, non-productive cough, sore throat, and rhinitis. Most people recover within one to two weeks without requiring any medical treatment. In the very young, the elderly and people suffering from medical conditions such as lung diseases, diabetes, cancer, kidney or heart problems, Influenza poses a serious risk. In these people, the infection may lead to severe complications of underlying diseases, pneumonia and death.

The currently circulating influenza viruses that cause human disease are divided into two groups: A and B. Influenza A has 2 subtypes which are important for humans: A (H3N2) and A (H1N1), of which the former is currently associated with most deaths. Influenza viruses are defined by 2 different protein components, known as antigens, on the surface of the virus. They are spike-like features called haemagglutinin (H) and neuraminidase (N) components.

PRINCIPLE OF THE TEST

The CerTest Influenza A+B Blister Test is a qualitative immunochromatographic assay for the determination of Influenza type A and type B antigens. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, to recognize these antigens.

During testing, the sample is allowed to react with the coloured conjugate (anti-type A mouse monoclonal antibodies-red microspheres and anti-type B mouse monoclonal antibodies-blue microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured particles. Different coloured lines will be visible, depending upon the type of antigens contents of the sample. These lines are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- CerTest Influenza A+B Blister Test is intended for professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not use test devices if the aluminium pouch has been damaged during storage.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new swab must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Nasal swab specimen:

Collect specimen with a sterile swab from one nostril. Insert the swab approximately 3 cm into the nostril rotating against the nasal wall. Process the swab as soon as possible after collecting the specimen.

Nasal Wash or Aspirate specimen:

For Adult: Place the irrigator up to the nose. Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side. Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity. Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.

For children: use an aspiration bulb or bulb syringe to instill the saline water into one nostril leaning the children head. Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container. Repeat for the other nostril and transfer the fluid into the same specimen container.

Samples should be processed as soon as possible after collection. The samples can be stored in the refrigerator (2-4 °C) for 8 hours prior to testing.

MATERIALS PROVIDED

- Blister tests
- Reagent A (saline solution) and Reagent B (sample diluent)
- Testing tubes
- Sterile swabs
- Disposable pipettes
- Positive Control swab
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Disposable gloves
- Timer

TEST PROCEDURE


Allow the tests, swabs and controls to reach to room temperature (15-30°C) prior to testing. Do not open the package until ready to perform the assay. Only bring to room temperature the number of tests required to assay before opening it.

- Procedure A using nasal swab samples:

1. Add 10 drops (1) Reagent A and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 10 times. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
3. Add 3 drops (3) Reagent B. Mix the solution with the pipette at least 10 times.
4. Cut the blister to obtain a test single pack, hold the non sealed side and open it peeling off the upper foil.

4.a. Don't remove the test from the blister cavity and use it as soon as possible. Place the test on a flat surface. Use a separate pipette and device for each sample or control. Dispense exactly 5 drops or 150µL on the white end of the test (4a). Read the result at 10 minutes.

4.b. Place the test vertically with the white end submerged into the sample taking care of not surpassing the limit of immersion indicated with the arrows (4b). Read the result at 10 minutes.


- Procedure B using nasal wash or aspirate samples:


1. Add 6 drops (1) of the nasal wash or aspirate samples with a pipette and 3 drops (2) of Reagent B in a testing tube.
2. Mix the solution with the pipette at least 10 times. Best results are obtained when the specimen is vigorously extracted in the solution (3).
3. Cut the blister to obtain a test single pack, hold the non sealed side and open it peeling off the upper foil.

3.a. Don't remove the test from the blister cavity and use it as soon as possible. Place the test on a flat surface. Use a separate pipette and device for each sample or control. Dispense exactly 5 drops or 150µL on the white end of the test (3a). Read the result at 10 minutes.

3.b. Place the test vertically with the white end submerged into the sample taking care of not surpassing the limit of immersion indicated with the arrows (3b). Read the result at 10 minutes.

INTERPRETATION OF RESULTS (Please refer to the illustration below)

NEGATIVE: Only one GREEN band (control line) appears in the white central zone of the reaction test (control region).

Influenza A POSITIVE: In addition to the GREEN control band, a distinguishable RED band (type A test line) also appears in the white central zone of the reaction test (tests lines).

Influenza B POSITIVE: In addition to the GREEN control band, a BLUE band (type B test line) also appears in the result region (tests lines).

Influenza A+B POSITIVE: All the lines above described (a GREEN control band in the control region, a RED band and a BLUE band in the result region) could appear at the same time during the test performance due to a simultaneous infection of Influenza A+B.

INVALID: A total absence of the control coloured band (GREEN) regardless the appearance or not of the tests lines (RED/BLUE). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red or blue coloured band in the result line region will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal and external procedural controls are included in the kit:

- A green line appearing in the control region is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.
- Each kit also contains a positive control swab. Use this control to test that the extraction reagents and the test are working properly.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed pack.
2. This test provides a presumptive diagnosis for Influenza A+B infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
3. The CerTest Influenza A+B Test should be used only with nasal swabs, nasal wash and nasal aspirate samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasal specimens must be obtained.
4. A negative result may be obtained if the specimen is inadequate or antigen concentration is below the sensitivity of the test. Therefore, it is recommended that all negative CerTest Influenza A+B Test results undergo confirmatory testing using other method and/or by cell culture.

PERFORMANCE

Sensitivity

Different virus extract dilutions were tested directly in the sample diluent or spiked in a negative nasal specimen in accordance with the kit instructions.

The detection of Influenza type A and/or type B with CerTest Influenza A+B Test showed >99% of sensitivity compared with another commercial rapid test.

Specificity

The use of mouse monoclonal antibodies in the elaboration of CerTest Influenza A+B Test assures high degree of specificity for the detection of Influenza type A and type B antigens.

The detection of Influenza type A and type B with CerTest Influenza A+B Test showed >99% of specificity compared with another commercial rapid test.

REFERENCES

1. Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. Journal of Clinical Microbiology, Mar. 2000, p. 1161-1165 Vol. 38, No. 3.
2. Adriana Weinberg and Miranda L. Walker. Evaluation of Three Immunoassay Kits for Rapid Detection of Influenza Virus A and B. Clinical And Diagnostic Laboratory Immunology, Mar. 2005, p. 367-370 Vol. 12, No. 3.
3. Andreea C. Cazacu, Gail J. Demmler et al. Comparison of a New Lateral-Flow Chromatographic Membrane Immunoassay to Viral Culture for Rapid Detection and Differentiation of Influenza A and B Viruses in Respiratory Specimens. Journal of Clinical Microbiology, Aug. 2004, p. 3661-3664 Vol. 42, No. 8.
4. Kenneth E. Imen and James J. Kelleher. Use of Monoclonal Antibodies for Rapid Diagnosis of Respiratory Viruses in a Community Hospital. Clinical and Diagnostic Laboratory Immunology, May 2000, p. 396-403 Vol. 7, No. 3.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Batch code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		For IVD performance evaluation only

CerTest María de Luna 11, nave 16
 BIOTEC S.L. E-50018 Zaragoza (SPAIN)
 www.cerTest.es

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BIBLIOGRAPHIE

1. Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. Journal of Clinical Microbiology, Mar. 2000, p. 1161-1165 Vol. 38, No. 3
2. Adriana Weinberg and Miranda L. Walker. Evaluation of Three Immunoassay Kits for Rapid Detection of Influenza Virus A and B. Clinical And Diagnostic Laboratory Immunology, Mar. 2005, p. 367-370 Vol. 12, No. 3
3. Andreea C. Cazacu, Gail J. Demmler et al. Comparison of a New Lateral-Flow Chromatographic Membrane Immunoassay to Viral Culture for Rapid Detection and Differentiation of Influenza A and B Viruses in Respiratory Specimens. Journal of Clinical Microbiology, Aug. 2004, p. 3661-3664 Vol. 42, No. 8.
4. Kenneth E. Irmen and James J. Kelleher. Use of Monoclonal Antibodies for Rapid Diagnosis of Respiratory Viruses in a Community Hospital. Clinical and Diagnostic Laboratory Immunology, May 2000, p. 396-403 Vol. 7, No. 3.

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