RAPIDIAG 40/41 Adeno

REF

COR C-1003

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
BIBLIOGRAPHIE

- 1 Improvement of the specificity of enzyme immunoassays for the detection of Rotavirus and Adenovirus in fecal specimens.
Rabeneau, H., Knoll, B., Allwin, R., Doerr, HW. and Weber, B.
Intervirolgy, 1998; 41(2-3):55-62
- 2 Comparison of detection methods for Adenovirus from enteric clinical specimens.
Ahluwalia, GS., Scott-Taylor, TH., Klisko, B. and Hammond, GW.
Diagn Microbiol. Infect Dis., 1994; 18(3):161-166
- 3 Evaluation of rapid culture centrifugation method for Adenovirus detection in stools.
Durepaire, N., Ranger-Rogez, S. and Denis, F.
Diagn Microbiol. Infect Dis., 1996; 24(1):25-29
- 4 Importance of Rotavirus and Adenovirus types 40 and 41 in acute gastroenteritis in Korean children.
Kim, KH, Yang, JM, Joe, SI, Cho, YG, Glass, RI and Cho, YJ.
J. Clin. Microbiol. 1990; 28(10):2279-2284
- 5 Gastroenteritis caused by Adenoviruses 40/41: epidemiological and clinical aspects.
Pena, MJ., Elcuaz, R., Suarez, J. and Lafarga, B.
Enferm.Infecc. Microbiol.Clin. 1992; 10(8):481-485
- 6 Prevalence of group A Rotavirus, human calcivirus, astrovirus and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France.
Bon, F., Facsia, P., Dauvergne, M., Tenebaum, D., Planson, H., Petion, AM., Pothier, P. and Kohli, E.
J. Clin. Microbiol. 1999; 37(9):3055-3058

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RAPIDIAG 40/41 Adeno

COR C-1003

In vitro rapid diagnostic test for Adenovirus gastro-enteritis (serotypes 40/41) in stool samples

DEFINITION

Diarrhoea and gastro-enteritis in human beings can be caused by viruses (Rotavirus, Adenovirus, Astrovirus, Norwalk virus, etc), bacteria such as Salmonella and E. coli, and protozoa such as Cryptosporidium and Giardia.

Viruses cause 45% of the diarrhoea in children under 1 year of age and 40% of the diarrhoea in children under 4.

The prevalence of Adenovirus is 4-12%. This makes it the second leading cause of viral enteritis in children under two years of age.

Infection occurs via the faecal-oral route, but can result from the inhalation of aerosols as well. The incubation period lasts from 5 to 8 days and the symptoms of the stomach and intestinal inflammation are watery diarrhoea, vomiting, fever, and abdominal cramps.

The Adenoviruses are divided into six subgroups labelled A to F. Subgroup F is the most frequently involved in paediatric gastroenteritis. The 40/41 Adeno-Strip detects subgroup F.

ASSAY PRINCIPLE

The **RAPIDIAG 40/41 Adeno** test is ready to use and is based on the use of a homogeneous membrane system technology with colloidal gold particles. The faecal sample must be diluted in the dilution buffer that is supplied with the test. A nitrocellulose membrane is sensitized with antibody to Adenovirus.

The test's specificity is ensured by a monoclonal antibody directed against specific proteins of human Adenovirus 40/41 that is coupled to the colloidal gold. This conjugate is insolubilized on a polyester membrane.

When the strip is dipped into the liquid phase of the faecal suspension, the resolubilized conjugate migrates with the sample by passive diffusion and the conjugate and sample material come into contact with the anti-Adenovirus monoclonal antibody adsorbed to the nitrocellulose. If the sample contains Adenovirus 40/41, the conjugate-Adenovirus 40/41 complex remains bound to the anti-Adenovirus monoclonal antibody. The result – in the form of a red line that develops on the strip – is visible within ten minutes.

The solution continues to migrate to encounter a control reagent (an anti-chicken IgY antibody) that binds a control conjugate, there by producing a second dark red line and confirming that the test is working properly.

REAGENTS

RAPIDIAG 40/41 Adeno Strip  <p>Each strip is sensitized with a mouse monoclonal anti-Adenovirus antibody directed against the Hexon antigens of Adenovirus and with a control reagent. The purified reagents are adsorbed to the nitrocellulose. The anti-Adenovirus conjugate is produced with mouse monoclonal antibody directed against the group F antigens (serotypes 40/41). This antibody is purified and coupled to colloidal gold particles. These strips come in a bottle or a pouch with a desiccant.</p>	25
Dilution buffer   <p>Saline solution buffered to pH 7.5 with Tris and containing EDTA, NaN₃ (<0.1%), a detergent, and charged proteins. Ready to use</p>	15 mL

MATERIAL REQUIRED BUT NOT PROVIDED

- 3 or 5 mL test tubes
- inoculating loops for taking the faecal samples
- cardboard rack
- **RAPIDIAG CONTROL** 40/41 Adenovirus (Ref.: COR C-1083)

PRECAUTIONS

All operations linked to the use of the test must be performed in accordance with Good Laboratory Practices (GLP).

The **RAPIDIAG 40/41 Adeno** tests are for *in vitro* diagnostic only.

Avoid touching the nitrocellulose with your fingers.

Wear gloves when handling the samples.

Dispose of gloves, swabs, test tubes, and sensitized strips in accordance with GLP.

Never use reagents from another kit.

If strips are stored in container, the container must be recapped resealed as soon as the necessary number of strips for the operation has been removed, since the strips are sensitive to humidity. Make sure that the desiccant is present.

If strips are stored in individual pouches, pouch must be opened with care to avoid damaging the strip.

Two green lines indicate the antibody adsorption sites. They disappear in the course of the test.

Discard the buffer solution if it is contaminated with bacteria or mould.

The reagent's quality cannot be guaranteed beyond their shelf-life dates or if the reagents are stored under inappropriate conditions.

To avoid diluting the colloidal gold conjugate in the solution, take care not to immerse the strip above the line placed under the green arrow.

SAMPLES COLLECTION AND HANDLING

- The stool specimens must be tested as soon after they are collected as possible.
- The samples can be stored in the refrigerator at +2°C/+8 °C for 24 hours. For longer storage they must be kept frozen at -20°C.
- Make sure that the specimens are not treated with solutions containing formaldehyde or its derivatives.

STABILITY AND STORAGE

- This test kit must be stored at +4°C to +37°C upon receipt.
- Do not use kits beyond the expiration date indicated on the label.
- After opening the bottle containing the test strips, the strips remain stable for 15 weeks if they remain in their original bottle and stored between +4°C and +37°C in a dry environment
- **DO NOT FREEZE.**

SET UP

Preparation of the test:

If this test was kept at +4°C, let all the reagents warm up to room temperature before proceeding with the test.

Write the patient's name or specimen number on the test tube (foresee one test tube per sample).

Place the marked test tubes in a rack.

Procedure:

- ① Add 0.5 ml or 15 drops of the dilution buffer solution to each tube.
- ② Dip the inoculating loop containing the stool sample into the tube. The dilution ratio must be at most 4% w/v. For liquid samples, take 2 loops of 10 µL, for solid samples, take 1 loop.
- ③ Stir to homogenized the solution and let stand for 1-2 minutes.
- ④ Discard the inoculating loop and immerse the sensitized strip in the direction indicated by the arrow.
- ⑤ Let react for 10 minutes, Results must be read on wet strips after 10 minutes incubation

The strong positive samples are detected after 1 or 3 minutes.

INTERPRETATION OF RESULTS

The results are to be interpreted as follows:

- 1 upper line = negative**
- 2 lines = positive**
- 0 line = invalid***

* The absence of the control line, which is the upper line, makes the result invalid. In this case, the sample must be retested.

The intensity of the test line may vary according to the quantity of antigens found in the sample. Any signal, even weak, on the test line must be regarded as a positive result.

Nevertheless, the test is qualitative and cannot predict the quantity of antigens present in the sample. The clinical presentation and other test results must be taken into consideration to establish diagnosis. During the drying process, a very faint shadow may appear at the test line. It should not be regarded as a positive result.

To store the results, let the strip dry after removing the absorbent material at its base. During the drying process, a very faint shadow may appear at the test line.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

Limit of detection

The **RAPIDIAG 40/41 Adeno** strip's limit of detection is the dilution detection threshold below which it no longer detects Adenovirus 40/41. This limit was determined by comparing it with ELISA on four stool specimens positive for Adenovirus 40 and 41. The results showed that the **RAPIDIAG 40/41 Adeno** kit's limit of detection for two of the specimens tested was identical to that of the ELISA technique. It was below that of the ELISA technique for one of the stool specimens and above that of the ELISA technique for the last stool specimen. It should be pointed out that the **RAPIDIAG 40/41 Adeno** test is not sensitive to the background noise problems that plaque ELISA.

Sensitivity, specific affinity and agreement to other methods

The kit was validated by comparing **RAPIDIAG 40/41 Adeno** kit's results with those of an ELISA test.

The **RAPIDIAG 40/41 Adeno** kit's sensitivity and specificity were tested on 153 stool samples. The following results were obtained:

ELISA	Positive	Negative	Total
RAPIDIAG 40/41 Adeno			
Positive	9	1	10
Negative	0	143	143
Total	9	144	153

Sensitivity = 100 % (9/9)

Specificity = 99.30 % (143/144)

Reliability (Concordance) = 99.35 % (152/153)
(N = 153)

Accuracy

Intra-lot

Same positive sample and a buffer solution have been processed 15 times on sticks of the same **RAPIDIAG 40/41 Adeno** production lot in the same experimental conditions.

The dilution buffer was tested 15 times in parallel with the positive faecal samples.

The 15 tests made on the 40/41 Adenovirus positive faecal sample - dilution 1/48 – were found positive with appearance of two colored lines.

All observed results were correct as expected.

The 15 tests made with the dilution buffer were found all negative with appearance of a single colored line (control line).

Inter-lot

Some samples (positive and buffer) were processed on six different **RAPIDIAG 40/41 Adeno** production lots.

The dilution buffer was tested 3 times in parallel with the positive faecal samples.

All results were correct as expected.

In every test, 6 tested lots gave positive results with the positive faecal samples in 40/41 Adenovirus– at the dilution 1/48 -and gave negative results with the dilution buffer.

Interferences

Cross reactivity with samples positive for the following pathogens was tested and found to be negative:

- *Cryptosporidium parvum* (n = 9)
- *Campylobacter jejuni* (n= 10)
- *Giardia lamblia* (n = 10)
- *Rotavirus* (n = 25)
- *E. coli* 0157 : H7 (n = 2)
- *Salmonella thyphimurium* (n = 1)
- *Salmonella enteritidis* (n = 1)
- *Yersinia enterocolitica* (n = 3)
- *Helicobacter pylori* (n = 1)
- *Aeromonas hydrophila* (n = 1)

LIMITS

RAPIDIAG 40/41 Adeno kit results must be compared with all other available clinical and laboratory information.

A positive test does not rule out the possibility that other pathogens may be present.

The **RAPIDIAG 40/41 Adeno** is an acute-phase screening test. Stool specimens that are collected after this phase may contain antigen titres below the reagent's sensitivity threshold.

QUALITY CONTROL

[bmd](#); in accordance with Good Laboratory Practices, recommend checking the test's performance regularly in line with the laboratory's requirements.

To do this, the reagent **RAPIDIAG CONTROL 40/41 Adenovirus** (Ref.No COR C-1083), in which the strip is immersed, may be used. (for more information, refer to the COR C-1083 package insert).

TECHNICAL PROBLEMS / COMPLAINTS


If you encounter a technical problem, or if performances do not match with those indicated in this package insert:

- 1- Record the lot No of the kit concerned.
- 2- If necessary, store the suspicious sample in the freezer as soon as possible.
- 3- Contact the [bmd](#) technical support.

REFERENCES

- 1 Improvement of the specificity of enzyme immunoassays for the detection of Rotavirus and Adenovirus in fecal specimens.
Rabeneau, H., Knoll, B., Allwin, R., Doerr, HW. and Weber, B.
Intervirology, 1998; 41(2-3):55-62
- 2 Comparison of detection methods for Adenovirus from enteric clinical specimens.
Ahluwalia, GS., Scott-Taylor, TH., Klisko, B. and Hammond, GW.
Diagn Microbiol. Infect Dis., 1994; 18(3):161-166
- 3 Evaluation of rapid culture centrifugation method for Adenovirus detection in stools.
Durepaire, N., Ranger-Rogez, S. and Denis, F.
Diagn Microbiol. Infect Dis., 1996; 24(1):25-29
- 4 Importance of Rotavirus and Adenovirus types 40 and 41 in acute gastroenteritis in Korean children.
Kim, KH, Yang, JM, Joe, SI, Cho, YG, Glass, RI and Cho, YJ.
J. Clin. Microbiol. 1990; 28(10):2279-2284
- 5 Gastroenteritis caused by Adenoviruses 40/41: epidemiological and clinical aspects.
Pena, MJ., Elcuaz, R., Suarez, J. and Lafarga, B.
Enferm.Infecc. Microbiol.Clin. 1992; 10(8):481-485
- 6 Prevalence of group A Rotavirus, human calcivirus, astrovirus and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France.
Bon, F., Facsia, P., Dauvergne, M., Tenebaum, D., Planson, H., Petion, AM., Pothier, P. and Kohli, E.
J. Clin. Microbiol. 1999; 37(9):3055-3058

SUMMARY OF METHOD

SYMBOLS USED

	Read instructions for use		Temperature limitation	REF	Catalog Number
	Number of tests	IVD	In Vitro Diagnostic Use	LOT	Lot Number
DIL SPE	Diluent assay		Use by		Do not reuse
STRIP	Test Strip		Keep dry	CE	EC Declaration of Conformity
CONT NaN₃	Contains Sodium Azide				

biomédical diagnostics

Office

Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com