RAPIDIAG CONTROLS

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DEFINITION

RAPIDIAG CONTROL reagents have been designed for the use with immunochromatographic **RAPIDIAG** tests.

These controls are designed to produce a positive reaction with **RAPIDIAG** strips.

These controls allow to verify the presence (not the specificity) of antibodies directed against the antigens specified.

RAPIDIAG CONTROL		Antibody
COR C-1083	RAPIDIAG CONTROL ADENO 40/41	Anti-Adenovirus serotype 40/41 (COR C-1003)
COR C-1092	RAPIDIAG CONTROL GRIPPE A&B	Anti-Grippe A&B (COR C-1012)

PRINCIPLE

When the strip is dipped in the Control diluted with the dilution buffer of each kit, the solution migrates by passive diffusion and the control reagent binds first to a conjugated antibody then to the antibody coated on the test line, thereby causing a coloured line to appear on the specific test line and allows to valid the test.

REAGENTS

RAPIDIAG CONTROL	CONTROL TEST	1mL
-------------------------	---------------------	-----

STABILITY AND STORAGE

- Store the reagents at +2°C to +8°C when not in use.
- Do not use reagents beyond the expiration date.
- The reagent can be aliquoted and frozen at -20°C for several months. They can be thawed only once.

PROTOCOL FOR USE

RAPIDIAG CONTROL reagents should be diluted ½ with the dilution buffer of the kit.

RAPIDIAG CONTROL reagents should be analyzed in the same manner as patient samples, according to the instructions supplied.

RESULTS INTERPRETATION

The results are to be interpreted as follows

1 upper line = negative

2 lines = positive

0 line = invalid*

* *The absence of the control line, which is the upper line, makes the result invalid. In this case, the sample must be retested.*

SYMBOLS USED

	Temperature limitation
IVD	In Vitro Diagnostic Use
CE	EC Declaration of Conformity
REF	Catalog Number
LOT	Lot Number
	Use by

biomédical diagnostics

Office

Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com