MULTISCREEN 6 (AMP/COC/THC/MET/MDMA/MOP)

REF DT 002

REF DT 003

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REF DT 002

REF DT 003

A rapid, one step screening test for qualitative detection of 6 drugs in human urine

MULTISCREEN 10 MTD allows the detection of multiple drugs and metabolites in human urine: **Amphetamine, Cocaine, Marijuana, Methamphetamine, Methylenedioxymethamphetamine (Ecstasy) and Morphine 300.**

INTENDED USE & SUMMARY

Urine based screen tests for multiple drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method to screen urine for multiple drugs of abuse.

The **MULTISCREEN 6** Test Panel (urine) is a lateral flow chromatographic immunoassay for the qualitative detection of following drugs without the need of instruments.

| Test* | Calibrator | Cut-off (ng/mL) |
|---|------------------------------------|-----------------|
| Amphetamine (AMP) | d-Amphetamine | 1,000 |
| Cocaine (COC) | Benzoylcegonine | 300 |
| Marijuana (THC) | 11-nor- Δ^9 -THC-9-COOH | 50 |
| Methamphetamine (MET) | d-Methamphetamine | 1,000 |
| Methylenedioxy-methamphetamine {Ecstasy} (MDMA) | d,l Methylenedioxy methamphetamine | 500 |
| Morphine (MOP 300) | Morphine | 300 |

* This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

MULTISCREEN 6 Test Panel (Urine) provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

PRINCIPLE

The **MULTISCREEN 6** Test Panel (Urine) is an immunoassay based on the principle of competitive binding.

Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action.

A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles.

The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will show up in the test line region of the specific drug strip.

The colored line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles.

A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region.

To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

Each test in the test panel contains mouse monoclonal antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in each control line.

MATERIALS

Materials Provided

- **MULTISCREEN 6** Test panels
- Package insert

TEST DEVICE

Materials Required But Not Provided

- Specimen collection container
- Timer

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- The test panel should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test panel should be discarded according to local regulations.

STORAGE AND STABILITY

- Store as packaged in the sealed pouch either at room temperature or refrigerated (+2°C/+30°C).
- The test panel is stable through the expiration date printed on the sealed pouch.
- The test panel must remain in the sealed pouch until use.
- **Do not freeze.**
- Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.


Specimen Storage

Urine specimens may be stored at +2°C/+8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

DIRECTIONS FOR USE

Allow the test panel, urine specimen, and/or controls to equilibrate to room temperature (+15°C/+30°C) prior to testing.

1. Remove the test card from the sealed pouch and use it as soon as possible. Remove the cap from the end of the test card. With arrows pointing toward the urine specimen, immerse the strip(s) of the test card vertically in the urine specimen for at least 10-15 seconds. **Immerse the strip(s) to at least the level of the wavy lines, but not above the arrow(s) on the test card.**
2. Replace the cap and place the test card on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear.
3. **Read the drug strip results at 5 minutes.** Do not read results after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE*: A colored line in the control line region (C) and a colored line in the test line region (T) for a specific drug indicate a negative result.

This indicates that the drug concentration in the urine specimen is below the designated cut-off level for that specific drug.

NOTE: The shade of color in the test region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive result.

This indicates that the drug concentration in the urine specimen exceeds the designated cut-off for that specific drug.

INVALID: Control line fails to appear.

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact imd@bmd.com.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The **MULTISCREEN 6** Test Panel (Urine) provides only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. The test does not distinguish between drugs of abuse and certain medications.
7. A positive result might be obtained from certain foods or food supplements.

BIBLIOGRAPHY

Tietz NW. - Textbook of Clinical Chemistry
W.B. Saunders Company. 1986; 1735

Baselt RC. - Disposition of Toxic Multi-Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

Hawks RL, CN Chiang. - Urine testing for Drugs of Abuse
National Institute for Drug Abuse (NIDA), Research Monograph 73, 198

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the **MULTISCREEN 6** Test Panel (Urine) and commercially available drug rapid tests. Testing was performed on approximately 300 specimens previously collected from subjects present for drug screen testing. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

% Agreement with Commercial Kit

| | Positive (%) | Negative (%) | Total |
|---------|--------------|--------------|-------|
| AMP | 97 % | > 99 % | 98 % |
| COC | 95 % | > 99 % | 98 % |
| THC | 98 % | >99 % | 99 % |
| MET | 98 % | >99 % | 99 % |
| MDMA | >99 % | 99 % | 99 % |
| MOP 300 | >99 % | >99 % | >99 % |

% Agreement with GC/MS

| | Positive (%) | Negative (%) | Total |
|---------|--------------|--------------|-------|
| AMP | 97 % | 95 % | 96 % |
| COC | 96 % | 90 % | 93 % |
| THC | 97 % | 88 % | 91 % |
| MET | 99 % | 94 % | 96 % |
| MDMA | >99 % | 98 % | 99 % |
| MOP 300 | >99 % | 94 % | 97 % |

Analytical Specificity

The following tables lists the concentration of compounds (ng/mL) that are detected positive in urine by the **MULTISCREEN 6** Test Panel (Urine) at 5 minutes.

| AMPHETAMINE | | METHAMPHETAMINE | |
|--------------------------------------|--------|---|--------|
| d-Amphetamine | 1,000 | d-Methamphetamine | 1,000 |
| d,l-Amphetamine | 3,000 | p-Hydroxymethamphetamine | 30,000 |
| l-Amphetamine | 50,000 | l-Methamphetamine | 8,000 |
| Phentermine | 3,000 | Mephentermine | 50,000 |
| 3,4-Methylenedioxy-amphetamine (MDA) | 2,000 | 3,4-Methylenedioxy-methamphetamine (MDMA) | 2,000 |
| COCAINE | | METHYLENEDIOXY-METHAMPHETAMINE | |
| Benzoylcegonine | 300 | 3,4-Methylenedioxy-méthamphétamine (MDMA) | 500 |
| Cocaine | 780 | 3,4-Methylenedioxy-methamphetamine (MDA) | 3,000 |
| Cocethylene | 12,500 | 3,4-Methylenedioxy-methamphetamine (MDEA) | 300 |
| Ecgonine | 32,000 | | |

| MARIJUANA | | | |
|-------------------------------|--------|-----------------|--------|
| 11-nor- Δ^9 -THC-9COOH | 50 | Δ^8 -THC | 15,000 |
| Cannabinol | 20,000 | Δ^9 -THC | 15,000 |
| 11-nor- Δ^8 -THC-9COOH | 30 | | |

| MORPHINE 300 | | | |
|------------------------------|--------|------------------------------------|---------|
| Morphine | 300 | Morphine 3- β -D-glucuronide | 1,000 |
| Codeine | 300 | Norcodeine | 6,250 |
| Ethylmorphine | 6,250 | Normorphine | 100,000 |
| Hydrocodone | 50,000 | Oxycodone | 30,000 |
| Hydromorphone | 3,125 | Oxymorphone | 100,000 |
| Levorphanol | 1,500 | Procaine | 15,000 |
| 6-Monoacetylmorphine (6-MAM) | 400 | Thebaine | 6,250 |

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or **Amphetamines, Cocaine, Marijuana, Methamphetamine, Methylenedioxy-methamphetamine, Morphine 300** positive urine. The following compounds show no cross-reactivity when tested with the **MULTISCREEN 6** Test Panel (Urine) at a concentration of 100 μ g/mL.

Non Cross-Reacting Compounds

| | | |
|-----------------------|------------------------|--------------------------------|
| Acetophenetidin | Diphenhydramine | Oxolinic acid |
| N-Acetylprocainamide | Ethyl-p-aminobenzoate | Oxymetazoline |
| Acetylsalicylic acid | β -Estradiol | Papaverine |
| Aminopyrine | Estrone- 3- sulfate | Penicillin-G |
| Amoxicillin | Erythromycin | Perphenazine |
| Ampicillin | Fenoprofen | Phenelzine |
| L-Ascorbic acid | Furosemide | Prednisone |
| Apomorphine | Gentisic acid | d, l-Propranolol |
| Aspartame | Hemoglobin | d-Pseudoephedrine |
| Atropine | Hydralazine | Quinidine |
| Benzilic acid | Hydrochlorothiazide | Quinine |
| Benzoic acid | Hydrocortisone | Salicylic acid |
| Bilirubin | o-Hydroxyhippuric acid | Serotonin |
| d, l-Brompheniramine | 3-Hydroxytyramine | Sulfamethazine |
| Caffeine | d, l-Isoproterenol | Sulindac |
| Cannabidiol | Isoxsuprine | Tetracycline |
| Chloralhydrate | Ketamine | Tetrahydrocortisone, 3-Acetate |
| Chloramphenicol | Ketoprofen | Tetrahydrocortisone |
| Chlorothiazide | Labetalol | Tetrahydrozoline |
| d, l-Chlorpheniramine | Loperamide | Thiamine |
| Chlorpromazine | Meprobamate | Thioridazine |
| Cholesterol | Methoxyphenamine | d, l-Tyrosine |
| Clonidine | Methylphenidate | Tolbutamide |
| Cortisone | Nalidixic acid | Triamterene |
| l-Cotinine | Naproxen | Trifluoperazine |
| Creatinine | Niacinamide | Trimethoprim |
| Deoxycorticosterone | Nifedipine | d, l-Tryptophan |
| Dextromethorphan | Norethindrone | Uric acid |
| Diclofenac | Noscapine | Verapamil |
| Diflunisal | d, l-Octopamine | Zomepirac |
| Digoxin | Oxalic acid | |

SYMBOLS USED

| | | | |
|---|----------------------------------|---|--------------------------------|
|  | Catalog number |  | Storage temperature limitation |
|  | Lot number |  | Rapid test |
|  | Use by |  | Do not reuse |
|  | Number of tests |  | EC Declaration of Conformity |
|  | For in vitro diagnostic use only |  | Read instructions for use |

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tél : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com