ImmunoDOT™ TPO-TG

REF GB 5885

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BIBLIOGRAPHIE

Czarnocka B, et.al. "Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid diseases" *FEBS* 190:147 (1985)

Yoshida H, et.al. "Association of serum antithyroid antibodies with lymphocytic infiltration of the thyroid gland: study of seventy autopsied cases" *J Clin Endoc Metab* 46:859 (1978)

Hooper B, et.al. "Autoimmunity in a rural community" *Clin Exp Immunol* 12:79 (1972)

Doniach D. "Humoral and genetic aspects of thyroid autoimmunity" *Clin Endocrinol Metab* 4:267 (1975)

Centers for Disease Control/National Institutes of Health Manual *Biosafety in Microbiological and Biomedical Laboratories* (1983)

Hijmans W, et.al. "Serological overlap between lupus erythematosus, rheumatoid arthritis and thyroid autoimmune disease" *Brit Med J* 2:909 (1961)

Evered DC, et.al. "Grades of Hypothyroidism" *Brit Med J* 1:675 (1973)

Volpe R. "The role of autoimmunity in hypoendocrine and hyperendocrine function: with special emphasis in autoimmune thyroid disease" *Ann Int Med* 87:86 (1977)

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GB 5885

50

Qualitative detection of autoantibodies against human thyroglobulin and human thyroid peroxidase (microsome).

INTENDED USE

The **ImmunoDOT™ TPO/TG** is an enzyme immunoassay (EIA) test for screening and detection of autoantibodies against human thyroglobulin and human thyroid peroxidase (microsome) in serum and heparinized whole blood and is used as an aid in the diagnosis of thyroid disorders. This product is intended for use in physician offices and clinical laboratories.

SUMMARY AND EXPLANATION

Autoimmune thyroid gland disorders are characterized by detection of anti-thyroid antibodies, primarily against thyroglobulin and/or microsomal thyroid antigens. Recently it has been shown that thyroid peroxidase (TPO) is the protein responsible for microsomal antigenicity. In addition to chronic thyroiditis, thyroid autoantibodies may be found in other thyroid disorders. These autoantibodies may also occur in apparently normal subjects. Thyroid microsomal (TPO) autoantibodies occur in sera of most autoimmune thyroid disease patients and predict raised serum TSH levels in random populations. The presence of autoantibody does not imply active tissue destruction. Microsomal (TPO) antibody level correlates with the degree of lymphoid infiltration of the thyroid gland.

Thyroid autoantibodies are detected by a variety of immunoassays. The common methods are indirect hemagglutination (IHA), indirect fluorescent antibody (IFA), EIA And RIA techniques.

ImmunoDOT™ TPO/TG Test uses highly purified human thyroglobulin devoid of microsomal antigen and recombinant human thyroid peroxidase (microsomal antigen) which does not contain contaminating thyroglobulin and/or mitochondria found in other microsome antigen preparations. These purified antigens are used in a dipstick format enzyme immunoassay technique.

ASSAY PRINCIPLE

The **ImmunoDOT™ TPO/TG** Test utilizes an EIA dot technique for the detection of antibodies. The antigens are dispensed as discrete dots onto a solid membrane.

After adding specimen to a reaction vessel, an assay strip is inserted, allowing patient antibodies reactive with the test antigen to bind to the strip's solid support membrane.

In the second stage, the reaction is enhanced by removal of non-specifically bound materials.


During the third stage, alkaline phosphatase-conjugated anti-human antibodies are allowed to react with bound patient antibodies.

Finally, the strip is transferred to enzyme substrate reagent which reacts with bound alkaline phosphatase to produce an easily seen, distinct dot.

STORAGE CONDITIONS

- Store reagents and assay strips between +2°C and +8°C.
- Reagents must be at room temperature (between +15°C and +30°C) before use. They should be used within one hour following placement in the heated workstation (incubator)
- Avoid contamination of reagents as that can give erroneous results.

MATERIAL PROVIDED

ImmunoDOT™ TPO/TG Assay strip		50
 <p>STRIP</p> <ul style="list-style-type: none"> ⊕ Positive control ⊕ Anti-TPO (cut off value) ⊕ Anti-TPO (positive value) ⊕ Anti-TG (cut off value) ⊕ Anti-TG (positive value) ○ Negative control 		
Diluent <i>buffered diluent (pH 6.2-7.6), protein stabilizers, and <0.1% NaN₃.</i>	DIL SPE	2 x 50mL
Enhancer <i>sodium chloride and <0.1% NaN₃.</i>	SOLN ENH	2 x 50mL
Conjugate <i>alkaline phosphatase conjugated goat anti-human antibodies in buffered diluent (pH 6.2-8.5) and stabilizers.</i>	CONJ IgG	2 x 50mL
Developer <i>5-bromo-4-chloro-3-indolyl phosphate and p-nitro blue tetrazolium chloride in buffered diluent (pH 9.0-11.0), 0.8% N, N-Dimethylformamide, and <0.1% NaN₃</i>	SUBS BCIP-NBT	2 x 50mL
Reaction Vessels		200
Package Insert		1

MATERIALS REQUIRED BUT NOT PROVIDED

- Workstation (incubator)
- Specimen collection apparatus
- Timer
- Analytical quality distilled or deionized water to be used for the washes.
- Pipets
- Absorbent toweling to blot dry the assay strip
- Positive control serum

SAMPLE COLLECTION AND HANDLING

- The **ImmunoDOT™ TPO/TG Test** is performed on serum or heparinized whole blood. The test requires approximately 10µL of serum or 20 µL of whole blood.
- Serum and heparinized whole blood are collected according to standard practices
- Serum and heparinized whole blood may be stored at +2°C and +8°C for up to five days.
- Serum may be frozen below -20°C for extended periods. Freezing whole blood samples is not advised.

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use Only.

The **ImmunoDOT™ TPO/TG** reagents have been optimized for use as a system. Do not substitute reagents or strips from any other source including those from other **ImmunoDOT™** systems.

Dilution or modification of these reagents may also affect the performance of the test.

Do not use kits after the stated expiration date.

Analytical quality deionized or distilled water must be used for the washes.

Strictly adhere to the test procedures for optimal results. Do not shorten or prolong stated incubation times as this may result in poor assay performance.

Some assay components contain sodium azide (NaN₃) that may react with lead and copper plumbing to form highly explosive metal azides. After disposing of reagents, flush sink with a large volume of water to prevent azide buildup.

Human sera used in the preparation of this product were tested and found non-reactive for hepatitis B surface antigen and for antibodies to HIV-1, HIV-2, HTLV-1, and hepatitis C virus. Because no test method can offer complete assurance that infectious agents are absent, handle reagents and patient samples as if capable of transmitting infectious disease.

SETUP

1. Turn on Workstation and adjust temperature according to the instructions for use of your instrument.

2. Remove 4 Reaction Vessels per test from the product box and insert into appropriate slots in Workstation.

For the large Workstation, add water up to the fill line of the Wash Vessel provided.

For the small Workstation, use an appropriate container and sufficient water to cover all reactive windows of the assay strip.

3. Appropriately label the Assay Strips.

4. Place:

- 2 mL Diluent (#1) in Reaction Vessel #1
- 2 mL Enhancer (#2) in Reaction Vessel #2
- 2 mL Conjugate (#3) in Reaction Vessel #3

5. Incubate for 10 minutes.

6. Verify that the temperature of the heating block is in accordance with the instructions of the manual corresponding to the model of the workstation.

7. If the large Workstation is used, insert the label end of the Assay Strip into the Strip Holder, one per groove, taking care not to touch the assay windows.

8. Insert the strip in the distilled water bath (Wash) for 30 seconds, just before starting the assay procedure. Avoid the formation of bubbles on the surface.

ASSAY PROCEDURE

1. SAMPLE INCUBATION

- Add specimen (10µL of serum or 20µL of whole blood) to Reaction Vessel #1.
- Plunge assay strip into Reaction Vessel #1
- Using 10 quick up and down motions with the Assay Strip, mix reagents and specimen thoroughly.
- **Let stand for 5 minutes.**

2. WASH STEP

- Remove Assay Strip from Reaction Vessel and swish in the Wash (distilled water). Use a swift back and forth motion for 10 seconds allowing for optimal washing of the Assay Strip.
- Blot the strip on an absorbent towel.

3. ENHANCER

- Place Assay Strip into Reaction Vessel #2.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 5 minutes.**

4. WASH STEP

- Remove Assay Strip from Reaction Vessel #2 and swish in the Wash for 10 seconds as described (step #2). The wash should diffuse completely through the strip.
- Blot the strip on an absorbent towel.

5. CONJUGATE

- Place Assay Strip into Reaction Vessel #3.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 15 minutes.**
- **At the beginning of this incubation, place 2mL of reagent #4 into the Reaction Vessel #4.**

6. WASH STEP

- Remove Assay Strip from Reaction Vessel #3 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip
- DO NOT remove the Assay Strip from the wash.
- **Allow the Assay Strip to stand in the wash for 5 minutes.**
- Blot the strip on an absorbent towel.

7. DEVELOPER

- Remove Assay Strip from wash and place into Reaction Vessel #4.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 5 minutes.**


8. WASH STEP

- Remove Assay Strip from Reaction Vessel #4 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.

9. DRY

- Blot and allow Assay Strip to dry. It is imperative that tests of borderline specimens be interpreted after the Assay Strip has been allowed to dry. **A false positive dot may be identified if the assay strip is not dry when interpreted.**

10. READING THE ASSAY STRIP

Positive cut off: a dot with an **EASILY SEEN**, distinct border.

Positive: a dot with an **EASILY SEEN**, distinct border is visible in the center of the window. The outer perimeter of the window must be white to pale gray.

Negative: no dot - If no dot is seen or a dot is difficult to see, interpret it as negative.

Each dot of the assay strip must be read independently.

TEST VALIDATION

The top and bottom membrane windows of the Assay Strip are reagent controls.

The top window is a positive reagent control and must be positive for further interpretation.

The bottom window is the reagent negative control and must be negative for further interpretation.

Reagent controls assure that reagents are active and that the test has been performed properly. If either reagent control is invalid, the test must be repeated.

The intensity of the positive control dot must not be used as a calibrator.

Positive reactions in the antigen windows of the strip may be either darker or lighter than the positive control depending on the antibody titer.

QUALITY CONTROL

It is recommended to use internally or externally sourced control material for the different specificities. Multiparametric control (bmd, Cat: HM020) contains antibodies directed against one of the specificities; handle in the same way as samples.


The assay is performed with reagents warmed at +44°C/+ 48°C. The temperature of the Workstation itself should be higher than the assay temperature in the reaction vessel. The Workstation should be monitored to assure the appropriate temperature is maintained.

LIMITATIONS

- ⇒ Epidemiologic factors, clinical findings, and other laboratory results should be considered in addition to autoantibody laboratory results for diagnosis of the patient.
- ⇒ Thyroid autoantibodies may be present in non-thyroid disorders such as Sjögren's syndrome, pernicious anemia, Addison's disease, myasthenia gravis and diabetes mellitus and in apparently healthy subjects.

INTERPRETATION

There are two reaction levels for each analyte (thyroglobulin and TPO). Interpret each analyte reaction separately and according to the following criteria:

No positive Dots

Results: NEGATIVE	Anti-TPO IU/ml	Anti-TG IU/ml
Nonreactive for autoantibody	<40	<100

One positive Dot (Only cut off Dot)

Results: TO CONFIRM	Anti-TPO IU/ml	Anti-TG IU/ml
Equivocal but to confirm on a second sample	40-65	100-165

Two positive Dots (Cut off and positive)

Results: POSITIVE	Anti-TPO IU/ml	Anti-TG IU/ml
+ (dot with distinct border)	66-130	166-330
++ (weak uniform dot)	131-260	331-660
+++ (intense uniform dot)	261-520	661-1320
++++ (intense dot)	>520	>1320

When the two Dots are positive, only the spot "Value positive" is interpreted compared to the reference colored strip to have a result semi-quantitative. These values are indicative, but were obtained from a correlation with bmd ELISA kits: TG-LISA and TPO-LISA.

PERFORMANCE AND CHARACTERISTICS OF THE TEST

ImmunoDOT™ TPO/TG Test was tested with presumptive normal samples and microsome and thyroglobulin positive (IHA) samples. The results in normal samples (n=125) predict assay specificity in a random, normal population and are shown in **Erreur ! Source du renvoi introuvable.** These data predict the TPO assay to be 91% specific and the thyroglobulin to be 96% specific. These data are consistent with previous reports.

Table 1 – Normal samples

Specificity	Non reactive	To confirm	Positive
TPO	91.2% (114/125)	6.4% (8/125)	2.4% (3/125)
TG	96% (120/125)	4.0% (5/125)	0% (0/125)

Sixty-six (66) indirect hemagglutination (IHA) positive samples were evaluated in the **ImmunoDOT™ TPO/TG** Test for anti-TPO (microsome) to assess test sensitivity compared to the classic test method. These results are shown in Table 1.

Table 1: Thyroid Peroxidase (microsome)

IHA Titer	Non reactive	To confirm	Positive
100	0	4	1
400	0	6	12
1600	0	0	28
6400	0	0	9
25600	0	0	6

Assuming that an IHA titer greater than or equal to 1:1600 is clinically relevant, the **ImmunoDOT™ TPO/TG** Test for anti-TPO (microsome) is 100% (43/43) sensitive when compared to IHA. Unlike the IHA method, **ImmunoDOT™ TPO/TG** identified 13 of the 23 borderline (1:100-400) IHA results as clearly reactive samples.

Twenty-one thyroglobulin (Tg) IHA positive samples were tested in the **ImmunoDOT™ TPO/TG** Test for anti-thyroglobulin. These results are shown in Table 2.

Table 2: Thyroglobulin

IHA Titer	Non reactive	To confirm	Positive
80	1	0	2
160	1	1	0
320	0	2	2
640	0	0	4
>1280	0	0	8

Assuming that an IHA titer greater than or equal to 1:160 is clinically relevant, the **ImmunoDOT™ TPO/TG** Test for anti-thyroglobulin is 94% (17/18) sensitive when compared to IHA. In addition, **ImmunoDOT™ TPO/TG** identified two out of three borderline (1:80) IHA positives as clearly reactive samples.

REFERENCES

Czarnocka B, et.al. "Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid diseases" *FEBS* 190:147 (1985)

Yoshida H, et.al. "Association of serum antithyroid antibodies with lymphocytic infiltration of the thyroid gland: study of seventy autopsied cases" *J Clin Endoc Metab* 46:859 (1978)

Hooper B, et.al. "Autoimmunity in a rural community" *Clin Exp Immunol* 12:79 (1972)

Doniach D. "Humoral and genetic aspects of thyroid autoimmunity" *Clin Endocrinol Metab* 4:267 (1975)


Centers for Disease Control/National Institutes of Health Manual *Biosafety in Microbiological and Biomedical Laboratories* (1983)

Hijmans W, et.al. "Serological overlap between lupus erythematosus, rheumatoid arthritis and thyroid autoimmune disease" *Brit Med J* 2:909 (1961)

Evered DC, et.al. "Grades of Hypothyroidism" *Brit Med J* 1:675 (1973)

Volpe R. "The role of autoimmunity in hypoenocrine and hyperendocrine function: with special emphasis in autoimmune thyroid disease" *Ann Int Med* 87:86 (1977)

SYMBOLS USED

	EC Declaration of conformity
	Rapid test
	<i>In Vitro</i> Diagnostic Device
	Catalogue number
	Lot Number
	Expiry Date
	Number of test
	Consult Instructions
	Temperature limitation
	Biological risk
	Contains sodium azide
	Reconstitute with

SUMMARY OF METHOD

SETUP	<ul style="list-style-type: none"> • Turn on and prepare the workstation (according to the model). The temperature in vessels must lie within +44°C and +48°C. • Insert the strip in the distilled water bath (Wash) for 30 seconds. • Pipet: <ul style="list-style-type: none"> - 2 mL of diluent (#1) into reaction vessel #1 - 2 mL of enhancer (#2) into reaction vessel #2 - 2 mL of conjugate (#3) into reaction vessel #3 • Incubate for 10 min. 	INCUBATION 10 minutes
1. INCUBATION OF SERA	<ul style="list-style-type: none"> • Add 10 µL of sample sera or positive control to vessel #1. • Immerse strip in reaction vessel #1. • Mix using a swift up and down motion 10 times. • Incubate for 60 min. 	INCUBATION 5 minutes
2. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Blot the assay strip on an absorbent towel. 	10 seconds
3. ENHANCER	<ul style="list-style-type: none"> • Place strip into vessel #2. • Mix using a swift up and down motion 10 times. • Incubate for 5 min. 	INCUBATION 5 minutes
4. WASH STEP	<ul style="list-style-type: none"> • Wash the in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Blot the assay strip on an absorbent towel. 	10 seconds
5. CONJUGATE	<ul style="list-style-type: none"> • Place strip into vessel #3. • Mix using a swift up and down motion 10 times. • Incubate for 30 min. • At the start of this incubation, place 2 mL of reagent #4 into vessel #4. 	INCUBATION 15 minutes
6. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Immerse the strip in the Wash (distilled water). Let stand 5 minutes. • Blot the assay strip on an absorbent towel. 	10 seconds then INCUBATION 5 minutes
7. DEVELOPER	<ul style="list-style-type: none"> • Place strip into vessel #4. • Mix using a swift up and down motion 10 times. • Incubate for 5 min. 	INCUBATION 5 minutes
8. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. 	10 seconds
9. DRY STEP	<ul style="list-style-type: none"> • Remove the residual water with absorbent paper, and let dry before interpreting the result. 	WAIT 5 minutes

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com