# ImmunoDOT™ EBV MONO - G

REF GB 6085

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# ImmunoDOT™ EBV MONO-G

REF GB 6085

Qualitative detection of IgG antibodies to Epstein-Barr viral capsid antigen (EBV-VCA) and Epstein-Barr early nuclear antigen (EBV-EBNA) in human serum.

## INTENDED USE

The **ImmunoDOT™ EBV MONO-G Test** is a qualitative enzyme immunoassay (EIA) for the rapid detection of IgG antibodies to Epstein-Barr viral capsid antigen (EBV-VCA) and Epstein-Barr early nuclear antigen (EBV-EBNA) in human serum. This test is intended to be used for the serodiagnosis of any prior infections by EBV. When used in conjunction with the **ImmunoDOT™ EBV MONO-M Test** it is an aid in the serodiagnosis of infectious (EBV) mononucleosis.

## SUMMARY AND EXPLANATION

Infectious mononucleosis (IM) associated with Epstein-Barr Virus (EBV), is a transitory, acute and febrile lymphoproliferative syndrome with a benign course. It is widespread with over 95% of adults having been infected.

Transmission between people is most often through the saliva, and infection is seen earlier with poor sanitation conditions. After penetration into the oral pharynx, the virus multiplies in situ, where the B lymphocytes are secondarily infected during passage into the pharyngeal lymphoid tissue. Venereal transmission would appear to be possible.

EBV Infection may be expressed as:

- an asymptomatic seroconversion;
- a moderate mononucleosidic syndrome with slight fever, sore throat, and adenopathies, particularly in young children;
- an infectious mononucleosis (IM) with sore throat, adenopathies, asthenia (weakness), sensitivity to certain antibiotics (Ampicillin) particularly in adolescents;
- a mononucleosidic syndrome with chronic fatigue ;
- occasional malignant pathologies related to the persistence of the virus in the host such as lymphoma in the immunocompromised patient, undifferentiated nasopharyngeal carcinoma, oral hairy cell leukemia, and Burkitt's lymphoma.

A well defined kinetic assay shows that the antibodies directed against the various EBV antigens appear during the course of the primary infection.

Anti-VCA IgG and IgM antibodies are against capsid antigens. Anti-VCA IgM antibodies appear only during acute infection or rarely during reactivation and disappear within 4 to 8 weeks. Anti-VCA IgG antibodies appear in all primary infections with EBV. After infection, these antibodies decline and persist for life.

Anti-EBNA IgG are antibodies to nuclear antigens. These antibodies appear late after EBV infection. They are always absent during the acute phase of primary infection.

The **ImmunoDOT™ EBV MONO-G test** detects anti-VCA IgG and anti-EBNA IgG antibodies that may appear during or after EBV infection.

EBV infection is serologically diagnosed by concurrent detection of both IgM and IgG antibodies (EBNA IgG, VCA IgG, VCA IgM).

Cytomegalovirus (CMV) and *Toxoplasma gondii* controls have not been validated. Ignore the results obtained for these two parameters.

Testing for anti-CMV and anti-*Toxoplasma* IgG antibodies will have to be performed using validated serological tests.

## ASSAY PRINCIPLE

The **ImmunoDOT™ EBV MONO-G** uses an enzyme-linked immunoassay (EIA) dot technique for the detection of the anti-VCA IgG and anti-EBNA IgG antibodies. The different antigens are present as discrete dots on the assay strip.

First, an assay strip is incubated in a vessel with dilute patient serum, thus allowing any antibodies present in the sample to react with the test antigens immobilized on the assay strip.

In the second step, the specificity is enhanced by removal of non-specifically bound materials.

In the third step the strip is exposed to alkaline phosphatase conjugated to anti-human IgG antibodies so that the enzyme conjugate may bind to any previously bound patient antibodies.

Finally, the strip is transferred to the vessel containing the enzyme substrate reagent which then reacts with the conjugated alkaline phosphatase, producing distinct, clearly visible dots.

## STORAGE CONDITIONS

- Store reagents and assay strips between +2°C and +8°C.
- Reagents must be at room temperature (between +15°C and +30°C) before use. They should be used within an hour following placement in the heated workstation (incubator)
- Avoid contamination of reagents as that can give erroneous results.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Workstation (incubator)
- Specimen collection apparatus
- Timer
- Analytical quality distilled or deionized water to be used for the washes.
- Pipets
- Absorbent toweling to blot dry the assay strip

## SAMPLE COLLECTION AND HANDLING

- The **ImmunoDOT™ EBV MONO-G Test** is performed on serum. The test requires 10 µL of serum.
- Avoid Lipemic or hemolyzed serum specimens has not been shown to be an acceptable specimen.
- Collect the samples of serum aseptically and store them between +2°C and +8°C for a maximum of 48 hours.
- If the assay of the samples will not be completed within 48 hours, or for shipment of the specimens, freeze at -20°C

## MATERIAL PROVIDED

		GB 6085
<b>ImmunoDOT™ EBV MONO-G Assay Strip</b>		50
(*) partially purified EBV lysate from cell strain P3HR1		
<b>Diluent</b>	<input type="checkbox"/> DIL <input type="checkbox"/> SPE	2 x 50mL
Consists of buffer (pH 6.2-7.6) containing IgG absorbent (goat anti-human IgG), protein stabilizers, and <0.1% NaN <sub>3</sub> .		
<b>Enhancer</b>	<input type="checkbox"/> SOLN <input type="checkbox"/> ENH	2 x 50mL
Consists of sodium chloride and <0.1% NaN <sub>3</sub> .		
<b>Conjugate</b>	<input type="checkbox"/> CONJ <input type="checkbox"/> IgG	2 x 50mL
Consists of alkaline phosphatase conjugated goat anti-human IgG (heavy chain specific) in buffered diluent (pH 6.2-7.6), with protein stabilizers.		
<b>Developer [substrat]</b>	<input type="checkbox"/> SUBS <input type="checkbox"/> BCIP-NBT	2 x 50mL
Consists of 5-bromo-4-chloro-3-indolyl phosphate and p-nitro blue tetrazolium chloride in buffered diluent (pH 8.0-11.0), 0.8% dimethylformamide and <0.1% NaN <sub>3</sub>		
<b>Positive control</b>	<input type="checkbox"/> CONTROL <input type="checkbox"/> +	1 x 100μL
<b>Reaction Vessels</b>	200	
<b>Package Insert</b>	1	

## WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use Only.

The **ImmunoDOT™ EBV MONO-G** reagents have been optimized for use as a system. Do not substitute reagents or strips from any other source including those from other **ImmunoDOT** systems.

Dilution or modification of these reagents may also affect the performance of the test.

Do not use kits after the stated expiration date.

Analytical quality deionized or distilled water must be used for the washes.

Strictly adhere to the test procedures for optimal results. Do not shorten or prolong stated incubation times as this may result in poor assay performance.

Some assay components contain sodium azide (NaN<sub>3</sub>) that may react with lead and copper plumbing to form highly explosive metal azides. After disposing of reagents, flush sink with a large volume of water to prevent azide buildup.

Human sera used in the preparation of this product were tested and found non-reactive for hepatitis B surface antigen and for antibodies to HIV-1, HIV-2, HTLV-1, and hepatitis C virus. Because no test method can offer complete assurance that infectious agents are absent, handle reagents and patient samples as if capable of transmitting infectious disease.

## SETUP

- Turn on Workstation and adjust temperature according to the instructions for use of your instrument.
- Remove 4 Reaction Vessels per test from the product box and insert into appropriate slots in Workstation.  
*For the large Workstation*, add water up to the fill line of the Wash Vessel provided.  
*For the small Workstation*, use an appropriate container and sufficient water to cover all reactive areas/dots of the assay strip.
- Appropriately label the Assay Strips.
- Place:
  - 2 mL Diluent (#1) in Reaction Vessel #1
  - 2 mL Enhancer (#2) in Reaction Vessel #2
  - 2 mL Conjugate (#3) in Reaction Vessel #3
- Incubate for 10 minutes.
- Verify that the temperature of the heating block is in accordance with the instructions of the manual corresponding to the model of the workstation.
- Insert the strip in the distilled water bath (Wash) for 30 seconds, just before starting the assay procedure. Avoid the formation of bubbles on the surface.

## ASSAY PROCEDURE

### 1. SAMPLE INCUBATION

- Add 10 μL serum sample or positive control to Reaction Vessel #1.
- Plunge assay strip into Reaction Vessel #1
- Using 10 quick up and down motions with the Assay Strip, mix contents thoroughly.
- Let stand for 60 minutes.

### 2. WASH STEP

- Remove Assay Strip from Reaction Vessel and swish in the Wash (distilled water). Use a swift back and forth motion for 10 seconds allowing for optimal washing of the Assay Strip.
- Blot the strip on an absorbent towel.

### 3. ENHANCER

- Place Assay Strip into Reaction Vessel #2.
- Mix thoroughly with 10 quick up and down motions.
- Let stand for 5 minutes.

### 4. WASH STEP

- Remove Assay Strip from Reaction Vessel #2 and swish in the Wash for 10 seconds as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.
- Blot the strip on an absorbent towel.

### 5. CONJUGATE

- Place Assay Strip into Reaction Vessel #3.
- Mix thoroughly with 10 quick up and down motions.
- Let stand for 30 minutes.
- At the beginning of this incubation, place 2mL of reagent #4 into the Reaction Vessel #4.

### 6. WASH STEP

- Remove Assay Strip from Reaction Vessel #3 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip
- DO NOT remove the Assay Strip from the wash.
- Allow the Assay Strip to stand in the wash for 5 minutes.
- Blot the strip on an absorbent towel.

## 7. DEVELOPER

- Remove Assay Strip from wash and place into Reaction Vessel #4.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 5 minutes.**

## 8. WASH STEP

- Remove Assay Strip from Reaction Vessel #4 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.

## 9. DRY

- Blot and allow Assay Strip to dry. It is imperative that tests of borderline specimens be interpreted after the Assay Strip has been allowed to dry. **A false positive dot may be identified if the assay strip is not dry when interpreted.**

## 10. READING THE ASSAY STRIP

**Positive:** A dot with an **EASILY SEEN**, distinct border is visible in the center of the window. The outer perimeter of the window must be white to pale gray.

**Negative:** If no dot is seen or a dot is difficult to see, interpret it as negative.

Each dot of the assay strip is read independently.

## QUALITY CONTROL

The dots of positive and negative controls must appear correctly; otherwise the test needs to be repeated.

The color intensity of the positive control must not be used as a calibrator. The reactions in dots containing antigens could have a different intensity than the positive control, because of a difference in the concentration of antibodies.

**The positive control serum** included in this kit, is moderately positive for all antibodies. The performance of every kit can be confirmed upon receipt by using this control which is expected to give a positive result.

The analysis is performed with reagents warmed at +44°C/+ 48°C. The temperature of the workstation has to exceed the temperature in the reaction vessels. The temperature of the workstation must be watched.

## LIMITATIONS

- ⇒ This is a qualitative test, and therefore it cannot be used to detect changes in antibody titers.
- ⇒ Primarily maternal antibodies will be detected in children less than one year of age.
- ⇒ Lipemic or hemolysed samples or samples contaminated by bacteria could provide erroneous results.
- ⇒ The **ImmunoDOT™ EBV MONO G** test used alone is not in any case diagnostic for mononucleosis infectious due to the Epstein Barr virus.
- ⇒ The complete serologic profile of an infection to EBV should include determinations for the presence of antibodies directed against the different EBV-antigens appearing during the course of infection: EBV VCA IgM, EBV VCA IgG and EBV EBNA IgG.
- ⇒ If one wants to determine the state of reactivation in an immunocompromised patient or show that the EBV is associated with a malignant disease, it is essential to determine titers of the anti-EA (early antigen) antibodies.
- ⇒ **The CMV controls and toxoplasma gondii have not been validated. Ignore any results associated with these two parameters. Any testing for anti- CMV and anti-Toxoplasma antibodies must be performed using validated serologic tests.**

⇒ The interpretation of the results of the **ImmunoDOT™ EBV MONO G** test should be done in parallel with those obtained with the **ImmunoDOT™ EBV MONO M** test which provides for the detection of anti-VCA isotype IgM antibodies.

## REFERENCES

Centers for Disease Control/National Institutes of Health Manual Biosafety in Microbiological and Biomedical Laboratories (1999).

Lennette, ET, Manual of Clinical Microbiology, 4th Ed.(Lennette, Balows, Housler, and Shadomy, eds) pp 728, ASM Press, Washington DC (1985).

Evans, AS, Medcal Virology VII (de la Maza and Peterson, eds) pp 57, Elsevier New-York, NY (1988).

Naraqui,S, Textbook of Human Virology (Belshe,ed), pp 892, PSG Publishing Littleton, MA (1984).

Perham TGM,EO Caul, RJ Conway,et al, Brit.J.Haematol 20:307 (1971).

Kane, RC, WE Rousseau, GR Noble, et al., Infect Immun 11:719 (1975).

Klemola, D,R. Von Essen, J. Paloheimo, et al., Scand.J. Infect.Dis 1:137 (1969).

McCabe RE and JS Remington, Principles and Practice of Infectious Diseases (Mandel, Douglas and Bennett, eds) pp 1542, John Wiley and Sons, New -York NY (1985).

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
Lennette EH, Lennette DA, and Lennette ET. Epstein-Barr Virus, in Diagnostic Procedures for Viral, Rickettsial, and Chlamydial Infections ( 7 th Ed ), Chapter 18, American Public Health Association, Washington DC, Lennette, ET, 1995.

## SYMBOLS USED

	EC Declaration of conformity
	Rapid test
	In Vitro Diagnostic Device
	Catalogue number
	Lot Number
	Expiry Date
	Number of test
	Consult Instructions
	Temperature limitation
	Biological risk
	Contains sodium azide
	Reconstitute with

## INTERPRETATION

Reading and interpreting the VCA and EBNA antibodies dots. The **ImmunoDOT™ EBV MONO G** kit is a test for the rapid, qualitative detection of EBV VCA IgG and EBV EBNA IgG antibodies.

ImmunoDOT™ MONO M		ImmunoDOT™ MONO G		Hematological and clinical symptoms	INTERPRETATION EBV infection steps	Frequency of the pattern	REMARKS
ACs* MNI**	VCA IgM	EBNA IgG	VCA IgG				
-	-	-	-	• negative	<b>Seronegative patient (no infected)</b>	adult: rarely child: frequent	
-/+	+	-	-	• negative • indeterminate • characteristic	<b>Primary infection (early step)</b>	adult: rarely child and teenager: frequent	Serological Testing in 3 weeks
- +	+ +	- -	+ +	• negative • indeterminate • characteristic	<b>Primary infection</b>	adult: rarely child and teenager: frequent	Serological Testing in 3 weeks
- -	+ +	+ -/+	+ +	• negative • indeterminate	<b>End of primary infection (passing step)</b>	Frequent patterns for child and teenager	Serological Testing in 3 weeks
-	-	+	+	• negative	<b>Past infection</b>	Frequent pattern (95% of the adult population)	2 % of patients never have anti-EBNA antibodies

## OTHERS

ImmunoDOT™ MONO M		ImmunoDOT™ MONO G		Hematological and clinical symptoms	INTERPRETATION EBV infection	Pattern and frequency	REMARKS
ACs* MNI**	VCA IgM	EBNA IgG	VCA IgG				
+/-	+	+	+	• negative • indeterminate • characteristic	<ul style="list-style-type: none"> <li>• acute chronic infection?</li> <li>• Reactivation?</li> <li>• Recent infection with evolution?</li> </ul>	Frequent patterns in immunodeficient patient with transplantation Rarely in other patients	Test following serological parameters : 1. IgM VCA 2. IgG VCA 3. IgG EA 4. IgG and IgM CMV 5. IgG and IgM Toxo

(\* ACs = heterophile antibodies (\*\* MNI = Infectious MonoNucleosis)

### Remarks:

- The presence of anti-VCA IgG and anti-EBNA IgG antibodies is confirmatory for past infection
- CMV Control not validated:** The IgM anti-CMV detection does not allow the diagnosis of the infection by CMV. The CMV well is a control well; results should be confirmed by other validated IgG and IgM CMV tests. CMV infections are frequent and ubiquitous. The prevalence of the infection is higher in populations living in conditions of poor sanitation or in densely populated areas (4). In France, approximately 50% of the adults from

25 to 30 years have anti-CMV antibodies. The majority of cases of primary infection with CMV are without clinical symptoms or distinct symptoms. Often mononucleosis syndromes lack heterophile antibodies. CMV is the source of 70% post transfusion mononucleosis syndromes (5,6,7). CMV is the primary cause of infections in transplant patients and during malignant hemopathy. CMV infections are also very frequent during the course of AIDS.

- Toxoplasma gondii Control not validated:** The presence of anti-Toxoplasma IgG can indicate past infection by Toxoplasma. It is

imperative to confirm this state of immunity by detecting with validated tests anti-Toxoplasma gondii IgG and IgM antibodies. Toxoplasmosis exists in two forms: acquired toxoplasmosis and congenital toxoplasmosis. The acquired form of toxoplasmosis (usually benign and often not apparent) is responsible for extremely rare, serious mononucleosis syndromes that occur mainly in immunocompromised subjects. The frequency of human toxoplasmosis in adults in France is 50 to 80 %.

## PERFORMANCE AND CHARACTERISTICS OF THE TEST

**ImmunoDOT™ EBV MONO G** test allows the detection of IgG anti-EBV VCA and IgG anti-EBV EBNA antibodies. **The presence of CMV and toxoplasma antigens is used only for the controls.**

### 1. IgG anti-VCA antibodies

A prospective study was performed in a reference center on 51 patient samples. The test has been compared to immunofluorescence (IF).

Dot \ IF	Positive	Negative
Positive	35	3
Negative	2	11

**Sensitivity: 94.6 %**  
**Specificity: 78.6 %**  
**Agreement: 90.2%**

### 2. IgG EBNA antibodies

Study carried out by a center of reference for the EBV N = 51: The test has been compared to immunofluorescence (IF).

Dot \ IF	Positive	Negative
Positive	34	1
Negative	1	15

**Sensitivity : 97.1 %**  
**Specificity : 93.8 %**  
**Agreement : 96.1%**

### 3. Reproducibility IgG anti-VCA spots

Reproducibility was performed on 3 sample ranges: strongly positive, positive and negative for IgG anti-VCA dots:

INTRA-ASSAY			INTER-ASSAYS		
ELISA level	N	VCA appearance	N	VCA appearance	
480 U/mL	+++	5	Dark circle	5	Dark circle
210 U/mL	+	5	Light circle	5	Pale circle
105 U/mL	-	5	No spot	5	No spot

### 4. Reproducibility IgG anti-EBNA dots

Reproducibility was performed on 3 samples ranges: strongly positive, positive and negative for IgG anti-EBNA dots:

INTRA-ASSAY			INTER-ASSAYS		
ELISA level	N	EBNA appearance	N	EBNA appearance	
600 U/mL	+++	5	Dark circle	5	Dark circle
330 U/mL	+	5	Pale circle	5	Pale circle
155 U/mL	-	5	No spot	5	No spot

**SUMMARY OF METHOD**

<b>SETUP</b>	<ul style="list-style-type: none"> <li>• Turn on and prepare the workstation (according to the model). <b>The temperature in vessels must lie within +44°C and +48°C.</b></li> <li>• Insert the strip in the distilled water bath (Wash) for 30 seconds.</li> <li>• Pipet:             <ul style="list-style-type: none"> <li>- 2 mL of diluent (#1) into reaction vessel #1</li> <li>- 2 mL of enhancer (#2) into reaction vessel #2</li> <li>- 2 mL of conjugate (#3) into reaction vessel #3</li> </ul> </li> <li>• Incubate for 10 min.</li> </ul>	<b>INCUBATION 10 minutes</b>
<b>1. INCUBATION OF SERA</b>	<ul style="list-style-type: none"> <li>• Add 10 µL of sample sera or positive control to vessel #1.</li> <li>• Immerse strip in reaction vessel #1.</li> <li>• Mix using a swift up and down motion 10 times.</li> <li>• Incubate for 60 min.</li> </ul>	<b>INCUBATION 60 minutes</b>
<b>2. WASH STEP</b>	<ul style="list-style-type: none"> <li>• Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane.</li> <li>• Blot the assay strip on an absorbent towel.</li> </ul>	<b>10 seconds</b>
<b>3. ENHANCER</b>	<ul style="list-style-type: none"> <li>• Place strip into vessel #2.</li> <li>• Mix using a swift up and down motion 10 times.</li> <li>• Incubate for 5 min.</li> </ul>	<b>INCUBATION 5 minutes</b>
<b>4. WASH STEP</b>	<ul style="list-style-type: none"> <li>• Wash the in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane.</li> <li>• Blot the assay strip on an absorbent towel.</li> </ul>	<b>10 seconds</b>
<b>5. CONJUGATE</b>	<ul style="list-style-type: none"> <li>• Place strip into vessel #3.</li> <li>• Mix using a swift up and down motion 10 times.</li> <li>• Incubate for 30 min.</li> <li>• At the start of this incubation, <b>place 2 mL of reagent #4 into vessel #4.</b></li> </ul>	<b>INCUBATION 30 minutes</b>
<b>6. WASH STEP</b>	<ul style="list-style-type: none"> <li>• Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane.</li> <li>• Immerse the strip in the Wash (distilled water). Let stand 5 minutes.</li> <li>• Blot the assay strip on an absorbent towel.</li> </ul>	<b>10 seconds then INCUBATION 5 minutes</b>
<b>7. DEVELOPER</b>	<ul style="list-style-type: none"> <li>• Place strip into vessel #4.</li> <li>• Mix using a swift up and down motion 10 times.</li> <li>• Incubate for 5 min.</li> </ul>	<b>INCUBATION 5 minutes</b>
<b>8. WASH STEP</b>	<ul style="list-style-type: none"> <li>• Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane.</li> </ul>	<b>10 seconds</b>
<b>9. DRY STEP</b>	<ul style="list-style-type: none"> <li>• Remove the residual water with absorbent paper, and <b>let dry before interpreting the result.</b></li> </ul>	<b>WAIT 5 minutes</b>

**BioMédical Diagnostics SA**

**Office**

Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)