ImmunoDOT™ EBV MONO-M

REF GB 6585

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US Patent 247,934 (MA Fletcher)

Biggar RJ, W Henle, G Fleisher, J Bocker, ET Lennette , G Henle 1978. Primary Epstein-Barr virus infections and African infants. I. Decline of maternal antibodies and time of infection Int J Cancer 22:239.

Biggar RJ, G Henle, J Bocker, ET Lennette, G Fleisher, W Henle 1978. Primary Epstein-Barr virus infections in African infants. II. Clinical and serological observations during seroconversion. Int J Cancer 22: 244.

Fleisher G, W Henle, G Henle, ET Lennette, RJ Biggar 1979. Primary Epstein-Barr virus infection in American infants: Clinical and serological observations. J Infect Dis 139: 553.

Sumaya CV, Y Ench 1985. Epstein-Barr Virus Infectious Mononucleosis in Children II. Heterophil Antibody and Viral-Specific Responses. Pediatrics 75:1011.

Lennette ET 1995. Epstein-Barr Virus (EBV). In Diagnostic Procedures for Viral, Rickettsial, and Chlamydial Infections, 7th Ed.(EH Lennette, DA Lennette, ET Lennette) American Public Health Association, Washington D.C.

Lennette ET, L Rymo, Y Mohan, G Masucci, K Merk, L Timar, G Klein 1993. Disease-related Differences in Antibody Patterns Against EBV-encoded Nuclear Antigens EBNA 1, EBNA 2 and EBNA 6. Eur J Cancer 29A: 1584.

Friedman HM 1981. Cytomegalovirus: subclinical infection or disease? Am J Med 70:215.

Lajo A, C Borque, F Del Castillo, A Martin-Ancel 1994. Mononucleosis caused by Epstein-Barr virus and cytomegalovirus in children: a comparative study of 124 cases. Pediat Infect Dis J 13:56.

Evans AS 1969. Infectious mononucleosis and other mono-like syndromes. N Eng J Med 286:836.

Hofgartner WT, SR Swanzy, RM Bacina, J Condon, M Gupta, PE Matlock, DL Bergeron, JJ Plorde, TR Fritsche 1997. Detection of Immunoglobulin G (IgG) and IgM Antibodies to Toxoplasma gondii: Evaluation of Four Commercial Immunoassay Systems. J Clin Micro 35:3313.

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ImmunoDOT™ EBV MONO-M

REF GB 6585

Qualitative detection of anti-EBV VCA-IgM and heterophile antibodies that may appear during or after EBV infection.

INTENDED USE

The **ImmunoDOT™ EBV MONO-M Test** is a qualitative enzyme immunoassay (EIA) for the rapid detection of IgM antibodies to Paul-Bunnell heterophile, Epstein-Barr virus capsid antigen (EBV-VCA), and cytomegalovirus (CMV) in human serum. This test is intended to be used for the serodiagnosis of any prior infections by EBV. When used in conjunction with the **ImmunoDOT™ EBV MONO-G Test** it is an aid in the serodiagnosis of infectious (EBV) mononucleosis.

SUMMARY AND EXPLANATION

Infectious mononucleosis (IM) associated with Epstein-Barr Virus (EBV), is a transitory, acute and febrile lymphoproliferative syndrome with a benign course. It is widespread with over 95% of adults having been infected.

Transmission between people is most often through the saliva, and infection is seen earlier with poor sanitation conditions.. After penetration into the oral pharynx, the virus multiplies in situ, where the B lymphocytes are secondarily infected during passage into the pharyngeal lymphoid tissue. Venereal transmission would appear to be possible.

EBV Infection may be expressed as:

- an asymptomatic seroconversion;
- a moderate mononucleosidic syndrome with slight fever, sore throat, and adenopathies, particularly in young children;
- an infectious mononucleosis (IM) with sore throat, adenopathies, asthenia (weakness), sensitivity to certain antibiotics (Ampicillin) particularly in adolescents;
- a mononucleosidic syndrome with chronic fatigue ;
- occasional malignant pathologies related to the persistence of the virus in the host such as lymphoma in the immunocompromised patient, undifferentiated nasopharyngeal carcinoma, oral hairy cell leukemia,, and Burkitt's lymphoma.

A well defined kinetic assay shows that the antibodies directed against the various EBV antigens appear during the course of the primary infection.

Anti-VCA IgG and IgM antibodies are against capsid antigens. Anti-VCA IgM antibodies appear only during acute infection or rarely during reactivation and disappear within 4 to 8 weeks. Anti-VCA IgG antibodies appear in all primary infections with EBV. After infection, these antibodies decline and persist for life.

Anti-EBNA IgG are antibodies to nuclear antigens. These antibodies appear late after EBV infection. They are always absent during the acute phase of primary infection.

The **ImmunoDOT™ EBV MONO-M** test detects, anti-EBV VCA-IgM and heterophile antibodies that may appear during or after EBV infection.

Anti-EA antibodies to Early Antigens (EA) appear in 60 % of infectious mononucleosis and disappear quickly. However, an anti-EA antibody level is demonstrated during a viral reactivation among immunodeficient patients.

Asymptomatic forms are frequent for children. In no specific clinical forms, the diagnosis can be difficult. Indeed, other infectious agents such as cytomegalovirus (CMV), Toxoplasma gondii, Herpes Virus (HSV), Rubella, or Human Immunodeficiency Virus (HIV) may cause mononucleosis syndromes.

Although the first viral agent is Epstein-Barr Virus and the infectious mononucleosis incidence due to each one of these agents varies according to the age and to the socio-economic class.

EBV infection is serologically diagnosed by concurrent detection of both IgM and IgG antibodies (EBNA IgG, VCA IgG, VCA IgM).

Cytomegalovirus (CMV) control has not been validated. Ignore the results obtained for this parameter. Testing for anti-CMV antibodies will have to be performed using validated serological tests.

ASSAY PRINCIPLE

The **ImmunoDOT™ EBV MONO-M** uses an enzyme-linked immunoassay (EIA) dot technique for the detection of the anti-VCA IgM and heterophile antibodies. The different antigens are present as discrete dots on the assay strip.

First, an assay strip is incubated in a vessel with dilute patient serum, thus allowing any antibodies present in the sample to react with the test antigens immobilized on the assay strip.

In the second step, the specificity is enhanced by removal of non-specifically bound materials.

In the third step the strip is exposed to alkaline phosphatase conjugated to anti-human IgM antibodies so that the enzyme conjugate may bind to any previously bound patient antibodies.

Finally, the strip is transferred to the vessel containing the enzyme substrate reagent which then reacts with the conjugated alkaline phosphatase, producing distinct, clearly visible dots.


STORAGE CONDITIONS

- Store reagents and assay strips between +2°C and +8°C.
- Reagents must be at room temperature (between +15°C and +30°C) before use. They should be used within an hour following placement in the heated workstation (incubator)
- Avoid contamination of reagents as that can give erroneous results.

MATERIALS REQUIRED BUT NOT PROVIDED

- Workstation (incubator)
- Specimen collection apparatus
- Timer
- Analytical quality distilled or deionized water to be used for the washes.
- Pipets
- Absorbent toweling to blot dry the assay strip

MATERIAL PROVIDED

ImmunoDOT™ EBV MONO-M Assay Strip	50
 <p>STRIP</p> <ul style="list-style-type: none"> ⊗ Positive control ⊗ Level 1: Heterophile Antigen (purified from bovine red blood cells) (more sensitive level) ⊗ Level 2: Heterophile Antigen (less sensitive level) ⊗ EBV-VCA (gp 125, purified by affinity chromatography) ⊗ CMV antigen control (glycine lysate of strain AD-169) ○ Negative control 	
Diluent DIL SPE <i>Consists of buffer (pH 6.2-7.6) containing IgG absorbent (goat anti-human IgG, against rheumatoid factor or Immunoglobulin G), protein stabilizers, and <0.1% NaN₃</i>	2 x 50mL
Enhancer SOLN ENH <i>Consists of sodium chloride and <0.1% NaN₃</i>	2 x 50mL
Conjugate CONJ IgM <i>Consists of alkaline phosphatase conjugated goat anti-human IgM (heavy chain specific) in buffered diluent (pH 6.2-7.6), with protein stabilizers.</i>	2 x 50mL
Developer SUBS BCIP-NBT <i>Consists of 5-bromo-4-chloro-3-indolyl phosphate and p-nitro blue tetrazolium chloride in buffered diluent (pH 8.0-11.0), 0.8% dimethylformamide and <0.1% NaN₃.</i>	2 x 50mL
Positive EBV control CONTROL + <i>consist of human sera containing IgM anti-VCA antibodies and heterophile antibodies and <0.1 % d'Azide de Sodium.</i>	1 x 100µL
Positive CMV control CONTROL + <i>consist of human sera containing anti-CMV antibodies and <0.1 % d'Azide de Sodium.</i>	1 x 100µL
Reaction Vessels	200
Package Insert	1

SAMPLE COLLECTION AND HANDLING

- The **ImmunoDOT™ EBV MONO-M Test** is performed on serum. The test requires 10 µL of serum.
- Avoid lipemic or hemolyzed serum specimens has not been shown to be an acceptable specimen.
- Collect the samples of serum aseptically and store them between +2°C and +8°C for a maximum of 48 hours.
- If the assay of the samples will not be completed within 48 hours, or for shipment of the specimens, freeze at -20°C.

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use Only.

The **ImmunoDOT™ EBV MONO-M** reagents have been optimized for use as a system. Do not substitute reagents or strips from any other source including those from other **ImmunoDOT™** systems.

Dilution or modification of these reagents may also affect the performance of the test.

Do not use kits after the stated expiration date.

Analytical quality deionized or distilled water must be used for the washes.

Strictly adhere to the test procedures for optimal results. Do not shorten or prolong stated incubation times as this may result in poor assay performance.

Some assay components contain sodium azide (NaN₃) that may react with lead and copper plumbing to form highly explosive metal azides. After disposing of reagents, flush sink with a large volume of water to prevent azide buildup.

Human sera used in the preparation of this product were tested and found non-reactive for hepatitis B surface antigen and for antibodies to HIV-1, HIV-2, HTLV-1, and hepatitis C virus. Because no test method can offer complete assurance that infectious agents are absent, handle reagents and patient samples as if capable of transmitting infectious disease.

SETUP

1. Turn on Workstation and adjust temperature according to the instructions for use of your instrument.
2. Remove 4 Reaction Vessels per test from the product box and insert into appropriate slots in Workstation.
For the large Workstation, add water up to the fill line of the Wash Vessel provided.
For the small Workstation, use an appropriate container and sufficient water to cover all reactive areas/dots of the assay strip.
3. Appropriately label the Assay Strips.
4. Place:
 - 2 mL Diluent (#1) in Reaction Vessel #1
 - 2 mL Enhancer (#2) in Reaction Vessel #2
 - 2 mL Conjugate (#3) in Reaction Vessel #3
5. Incubate for 10 minutes.
6. Verify that the temperature of the heating block is in accordance with the instructions of the manual corresponding to the model of the workstation.
7. Insert the strip in the distilled water bath (Wash) for 30 seconds, just before starting the assay procedure. Avoid the formation of bubbles on the surface.

ASSAY PROCEDURE

1. SAMPLE INCUBATION

- Add 10 µL serum sample or positive control to Reaction Vessel #1.
- Plunge assay strip into Reaction Vessel #1
- Using 10 quick up and down motions with the Assay Strip, mix contents thoroughly.
- **Let stand for 60 minutes.**

2. WASH STEP

- Remove Assay Strip from Reaction Vessel and swish in the Wash (distilled water). Use a swift back and forth motion for 10 seconds allowing for optimal washing of the Assay Strip.
- Blot the strip on an absorbent towel.

3. ENHANCER

- Place Assay Strip into Reaction Vessel #2.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 5 minutes.**

4. WASH STEP

- Remove Assay Strip from Reaction Vessel #2 and swish in the Wash for 10 seconds as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.
- Blot the strip on an absorbent towel.

5. CONJUGATE

- Place Assay Strip into Reaction Vessel #3.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 30 minutes.**
- **At the beginning of this incubation, place 2mL of reagent #4 into the Reaction Vessel #4.**

6. WASH STEP

- Remove Assay Strip from Reaction Vessel #3 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip
- DO NOT remove the Assay Strip from the wash.
- **Allow the Assay Strip to stand in the wash for 5 minutes.**
- Blot the strip on an absorbent towel.

7. DEVELOPER

- Remove Assay Strip from wash and place into Reaction Vessel #4.
- Mix thoroughly with 10 quick up and down motions.
- **Let stand for 5 minutes.**

8. WASH STEP

- Remove Assay Strip from Reaction Vessel #4 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.

9. DRY

- Blot and allow Assay Strip to dry. It is imperative that tests of borderline specimens be interpreted after the Assay Strip has been allowed to dry. **A false positive dot may be identified if the assay strip is not dry when interpreted.**

10. READING THE ASSAY STRIP

Positive: A dot with an **EASILY SEEN**, distinct border is visible in the center of the window. The outer perimeter of the window must be white to pale gray.

Negative: If no dot is seen or a dot is difficult to see, interpret it as negative.

Each dot of the assay strip is read independently.

QUALITY CONTROL

The dots of positive and negative controls must appear correctly; otherwise the test needs to be repeated.

The color intensity of the positive control must not be used as a calibrator. The reactions in dots containing antigens could have a different intensity than the positive control, because of a difference in the concentration of antibodies.

The positive control serums included in this kit, are moderately positive for all antibodies. The performance of every kit can be confirmed upon receipt by using these controls which are expected to give a positive result.

The analysis is performed with reagents warmed at **+44°C/+ 48°C**. The temperature of the workstation has to exceed the temperature in the reaction vessels. The temperature of the workstation must be watched.

LIMITATIONS

- ⇒ This is a qualitative test and therefore it cannot be used to detect changes in antibody titers.
- ⇒ Primarily maternal antibodies will be detected in children less than one year of age.
- ⇒ Lipemic or hemolysed samples or samples contaminated by bacteria could provide erroneous results.
- ⇒ The **ImmunoDOT™ EBV MONO G** test used alone is not in any case diagnostic for a mononucleosis infectious due to the Epstein Barr virus.
- ⇒ The complete serologic profile of an infection to EBV should include determinations for the presence of antibodies directed against the different EBV-antigens appearing during the course of infection: EBV VCA IgM, EBV VCA IgG and EBV EBNA IgG.
- ⇒ If one wants to determine the state of reactivation in an immuno-compromised patient or show that the EBV is associated with a malignant disease, it is essential to determine titers of the anti-EA (early antigen) antibodies.
- ⇒ **The CMV control has not been validated. Ignore any results associated with this parameter. Any testing for anti-CMV antibodies must be performed using validated serologic tests.**
- ⇒ The interpretation of the results of the **ImmunoDOT™ EBV MONO M** test should be done in parallel with those obtained with the **ImmunoDOT™ EBV MONO G** test which provides for the detection of anti-VCA isotype IgM antibodies.
- ⇒ The performance of this device has not been established for semi-quantitative utility of H1 and H2.
- ⇒ Performance with this device has not been established for either prenatal screening or newborn testing, in a non-clinical laboratory environment (e.g., point of care testing).

INTERPRETATION

1. Reading and interpreting the VCA antibodies dots. The **ImmunoDOT™ EBV MONO M** kit is a test for the rapid, qualitative detection of EBV VCA IgM antibodies.
2. Reading and interpretation of the heterophil antibody wells : The two heterophil antigen wells (1 and 2 levels) are an help for interpretation. The level 1 is a weak heterophil antibody level and the level 2 is an higher level.
3. CMV Control well allows the detection of IgM anti-CMV antibodies. It's an indicative control.

ImmunoDOT™ MONO M		ImmunoDOT™ MONO G		Hematological and clinical symptoms	INTERPRETATION EBV infection steps	Frequency of the pattern	REMARKS
ACs* MNI**	VCA IgM	EBNA IgG	VCA IgG				
-	-	-	-	<ul style="list-style-type: none"> negative 	Seronegative patient (no infected)	adult: rarely child: frequent	
-/+	+	-	-	<ul style="list-style-type: none"> negative indeterminate characteristic 	Primary infection (early step)	adult: rarely child and teenager: frequent	Serological Testing in 3 weeks
- +	+ +	- -	+ +	<ul style="list-style-type: none"> negative indeterminate characteristic 	Primary infection	adult: rarely child and teenager: frequent	Serological Testing in 3 weeks
- -	+ +	+ -/+	+ +	<ul style="list-style-type: none"> negative indeterminate 	End of primary infection (passing step)	Frequent patterns for child and teenager	Serological Testing in 3 weeks
-	-	+	+	<ul style="list-style-type: none"> negative 	Past infection	Frequent pattern (95% of the adult population)	2 % of patients never have anti-EBNA antibodies

OTHERS

ImmunoDOT™ MONO M		ImmunoDOT™ MONO G		Hematological and clinical symptoms	INTERPRETATION EBV infection	Pattern and frequency	REMARKS
ACs* MNI**	VCA IgM	EBNA IgG	VCA IgG				
+/-	+	+	+	<ul style="list-style-type: none"> negative indeterminate characteristic 	<ul style="list-style-type: none"> acute chronic infection? Reactivation ? Recent infection with evolution? 	Frequent patterns in immunodeficient patient with transplantation Rarely in other patients	Test following serological parameters : 1. IgM VCA 2. IgG VCA 3. IgG EA 4. IgG and IgM CMV 5. IgG and IgM Toxo

(*) ACs = heterophile antibodies (**) MNI = Infectious MonoNucleosis

Remarks :

- IgG anti-VCA and IgM antibodies combined with heterophile antibodies or no confirm a current infection. For children, the primary infection is often with indeterminate symptoms.
- Heterophile antibodies negative results for do not rule out the diagnosis of an accurate mononucleosis infection.
- Heterophile antibodies positive results may occur a recent or current mononucleosis infection. They appear into 2 to 3 weeks and disappear into 1 to 3 months. In addition, CMV, HSV, toxoplasmosis, hepatitis B, rubella infections, also induce heterophile antibodies. It should however be noted that in 20 to 30 % of the cases, the heterophile antibodies do not appear (or very tardily) especially for child.
- **CMV Control:** The IgM anti-CMV detection does not allow the diagnosis of the infection by CMV. **The CMV well is a control well; results should be confirmed by other validated IgG and IgM CMV tests.** CMV infections are frequent and ubiquitous. The prevalence of the infection is higher in populations living in conditions of poor sanitation or in densely populated areas (4). In France, approximately 50% of the adults from 25 to 30 years have anti-CMV antibodies. The majority of cases of primary infection with CMV are without clinical symptoms or distinct symptoms. Often mononucleosis syndromes lack heterophile antibodies. CMV is the source of 70% post transfusion mononucleosis syndromes (5,6,7). CMV is the primary cause of infections in transplant patients and during malignant hemopathy. CMV infections are also very frequent during the course of AIDS.

PERFORMANCE AND CHARACTERISTICS OF THE TEST

ImmunoDOT™ EBV MONO M test allows the detection of IgM anti-EBV VCA and heterophile antibodies. **The presence of anti-CMV IgM antibodies is used only for the controls.**

1. IgM anti-VCA antibodies

A prospective study was performed in a reference center on 95 patient samples. The test has been compared to immunofluorescence (IF).

Dot \ IF	Positive	Negative
Positive	40	2
Negative	1	52

Sensitivity 40/41: 97.6 %
Specificity 52/54: 96.3 %
Agreement 92/95: 96.8 %

2. Heterophile antibodies

A prospective study was performed in a reference center on 61 patient samples.

The specificity was estimated on a seronegative population for EBV and a population with post infection serology.

The test sensibility was estimated on a positive population for Paul-Bunnell and Davidsohn (PBD) technique and with a confirmed primary infection.

Dot \ IF	Positive	Negative
Positive	26	1
Negative	1	33

Sensitivity 26/27: 96 %
Specificity 33/34: 97 %
Agreement 59/61: 97 %

3. Reproducibility IgM anti-VCA spots

Positive and negative samples in ELISA were tested 5 times with **ImmunoDOT™ EBV MONO M** Test. Reproducibility was performed on 3 samples: high positive, positive and negative :

Results for

Elisa level	INTRA-ASSAY		INTER-ASSAYS	
	N	VCA appearance	N	VCA appearance
500 U/mL +++	5	Dark circle	5	Dark circle
280 U/mL +	5	Light circle	5	Pale circle
30 U/mL -	5	No spot	5	No spot

4. Reproducibility Heterophile antibody dots

Reproducibility was performed on 3 samples ranges: strongly positive, positive and negative for IgG anti-EBNA dots:

PBD level	INTRA-ASSAY		INTER-ASSAYS	
	N	Heterophile appearance	N	Heterophile appearance
+++	5	Dark circle: level 1 + 2	5	Dark circle: level 1 + 2
+	5	Pale circle: level 1	5	Pale circle: level 1
-	5	No spot:	5	No spot

REFERENCES

US Patent 247,934 (MA Fletcher)

Biggar RJ, W Henle, G Fleisher, J Bocker, ET Lennette, G Henle 1978. Primary Epstein-Barr virus infections and African infants. I. Decline of maternal antibodies and time of infection Int J Cancer 22:239.

Biggar RJ, G Henle, J Bocker, ET Lennette, G Fleisher, W Henle 1978. Primary Epstein-Barr virus infections in African infants. II. Clinical and serological observations during seroconversion. Int J Cancer 22: 244.

Fleisher G, W Henle, G Henle, ET Lennette, RJ Biggar 1979. Primary Epstein-Barr virus infection in American infants: Clinical and serological observations. J Infect Dis 139: 553.

Sumaya CV, Y Ench 1985. Epstein-Barr Virus Infectious Mononucleosis in Children II. Heterophil Antibody and Viral-Specific Responses. Pediatrics 75:1011.

Lennette ET 1995. Epstein-Barr Virus (EBV). In Diagnostic Procedures for Viral, Rickettsial, and Chlamydial Infections, 7th Ed.(EH Lennette, DA Lennette, ET Lennette) American Public Health Association, Washington D.C.

Lennette ET, L Rymo, Y Mohan, G Masucci, K Merk, L Timar, G Klein 1993. Disease-related Differences in Antibody Patterns Against EBV-encoded Nuclear Antigens EBNA 1, EBNA 2 and EBNA 6. Eur J Cancer 29A: 1584.

Friedman HM 1981. Cytomegalovirus: subclinical infection or disease? Am J Med 70:215.

Lajo A, C Borque, F Del Castillo, A Martin-Ancel 1994. Mononucleosis caused by Epstein-Barr virus and cytomegalovirus in children: a comparative study of 124 cases. Pediat Infect Dis J 13:56.

Evans AS 1969. Infectious mononucleosis and other mono-like syndromes. N Eng J Med 286:836.

Hofgartner WT, SR Swanzy, RM Bacina, J Condon, M Gupta, PE Matlock, DL Bergeron, JJ Plorde, TR Fritsche 1997. Detection of Immunoglobulin G (IgG) and IgM Antibodies to Toxoplasma gondii: Evaluation of Four Commercial Immunoassay Systems. J Clin Micro 35:3313.

SYMBOLS USED

EC Declaration of conformity

Rapid test

In Vitro Diagnostic Device

Catalogue number

Lot Number

Expiry Date

Number of test

Consult Instructions

Temperature limitation

Biological risk

Contains sodium azide

Reconstitute with

SUMMARY OF METHOD

SETUP	<ul style="list-style-type: none"> • Turn on and prepare the workstation (according to the model). The temperature in vessels must lie within +44°C and +48°C. • Insert the strip in the distilled water bath (Wash) for 30 seconds. • Pipet: <ul style="list-style-type: none"> - 2 mL of diluent (#1) into reaction vessel #1 - 2 mL of enhancer (#2) into reaction vessel #2 - 2 mL of conjugate (#3) into reaction vessel #3 • Incubate for 10 min. 	INCUBATION 10 minutes
1. INCUBATION OF SERA	<ul style="list-style-type: none"> • Add 10 µL of sample sera or positive control to vessel #1. • Immerse strip in reaction vessel #1. • Mix using a swift up and down motion 10 times. • Incubate for 60 min. 	INCUBATION 60 minutes
2. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Blot the assay strip on an absorbent towel. 	10 seconds
3. ENHANCER	<ul style="list-style-type: none"> • Place strip into vessel #2. • Mix using a swift up and down motion 10 times. • Incubate for 5 min. 	INCUBATION 5 minutes
4. WASH STEP	<ul style="list-style-type: none"> • Wash the in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Blot the assay strip on an absorbent towel. 	10 seconds
5. CONJUGATE	<ul style="list-style-type: none"> • Place strip into vessel #3. • Mix using a swift up and down motion 10 times. • Incubate for 30 min. • At the start of this incubation, place 2 mL of reagent #4 into vessel #4. 	INCUBATION 30 minutes
6. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Immerse the strip in the Wash (distilled water). Let stand 5 minutes. • Blot the assay strip on an absorbent towel. 	10 seconds then INCUBATION 5 minutes
7. DEVELOPER	<ul style="list-style-type: none"> • Place strip into vessel #4. • Mix using a swift up and down motion 10 times. • Incubate for 5 min. 	INCUBATION 5 minutes
8. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Water (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. 	10 seconds
9. DRY STEP	<ul style="list-style-type: none"> • Remove the residual water with absorbent paper, and let dry before interpreting the result. 	WAIT 5 minutes

BioMédical Diagnostics SA

Office

Actipole 25
4 blvd de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com