ImmunoDOT™ Typage du Virus Herpès Simplex IgG

REF GB 7125

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BIBLIOGRAPHIE

Nahmias, A., and W. Dowdle. 1968 Antigenic and biologic differences in herpesvirus hominis. *Prog. Med. Virol.* 10:110-159.

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McGeoch, D.J., H.W. Moss, D. McNab, and M.C. Frame. 1987. DNA sequence and genetic content of the HindIII region in the short unique component of herpes simplex virus type 2 genome: identification of the gene encoding glycoprotein G, and evolutionary comparisons. *J. Gen. virol.* 68:19-38.

Rapoport, T.A. 1986. Protein translocation across and integration into membranes. *Crit Rev biochem.* 20:73-137.

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
ImmunoDOT™ Herpes Simplex Virus Typing IgG

GB 7125

Enzyme immunoassay (EIA) detecting HSV or glycoprotein G (gG) type specific IgG or IgM antibodies.

INTENDED USE

ImmunoDOT™ HSV IgG test is an enzyme immunoassay (EIA) detecting HSV or glycoprotein G (gG) type specific IgG antibodies. The test detects the presence or absence of past HSV exposure and specifically determines whether past infection(s) is due to HSV Type 1, Type 2 or both Type 1 and 2.

SUMMARY AND EXPLANATION

There are two herpes antigenic types (1,2). "Definitive diagnosis of genital herpes infections is fundamental to the management of patients and the development of strategies to prevent transmission to partners and neonates" (3). Such diagnosis has proven inaccurate when based solely on clinical history and impression (4). Instead, virus, antigen or nucleic acid detection and classification is used for patients presenting with lesions or type-specific serological tests may be used when lesions are absent.

For type specific serology, either western blot (5,6,7) or determination of type specific protein by molecular biology (8,9) is used. Acceptable type specific classification is not possible using whole virus lysate, the commonly used antigen of early HSV serology kits. The most commonly used type specific protein is glycoprotein G. ImmunoDOT HSV Typing test uses HSV gG type 1 and type 2 specific proteins.

ASSAY PRINCIPLE

The **ImmunoDOT™ HSV IgG** assay uses an enzyme-linked immunoassay (EIA) dot technique for the detection of antibodies to HSV-1 gG, HSV-1 virus lysate proteins, HSV-2 gG or HSV-2 virus lysate proteins.

An assay strip is incubated with dilute patient serum, allowing patient antibodies reactive with the test antigens to bind to the solid support membrane.

In the second stage, the reaction is enhanced by removal of non-specifically bound materials.

During the third stage, alkaline phosphatase-conjugated antihuman antibodies are allowed to react with bound patient antibodies.

Finally, the strip is transferred to enzyme substrate reagent which reacts with bound alkaline phosphatase to produce an easily seen, distinct dot.

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use Only.

These reagents have been optimized for use as a system. Do not interchange with other reagents or other ImmunoDOT™ assay system.

Dilution or adulteration of these reagents may also affect the performance of the test.

Do not use kit if evidence of microbial contamination is present. Do not use any kits beyond the stated expiration date.

Analytic quality deionized or distilled water must be used as Clarifier.

Close adherence to the test procedure will assure optimal performance. Do not shorten or lengthen stated incubation times since this may result in poor assay performance.

Some assay components contain sodium azide (NaN₃) that may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide buildup.

Human sera used in the preparation of this product were tested and found non-reactive for hepatitis B surface antigen and for antibodies to HIV-1, HIV-2, HTLV-1, and hepatitis C virus. Because no test method can offer complete assurance that infectious agents are absent, handle reagents and patient samples as if capable of transmitting infectious disease.

STORAGE CONDITIONS

- Store reagents and assay strips between +2°C and +8°C.
- Reagents must be at room temperature (between +15°C and +30°C) before use. They must be used within hour which follows their implementation in the warmed workstation.
- Avoid contamination of reagents as that can give erroneous results.


SAMPLE COLLECTION AND HANDLING

- The **ImmunoDOT™ HSV IgG Test** is performed on serum. The test requires 50µL of serum.
- Avoid lipemic or hemolyzed serum specimens has not been shown to be an acceptable specimen.
- Store samples at room temperature for no longer than eight hours. If the assay will not be completed within eight hours, store the sample between +2°C and +8°C.
- If the assay of the samples will not be completed within 48 hours, freeze at -20°C.

MATERIALS REQUIRED BUT NOT PROVIDED

- Workstation (incubator)
- Specimen collection apparatus
- Timer
- Analytical quality distilled or deionized water to be used for the washes.
- Pipets
- Absorbent toweling to blot dry the assay strip

MATERIAL PROVIDED

ImmunoDOT™ HSV IgG Assay Strips		25
		
<i>Antigens in order beginning with the window next to the label: Positive reagent control (human serum); [Total 1] HSV-1, [Total 2] HSV-2 lysate, [Specific 1] HSV-1 gG, [Specific 2] HSV-2 gG and Negative reagent control.</i>		
Diluent	<input type="checkbox"/> DIL <input type="checkbox"/> SPE	1 x 50mL
Consists of buffered diluent containing protein stabilizers, and <0.1% NaN ₃ .		
Enhancer	<input type="checkbox"/> SOLN <input type="checkbox"/> ENH	1 x 50mL
Consists of sodium chloride and <0.1% NaN ₃ .		
Conjugate	<input type="checkbox"/> CONJ <input type="checkbox"/> IgG	1 x 50mL
Consists of alkaline phosphatase conjugated goat antihuman IgM (heavy chain specific) in buffered diluent, protein stabilizers.		
Developer	<input type="checkbox"/> SUBS <input type="checkbox"/> BCIP-NBT	1 x 50mL
Consists of 5-bromo-4-chloro-3-indolyl phosphate and p-nitro blue tetrazolium chloride in buffered diluent and <0.1% NaN ₃ .		
HSV-1 Positive Control	<input type="checkbox"/> CONTROL <input type="checkbox"/> +	1 x 300μL
Consists of HSV-1 IgG positive human serum that is HSV-2 negative containing protein stabilizers and <0.1% NaN ₃ .		
HSV-2 Positive Control	<input type="checkbox"/> CONTROL <input type="checkbox"/> +	1 x 300μL
Consists of HSV-2 IgG positive human serum that is HSV-1 negative containing protein stabilizers and <0.1% NaN ₃ .		
Reaction Vessels		100
Package Insert		1

SETUP

- Turn on Workstation and adjust temperature according to the instructions for use of your instrument.
- Remove 4 Reaction Vessels per test from the product box and insert into appropriate slots in Workstation.
For the large Workstation, add water up to the fill line of the Wash Vessel provided.
For the small Workstation, use an appropriate container and sufficient water to cover all reactive areas/dots of the assay strip.
- Appropriately label the Assay Strips.
- Place:
 - 2 mL Diluent (#1) in Reaction Vessel #1
 - 2 mL Enhancer (#2) in Reaction Vessel #2
 - 2 mL Conjugate (#3) in Reaction Vessel #3

- Incubate for 10 minutes.
- Verify that the temperature of the heating block is in accordance with the instructions of the manual corresponding to the model of the workstation.
- Insert the strip in the distilled water bath (Wash) for 30 seconds, just before starting the assay procedure. Avoid the formation of bubbles on the surface.
- If the large Workstation is used, insert the label end of the assay strip into the Strip Holder, one per groove, taking care not to touch the assay windows.

ASSAY PROCEDURE

1. SAMPLE INCUBATION

- Add 50 μL serum to Reaction Vessel #1.
- Plunge assay strip into Reaction Vessel #1
- Using 10 quick up and down motions with the Assay Strip, mix contents thoroughly.
- Let stand for 60 minutes.**

2. WASH STEP

- Remove Assay Strip from Reaction Vessel and swish in the Wash (distilled water). Use a swift back and forth motion for 10 seconds allowing for optimal washing of the Assay Strip.
- Blot the strip on an absorbent towel.

3. ENHANCER

- Place Assay Strip into Reaction Vessel #2.
- Mix thoroughly with 10 quick up and down motions.
- Let stand for 5 minutes.**

4. WASH STEP

- Remove Assay Strip from Reaction Vessel #2 and swish in the Wash for 10 seconds as described (step #2). The wash should diffuse completely through the strip.
- Blot the strip on an absorbent towel.

5. CONJUGATE

- Place Assay Strip into Reaction Vessel #3.
- Mix thoroughly with 10 quick up and down motions.
- Let stand for 30 minutes.**
- At the beginning of this incubation, place 2mL of reagent #4 into the Reaction Vessel #4.**

6. WASH STEP

- Remove Assay Strip from Reaction Vessel #3 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.
- DO NOT remove the Assay Strip from the wash.
- Allow the Assay Strip to stand in the wash for 5 minutes.**
- Blot the strip on an absorbent towel.

7. DEVELOPER

- Remove Assay Strip from wash and place into Reaction Vessel #4.
- Mix thoroughly with 10 quick up and down motions.
- Let stand for 5 minutes.**

8. WASH STEP

- Remove Assay Strip from Reaction Vessel #4 and swish in the Wash as described (step #2) for 10 seconds. The wash should diffuse completely through the strip.

9. DRY

- Blot and allow Assay Strip to dry. It is imperative that tests of borderline specimens be interpreted after the Assay Strip has been allowed to dry. **A false positive dot may be identified if the assay strip is not dry when interpreted.**

10. READING THE ASSAY STRIP

Positive: A dot with an **EASILY SEEN**, distinct border is visible in the center of the window. The outer perimeter of the window must be white to pale gray.

Negative: If no dot is seen or a dot is difficult to see, interpret it as negative.

Each dot of the assay strip is read independently.

In order to minimize the possibility of "over interpreting" positive test results it is recommended that during initial validation of the assay (as may be required by the laboratory by regulation), the laboratory test a series of presumptive negative samples and each technician interpret the assay strips in a blinded fashion. Please call bmd Technical Service for further clarification (To report results, refer to Interpretation Section).

QUALITY CONTROL

The assay is performed with reagents at **+44°C/+48°C**. The temperature of the Workstation itself should be higher than the assay temperature in the reaction vessel. The Workstation should be monitored to assure the appropriate temperature is maintained.

NCCLS C24-A may be consulted for guidance on appropriate quality control practices. These should be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Unless otherwise required, it is recommended that control sera be tested upon receipt of a kit. If the control is not reactive, results should not be reported and bmd Technical Service should be contacted before the kit is used again.

The kit uses reagent controls to assure performance each time a test is performed. The Positive Control window (well #1) contains human serum and tests reagent reactivity. It must be reactive but the intensity must not be used as a calibrator. As a negative reagent control check, the backgrounds around dots and the bottom window (#6) must be white. If the positive is not reactive or the negative is reactive, do not interpret the assay strip.

INTERPRETATION

Interpretation	Total		Specific	
	1	2	1	2
Negative	-	-	-	-
HSV Type 1	+	+	+	-
HSV Type 2	+	+	-	+
HSV Type 1&2	+	+	+	+
HSV Positive Type Unknown	+	+	-	-
Non Interpretable	-	-	+	-
Non Interpretable	-	-	-	+
Non Interpretable	-	-	+	+

- All results from this and other serologies must be correlated with clinical history, epidemiological data and other data available to the attending physician in evaluating the patient.
- It is expected that some patient samples have both HSV Total 1 and 2 positive and that both Specific 1 and 2 negative. This indicates either a recent infection or that the patient never formed type specific antibody. If a recent infection, isolation, antigen or nucleic acid detection methods are indicated. Otherwise, another method such as western blot may be used. However, it is also possible that no type specific antibody is present and therefore typing is not possible.
- As with other serological tests, negative results do not rule out the diagnosis of herpes simplex disease.

- A single positive result only indicates previous immunologic exposure. Other methods such as nucleic acid, antigen or viral culture are required to establish current infection.
- It is unlikely that patterns other than presented will occur. Repeat the testing and if the pattern repeats, contact the manufacturer.

LIMITATIONS

The performance of this assay has not been established for ruling out diseases with similar symptoms, e.g., *Candida albicans*, *Bacteriodes* species, *G. vaginalis*, *Mobiluncus* species. Instead, also use culture or other appropriate methods.

EXPECTED RESULTS

It has been shown that gG specific antibody typically develops over the course of six months. Appearance is slower and later than antibodies to other, non-type specific proteins. Therefore, it is expected that some patients will be HSV antibody positive yet anti-gG negative.

Eighty-four sera collected from U.S. blood donors were tested using **ImmunoDOT™ HSV IgG** kit. These results are shown in Table 1.

Table 1

Interpretation	Percent
Negative	15.5%
HSV Type 1	36.9%
HSV Type 2	11.9%
HSV Type 1&2	32.1%
HSV Positive Type Unknown	3.6%

PERFORMANCE AND CHARACTERISTICS OF THE TEST

Eighty-four sera collected from U.S. blood donors were tested using: **ImmunoDOT™ HSV IgG** kit, an alternate commercial microtiter kit for specific detection of either type 1 or type 2 gG antibodies [Kit A], and two alternate commercial microtiter kits for detection of HSV lysate antibodies (Not able to specify type) [Kits B and C].

Total (Traditional Method)

Comparison between **ImmunoDOT™** Total 1 and 2 results and two commercial microtiter kits (Kits B and C) is shown in Table 2. Agreement is 100% and therefore relative sensitivity and specificity are 100%.

Table 2

Alternate Kit	ImmunoDOT™	
	Negative	Positive
Negative	13	0
Positive	0	71

Specific (HSV Serology Type)

ImmunoDOT™ Specific results and the alternate commercial type specific microtiter kit comparison, using the eighty-four blood donor samples, is shown in Table 3. One sample is HSV antibody negative using kits B and C, reported negative by **ImmunoDOT™** but reported HSV Type 2 by the Alternate gG assay. This is considered a false positive and not shown in Table 3. Agreement between the two gG assays is 90% (75/83).

Table 3

ImmunoDOT™	Alternate	Number
Negative	Negative	15 (18%)
Type 1	Type 1	32 (39%)
Type 2	Type 2	9 (11%)
Type 1&2	Type 1&2	22 (27%)
Type 1	Negative	0
Type 1	Type 2	0
Type 1	Type 1&2	0
Type 2	Negative	0
Type 2	Type 1	0
Type 2	Type 1&2	2 (2%)
Type 1&2	Negative	0
Type 1&2	Type 1	2 (2%)
Type 1&2	Type 2	0
Negative	Type 1	1 (1%)
Negative	Type 2	0
Negative	Type 1&2	0

There is 94% (78/84) agreement between the two methods. Assuming alternate test is correct (relative performance) is listed in Table 4.

Table 4

Interpretation	Relative Sensitivity	Relative Specificity
Negative	100% (15/15)	100% (15/15)
Type 1	95% (56/59)	100% (15/15)
Type 2	94% (31/33)	100% (15/15)

Precision

Like any visual interpreted test, dot intensity is directly related to precision. The darkest dots are most reliable while weaker reactions (less intensity) are proportionately less reliable (equivocal or borderline).

REFERENCES

- Nahmias, A., and W. Dowdle. 1968 Antigenic and biologic differences in herpesvirus hominis. *Prog. Med. Virol.* 10:110-159.
- Schnewels, K.E. 1962. Serologische untersuchungen zur typendifferenzierung des herpesvirus hominis. *Z. Immunitaetsforsch.* 124:24.48.
- Ashley, R.L., and A Wald. 1999. Genital herpes: review of the epidemic and potential use of type-specific serology. *Clin. Microbiol. Rev.* 12(1): 1-8.
- Brown, Z. A., J.K. Benedetti, D.H. Watts, S. Selke, S. Berry, R. Ashley, and L. Corey. 1995. A comparison between detailed and simple histories in the diagnosis of genital herpes complicating pregnancy. *Am. J. Obstet. Gynecol.* 172:1299-1303.
- Bernstein, D.I., M.A. Lovett, and Y.J. Bryson. 1984. Serologic analysis of first-episode nonprimary genital herpes simplex virus infection. *Am. J. Med.* 77:1055-1060.
- Ho, d.W.T., P.R. Fireld, E. Sjogren-Jansson, S. Jeansson, and A.L. Cunningham. 1992. Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). *J. Virol. Methods* 36:249-264.
- Holmberg, S.D., J.A. Stewart, G. russel, R.H. Byers, F.K. Lee, P.M. O'Malley, and A. J. Nahmias. 1988. Prior herpes simplex virus type 2 infection as a risk factor for HIV infection. *JAMA* 259:1048-1050.
- McGeoch, D.J., H.W. Moss, D. McNab, and M.C. Frame. 1987. DNA sequence and genetic content of the HindIII region in the short unique component of herpes simplex virus type 2 genome: identification of the gene encoding glycoprotein G, and evolutionary comparisons. *J. Gen. virol.* 68:19-38.
- Rapoport, T.A. 1986. Protein translocation across and integration into membranes. *Crit Rev biochem.* 20:73-137.

SYMBOLS USED

	EC Declaration of conformity
	Rapid test
	<i>In Vitro</i> Diagnostic Device
	Catalogue number
	Lot Number
	Expiry Date
	Number of test
	Consult Instructions
	Temperature limitation
	Biological risk
	Contains sodium azide

SUMMARY OF METHOD

SETUP	<ul style="list-style-type: none"> • Turn on and prepare the workstation (according to the model). The temperature in vessels must lie within +44°C and +48°C. • Place reaction Vessels into slots in Workstation and add water to the Clarifier Vessel. • Insert the strip in the distilled water bath (Wash) for 30 seconds. • Pipet: <ul style="list-style-type: none"> - 2 mL of diluent (#1) into reaction vessel #1 - 2 mL of enhancer (#2) into reaction vessel #2 - 2 mL of conjugate (#3) into reaction vessel #3 • Incubate for 10 min. 	INCUBATION 10 minutes
1. INCUBATION OF SERA	<ul style="list-style-type: none"> • Add 50 µL of sample sera or positive control to vessel #1. • Immerse strip in reaction vessel #1. • Mix using a swift up and down motion 10 times. • Incubate for 60 min. 	INCUBATION 60 minutes
2. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Blot the assay strip on an absorbent towel. 	10 seconds
3. ENHANCER	<ul style="list-style-type: none"> • Place strip into vessel #2. • Mix using a swift up and down motion 10 times. • Incubate for 5 min. 	INCUBATION 5 minutes
4. WASH STEP	<ul style="list-style-type: none"> • Wash the in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Blot the assay strip on an absorbent towel. 	10 seconds
5. CONJUGATE	<ul style="list-style-type: none"> • Place strip into vessel #3. • Mix using a swift up and down motion 10 times. • Incubate for 30 min. • At the start of this incubation, place 2 mL of developer into vessel #4. 	INCUBATION 30 minutes
6. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Wash (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. • Immerse the strip in the Wash (distilled water). Let stand 5 minutes. • Blot the assay strip on an absorbent towel. 	10 seconds then INCUBATION 5 minutes
7. DEVELOPER	<ul style="list-style-type: none"> • Place strip into vessel #4. • Mix using a swift up and down motion 10 times. • Incubate for 5 min. 	INCUBATION 5 minutes
8. WASH STEP	<ul style="list-style-type: none"> • Wash the strip in the Water (distilled water) using a swift back and forth motion from left to right so that the water diffuses completely through the membrane. 	10 seconds
9. DRY STEP	<ul style="list-style-type: none"> • Remove the residual water with absorbent paper, and let dry before interpreting the result. 	WAIT 5 minutes

BioMédical Diagnostics SA

Office
Actipole 25
4 blvd de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com