CARDIO-LISA IgG

REF HM 001

DÉFINITION

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FONT J.
Anticardiolipin antibodies in patients with autoimmune diseases : isotype distribution and clinical associations.
Clin. Rheum. 1989, 8 : 475-785

HESS D.C.
Anticardiolipin antibodies : a study of frequency in TIA and stroke.
Neurology 1991, 23 (suppl) : 129-132

INFANTE RIVARD C. et al.
Lupus anticoagulant, anticardiolipin antibodies and fetal loss. A case control study.
N. Engl. J. Med. 1991, 325 : 1063-1066

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
96

DEFINITION

CARDIO-LISA IgG (bmd) is an enzyme linked immunoassay (ELISA) for the quantitative determination of anti-cardiolipin antibodies in human serum or plasma samples.

DIAGNOSTIC VALUE

Anti-cardiolipin autoantibodies are present in clinical contexts comparable to those in which lupus anticoagulants (anti-prothrombinase) are found. Moreover, they are frequently associated. In such patients, it is not rare to observe a dissociated syphilitic serology with positive VDRL and negative TPHA. The most frequent associated diseases are Disseminated Lupus Erythematosus, drug-induced lupus, thromboembolic disease and thrombosis in general, recurrent abortion, autoimmune peripheral cytopenia, primary biliary cirrhosis, certain acute infections and HIV and HTLV retrovirus infections. With the exception of infections, in which anti-cardiolipin antibodies seem only transiently present, the antibodies would seem to have a predictive value with respect to thrombosis. They exacerbate the life threat in myocardial infarction and enhance the risk of graft rejection following bypass procedures. The term, "anti-cardiolipin antibody syndrome" has been suggested to designate the morbid signs and symptoms that are associated with the presence of the antibodies. Anti-cardiolipin antibodies are generally of classes IgG, IgG and IgM, or IgG and IgA.

Indications for screening of anti-cardiolipin antibodies

The indications are mainly diagnostic and prognostic. The main clinical indications are:

- Thrombosis
- Disseminated Lupus Erythematosus
- Drug-induced lupus
- Hematological lupus with anti-mitochondrial antibodies, type 5
- Soulier-Boffa syndrome
- Recurrent abortion

At biological level:

- A lupus anticoagulant (anti-prothrombinase) confirmed or expected in the context of an increased activated partial thromboplastin time (APTT)
- A dissociated syphilitic serology: VDRL positive and TPHA negative
- Anti-mitochondrial antibodies of types 1 or 5
- Peripheral cytopenia, particularly that associated with autoimmune disease.

Anti-cardiolipin antibodies, lupus anticoagulant, dissociated syphilitic serology, type 5 anti-mitochondrial antibodies and a positive Coombs' test result may be associated in the same patient, but that is rather rare.

In contrast, cases in which only one or two antibodies are present are relatively frequent. For that reason, screening for antibodies must be conducted individually, since the absence of one does not necessarily exclude the presence of the other. An attempt has been made to explain their association as immunization against one or several phospholipids common to the screening tests. An alternative hypothesis is that antibodies with cross-reactivity to several phospholipids are involved. Human monoclonal antibodies yielding cross-reactions with DNA, cardiolipin and other phospholipids have been reported (Schoenfeld et al., 1983). That type of antibody is reported to be directed against the phosphodiester bond.

Anti-cardiolipin antibodies may be present in normal subjects. If that is the case, they are present at low levels and, frequently, transiently. Their significance is unknown.

ASSAY PRINCIPLE

Cardiolipin (diphosphatidylglycerol) is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the antigen coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Binding of anti-cardiolipin antibodies is detected using an anti-human IgG alkaline phosphatase labeled conjugate. After incubation, the wells are washed again to eliminate any excess of conjugate.
- The bound enzyme is revealed by addition of substrate (Para Nitro Phenyl Phosphate, producing a yellow colour. The colour intensity is proportional to anti-cardiolipin antibodies.
- After stopping the reaction by 1N sodium hydroxide, the optical density is read by a spectrophotometer at 405nm

Anti-cardiolipin levels of samples patients are expressed in GPL international units (defined by Harris et al.), using a standard curve.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

REAGENTS

96 wells microplate coated with Cardioliipin (12x8 individual breakaway wells with holder) MP	12 strips
4 Standards titered in GPL international units <u>Ready to use</u> <i>The values are shown on the vial label.</i> CAL n	4 x 1.5mL
1 Positive control titered in GPL international units <u>To dilute</u> <i>The values expected are shown on the vial label.</i> CONTROL +	1 x 350µL
1 vial of Negative control <u>To dilute</u> CONTROL -	1 x 350µL
1 vial of Goat anti-human IgG conjugate PAL <u>Ready to use</u> CONJ IgG	1 x 15mL
1 vial of Substrate containing 6x5mg PNPP tablets <u>To dilute</u> SUBS PNPP	1
1 vial of Phosphate buffer pH 7,2 (concentrated 10x) BUF WASH 10x <u>To dilute in distilled water</u>	1 x 100mL
1 vial of Substrate dilution buffer <u>Ready to use</u> BUF SUBS	1 x 35mL
1 vial of Solution A: dilution and washing buffer <u>To dilute</u> BUF Sol.A	1 x 15mL
1 vial of Stop solution (Sodium hydroxide 1N) <u>Ready to use</u> SOLN STOP	1 x 9mL

MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- timer
- incubator adjusted at +37°C
- physiological serum (NaCl 0.9%)
- microplate spectrophotometer with 405 nm filter
- 8 channel pipettes

SETUP

All the reagents should be prepared as required:

1. Preparation of Phosphate buffer solution

- Dilute concentrated Phosphate buffer 1/10 in distilled water.
- Storage period: 3 month at +2°C/+8°C.
NB. If there are crystals in the concentrated solution, warm the bottle to 37°C for 15 minutes before use.

2. Preparation of washing buffer

- Dilute solution A 1/100 in Phosphate buffer solution.

3. Preparation of dilution buffer (solution D)

- Dilute solution A 1/10 in Phosphate buffer.

4. Preparation of samples and controls.

- Dilute to 1/50 in dilution buffer (solution D) :
Ex: 10µL sample qsp 0,5mL solution D and shake vigorously at the vortex.

5. Preparation of ready-to-use conjugate.

- Estimate amount required for handling and transfer to a tube.

6. Preparation of the substrate.

- Dissolve the PNPP tablets in the substrate dilution buffer at a ratio of 1 tablet per 5 .
- Protect from light until use.

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C - +25°C) before handling. Check that all plates are well drained after each wash. Do not put plates on top of each other during incubation at +37°C.

Check that no condensed water is left at the bottom of the wells before reading.

Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

Reagents in solution (except for substrate buffer and stop solution) contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

CARDIO-LISA IgG has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

CARDIO-LISA IgG has been optimized for the use as describe in this procedure. Do not substitute reagents from other manufacturers. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Set out:

- one blank reagent well (made up with solution D: dilution buffer)
- in duplicate:
 - 4 standard wells
 - 1 well for each control
 - to obtain more precise results, it is also advisable to run samples in duplicate.
- then detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

2. Sample incubation

Add 100µL of standards, diluted controls or samples.

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 200µL of washing buffer

3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature.

Wash step:

Repeat the procedure described previously

Make a last wash with 200µL of 0,9% sodium chloride.

4. Incubation of substrate

Add 100µL substrate into each well

Incubate for 30 minutes in the dark at 37°C.

5. Stop solution

Add 50µL of 1N sodium hydroxide to each well.

6. Reading

Read the optical density of each well at 405 nm against the reagent blank set at zero absorbance within 30 minutes.

RESULTS

Standard curve

Determine the mean optical density for all duplicate readings.

Trace the standard curve on log-log paper, plotting the units of the 4 standard points along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).

Comments

- The mean OD of the standard 4 should be at least 0.8.
- The samples with values greater than that of standard 4 may be diluted to 1/200 to obtain a more precise result. The number of units should be multiplied by 4.

INTERPRETATION OF RESULTS

GPL units	< 10	10-15	> 15
Results	Negative	Weakly positive*	Positive

*Weakly positive results should be controlled on a second specimen obtained a few weeks later and interpreted in the light of the results of additional investigations and the clinical context.

QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive control concentration should fall within the range printed on the positive control label. If either reagent control is invalid, the test should be repeated.

bmd offers an external multiparameter quality control (Immunotrol I, Cat: HM036). It contains antibodies IgG directed against cardiolipin antigen. This control should be handled as samples.

REFERENCES

- ABUAF N. et al.
Vers une standardisation du dosage des anticorps anti-cardiolipine.
STV 1993, 5 : 663-671
- BEVER E.M.
Cofacteurs involved in the antiphospholipid syndrome.
Editorial.
Lupus 1992, 1 : 51-53
- FONT J.
Anticardiolipin antibodies in patients with autoimmune diseases : isotype distribution and clinical associations.
Clin. Rheum. 1989, 8 : 475-785
- HESS D.C.
Anticardiolipin antibodies : a study of frequency in TIA and stroke.
Neurology 1991, 23 (suppl) : 129-132
- INFANTE RIVARD C. et al.
Lupus anticoagulant, anticardiolipin antibodies and fetal loss. A case control study.
N. Engl. J. Med. 1991, 325 : 1063-1066
- LARAKI R.
Le syndrome des antiphospholipides
STV 1991, 325 : 363-369
- MEYER O., ROUQUETTE A-M et YOUINOU P.
Autoanticorps, Marqueurs des Maladies autoimmunes. BMD Editions (1999).

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	100µL	Solution D: blank reagent	30 min at room temperature
	100µL	Standards (double)	
	100µL	Diluted controls (double)	
	100µL	Diluted samples	
WASHING	Wash 3 times with washing buffer (200µL/ well)		
INCUBATION OF CONJUGATE	100µL	Anti-human IgG PAL conjugate	30 min at room temperature
WASHING	Wash 3 times with washing buffer (200 µL/ well) Wash once with NaCl 0.9% (200 µL/ well)		
ENZYME REACTION	100µL	Substrate solution	30 min at +37°C
STOP REACTION	50µL	NaOH 1N	

SYMBOLS USED

EC Declaration of conformity

ELISA Test

In Vitro Diagnostic Device

Catalogue number

Lot Number

Expiry Date

Number of test

Consult Instructions

Manufactured by

Temperature limitation

Biological risk

Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com