DNA-LISA

REF **HM 002**

DÉFINITION

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BIBLIOGRAPHIE

- AHMAN A et al. Molecular expression systems for anti-DNA antibodies - 1 - 2. Lupus 2002, 11, 824-832, 833-842.
- CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance. Clin. and Exp. Rheumatol. 1997 ;15, 457-465
- HOLLINGSWORTH P.N. et al. Antinuclear antibodies In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.74-90.
- HOMBURGER H.A. Cascade testing for autoantibodies in connective tissue disease. Mayo Clin. Proc., 1995, 70, 183-184.
- ISENBERG et al. Clinical laboratory assays for measuring anti-dsDNA antibodies. Where are we now? Lupus 2002, 11, 797-800.
- RAHMAN A et al. Anti-DNA antibodies – structure and function. Lupus 2002, 11, 776-779.
- RAVIRAJAN CT et al. An analysis of clinical disease activity and nephritis associated serum autoantibody profiles in patients with systemic lupus erythematosus: a cross-sectional study. Rheumatology 2001, 40, 1405-1412.
- SMEENK RJ et al. dsDNA autoantibodies. In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234.

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REF **HM 002**

DEFINITION

The **DNA-LISA** kit is an enzyme linked immunoassay for the quantitative determination of anti-dsDNA antibodies in human serum.

DIAGNOSTIC VALUE

Connective diseases are systemic auto-immune diseases. Their classification is based upon clinical and biological data. Biological criteria is based on the detection of anti-nuclear antibodies.

Among them, anti-dsDNA (ds: double strength) antibodies are recognized as the major serologic marker of Systemic Lupus Erythematosus (SLE). Their specificity and their sensitivity give them a high diagnostic value. Therefore they are a part of the clinical and biological criteria established in 1982 by the American Rheumatism Association (ARA) for the diagnosis of the SLE.

The detection of anti-dsDNA antibodies is particularly useful in two different ways : as an aid to the diagnosis of SLE and as a tool to monitor the course of the disease.

For the second purpose repeated serum sampling of individual patients can be very informative about the clinical course of the disease because a relationship exists between anti-DNA and diseases activity : flares of SLE are generally preceded by a rise in anti-dsDNA levels, followed by a steep drop during the exacerbation (particularly in nephritis).

Furthermore different treatments of patients have varying influences on anti-dsDNA levels and can be adapted by a regular follow-up of these antibodies.

SAMPLES COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic or hemolyzed sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 7 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

ASSAY PRINCIPLE

dsDNA antigen is coated on to wells of a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well for binding. After incubation, unbound antibodies are removed by washing.
- Then, a murine monoclonal anti-human IgG conjugated to horseradish peroxidase is added to each of the wells. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.

- The chromogenic step is performed by the addition of substrate ABTS (2,2 Azino-di-[3-sulfonate déthyl benzthiazoline]). During this step, a colour will develop in proportion to the amount of antibody in the sample.

- Addition of Stop Solution H₂SO₄ (0.25 M) serves to stop the enzymatic reaction.

- After stopping the reaction by H₂SO₄ (0.25M), the optical density is read by a spectrophotometer at 405nm

- Anti-dsDNA levels of sample patients are calculated using a standard curve of 4 standards.

Titers of anti-dsDNA standards are expressed in IU/ml (International Standard for dsDNA Wo/80).

REAGENTS

dsDNA antigen coated microplate with individual breakaway wells (12x8), with holder. MP	12 strips
4 standards titered in international units (WHO) <u>Ready to use</u> Each titer is printed on the vial label. CAL n	4 x 1,5mL
1 vial of positive control titered in international units (WHO) <u>To dilute</u> Expected values after dilution are printed on the vial label. CONTROL +	1 x 350µL
1 vial of negative control <u>To dilute</u> CONTROL -	1 x 350µL
1 vial of phosphate-Tween buffer pH 7.2: to be reconstituted with distilled water. <u>10X concentrated.</u> BUF WASH 10x	1 x 100mL
1 vial of anti-human IgG conjugated to horseradish peroxidase <u>Ready to use.</u> CONJ IgG	1 x 15mL
1 vial of substrat solution (ABTS) <u>Ready to use</u> SUBS ABTS	1 x 15mL
1 vial of stop solution H ₂ SO ₄ (0.25M) <u>Ready to use</u> SOLN STOP	1 x 15mL

MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- microplate spectrophotometer with 405nm filter
- 8 channel pipettes

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

SETUP

All the reagents should be prepared as required:

1. Dilution and washing buffer.

- Dilute the concentrated Phosphate-Tween at 1/10 in distilled water.
- Storage period: 3 months at +2°C/+8°C.

NB. If there are crystals in the concentrated solution, warm the bottle to +37°C for 15 minutes before use.

2. Preparation of samples and controls

- Dilute samples and controls at 1/101 in "dilution and washing buffer".
(10µL of sample + 1mL of "dilution and washing buffer")
- Vortex vigorously.

3. Use of ready-to-use conjugate.

- Estimate amount required for handling and transfer to a tube.

4. Use of ready-to-use substrate.

- Estimate amount required for handling and transfer to a tube.

PRECAUTIONS

Place all reagents outside of their packaging location. They must be necessarily returned to ambient temperature (+18°C/+25°C) at least one half an hour before beginning the test.

⚠ The temperature of the reagents could disturb the final result.

Check that all plates are well drained after each wash.

Avoid to use reagents if signs of contamination or other visible changes occur.

Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

Reagents in solution (except for substrate buffer and stop solution) contain less than 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

DNA-LISA has been developed according to E.C Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

DNA-LISA has been optimized for the use as described in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations, set out:

- 1 "blank reagent" well (made up with TDL)
- 4 "standards" wells (in duplicate)
- 1 well for each control (in duplicate)
- To improve the results reproducibility it is recommended to run the samples in duplicate.
- Remove the exact number of wells needed for the series and place them on the plastic support.

2. Samples incubation

Add 100µL of standards or calibrator, diluted controls or samples.

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with dilution and washing buffer (300µL/well).

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.

3. Conjugate incubation

Add 100µL of conjugate.

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with dilution and washing buffer (300µL/wells).

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.

4. Substrate incubation

Add 100µL substrate into each well.

Incubate for 30 minutes at room temperature.

5. Stop reaction

Add 100µL of H₂SO₄ (0,25M) to each well.

6. Reading

Read the optical density of each well at 405 nm against the reagent blank set at zero absorbance within 30 minutes.

RESULTS

Standard curve

- Determine the mean optical density for all duplicate readings.
- Trace the standard curve on log-log paper, plotting the units of the 4 standard points along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).

Comments:

- The mean OD of the standard 4 should be at least 0.8.
- The linear coefficient of regression should be higher than 0.990 for the validation of tests.
- The samples with values greater than that of standard 4 may be diluted to 1/400 to obtain a more precise result. The number of units should be multiplied by 4.

INTERPRETATION OF RESULTS

International units (WHO)	< 35 IU/mL	35-55 IU/mL	> 55 IU/mL
Results	Negative	Borderline*	Positive

*: Borderline results should be controlled on a second sample and the interpretation of the results should be done in the frame of additional testing and taking into account the clinical status of the patient.

LIMITS

Because of differences which may exist between laboratories and locations with respect to population, diet, laboratory technique and selection of reference groups, it is important for each laboratory to establish the validity of the above suggested range.

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative controls should be included with each run of the test. Positive control concentration should fall within the range printed on the Positive control label. If either reagent control is invalid, the test should be repeated.

bmd offers a multiparametric quality control (Immunotrol I, Cat. No: HM036 or Immunotrol IV, Cat. No: HM051) which contain antibodies directed against dsDNA antigens. They must be handled in the same way as samples.

CHARACTERISTICS AND PERFORMANCE OF THE TEST

1. SENSITIVITY, SPECIFIC AFFINITY AND AGREEMENT TO OTHER METHODS

The new peroxidase **DNA-LISA** kit has been compared to a marketed ELISA kit.

The study was performed on 338 samples of which:

⇒ 47 positive samples to anti-dsDNA antibodies with the reference method

⇒ 291 negative samples to anti-dsDNA antibodies with the reference method

Results

		ELISA marketed Kit		
		+	-	
bmd DNA-LISA	+	45	14	Relative sensibility : 97.8%
	-	1	278	Relative specificity : 95.2% Accuracy : 95.6%

2. PRECISION

Estimated on three positive samples : strongly, moderate and weakly positive.

Intra-assay (10 tests in the same assay)		Inter-assays (5 tests in 4 different assays)	
Titre	CV (%)	Titre	CV (%)
320	2.0	350	5.0
86	5.0	108	11.0
48	6.0	51	11.0


REFERENCES

- AHMAN A et al. Molecular expression systems for anti-DNA antibodies - 1 – 2. Lupus 2002, 11, 824-832, 833-842.
- CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance. Clin. and Exp. Rheumatol. 1997 ;15, 457-465
- HOLLINGSWORTH P.N. et al. Antinuclear antibodies In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.74-90.
- HOMBURGER H.A. Cascade testing for autoantibodies in connective tissue disease. Mayo Clin. Proc., 1995, 70, 183-184.
- ISENBERG et al. Clinical laboratory assays for measuring anti-dsDNA antibodies. Were are we now? Lupus 2002, 11, 797-800.
- RAHMAN A et al. Anti-DNA antibodies – structure and function. Lupus 2002, 11, 776-779.
- RAVIRAJAN CT et al. An analysis of clinical disease activity and nephritis associated serum autoantibody profiles in patients with systemic lupus erythematosus: a cross-sectional study. Rheumatology 2001, 40, 1405-1412.
- SMEENK RJ et al. dsDNA autoantibodies. In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234

SUMMARY OF METHOD

	<i>AMOUNT TO BE DISTRIBUTED</i>	<i>REAGENTS</i>	<i>INCUBATION CONDITIONS</i>
INCUBATION OF SAMPLES	100µL	Dilution buffer (blank reagent)	30 mn room temperature
	100µL	Standards (in duplicate)	
	100µL	Diluted controls (in duplicate)	
	100µL	Diluted samples	
WASHING	Wash 3 times with dilution and washing buffer (300µL/well)		
⚠ Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.			
INCUBATION OF CONJUGATE	100µL	Anti-IgG-PEROX Ready to use	30 mn room temperature
WASHING	Wash 3 times with dilution and washing buffer (300µL/well)		
⚠ Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.			
ENZYME REACTION	100µL	Substrate solution ABTS	30 mn room temperature.
STOP REACTION	100µL	Stop Solution H ₂ SO ₄ (0.25M)	

SYMBOLS USED

	EC Declaration of conformity		Number of test
	ELISA Test		Consult Instructions
	<i>In Vitro</i> Diagnostic Device		Manufactured by
	Catalogue number		Temperature limitation
	Lot Number		Biological risk
	Expiry Date		Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com