# PHOSPHO-LISA IgG et IgM

REF

HM 007

96

## DÉFINITION

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## BIBLIOGRAPHIE

- ALARCON-SEGOVIA D, CABRAL A.R.  
Functional and immunochemical heterogeneity of antiphospholipid antibodies : A classification. J. Rheum., 19, 1166-1169, 1992.
- ASHERSON R.A., KHAMASTA M.A., ROS J.  
The "primary" antiphospholipid syndrome : major clinical and serological features. Med., 68, 366-373, 1989.
- GALLI M., BARBUI T., CONFURIUS P. ZWALL R.F.A., BEVERS E.M. Anticoagulant mechanisms of antiphospholipids antibodies.  
Lupus Anticoagulants : Recent advances in diagnosis and understanding -ISTH Congress, Satellite Symposium- NEW-YORK 1993
- HARRIS E.N. Annotation - Antiphospholipid antibodies - Br. J. Haematol., 74, 1-9, 1990
- LAROCHE P, BERARD M, ROUQUETTE AM, DESGRUELLE, C. BOFFA MC.  
Advantage of using both anionic and switterionic phospholipid antigens for the detection of antiphospholipid antibodies. Am J Clin Pathol 106, 549-554, 1996
- LOIZOU S., MACKWORTH-YOUNG C.G., COFINER C., WALPORT M.J. Heterogeneity of binding reactivity to different phospholipids of antibodies from patients with systemic lupus erythematosus. Clin. Exp. Immunol., 80, 171-176, 1990
- Mc INTYRE J.A., WAGENKNECHT D.R., SUGI T.  
Role of lipoprotein cofactors in the detection of antiphospholipid antibodies. Lupus Anticoagulants : Recent Advances in diagnosis and understanding -ISTH Congress, Satellite Symposium- NEW-YORK 1993
- MEYER O., ROUQUETTE A-M et YOUINOUE P.  
Autoanticorps, Marqueurs des Maladies autoimmunes. BMD Editions (1999).
- TRIPLETT D.A., BRANDT J.T., MUSGRAVE K.A., ORR C.A.  
The relationship between lupus anticoagulants and antibodies to phospholipids. J. Am. Med. Ass., 259, 550-554, 1988.

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96

## INTENDED USE

**PHOSPHO-LISA IgG / IgM (BMD)** is an enzyme linked immunoassay (ELISA) for the quantitative determination of IgG and IgM anti-phospholipid antibodies in human serum or plasma samples.

A mixture of four phospholipids (three anionic and one zwitterionic) is used to coat the microplate in order to reveal the anti-phospholipid antibodies involved in different clinical conditions

## DIAGNOSTIC VALUE

Anti-phospholipids (APL) are IgG, IgM or IgA immunoglobulins, whose antigen targets are anionic phospholipids (PL) (cardiolipin, phosphatidylserine, phosphatidylinositol) and also zwitterionic phospholipids such as phosphatidylethanolamine and, to a lesser extent, phosphatidylcholine.

APL are inhibitors, interfering in the phospholipid-dependent stages of coagulation.

APL detection using the ELISA method calls for the presence of serum co-factors. The interaction between co-factors and PL is revealed by the formation of protein-PL complexes which induce the formation of neoantigens specifically recognised by the APL.

The co-factor for anionic PL has been identified as  $\beta_2$ -glycoprotein I or apolipoprotein H. Preliminary studies suggest that there is a different co-factor for phosphatidylethanolamine which is present in human serum.

The term "Anti-phospholipid syndrome" describes the association of one or several clinical criteria, venous and/or arterial thrombosis, repeated early abortion, thrombopenia with the presence of an APL with or without a circulating anticoagulant.

Apart from this syndrome, APLs have been implicated in a variety of clinical conditions:

- autoimmune disease
- cancer
- thrombopenia
- infection (syphilis, AIDS)
- consumption of some drug

## SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum or plasma.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

## ASSAY PRINCIPLE

A mixture of four phospholipids is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the antigen coated well in duplicate to research IgG and IgM antibodies. After incubation, unbound proteins are removed by washing.
- Binding of anti-phospholipid antibodies is detected using anti-human IgG or IgM alkaline phosphatase labelled conjugates. After incubation, the wells are washed again to eliminate any excess of conjugate.
- The bound enzyme is revealed by addition of substrate (Para Nitro Phenyl Phosphate), producing a yellow colour. The colour intensity is proportional to anti-phospholipid antibodies level in the sample.
- After stopping the reaction by 1N sodium hydroxide, the optical density is read by a spectrophotometer at 405nm

Anti-phospholipid levels of samples patients are expressed in GPL or MPL units, using a standard curve.

## STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- timer
- incubator adjusted at +37°C
- physiological serum (NaCl 0.9%)
- microplate spectrophotometer with 405 nm filter
- 8 channel pipettes

## REAGENTS

Microplate coated with Phospholipids (12x8 individual breakaway wells with holder) <span style="float: right;">MP</span>	12 strips
4 Standards in GPL international units <u>Ready to use</u> <span style="float: right;">CAL n</span> <i>The values are shown on the vial label.</i>	4 x 1mL
4 Standards in MPL international units <u>Ready to use</u> <span style="float: right;">CAL n</span> <i>The values are shown on the vial label.</i>	4 x 1mL
1 vial of Positive control in international GPL units <u>To dilute</u> <span style="float: right;">CONTROL +</span> <i>The values expected are shown on the vial label.</i>	1 x 350µL
1 vial of Positive control in international MPL units <u>To dilute</u> <span style="float: right;">CONTROL +</span> <i>The values expected are shown on the vial label.</i>	1 x 350µL
1 vial of Negative control <u>To dilute</u> <span style="float: right;">CONTROL -</span>	1 x 350µL
1 vial of Goat anti-human IgG conjugate PAL <u>Ready to use</u> <span style="float: right;">CONJ IgG</span>	1 x 7.5mL
1 vial of Goat anti-human IgM conjugate PAL <u>Ready to use</u> <span style="float: right;">CONJ IgM</span>	1 x 7.5mL
1 vial of Substrate containing 6x5mg PNPP tablets <u>To dilute</u> <span style="float: right;">SUBS PNPP</span>	1
1 vial of Phosphate buffer pH 7,2 (concentrated 10x) <span style="float: right;">BUF WASH 10x</span> <u>To dilute in distilled water</u>	1 x 100mL
1 vial of Substrate dilution buffer <u>Ready to use</u> <span style="float: right;">BUF SUBS</span>	1 x 35mL
1 vial of <b>Solution A</b> : dilution and washing buffer <u>To dilute</u> <span style="float: right;">BUF Sol.A</span>	1 x 15mL
1 vial of Stop solution (Sodium hydroxide 1N) <u>Ready to use</u> <span style="float: right;">SOLN STOP</span>	1 x 9mL

## SETUP

All the reagents should be prepared as required:

### 1. Preparation of Phosphate buffer solution

- Dilute concentrated Phosphate buffer 1/10 in distilled water.
- Storage period: 3 month at +2°C/+8°C.  
*NB. If there are crystals in the concentrated solution, warm the bottle to 37°C for 15 minutes before use.*

### 2. Preparation of washing buffer

- Dilute solution A 1/100 in Phosphate buffer solution.

### 3. Preparation of dilution buffer (solution D)

- Dilute solution A 1/10 in Phosphate buffer.

### 4. Preparation of samples and controls.

- Dilute to 1/50 in dilution buffer (solution D) :  
Ex: 10µL sample qsp 0.5mL solution D and shake vigorously at the vortex.

### 5. Preparation of ready-to-use conjugate.

- Estimate amount required for handling and transfer to a tube.

### 6. Preparation of the substrate.

- Dissolve the PNPP tablets in the substrate dilution buffer at a ratio of 1 tablet per 5 .
- Protect from light until use.

## PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C-+25°C) before handling. Check that all plates are well drained after each wash. Do not put plates on top of each other during incubation at 37°C.

Check that no condensed water is left at the bottom of the wells before reading.

*Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.*

Reagents in solution (except for substrate buffer and stop solution) contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**PHOSPHO-LISA IgG and IgM** has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

**PHOSPHO-LISA IgG and IgM** has been optimized for the use as describe in this procedure. Do not substitute reagents from other manufacturers. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

## METHOD

### 1. Preparing the test

Use the work sheet to note the sample locations.

- For each isotype, set out:
  - One blank reagent well (made up with dilution buffer)
  - 4 standard wells for the detection of IgG antibodies (in duplicate)
  - 4 standard wells for the detection of IgM antibodies (in duplicate)
  - 1 well for IgG and IgM positive control (in duplicate)
  - 1 well for the negative control (in duplicate)
  - to obtain more precise results, it is also advisable to run samples in duplicate.
- then detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

### 2. Sample incubation

Add 100µL of standards, diluted controls or samples.

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 200µL of washing buffer

### 3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature

Wash step:

Repeat the procedure described previously

Make a last wash with 200µl of 0,9% sodium chloride.

### 4. Incubation of substrate

Add 100µL substrate into each well

Incubate for 30 minutes in the dark at +37°C.

### 5. Stop solution

Add 50µL of 1N sodium hydroxide to each well.

### 6. Reading

Read the optical density of each well at 405 nm against the reagent blank set at zero absorbance within 30 minutes.

## RESULTS

### Standard curve

Determine the mean optical density for all duplicate readings.

Trace the standard curve on log-log paper, plotting the units of the 4 standard points along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).

### Comments

- The mean OD of the standard 4 should be at least 0.8.
- The samples with values greater than that of standard 4 may be diluted to 1/200 to obtain a more precise result. The number of units should be multiplied by 4.

## INTERPRETATION OF RESULTS

GPL or MPL units	< 10	10-15	> 15
Results	Negative	Weakly positive*	Positive

\*Weakly positive results should be controlled on a second specimen obtained a few weeks later and interpreted in the light of the results of additional investigations and the clinical context.

### QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive control concentration should fall within the range printed on the positive control label. If either reagent control is invalid, the test should be repeated.

bmd offers external multiparameter quality controls. Immunotrol I, cat. No: **HM036** which contains antibodies **IgG directed against phospholipids** and Immunotrol II, cat. No: **HM037**, which contains antibodies **IgM directed against phospholipids**. These controls should be handled as samples.

### CHARACTERISTICS AND PERFORMANCE OF THE TEST

#### 1. MEASURE RANGE (FOR 1:50 SAMPLES DILUTION)

Linear range : 0 to 100 units

Security range: > 1000 units

#### 2. THRESHOLD VALUES

Selected population : 100 sera from healthy individuals

- the positive threshold corresponds to the mean of the OD + 3 standard deviations, obtained for the population tested (i.e, 97th percentile OD distribution)

-between 3 and 5 standard deviations, samples are considered as weakly positive.

#### 3. ACCURACY

Estimated by recovery test.

Principle: negative sample is supercharged with increasing quantities of a positive sample.

The overlap percentage equals (titre obtained/theoretical titre) x 100.

	Mean percentage
IgG Isotype	96%
IgM Isotype	99%

Linear recovery : from 5 to 100 units

#### 4. REPRODUCIBILITY

Estimated on the basis of 3 samples : positive, weakly positive and negative.

	INTRA-ASSAY			INTER-ASSAY	
	Titre	n	CV (%)	n	CV (%)
IgG Isotype	11 GPL	5	4.7	5	6.0
	26.5 GPL	5	3.6	5	6.0
	43 GPL	5	5.2	5	7.4

## REFERENCES

- ALARCON-SEGOVIA D, CABRAL A.R.  
Functional and immunochemical heterogeneity of antiphospholipid antibodies : A classification. J. Rheum., 19, 1166-1169, 1992.
- ASHERSON R.A., KHAMASTA M.A., ROS J.  
The "primary" antiphospholipid syndrome : major clinical and serological features. Med., 68, 366-373, 1989.
- GALLI M., BARBUI T., CONFURIUS P. ZWALL R.F.A., BEVERS E.M. Anticoagulant mechanisms of antiphospholipids antibodies. Lupus Anticoagulants : Recent advances in diagnosis and understanding -ISTH Congress, Satellite Symposium- NEW-YORK 1993
- HARRIS E.N.  
Annotation - Antiphospholipid antibodies. Br. J. Haematol., 74, 1-9, 1990
- LAROCHE P, BERARD M, ROUQUETTE AM, DESGRUELLE, C. BOFFA MC. Advantage of using both anionic and switterionic phospholipid antigens for the detection of antiphospholipid antibodies. Am J Clin Pathol 106, 549-554, 1996
- LOIZOU S., MACKWORTH-YOUNG C.G., COFINER C., WALPORT M.J. Heterogeneity of binding reactivity to different phospholipids of antibodies from patients with systemic lupus erythematosus. Clin. Exp. Immunol., 80, 171-176, 1990
- Mc INTYRE J.A., WAGENKNECHT D.R., SUGI T.  
Role of lipoprotein cofactors in the detection of antiphospholipid antibodies. Lupus Anticoagulants : Recent Advances in diagnosis and understanding -ISTH Congress, Satellite Symposium- NEW-YORK 1993
- MEYER O., ROUQUETTE A-M et YOUINOU P.  
Autoanticorps, Marqueurs des Maladies autoimmunes. BMD Editions (1999).
- TRIPLETT D.A., BRANDT J.T., MUSGRAVE K.A., ORR C.A.  
The relationship between lupus anticoagulants and antibodies to phospholipids. J. Am. Med. Ass., 259, 550-554, 1988.

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#### 5. SENSITIVITY, SPECIFIC AFFINITY AND AGREEMENT WITH ANOTHER METHOD

Sera from 141 patients have been selected for potential AP related disorders :

Autoimmune diseases (72 including 56 Systemic Lupus Erythematosus ), arterial or venous thrombosis (18), recurrent foetal losses (11), vasculitides (23), neurologic syndrome (4), neoplasia (3), HIV (3), miscellaneous (7).

#### Measurement of APL using individual phospholipids as antigens

Anti-Cardiolipine (CL) antibodies were measured using commercial ELISA kits (CARDIO-LISA IgG and IgM from bmd). Antibodies directed against anti-phosphatidylinositol (PI), anti-phosphatidylserine (PS) et anti-phosphatidylethanolamine (PE) were measured following the same procedure.

79/141 (56%) were positive for at least one of the four antigenic specificities :

#### ► Isotypes distribution pattern

The most common APL isotype was IgG : alone (67%) or associated with IgM (28%). Only 5% of APL were exclusively of the IgM isotype

#### Antigenic specificities distribution pattern (in percentage)

	CL	PI	PS	PE
IgG	92	96	89	25
IgM	77	73	81	38

#### ► Distribution pattern regarding the antigenic target (multiple or only one target)

Most of positive samples were directed against the three anionic phospholipids :

	IgG (n)	IgM (n)
Multiple targets : APL directed against anionic phospholipids	49	15
Multiple targets : APL directed against anionic phospholipids and PEP	18	8
Only one target :	8 with:	3 with:
-anti-CL	2	2
-anti-PI	5	0
-anti-PS	0	0
-anti-PE	1	1

#### Simultaneous measurement of the four specificities using the PHOSPHO-LISA kit

All 141 samples tested previously were evaluated using PHOSPHOLISA. Sensitivity and specificity of the mixed phospholipids assay were compared to those obtained with each individual phospholipids test.

Sensitivity : 100% for both IgG and IgM isotypes : all the samples positive for at least one of the four specificities were also positive with PHOSPHO-LISA.

Specificity : 100% for IgM antibodies et 98,5% for the IgG isotype : Among the 66 samples negative for each of the four phospholipids, 65 were negative with PHOSPHO-LISA. The only one sample, negative with each individual assay , was weakly positive with PHOSPHO-LISA (13 units).

## SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	IgG IsoTYPE	IgG IsoTYPE	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	100µL	Solution D: blank reagent	Solution D: blank reagent	30 min at room temperature
	100µL	Standards (double)	Standards (double)	
	100µL	Diluted controls and samples	Diluted controls and samples	
WASHING	Wash 3 times with washing buffer (200µL/ well)			
INCUBATION OF CONJUGATE	100µL	Anti-IgG PAL conjugate	Anti-IgM PAL conjugate	30 min at room temperature
WASHING	Wash 3 times with washing buffer (200µL/ well) Wash once with NaCl 0.9% (200µL/ well)			
ENZYME REACTION	100µL	Substrate solution		30 min at +37°C
STOP REACTION	50µL	NaOH 1N		

## SYMBOLS USED

	EC Declaration of conformity		Number of test
	ELISA Test		Consult Instructions
	<i>In Vitro</i> Diagnostic Device		Manufactured by
	Catalogue number		Temperature limitation
	Lot Number		Biological risk
	Expiry Date		Reconstitute with

***BioMédical Diagnostics SA***

**Office**  
Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)