PHOSPHO-LISA IgG

REF HM 007G

DÉFINITION

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BIBLIOGRAPHIE

ALARCON-SEGOVIA D, CABRAL A.R.
Functional and immunochemical heterogeneity of antiphospholipid antibodies : A classification.
J. Rheum., 19, 1166-1169, 1992.

ASHERSON R.A., KHAMASTA M.A., ROS J.
The "primary" antiphospholipid syndrome : major clinical and serological features.
Med., 68, 366-373, 1989.

GALLI M., BARBUI T., CONFURIUS P. ZWALL R.F.A., BEVERS E.M.
Anticoagulant mechanisms of antiphospholipid antibodies.
Lupus Anticoagulants : Recent advances in diagnosis and understanding -ISTH Congress, Satellite Symposium-NEW-YORK 1993

HARRIS E.N.
Annotation - Antiphospholipid antibodies
Br. J. Haematol., 74, 1-9, 1990

LAROCHE P, BERARD M, ROUQUETTE AM, DESGRUELLE, C. BOFFA MC.
Advantage of using both anionic and switterionic phospholipid antigens for the detection of antiphospholipid antibodies
Am J Clin Pathol 106, 549-554, 1996

LOIZOU S., MACKWORTH-YOUNG C.G., COFINER C., WALPORT M.J.
Heterogeneity of binding reactivity to different phospholipids of antibodies from patients with systemic lupus erythematosus.
Clin. Exp. Immunol., 80, 171-176, 1990

Mc INTYRE J.A., WAGENKNECHT D.R., SUGI T.
Role of lipoprotein cofactors in the detection of antiphospholipid antibodies.
Lupus Anticoagulants : Recent Advances in diagnosis and understanding -ISTH Congress, Satellite Symposium-NEW-YORK 1993

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REF HM 007G

INTENDED USE

PHOSPHO-LISA IgG (bmd) is an enzyme linked immunoassay (ELISA) for the quantitative determination of IgG anti-phospholipid antibodies in human serum or plasma samples.

A mixture of four phospholipids (three anionic and one zwitterionic) is used to coat the microplate in order to reveal the anti-phospholipid antibodies involved in different clinical conditions.

DIAGNOSTIC VALUE

Anti-phospholipids (APL) are IgG, IgM or IgA immunoglobulins, whose antigen targets are anionic phospholipids (PL) (cardiolipin, phosphatidylserine, phosphatidylinositol) and also zwitterionic phospholipids such as phosphatidylethanolamine and, to a lesser extent, phosphatidylcholine.

APL are inhibitors, interfering in the phospholipid-dependent stages of coagulation.

APL detection using the ELISA method calls for the presence of serum co-factors. The interaction between co-factors and PL is revealed by the formation of protein-PL complexes which induce the formation of neoantigens specifically recognised by the APL.

The co-factor for anionic PL has been identified as β_2 -glycoprotein I or apolipoprotein H. Preliminary studies suggest that there is a different co-factor for phosphatidylethanolamine which is present in human serum.

The term "Anti-phospholipid syndrome" describes the association of one or several clinical criteria, venous and/or arterial thrombosis, repeated early abortion, thrombopenia with the presence of an APL with or without a circulating anticoagulant.

Apart from this syndrome, APLs have been implicated in a variety of clinical conditions:

- autoimmune disease
- cancer
- thrombopenia
- infection (syphilis, AIDS)
- consumption of some drug

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum or plasma.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

ASSAY PRINCIPLE

A mixture of four phospholipids is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the antigen coated well in duplicate to detect IgG antibodies. After incubation, unbound proteins are removed by washing.
- Binding of anti-phospholipid antibodies is detected using anti-human IgG alkaline phosphatase labelled conjugates. After incubation, the wells are washed again to eliminate any excess of conjugate.
- The bound enzyme is revealed by addition of substrate (Para Nitro Phenyl Phosphate), producing a yellow colour. The colour intensity is proportional to anti-phospholipid antibodies level in the sample.
- After stopping the reaction by 1N sodium hydroxide, the optical density is read by a spectrophotometer at 405nm

Anti-phospholipid levels of samples patients are expressed in GPL, using a standard curve.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- timer
- incubator adjusted at +37°C
- physiological serum (NaCl 0.9%)
- microplate spectrophotometer with 405 nm filter
- 8 channel pipettes

REAGENTS

Microplate coated with Phospholipids (12x8 individual breakaway wells with holder) MP	12 strips
4 Standards in GPL international units Ready to use CAL n <i>The values are shown on the vial label.</i>	4 x 1.5mL
1 vial of Positive control in international GPL units To dilute CONTROL + <i>The values expected are shown on the vial label.</i>	1 x 350µL
1 vial of Negative control To dilute CONTROL -	1 x 350µL
1 vial of Goat anti-human IgG conjugate PAL Ready to use CONJ IgG	1 x 15mL
1 vial of Substrate containing 6x5mg PNPP tablets To dilute SUBS PNPP	1
1 vial of Phosphate buffer pH 7,2 (concentrated 10x) BUF WASH 10x To dilute in distilled water	1 x 100mL
1 vial of Substrate dilution buffer Ready to use BUF SUBS	1 x 35mL
1 vial of Solution A : dilution and washing buffer To dilute BUF Sol.A	1 x 15mL
1 vial of Stop solution (Sodium hydroxide 1N) Ready to use SOLN STOP	1 x 9mL

SETUP

All the reagents should be prepared as required:

1. Preparation of Phosphate buffer solution

- Dilute concentrated Phosphate buffer 1/10 indistilled water.
- Storage period: 3 month at +2°C/+8°C.

NB. If there are crystals in the concentrated solution, warm the bottle to 37°C for 15 minutes before use.

2. Preparation of washing buffer

- Dilute solution A 1/100 in Phosphate buffer solution.

3. Preparation of dilution buffer (solution D)

- Dilute solution A 1/10 in Phosphate buffer.

4. Preparation of samples and controls.

- Dilute to 1/50 in dilution buffer (solution D) :
Ex: 10µL sample qsp 0,5mL solution D and shake vigorously at the vortex.

5. Preparation of ready-to-use conjugate.

- Estimate amount required for handling and transfer to a tube.

6. Preparation of the substrate.

- Dissolve the PNPP tablets in the substrate dilution buffer at a ratio of 1 tablet per 5.
- Protect from light until use.

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C--+25°C) before handling.

Check that all plates are well drained after each wash.

Do not put plates on top of each other during incubation at 37°C.

Check that no condensed water is left at the bottom of the wells before reading.

Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

Reagents in solution (except for substrate buffer and stop solution) contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

PHOSPHO-LISA IgG has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

PHOSPHO-LISA IgG has been optimized for the use as describe in this procedure. Do not substitute reagents from other manufacturers. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

- Set out :
 - one blank reagent well (made up with solution D: dilution buffer)
 - 4 standard wells for the detection of IgG antibodies (in duplicate)
 - 1 well for each positive control (in duplicate)
 - 1 well for the negative control(in duplicate)
 - to obtain more precise results, it is also advisable to run samples in duplicate.
- then detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

2. Sample incubation

Add 100µL of standards, diluted controls or samples.

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 200µL of washing buffer

3. Incubation of conjugate

Add 100µL of conjugate
 Incubate for 30 minutes at room temperature
 Wash step:
 Repeat the procedure described previously
 Make a last wash with 200µL of 0,9% sodium chloride.

4. Incubation of substrate

Add 100µL substrate into each well
 Incubate for 30 minutes in the dark at 37°C.

5. Stop solution

Add 50µL of 1N sodium hydroxide to each well.

6. Reading

Read the optical density of each well at 405 nm against the reagent blank set at zero absorbance within 30 minutes.

2. THRESHOLD VALUES

Selected population : 100 sera from healthy individuals
 - the positive threshold corresponds to the mean of the OD + 3 standard deviations, obtained for the population tested (i.e, 97th percentile OD distribution)
 - between 3 and 5 standard deviations, samples are considered as weakly positive.

3. ACCURACY

Estimated by recovery test.
 Principle: negative sample is supercharged with increasing quantities of a positive sample.
 The overlap percentage equals (titre obtained/ theoretical titre) x 100.

	Mean percentage
IgG Isotype	96%

Linear recovery : from 5 to 100 units

RESULTS

Standard curve

Determine the mean optical density for all duplicate readings.
 Trace the standard curve on log-log paper, plotting the units of the 4 standard points along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).

Comments

- The mean OD of the standard 4 should be at least 0.8.
- The samples with values greater than that of standard 4 may be diluted to 1/200 to obtain a more precise result. The number of units should be multiplied by 4.

4. REPRODUCIBILITY

Estimated on the basis of 3 samples : positive, weakly positive and negative.

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INTERPRETATION OF RESULTS

GPL or MPL units	< 10	10-15	> 15
Results	Negative	Weakly positive*	Positive

*Weakly positive results should be controlled on a second specimen obtained a few weeks later and interpreted in the light of the results of additional investigations and the clinical context.

QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive control concentration should fall within the range printed on the positive control label. If either reagent control is invalid, the test should be repeated.

bmd offers an external multiparameter quality control (Immunotrol I, Cat: HM036). It contains antibodies IgG directed against phospholipid and cardiolipin antigens. This control should be handled as samples.

CHARACTERISTICS AND PERFORMANCE OF THE TEST

1. MEASURE RANGE (FOR 1:50 SAMPLES DILUTION)

Linear range : 0 to 100 units
 Security range: > 1000 units

5. SENSITIVITY, SPECIFIC AFFINITY AND AGREEMENT WITH ANOTHER METHOD

Sera from 141 patients have been selected for potential AP related disorders :
 Autoimmune diseases (72 including 56 Sytemic Lupus Erythematosus), arterial or venous thrombosis (18), recurrent fetal losses (11), vasculitides (23), neurologic syndrome (4), neoplasia (3), HIV (3), miscellaneous (7).

Measurement of APL using individual phospholipids as antigens

Anti-Cardiolipine (CL) antibodies were measured using commercial ELISA kits (CARDIO-LISA IgG and IgM from BMD). Antibodies directed against anti-phosphatidylinositol (PI), anti-phosphatidylserine (PS) et anti-phosphatidylethanolamine (PE) were measured following the same procedure.

79/141 (56%) were positive for at least one of the four antigenic specificities :

➤ Isotypes distribution pattern

The most common APL isotype was IgG : alone (67%) or associated with IgM (28%).

Antigenic specificities distribution pattern (in percentage)

	CL	PI	PS	PE
IgG	92	96	89	25

➤ **Distribution pattern regarding the antigenic target (multiple or only one target)**

Most of positive samples were directed against the three anionic phospholipids :

	IgG (n)
Multiple targets : APL directed against anionic phospholipids	49
Multiple targets : APL directed against anionic phospholipids and PEP	18
Only one target :	8 with:
-anti-CL	2
-anti-PI	5
-anti-PS	0
-anti-PE	1

Simultaneous measurement of the four specificities using the PHOSPHO-LISA kit

All 141 samples tested previously were evaluated using PHOSPHOLISA. Sensitivity and specificity of the mixed phospholipids assay were compared to those obtained with each individual phospholipid test.

Sensitivity : 100% for both IgG and IgM isotypes : all the samples positive for at least one of the four specificities were also positive with PHOSPHO-LISA.

Specificity : 100% for IgM antibodies et 98,5% for the IgG isotype : Among the 66 samples negative for each of the four phospholipids, 65 were negative with PHOSPHO-LISA. The only one sample, negative with each individual assay , was weakly positive with PHOSPHO-LISA (13 units).

REFERENCES

ALARCON-SEGOVIA D, CABRAL A.R.
Functional and immunochemical heterogeneity of antiphospholipid antibodies : A classification.
J. Rheum., 19, 1166-1169, 1992.

ASHERSON R.A., KHAMASTA M.A., ROS J.
The "primary" antiphospholipid syndrome : major clinical and serological features.
Med., 68, 366-373, 1989.

GALLI M., BARBUI T., CONFURIUS P. ZWALL R.F.A., BEVERS E.M.
Anticoagulant mechanisms of antiphospholipids antibodies.
Lupus Anticoagulants : Recent advances in diagnosis and understanding -ISTH Congress, Satellite Symposium- NEW-YORK 1993

HARRIS E.N.
Annotation - Antiphospholipid antibodies
Br. J. Haematol., 74, 1-9, 1990

LAROCHE P, BERARD M, ROUQUETTE AM, DESGRUELLE, C. BOFFA MC.
Advantage of using both anionic and switterionic phospholipid antigens for the detection of antiphospholipid antibodies
Am J Clin Pathol 106, 549-554, 1996

LOIZOU S., MACKWORTH-YOUNG C.G., COFINER C., WALPORT M.J.
Heterogeneity of binding reactivity to different phospholipids of antibodies from patients with systemic lupus erythematosus.
Clin. Exp. Immunol., 80, 171-176, 1990

Mc INTYRE J.A., WAGENKNECHT D.R., SUGI T.
Role of lipoprotein cofactors in the detection of antiphospholipid antibodies.
Lupus Anticoagulants : Recent Advances in diagnosis and understanding -ISTH Congress, Satellite Symposium- NEW-YORK 1993

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TRIPLETT D.A., BRANDT J.T., MUSGRAVE K.A., ORR C.A.
The relationship between lupus anticoagulants and antibodies to phospholipids.
J. Am. Med. Ass., 259, 550-554, 1988.

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	IgG ISOTYPE	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	100µL	Solution D: blank reagent	30 min at room temperature
	100µL	Standards (double)	
	100µL	Diluted controls and samples	
WASHING	Wash 3 times with washing buffer (200µL/ well)		
INCUBATION OF CONJUGATE	100µL	Anti-IgG PAL conjugate	30 min at room temperature
WASHING	Wash 3 times with washing buffer (200µL/ well) Wash once with NaCl 0.9% (200µL/ well)		
ENZYME REACTION	100µL	Substrate solution	30 min at +37°C
STOP REACTION	50µL	NaOH 1N	

SYMBOLS USED

	EC Declaration of conformity		Number of test
	ELISA Test		Consult Instructions
	<i>In Vitro</i> Diagnostic Device		Manufactured by
	Catalogue number		Temperature limitation
	Lot Number		Biological risk
	Expiry Date		Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com