ENA-LISA

REF **HM 008**

12 x 8 tests

DÉFINITION

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HANS L.P. et al. A comparison of ELISA aSS-Ays as routine diagnostic test for detection of autoantibodies against extractable nuclear antigens. Clin. Biochem., 32, 3, 179-183, 1999.

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SONTHEIMER et al. Antinuclear Antibodies: clinical correlations and biological significance. Adv. Dermatol. 1992; 7, 3-52.

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12 x 8 tests

DEFINITION

E.N.A.-LISA is an enzyme linked immunoassay (ELISA) for specific detection of human antibodies directed against: SS-A (Ro), SS-B (La), Sm, Sm/RNP, Scl-70 and Jo-1.

DIAGNOSTIC VALUE

The detection and characterisation of anti-ENA (Extractable Nuclear Antigens) antibodies are an important element for the clinician in the diagnosis and classification of sub-groups of connective tissue diseases.

Antibodies associated with different autoimmune pathologies

- **Anti-SS-A and anti-SS-B antibodies**
Either one or both of these are observed, in the same pathological conditions Sjögren's syndrome and Systemic Lupus Erythematosus (SLE) (6). Anti-SSA are also found in mothers who are carrying a baby with neonatal lupus syndrome including heart block.
- **Anti-RNP antibody**
Found primarily in Mixed Connective Tissue Disease or Sharp's Syndrome, for which they are the biological marker. They are also observed in SLE.

Specific antibodies for a single connective tissue disease

- **Anti-Sm antibodies** characterize severe forms of SLE.
- **Anti-Scl70 antibodies** are specific for diffuse Scleroderma.
- **Anti-Jo1 antibodies** are found exclusively in sera from patients suffering from Polymyositis.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

ASSAY PRINCIPLE

Each different extractable nuclear antigen is coated on to specific coloured wells of a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well for each specific antigen, which allows to bind. After incubation, unbound proteins are removed by washing.

- Then, a murine monoclonal anti-human IgG conjugated to horseradish peroxidase is added to each of the wells. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of substrate ABTS (2,2 Azino-di-[3-sulfonate déthyl benzthiazoline]). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution H₂SO₄ (0.25 M) serves to stop the enzymatic reaction.
- After stopping the reaction by H₂SO₄ (0.25M), the optical density is read by a spectrophotometer at 405nm
- A calibration device allows to expressed in arbitrary units the titre of anti-ENA autoantibodies in samples patients.

REAGENTS

Microplate coated with Extractable Nuclear Antigens (12x8 individual breakaway wells). (*)	MP	12 strips
1 vial of Calibrator tittered in arbitrary units for 6 suspected specificities to be measured.	CAL n	1 x 1,5mL
Ready to use Each specificity value is shown on the vial label.		
1 vial of Anti-human IgG conjugated to peroxidase	CONJ IgG	1 x 15mL
Ready to use.		
1 vial of Phosphate-Tween buffer pH 7,2 (concentrated 10x) To dilute in distilled water	BUF WASH 10x	1 x 100mL
1 vial of Substrat solution (ABTS)	SUBS ABTS	1 x 15mL
Ready to use		
1 vial of Stop solution H ₂ SO ₄ (0.25M)	SOLN STOP	1 x 15mL
Ready to use		

MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- microplate spectrophotometer with 405 nm filter
- 8 channel pipettes

(*) Identification of the wells is as following:

Wells	Colour	Specificity	Sample identification
A	clear	Blank reagent	Dilution buffer
B	red	calibrator	Calibrator
C	blue	SS-A	Diluted sample
D	yellow	SS-B	Diluted sample
E	green	Sm	Diluted sample
F	orange	Sm/RNP*	Diluted sample
G	white	Scl-70	Diluted sample
H	black	Jo-1	Diluted sample

(*) To test in parallel with the antigen Sm for the detection of the anti-RNP antibodies.

Determination of the 6 specificities, for a same patient, is carried out using a complete strip (A1 well to H1 well for example).

Determination of several specificities, for one or more patients: extract the number of wells corresponding to each one and two wells, A and B. Return unused strips into their plastic bag.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

SETUP

All the reagents should be prepared as required:

1. Dilution and washing buffer.

- Dilute concentrated Phosphate-Tween Buffer 1/10 in distilled water.
- Storage period: 3 months at +2°C/+8°C (avoid to use if signs of contamination or other visible changes occur)

NB. If there are crystals in the concentrated solution, warm the bottle to +37°C for 15 minutes before use.

2. Preparation of samples

- Dilute the samples at 1/101 in dilution buffer (10µL sample + 1mL dilution buffer).
- Vortex vigorously.

3. Use of ready-to-use conjugate.

- Estimate amount required for handling and transfer to a tube.

4. Use of ready-to-use substrate.

- Estimate amount required for handling and transfer to a tube.

PRECAUTIONS

Place all reagents outside of their packaging location. They must be necessarily returned to ambient temperature (+18°C/+25°C) at least one half an hour before beginning the test.

⚠ The temperature of the reagents could disturb the final result.

Check that all plates are well drained after each wash.

Avoid to use reagents if signs of contamination or other visible changes occur.

Human sources for the preparation of standards have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

Reagents in solution (except for substrate buffer and stop solution) contain less than 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures with copper or lead piping. Rinse thoroughly after flushing.

ENA-LISA has been developed according to CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

ENA-LISA has been optimized for the use as described in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations, set out:

- one blank reagent well (made up with dilution buffer)
- 1 well for calibrator
- 1 well for each specificity and each patient

2. Calibrator and Sample incubation

Add 100µL of dilution buffer in the well A

Add 100µL of diluted calibrator in the well B

Add 100µL of diluted sample in the wells corresponding to required specificities: wells C to H

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing and dilution buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.

3. Incubation of conjugate

Add 100µL of conjugate into each well

Incubate for 30 minutes at room temperature.

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing and dilution buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.

4. Incubation of substrate

Add 100µL substrate into each well.
Incubate for 30 minutes at room temperature.

5. Stop reaction

Add 100µL of H₂SO₄ (0,25M) to each well.

6. Reading

Read the optical density of each well at 405 nm against the reagent blank set at zero absorbance within 30 minutes.

RESULTS

- OD for the calibrator should be at least 0.8.
- Calculate the conversion factor for each specificity: **Calibrator value (U-SS-A, U-SS-B, U-Sm...)/ calibrator OD**
- Multiply OD of each well sample by the corresponding factor and deduce the required anti-ENA autoantibody titre.

Sample calculation:

For a calibrator titrated to 60 U-SS-A, showing an OD of 1.25, the SS-A conversion factor equals 60/1.25, which is 48.

A sample with OD 0.5 (well SS-A) has a titre of 0.5 x 48, which is 24 U-SS-A.

INTERPRETATION OF RESULTS

Arbitrary Units	< 30	30 - 40	> 40
Results	Negative	Borderline *	Positive **

* Borderline results should be controled on a second sample taken several weeks later and interpreted in the light of complementary examination and the clinical context.

** Linearity of the anti-ENA antibody-antigen response is optimal for titres which do not exceed the value of the calibrator. Above this, accuracy and reproducibility of results can be affected particularly by reaction conditions (such as temperature). Therefore results should be expressed as "over 100 units".

LIMITS

Because of differences which may exist between laboratories and locations with respect to population, diet, laboratory technique and selection of reference groups, it is importante for each laboratory to establish the validity of the above suggested range.

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

QUALITY CONTROL

If 4 or more specificities are found with the same patient (example: SS-A – SS-B - Sm/RNP - Scl-70), it is recommended to control this sample with an other method.

bmd offers external multiparameter quality controls ; Immuno-Trol IV control (Cat. No: HM051) which contains antibodies directed against SS-A, SS-B, Sm and Sm/RNP and the Immuno-Trol V control (Cat. No: HM052) which contains antibodies directed against Scl-70 and Jo1. These controls should be handled as samples.

CHARACTERISTICS AND PERFORMANCE OF THE TEST

1. SENSITIVITY, SPECIFIC AFFINITY AND AGREEMENT WITH ANOTHER METHOD

The ENA-LISA test was compared with six reference methods.

The study involved 338 samples (for SS-A, SS-B, Sm/RNP, Scl-70 and Jo-1) and 368 samples (for Sm)

⇒ 212 positive samples (for SS-A, SS-B, Sm/RNP, Scl-70 and Jo-1) and 240 positive samples (for Sm) with at least one anti-ENA antibody

⇒ 126 negative samples (for SS-A, SS-B, Sm/RNP, Scl-70 and Jo-1) and 128 negative samples (for Sm) for which any anti-ENA antibody was detected by the reference method.

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HANS L.P. et al. A comparison of ELISA aSS-Ays as routine diagnostic test for detection of autoantibodies against extractable nuclear antigens.

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HOLLINGSWORTH PN et al. Antinuclear antibodies.

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SONTHEIMER et al. Antinuclear Antibodies: clinical correlations and biological significance.

Adv. Dermatol. 1992; 7, 3-52.

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	100µL	Dilution buffer: well A (blank reagent)	30 mn at room temperature
	100µL	Calibrator: well B	
	100µL	Diluted positive samples: wells C to H	
WASHING	Wash 3 times with dilution and washing buffer (300µL/well)		
⚠ Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.			
INCUBATION OF CONJUGATE	100µL	Anti-h IgG-PEROX Ready to use	30 mn at room temperature
WASHING	Wash 3 times with dilution and washing buffer (300µL/well)		
⚠ Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess liquid.			
ENZYME REACTION	100µL	Substrate solution ABTS	30 mn at room temperature
STOP REACTION	100µL	Stop Solution H ₂ SO ₄ (0.25M)	

SYMBOLS USED

EC Declaration of conformity

ELISA Test

In Vitro Diagnostic Device

Catalogue number

Lot Number

Expiry Date

Number of test

Consult Instructions

Manufactured by

Temperature limitation

Biological risk

Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com