# DNA-DOT

**REF** HM 015

Français

## DÉFINITION

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### BIBLIOGRAPHIE

AHMAN A et al.

Molecular expression systems for anti-DNA antibodies - 1 - 2. Lupus 2002, 11, 824-832, 833-842.

CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance.

Clin. and Exp. Rheumatol. 1997 ;15, 457-465

HOLLINGSWORTH P.N. et al. Antinuclear antibodies

In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.74-90.

HOMBURGER H.A. Cascade testing for autoantibodies in connective tissue disease.

Mayo Clin. Proc., 1995, 70, 183-184.

ISENBERG et al.

Clinical laboratory assays for measuring anti-dsDNA antibodies. Were are we now?

Lupus 2002, 11, 797-800.

RAHMAN A et al.

Anti-DNA antibodies – structure and function.

Lupus 2002, 11, 776-779.

RAVIRAJAN CT et al.

An analysis of clinical disease activity and nephritis associated serum autoantibody profiles in patients with systemic lupus erythematosus: a cross-sectional study.

Rheumatology 2001, 40, 1405-1412.




SMEENK RJ et al.

dsDNA autoantibodies.

In: Peter JB and Shoenfeld Y. Autoantibodies.

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English

## DEFINITION

The DNA-DOT kit is a qualitative enzyme immunoassay (EIA) for the detection of anti-dsDNA antibodies.

## DIAGNOSTIC VALUE

Connective diseases are systemic auto-immune diseases. Their classification is based upon clinical and biological data. Biological criteria is based on the detection of anti-nuclear antibodies.

Among them, anti-dsDNA (ds: double strength) antibodies are recognized as the major serologic marker of Systemic Lupus Erythematosus (SLE). Their specificity and their sensitivity give them a high diagnostic value. Therefore they are a part of the clinical and biological criteria established in 1982 by the American Rheumatism Association (ARA) for the diagnosis of the SLE.

The detection of anti-dsDNA antibodies is particularly useful in two different ways : as an aid to the diagnosis of SLE and as a tool to monitor the course of the disease.

For the second purpose repeated serum sampling of individual patients can be very informative about the clinical course of the disease because a relationship exists between anti-DNA and diseases activity : flares of SLE are generally preceded by a rise in anti-dsDNA levels, followed by a steep drop during the exacerbation (particularly in nephritis).

Furthermore different treatments of patients have varying influences on anti-dsDNA levels and can be adapted by a regular follow-up of these antibodies.

## ASSAY PRINCIPLE


The antigens are coated on distinct membranes on the assay strip.

- First, the assay strip is incubated in a reaction tube containing the diluted sample of the patient. If this sample contains the antibodies, these antibodies will recognize and bind to the corresponding antigen. After incubation, a first wash step removes all of the unbound proteins.
- An alkaline phosphatase labeled protein A conjugate binds to the captured antibodies. The excess of unbound conjugate is removed by a second wash step.
- The bound conjugate is visualised with an insoluble substrate for the phosphatase enzyme (BCIP/NBT). This will give rise to a blue colored circle, confirming the presence of auto-antibodies in the sample tested.

## SAMPLES COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

## REAGENTS

	HM015
Assay strips	25
 <p>control Ag DNA</p> <p>Antigen coated : - dsDNA purified from calf thymus</p> <p><b>STRIP</b></p>	
1 vial of PAL Conjugate (Phosphatase alkaline labeled protein A) To be diluted	650µL
CONJ IgG	
1 vial of Substrate (BCIP/NBT) Ready to use	18mL
SUBS BCIP-NBT	
1 vial of Phosphate-Tween buffer (concentrate 10x) To be diluted in distilled water	30mL
BUF WASH 10x	
Reaction tubes (1,5mL)	100
Self-stick labels for returned results	25

## MATERIAL REQUIRED BUT NOT PROVIDED

- Pipettes
- Timer
- Distilled water (to dilute the concentrated Phosphate-Tween buffer)
- Adsorbent towels (to blot dry the assay strip)
- Support for 96 reaction tubes (bmd, cat. HM501)

## STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## PRECAUTIONS

Be sure that the strips are dry after each wash.

Reagents in solution (except for substrate buffer) contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**DNA-DOT** has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labelling of dangerous preparations.

**DNA-DOT** has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

## SETUP

All reagents must be at room temperature before use.

### 1. Preparation of the wash and dilution buffer

- Dilute 10 times the concentrated Phosphate-Tween buffer in distilled water.
- Storage period: 3 months at +2°C/+8°C

### 2. Preparation of samples

- Mark one strip per sample.
- Mark 3 reaction tubes per sample and number from 1 to 3:
  - ≈ tube 1: sample incubation.
  - ≈ tube 2: conjugate incubation.
  - ≈ tube 3: substrate incubation.
- Dispense 600µL of wash and dilution buffer in tubes 1 and 2.
- Mark another microtube reserved for the two first washing steps and dispenses 600µL of wash and dilution buffer.

### Remarks and precautions:

- ❖ *The substrate is light sensitive. It requires to be dispense immediately at the time of his use in the tubes 3.*
- ❖ *Certain environmental conditions as the temperature may influence the final result. For temperatures > +30°C, it is recommended to be vigilant on the coherence of the results obtained.*
- ❖ *Certain operational conditions may influence the intensity of the blue color and should therefore be avoided:*
  - *the agitation time in each wash step is less than a minute*
  - *the membrane sticking to the wall of the vessel during incubation steps.*
  - *incomplete mixing of sample and conjugate*

## INTERPRETATION OF RESULTS

### 1. Validation criterion for the manipulation:

- The positive control should give a blue colored circle with a well defined outline.

### 2. Reading of the strips:

In order to interpret correctly make sure that the strips are completely dry.

APPEARANCE OF THE WELL	RESULTS
no circle	absence of auto-antibodies NEGATIVE
presence of a circle of which the staining is variable from bright to dark	presence of auto-antibodies POSITIVE

*Each membrane showing a difficult to distinguish circle should be controlled with a second blood sample, drawn a few weeks later. The interpretation of the results should be done in relation to other examinations and in the clinical context.*

## QUALITY CONTROL

bmd offers a multiparametric quality controls (Immunotrol I, Cat. No HM036 et Immunotrol IV, Cat.No HM051). It contains anti-native DNA antibodies; handle in the same way as samples.

## CHARACTERISTICS AND PERFORMANCE OF THE TEST

### 1. INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH AN OTHER METHOD

The **DNA-DOT** has been compared with the DNA-LISA (Cat. HM002) commercialised by bmd.

The study is based on 189 serum samples of patients suffering of pathologies associated with the presence of anti-native DNA antibodies, healthy individuals and patients suffering of other autoimmune diseases.

All the discordant results between the two methods are been controlled by indirect immunofluorescence on *Crithidia luciliae*.

	DNA-LISA	
	Positive	Negative
DNA-DOT	68	3
	1	117

### Controlled by immunofluorescence

- 1 sample negative on the DOT and positive by ELISA : negative on IF
- 3 samples positive on the dot and negative by ELISA : 2 were found negative on IF and 1 was found positive.

### Indication of the diagnostic value of the DNA-DOT :

in comparison with the DNA-LISA (cut off value of 10UA/mL) :

Relative sensibility	Relative specificity	Concordance	Positive predictive value	Negative predictive value
98.5%	97.5%	97.9%	95.8%	99.2%

### 2. ACCURACY

Estimated by a test of recovery :

Principle : negative sample with growing quantities of a positive sample.

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### 3. REPRODUCIBILITY

Estimated on 3 samples : positive, weak positive and negative.

DNA-LISA Titre	INTRA-ASSAY		INTER-ASSAYS	
	n	DNA-DOT appearance	n	DNA-DOT appearance
46U/mL	5	dark circle	5	dark circle
17U/mL	5	bright circle	5	bright circle
<7U/mL	5	no circle	5	no circle

### 4. PERFORMANCE COMPARISON OF DNA-DOT, NEW VERSION (SHORT PROCEDURE) TO DNA.-DOT, PREVIOUS VERSION (LONG PROCEDURE)

The study is based on 94 serum samples of patients suffering of pathologies associated with the presence of anti-native DNA antibodies, healthy individuals and patients suffering of other autoimmune diseases.

DNA-DOT (long procedure)	DNA-DOT (short procedure)		
		Positive	Negative
	Positive	51	1
Negative	1	41	

Indication of the diagnostic value of the DNA-DOT (short procedure):




in comparison with the DNA-DOT (long procedure) :

Relative sensibility	Relative specificity	Concordance
98.1%	97.6%	97.8%

### BIBLIOGRAPHY

- AHMAN A et al.  
Molecular expression systems for anti-DNA antibodies - 1 – 2.  
Lupus 2002, 11, 824-832, 833-842.
- CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance.  
Clin. and Exp. Rheumatol. 1997 ;15, 457-465
- HOLLINGSWORTH P.N. et al. Antinuclear antibodies  
In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.74-90.
- HOMBURGER H.A. Cascade testing for autoantibodies in connective tissue disease.  
Mayo Clin. Proc., 1995, 70, 183-184.
- ISENBERG et al.  
Clinical laboratory assays for measuring anti-dsDNA antibodies. Were are we now?  
Lupus 2002, 11, 797-800.
- RAHMAN A et al.  
Anti-DNA antibodies – structure and function.  
Lupus 2002, 11, 776-779.
- RAVIRAJAN CT et al.  
An analysis of clinical disease activity and nephritis associated serum autoantibody profiles in patients with systemic lupus erythematosus: a cross-sectional study.  
Rheumatology 2001, 40, 1405-1412.
- SMEENK RJ et al.  
dsDNA autoantibodies.  
In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234.

### SYMBOLS USED

	Catalog number
	Batch code
	Number of tests
	Rapid test
	In vitro diagnostic device
	Use by
	Storage temperature limitation
	Read instructions for use
	EC Declaration of Conformity
	Reconstitute with water
	Contains sodium azide
	Contains Proclin

**SUMMARY OF METHOD**

	<b>OPERATIONS STEPS</b>	<b>INCUBATION TIME</b>
<b>INCUBATION OF SAMPLES</b>	<ul style="list-style-type: none"> <li>- Add 20µL sample to tube 1 containing wash and dilution buffer.</li> <li>- Put in the identified strip.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	10 minutes at R.T
During sample incubation, add 20µL conjugate to tube 2 containing wash and dilution buffer		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>INCUBATION WITH CONJUGATE</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 2.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	10 minutes at R.T
At the end of the conjugate incubation, dispense 600µL of substrate in tube 3		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>INCUBATION OF SUBSTRATE</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 3.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	5 minutes at R.T
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Wash the strip under tap water.</li> </ul>	
<b>DRY</b>	<ul style="list-style-type: none"> <li>- Remove water by blotting the strip with an adsorbent towel and <b>let dry completely</b> before reading the result.</li> </ul>	about 30 minutes at R.T

R.T : Room Temperature between +18°C and + 25°C. \* : Up and down agitations should be slow.

**BioMédical Diagnostics SA**

**Office**  
 Actipole 25  
 4 bld de Beaubourg  
 77435 Marne la Vallée Cx2  
 France

Tel : 33 1 64 62 10 12  
 Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
 Internet : [www.bmd-net.com](http://www.bmd-net.com)