# ANCA-MBG-DOT

**REF** HM 025

## DEFINITION

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## BIBLIOGRAPHIE

- HAGEN EC et al.  
Development and standardization of solid-phase assays for the detection of anti-neutrophil cytoplasmic antibodies (ANCA) for clinical application : report of a large clinical evaluation study. Proceedings of the 6th International ANCA Workshop. Paris, 28 June - 1 July 1995. *Clin Exp Immunol*, 101/1: 29, 1995.
- HELLMARK T et al.  
Characterization of anti-GBM antibodies involved in Goodpasture's syndrome. *Kidney Int.* 1994, 46, 823-829.
- HERODY M et al.  
Anti-GBM disease: predictive value of clinical, histological and serological data. *Clin. Nephrol.* 1993, 40, 249-255.
- HUDSON BG et al.  
Molecular characteristics of the Goodpasture autoantigen. *Kidney Int.*, 1993, 43, 135-139.
- HUMBEL RL.  
Néphropathies d'origine immunologique. Autoanticorps et maladies autoimmunes. *Ed. Elsevier* 1994, 171-184.

JENNETTE JC et al.  
Pathogenic potential of anti-neutrophil cytoplasmic auto-antibodies. *lab Invest* 70: 135-137, 1994.

JOHANSSON C et al.  
The structural organization of type IV collagen. Identification of three NCI populations in the glomerular basement membrane. *J.Biol. Chem.* 1992, 267, 24533-24537.

KELLY D et al. Goodpasture syndrome: molecular and clinical advances. *Medicine (Baltimore)* 1994, 73, 171-185.

KLEPPEL MM et al. Human tissue distribution of novel basement membrane collagen. *Am. J. Pathol.* 1989, 134, 813-825.

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**REF** HM 025

## DEFINITION

The **ANCA-MBG-DOT** kit is a qualitative enzyme immunoassay (EIA) for the simultaneous detection of :

- **anti-neutrophil cytoplasm antibodies** (ANCA) directed against Myeloperoxidase (MPO) and Serine Proteinase 3 (PR3).
- **anti-glomerular basement membrane** (GBM) antibodies using a purified antigen based on type IV

## DIAGNOSTIC VALUE

### • ANCA

ANCA antibodies form a family of autoantibodies for which the diagnostic and predictive value in differentiating necrotising vasculitis and ‘overlap’ systemic vasculitis, associated with glomerulonephritis, is now well established.

The autoantibodies are directed against antigens present in the primary azurophilic granules (alpha granules) of polynuclear neutrophils and in the lysosomes of monocytes.

ANCA antibodies are first detected by means of indirect immunofluorescence using human ethanol-fixed neutrophil slides. The presence of ANCA is characterised in two distinct staining patterns: one cytoplasmic and the other perinuclear.

These two types called c-ANCA and p-ANCA respectively, react predominantly antibodies with proteinase 3 (PR3) and myeloperoxidase (MPO).

The detection of ANCA with Immuno Fluorescent Assay (IFA) only, may result in problems in the differential diagnosis because very similar patterns are created by other autoantibodies. This also holds true in Wegener’s granulomatosis and vasculitis, where the association of ANCA with other autoantibodies is described (rheumatoid factor, antinuclear antibodies, anti-smooth muscle antibodies, anti-thyroid, etc.) as in collagenosis, and autoimmune hepatitis where, in addition to specific autoantibodies, ANCA could also develop.

The diagnostic sensitivity of ANCA is of particular importance with Wegener’s granulomatosis (anti-PR3), Churg Strauss syndrome, microscopic periarteritis and idiopathic crescentic glomerulonephritis (anti-MPO).

These diseases represent a medical emergency and require a safe strategy to identify ANCA in two ways:

- ✦ *By indirect immunofluorescence where two known patterns indicate a general diagnosis.*
- ✦ *A reliable means to differentiate the specificity of the antibodies against the target antigens (MPO and PR3) using specific immunologic techniques and purified antigens.*

### • GBM

The measurement of anti-glomerular basement membrane (GBM) autoantibodies is intended as an aid in the diagnosis of Goodpasture’s syndrome, based on clinical signs of rapidly progressive renal failure and lung haemorrhage. It is an autoimmune disease caused by autoantibodies reacting with the components of the glomerular basement membrane (alpha chain 3 of type IV collagen).

This disease progresses rapidly and in a fatal manner in more than 75% of the cases without treatment. An early diagnosis, based on the presence of the specific marker, anti-GBM autoantibodies is indispensable in minimizing significant damage to the kidney and often the lung. In addition, the continuing presence of these autoantibodies indicates that there should be a delay in renal transplantation.


The detection of anti-GBM antibodies allows one to distinguish Goodpasture’s syndrome from other nephropathies of immunologic origin (glomerulonephritis associated with anti-neutrophil cytoplasmic antibodies) which requires a totally different treatment. A negative result that allows the clinician to rule Goodpasture’s syndrome as a possible diagnosis is very important.

## ASSAY PRINCIPLE

The antigens are coated on distinct membranes on the assay strip.

- First, the assay strip is incubated in a reaction tube containing the diluted sample of the patient. If this sample contains at least one of the antibodies, these antibodies will recognize and bind to the corresponding antigen.
- After incubation, a first wash step removes all of the unbound proteins.
- An alkaline phosphatase labeled protein A conjugate binds to the captured antibodies. The excess of unbound conjugate is removed by a second wash step.
- The bound conjugate is visualised with an insoluble substrate for the phosphatase enzyme (BCIP/NBT). This will give rise to a blue colored circle, confirming the presence of auto-antibodies in the sample tested.

## REAGENTS

|  |   | HM025 |
|--|---|-------|
| Assay strips   |   |       |
|        | <p><i>Antigens coated:</i></p> <ul style="list-style-type: none"> <li>- Human purified myeloperoxidase,</li> <li>- Human purified proteinase 3,</li> <li>- Bovine purified “Goodpasture” antigen</li> </ul> <p><b>STRIP</b></p> | 10    |
| 1 vial of PAL Conjugate (Phosphatase alkaline labeled protein A)<br>To dilute              | <b>CONJ</b> <b>IgG</b>  | 300µL |
| 1 vial of Substrate (BCIP/NBT)<br>Ready to use   | <b>SUBS</b> <b>BCIP-NBT</b>   | 12mL  |
| 1 vial of Phosphate-Tween buffer pH 7,2 (concentrated 10x)<br>To dilute in distilled water | <b>BUF</b> <b>WASH</b> <b>10x</b>   | 30mL  |
| Reaction tubes (1,5mL)   |   | 40    |
| Self-stick labels for returned results   |   | 10    |

## SAMPLES COLLECTION AND HANDLING

---

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

## MATERIAL REQUIRED BUT NOT PROVIDED

---

- Pipettes
- Timer
- Distilled water (to dilute the concentrated Phosphate-Tween buffer)
- Adsorbent towels (to blot dry the assay strip)
- Support for 96 microtubes (bmd, cat. HM501)

## PRECAUTIONS

---

Be sure that the strips are dry after each wash.

Reagents in solution (except for substrate buffer) contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**ANCA-MBG-DOT** has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labelling of dangerous preparations.

**ANCA-MBG-DOT** has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

## STABILITY AND STORAGE

---

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## SETUP

---

All reagents must be at room temperature before use.

### 1. Preparation of the wash and dilution buffer

- Dilute 10 times the concentrated Phosphate-Tween buffer in distilled water.
- Storage period: 3 months at +2°C/+8°C

### 2. Preparation of samples

- Mark one strip per sample.
- Mark 3 reaction tubes per sample and number from 1 to 3:
  - ∞ tube 1: sample incubation.
  - ∞ tube 2: conjugate incubation.
  - ∞ tube 3: substrate incubation.
- Dispense 1mL of wash and dilution buffer in tubes 1 and 2.
- Mark another microtube reserved for the two first washing steps and dispenses 1mL of wash and dilution buffer.

### Remarks and precautions:

- ❖ The substrate is light sensitive. It requires to be dispense immediately at the time of his use in the tubes 3.
- ❖ Certain environmental conditions as the temperature may influence the final result. For temperatures > +30°C, it is recommended to be vigilant on the coherence of the results obtained.
- ❖ Certain operational conditions may influence the intensity of the blue color and should therefore be avoided:
  - the agitation time in each wash step is less than a minute
  - the membrane sticking to the wall of the vessel during incubation step.
  - incomplete mixing of sample and conjugate

## INTERPRETATION OF RESULTS

---

### 1. Validation criterion for the manipulation:

- The positive control should give a blue colored circle with a well defined outline.

### 2. Reading of the strips:

In order to interpret correctly make sure that the strips are completely dry.

For each specificity, interpret accordint to this scheme:

| APPEARANCE OF THE WELL   | RESULTS                                 |
|--|---|
| no circle  | absence of auto-antibodies<br>NEGATIVE  |
| presence of a circle of which the staining is variable from bright to dark | presence of auto-antibodies<br>POSITIVE |

*Each membrane showing a difficult to distinguish circle should be controlled with a second blood sample, drawn a few weeks later. The interpretation of the results should be done in relation to other examinations and in the clinical context.*

## QUALITY CONTROL

---

bmd offers an external multiparametric quality control (Immuno-Trol I, Cat: HM036), it contains antibodies directed against ANCA (MPO and PR3 antigens) and GBM. This control should be handled as samples.

## CHARACTERISTICS AND PERFORMANCE OF THE TEST

---

### 1. INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH AN OTHER METHOD

- ANCA

The MPO and PR3 specificities have been compared with an ELISA on microtiterplates.

The study is based on 172 serum samples of patients suffering of pathologies associated with the presence of ANCA antibodies, healthy individuals and patients suffering of other autoimmune diseases.

All the discordant results between the two methods are been controlled by indirect immunofluorescence.

| Anti-MPO |          |              |          |
|----------|----------|--------------|----------|
|          |          | ELISA method |          |
|          |          | Positive     | Negative |
| ANCA-DOT | Positive | 32           | 1        |
|          | Negative | 0            | 139      |

| Anti-PR3 |          |              |          |
|----------|----------|--------------|----------|
|          |          | ELISA method |          |
|          |          | Positive     | Negative |
| ANCA-DOT | Positive | 20           | 8        |
|          | Negative | 2            | 142      |

#### Anti-MPO

#### Anti-PR3

#### Indication of the diagnostic value for selected samples

Relative sensitivity: 100%

Relative specific affinity: 99%

Concordance: 99%

Relative sensitivity: 91%

Relative specific affinity: 95%

Concordance: 94%

#### Controlled by immunofluorescence

↪ 1 sample positive on the DOT and negative by ELISA. This serum is positive for PR3 on the two methods. It gives p-ANCA and c-ANCA fluorescence pattern.

↪ 2 samples negative on the DOT and positive by ELISA, give a p-ANCA pattern.

↪ 8 samples positive on the DOT and negative by ELISA, 7 give a c-ANCA pattern.

#### • GBM

GBM has been compared with an ELISA on microtiterplates.

The study is based on 96 serum samples of patients suffering of Goodpasture's syndrome, healthy individuals and patients suffering of other autoimmune diseases.

| Anti-GBM |          |              |          |
|----------|----------|--------------|----------|
|          |          | ELISA method |          |
|          |          | Positive     | Negative |
| GBM-DOT  | Positive | 10           | 0        |
|          | Negative | 0            | 86*      |

\*: one sample was on the limit of 20 Units on ELISA

#### Indication of the diagnostic value for selected samples

Relative sensitivity: 100 %

Relative specific affinity: 100 %

Concordance: 100 %

#### 2. ACCURACY

Estimated for each specificity by a test of recoverment :

**Principle:** negative sample with growing quantities of a positive sample.

| Titre of the positive samples |                  |                   |
|-------------------------------|------------------|-------------------|
| anti-MPO: 84U/mL              | anti-PR3: 78U/mL | anti-GBM: 170U/mL |

| % negative sera | % positive sera | anti-MPO         | anti-PR3         | Anti-GBM         |
|-----------------|-----------------|------------------|------------------|------------------|
| 100             | 0               | no circle        | no circle        | no circle        |
| 90              | 10              | No circled halot | No circled halot | No circled halot |
| 80              | 20              | Pale circled     | Pale circled     | Pale circled     |
| 70              | 30              | Pale circled     | Pale circled     | Pale circled     |
| 60              | 40              | Pale circled     | Pale circled     | Pale circled     |
| 50              | 50              | Dark circled     | Dark circled     | Dark circled     |

#### 3. REPRODUCIBILITY

Estimated on 3 samples: positive, dubious or weak positive and negative.

|          | INTRA-ASSAY |   |                         | INTER-ASSAY |                         |
|----------|-------------|---|-------------------------|-------------|-------------------------|
|          | Titre       | n | ANCA-GBM-DOT appearance | n           | ANCA-GBM-DOT appearance |
| anti-MPO | 84U/ml      | 5 | dark circle             | 5           | dark circle             |
|          | 17U/ml      | 5 | pale circle             | 5           | pale circle             |
|          | 2U/ml       | 5 | no spot                 | 5           | no spot                 |
| anti-PR3 | 78U/ml      | 5 | dark circle             | 5           | dark circle             |
|          | 14U/ml      | 5 | pale circle             | 5           | pale circle             |
|          | 3U/ml       | 5 | no spot                 | 5           | no spot                 |
| Anti-GBM | 290U/ml     | 5 | dark circle             | 5           | dark circle             |
|          | 49U/ml      | 5 | pale circle             | 5           | pale circle             |
|          | <10U/ml     | 5 | no circle               | 5           | no circle               |

The comparison of the performance GBM-DOT, new version (little strip) to GBM-DOT, previous version (large strip) shows an agreement of 100% between the two versions.

#### SYMBOLS USED

|  |                                |
|--|--------------------------------|
|  | Catalog number                 |
|  | Batch code                     |
|  | Number of tests                |
|  | Rapid test                     |
|  | In vitro diagnostic device     |
|  | Use by                         |
|  | Storage temperature limitation |
|  | Read instructions for use      |
|  | EC Declaration of Conformity   |
|  | Reconstitute with water        |
|  | Contains sodium azide          |
|  | Contains Proclin               |

#### BIBLIOGRAPHY

HAGEN EC et al.

Development and standardization of solid-phase assays for the detection of anti-neutrophil cytoplasmic antibodies (ANCA) for clinical application : report of a large clinical evaluation study. Proceedings of the 6th International ANCA Workshop. Paris, 28 june - 1 july 1995. *Clin Exp Immunol*, 101/1: 29, 1995.

HELLMARK T et al.

Characterization of anti-GBM antibodies involved in Goodpasture's syndrome. *Kidney Int.* 1994, 46, 823-829.

HERODY M et al.

Anti-GBM disease: predictive value of clinical, histological and serological data. *Clin. Nephrol.* 1993, 40, 249-255.

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JENNETTE JC et al.  
Pathogenic potential of anti-neutrophil cytoplasmic auto-antibodies. *lab Invest* 70: 135-137, 1994.

JOHANSSON C et al.  
The structural organization of type IV collagen. Identification of three NC1 populations in the glomerular basement membrane. *J.Biol. Chem.* 1992, 267, 24533-24537.

KELLY D et al. Goodpasture syndrome: molecular and clinical advances. *Medicine (Baltimore)* 1994, 73, 171-185.

KLEPPEL MM et al. Human tissue distribution of novel basement membrane collagen. *Am. J. Pathol.* 1989, 134, 813-825.

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SAXENA R et al.  
Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome. *Intern. Med.* 1995, 237, 143-152.

## SUMMARY OF METHOD

|  | OPERATIONS STEPS  | INCUBATION TIME         |
|--|---|-------------------------|
| <b>INCUBATION OF SAMPLES</b>   | <ul style="list-style-type: none"> <li>- Add 20µL sample to tube 1 containing wash and dilution buffer.</li> <li>- Put in the identified strip.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul> | 45 minutes at R.T       |
| During sample incubation, add 20µL conjugate to tube 2 containing wash and dilution buffer |   |                         |
| <b>WASH</b>  | <ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>                                | 1 minute at R.T         |
| <b>INCUBATION OF CONJUGATE</b>   | <ul style="list-style-type: none"> <li>- Put the strip in tube 2.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>   | 20 minutes at R.T       |
| At the end of the conjugate incubation, dispense 1mL of substrate in tube 3                |   |                         |
| <b>WASH</b>  | <ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>                                | 1 minute at R.T         |
| <b>INCUBATION OF SUBSTRATE</b>   | <ul style="list-style-type: none"> <li>- Put the strip in tube 3.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>   | 10 minutes at R.T       |
| <b>WASH</b>  | <ul style="list-style-type: none"> <li>- Wash the strip under tap water.</li> </ul>   |                         |
| <b>DRY</b>   | <ul style="list-style-type: none"> <li>- Remove water by blotting the strip with an adsorbent towel and <b>let dry completely</b> before reading the result.</li> </ul>   | about 30 minutes at R.T |

R.T : Room Temperature between +18°C and + 25°C.

\* : Up and down agitations should be slow.

### BioMédical Diagnostics SA

#### Office

Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)