# FR-LISA IgM

REF

HM 028

2 X

## DÉFINITION

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2 X 

## DEFINITION

**FR-LISA IgM** is an enzyme linked immunoassay (ELISA) for the quantitative determination of IgM class rheumatoid factors directed against Fc determinants from human or animal G immunoglobulins in human serum samples.

## DIAGNOSTIC VALUE

The rheumatoid factors (RF) constitute a population of autoantibodies which react with different antigen sites located on the Fc fragment of class G immunoglobulins. These antigen determinants may be specific to human IgG or common to the IgG of various animal species. The rheumatoid factors may belong to all the immunoglobulin classes, but the majority is IgM.

The rheumatoid factors are essentially associated with inflammatory rheumatic pathologies. They occur with increased frequency in the serum of patients suffering from rheumatoid polyarthritis. The specific affinity of the RF for rheumatoid polyarthritis is increased by the finding that there is an associated positivity for human origin G immunoglobulin and for rabbit immunoglobulin G.

Widely used diagnostic methods are latex agglutination tests which identify the IgM class RF directed against the human antigen determinants and the Waaler-Rose haemagglutination reaction which identifies IgM class RF directed against rabbit antigen determinants.

However, immunoenzymatic methods, which are more sensitive and more specific, have the advantage of being able to measure the two types of specific affinity simultaneously.

Apart from rheumatic pathologies, RF may appear in the course of some infectious diseases of bacterial or viral origin or during inflammatory pathologies. Their prevalence then is in the region of 2%. They are small in number and their reactivity is directed solely against human origin IgG. In the normal population, RF prevalence is between 2 and 5%.

## STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in the ir own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## ASSAY PRINCIPLE

Each antigen (human or animal Fc fragment) is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the antigen coated well for each specific affinity, which allows to bind. After incubation, unbound proteins are removed by washing.
- Binding of rheumatoid factors is detected using an anti-human IgM alkaline phosphatase labeled conjugate. After incubation, the wells are washed again to eliminate any excess of conjugate.
- The bound enzyme is revealed by addition of substrate (Para Nitro Phenyl Phosphate, producing a yellow color. The color intensity is proportional to rheumatoid factors level in the sample.
- After stopping the reaction by 1N sodium hydroxide, the optical density is red by a spectrophotometer at 405nm
- Rheumatoid factors levels of sample patients are expressed in international units (First British Standard for Rheumatoid Arthritis 64/2) for each specificity, using a standard curve or a calibrator.

## SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

## MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- timer
- incubator adjusted at +37°C
- physiological serum (NaCl 0.9%)
- microplate spectrophotometer with 405nm filter
- 8 channel pipettes

## REAGENTS

Microplate coated with Human Fc fragment (12x8 individual breakaway wells with holder): <b>colorless wells</b> <span style="border: 1px solid black; padding: 2px;">MP</span>	12 strips
Microplate coated with Animal Fc fragment (12x8 individual breakaway wells with holder): <b>red wells</b> <span style="border: 1px solid black; padding: 2px;">MP</span>	12 strips
4 Standards for the human specificity tittered in international units <u>Ready to use</u> <span style="border: 1px solid black; padding: 2px;">CAL</span> <span style="border: 1px solid black; padding: 2px;">n</span> <i>The values are shown on the vial label.</i>	4 x 1.5mL
4 Standards for the animal specificity tittered in international units <u>Ready to use</u> <span style="border: 1px solid black; padding: 2px;">CAL</span> <span style="border: 1px solid black; padding: 2px;">n</span> <i>The values are shown on the vial label.</i>	4 x 1.5mL
1 vial Positive control for the human specificity tittered in international units <u>To dilute</u> <span style="border: 1px solid black; padding: 2px;">CONTROL</span> <span style="border: 1px solid black; padding: 2px;">+</span> <i>The values expected are shown on the vial label.</i>	1 x 350µL
1 vial Positive control for the animal specificity tittered in international units <u>To dilute</u> <span style="border: 1px solid black; padding: 2px;">CONTROL</span> <span style="border: 1px solid black; padding: 2px;">+</span> <i>The values expected are shown on the vial label.</i>	1 x 350µL
1 vial of Negative control common to both specificities <u>To dilute</u> <span style="border: 1px solid black; padding: 2px;">CONTROL</span> <span style="border: 1px solid black; padding: 2px;">-</span>	1 x 350µL
1 vial Goat anti-human IgM conjugate for the human specificity <u>Ready to use</u> <span style="border: 1px solid black; padding: 2px;">CONJ</span> <span style="border: 1px solid black; padding: 2px;">IgM</span>	1 x 15mL
1 vial of Goat anti-human IgM conjugate for the animal specificity <u>Ready to use</u> <span style="border: 1px solid black; padding: 2px;">CONJ</span> <span style="border: 1px solid black; padding: 2px;">IgM</span>	1 x 15mL
1 vial of Substrate containing 12x5mg PNPP tablets <u>To dilute</u> <span style="border: 1px solid black; padding: 2px;">SUBS</span> <span style="border: 1px solid black; padding: 2px;">PNPP</span>	1
1 vial of Phosphate-Tween buffer pH 7,2 (concentrated 10x) <u>To dilute in distilled water</u> <span style="border: 1px solid black; padding: 2px;">BUF</span> <span style="border: 1px solid black; padding: 2px;">WASH</span> <span style="border: 1px solid black; padding: 2px;">10x</span>	1 x 100mL
1 vial of Substrate dilution buffer <u>Ready to use</u> <span style="border: 1px solid black; padding: 2px;">BUF</span> <span style="border: 1px solid black; padding: 2px;">SUBS</span>	1 x 62mL
1 vial of Stop solution (Sodium hydroxide 1N) <u>Ready to use</u> <span style="border: 1px solid black; padding: 2px;">SOLN</span> <span style="border: 1px solid black; padding: 2px;">STOP</span>	1 x 15mL

## SETUP

All the reagents should be prepared as required:

### 1. Dilution and washing buffer.

- Dilute concentrated Phosphate-Tween Buffer 1/10 in distilled water.
- Storage period: 3 month at +2°C/+8°C.  
*NB. If there are crystals in the concentrated solution, warm the bottle to+ 37°C for 15 minutes before use.*

### 2. Preparation of samples and controls

- Dilute to 1/201 in dilution buffer (10µL sample + 2 mL dilution buffer).
- Vortex vigorously.

### 3. Preparation of ready-to-use conjugate.

- Estimate amount required for handling and transfer to a tube.

### 4. Preparation of the substrate.

- Dissolve the PNPP tablets in the substrate dilution buffer at a ratio of 1 tablet per 5 mL. Protect from light until use.

## PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C/+25°C) before handling.

Check that all plates are well drained after each wash.

Do not put plates on top of each other during incubation at +37°C.

Check that no condensed water is left at the bottom of the wells before reading.

*Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.*

Reagents in solution (except for substrate buffer and stop solution) contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**FR-LISA IgM** has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

**FR-LISA IgM** has been optimized for the use as describe in this procedure. Do not substitute reagents from other manufacturers. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

## METHOD

### 1. Preparing the test

Use the work sheet to note the sample locations.

- For each specific affinity test, set out:
  - one blank reagent well (made up with dilution buffer)
  - 4 standard wells or 1 "standard 4" well, depending on the chosen method of calibration (**in duplicate**) (cf. "Interpretation of results")
  - 1 well for each control (**in duplicate**)
  - to obtain more precise results, it is also advisable to run samples in duplicate.
- then detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

### 2. Sample incubation

Add 100 µl of standard(s), diluted controls or samples.

Incubate for 30 minutes at room temperature.

*Wash step:*

Remove the contents of the wells by rapid inversion.

Wash 3 times with 200µL of washing buffer (200µL/wells)

### 3. Incubation of conjugate

Add 100µL of each conjugate.

Incubate for 30 minutes at room temperature.

*Wash step:*

Remove the contents of the wells by rapid inversion.

Wash 3 times with 200µL of washing buffer (200µL/wells).

Make a last wash with 200µL of 0,9% sodium chloride.

### 4. Incubation of substrate

Add 100µL substrate into each well.

Incubate for 15 minutes in the dark at +37°C

### 5. Stop solution

Add 50µL of 1N sodium hydroxide to each well.

### 6. Reading

Read the optical density of each well at 405 nm against the reagent blank set at zero absorbance within 30 minutes.

## RESULTS

For each specificity, the results can be calculated by one of the following two methods:

### Standard curve

- Determine the mean optical density for all duplicate readings.
- Trace the standard curve on log-log paper, plotting the units of the 4 standard points along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).

### One point calibration

- Calculate the mean OD for standard 4.
- Calculate the conversion factor equal to: titre at standard 4/ OD for standard 4.
- Multiply the OD of each sample well by this factor. Deduce the titre in rheumatoid factor.

### Sample calculation:

- For a calibrator equal to 250 UI/ml, giving a mean OD of 2.5, the conversion factor is equal to 250/2.5, or 100.

A sample with OD 1.5, corresponds to a titre of 1.5 x 100, or 150 UI/ml.

*Examples of results obtained using the two methods of calculation*

	Standard curve	One point calibration
<b>Strongly positive sample</b>	238UI/mL	260UI/mL
<b>Positive sample</b>	93UI/mL	109UI/mL
<b>Weakly positive sample</b>	32UI/mL	39UI/mL

### Comments:

- The mean OD of the standard 4 should be at least 0.8.
- The samples with values greater than that of standard 4 may be diluted to 1/400 to obtain a more precise result. The number of units should be multiplied by 2.

## INTERPRETATION OF RESULTS

Specific affinity	negative	limit	positive
human	<15 UI/ml	15-20 UI/ml	> 20 UI/ml
animal	<25 UI/ml	25-30 UI/ml	> 30 UI/ml

## QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive control concentration should fall within the range printed on the positive control label. If either reagent control is invalid, the test should be repeated.

bmd offers an external multiparameter quality control (Immunotrol III, Cat: HM050) for determination of human and animal anti-Rheumatoid Factors (IgM isotype) autoantibodies. This control should be handled as samples.

## CHARACTERISTICS AND PERFORMANCE OF THE TEST

### 1. SENSITIVITY, SPECIFIC AFFINITY AND AGREEMENT WITH ANOTHER METHOD

The "**FR-LISA IgM**" test was compared with two methods, depending on the specific affinity investigated:

- *human specific affinity*: quantitative ELISA microplate method: EUROGENETICS RF (EUROGENETICS FRANCE -45072 Orléans La Source)
- *animal specific affinity*: quantitative haemagglutination microplate method based on the principle of the WAALER-ROSE reaction: POLYARTITRE (FUMOUCHE - 92600 Asnières).

The study involved 245 samples, taken from:

- ⇒ patients for whom the clinical diagnosis of rheumatoid polyarthritis had already been established
- ⇒ healthy individuals, or patients with other autoimmune diseases.

In the case of human specific affinity, any discordant result obtained from the two methods was verified by an agglutination method using latex particles and marketed by BEHRING.

Mean recovery: 103%.

<b>Human specific affinity</b>		Method of reference	
		Positives	Negatives
FR-LISA	Positives	80	1
	Negatives	3	161

<b>Animal specific affinity</b>		Method of reference	
		Positives	Negatives
FR-LISA	Positives	53	0
	Negatives	3	189

### Diagnostic value

#### Human specificity

Relative sensitivity: 96.3%

Relative specificity: 99.4%

Agreement: 98.4%

The 4 discordant samples were verified by the haemagglutination method:

- of the 3 results negative by FR-LISA and positive by the reference method, 2 were found negative and one at the threshold limit of positive (1/160).
- the result positive by FR-LISA and negative by the reference method, was found negative.

#### Animal specificity

Relative sensitivity: 94.6%

Relative specificity: 100 %

Agreement: 98.8%

The 3 discordant samples yielded values close to the threshold limit for the two techniques.

### 2. ACCURACY

(Estimated by recovery test)

*Principle:* negative sample is supercharged with increasing quantities of a positive sample.

The overlap percentage equals (titre obtained / theoretical titre) x 100.

The theoretical titre was calculated on the basis of the values obtained for the two samples tested separately.

% negative serum	% positive serum	HUMAN SPECIFICITY		
		negative sample (7 IU/ml) and positive sample (181 IU/ml)		
		Titre obtained	Theoretical titre	Overlap (%)
100	0	7	7	100
80	20	39	42	93
60	40	73	76	96
40	60	117	112	104
20	80	153	146	105
0	100	181	181	100

Mean recovery: 99.6%.

% negative serum	% positive serum	HUMAN SPECIFICITY		
		negative sample (7 IU/ml) and positive sample (181 IU/ml)		
		Titre obtained	Titre obtained	Titre obtained
100	0	8	8	100
80	20	31	29	107
60	40	54	52	104
40	60	76	73	104
20	80	98	95	103
0	100	117	117	100

### 3. REPRODUCIBILITY

(Estimated on the basis of 3 samples, strongly positive, weakly positive and negative)

	INTRA-ASSAY			INTER-ASSAYS		
	Titre	n	CV (%)	n	CV (%)	
Human specificity	219 UI/ml	5	2.5	5	3.5	
	44 UI/ml	5	3.7	5	4.4	
	14 UI/ml	5	5.6	5	6.5	
Animal specificity	117 UI/ml	5	3.9	5	4.5	
	38 UI/ml	5	4.8	5	6.4	
	22 UI/ml	5	4.5	5	4.9	

### 4. METHOD FOR DETERMINING THE THRESHOLD VALUES ACCORDING TO STANDARD 64/2

Selected population: 100 sera from healthy individuals previously tested by the method of reference. Those found positive were excluded from the calculations below.

The sera (93 for human specificity, 99 for animal specificity) and standard 64/2 were tested in the same run.

- the **negative threshold** corresponds to the mean of the OD + 3 standard deviations, obtained for the population tested: 15 and 25 IU/ml respectively for human and animal specificity.
- the **positive threshold** corresponds to the mean of the OD + 5 standard deviations, obtained for the population tested: 20 and 30 IU/ml respectively for human and animal specificity.
- Between these values, sample titres are considered **borderline**.

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
## SUMMARY OF METHOD

	<i>VOLUME TO BE DISTRIBUTED</i>	<i>HUMAN SPECIFICITY</i>	<i>ANIMAL SPECIFICITY</i>	<i>INCUBATION TIME</i>
<i>INCUBATION OF SAMPLES</i>	100µL	dilution buffer	dilution buffer	30 min at Room Temperature
	100µL	blank reagent standard(s)	blank reagent standard(s)	
	100µL	diluted controls and samples	diluted controls and samples	
<i>WASHING</i>	Wash 3 times with washing buffer (200µL/ well)			
<i>INCUBATION OF CONJUGATES</i>	100µL	Anti-IgM PAL conjugate	Anti-IgM PAL conjugate	30 min at Room Temperature
<i>WASHING</i>	Wash 3 times with washing buffer (200µL/ well) Wash once with NaCl 0.9% (200µL/ well)			
<i>ENZYME REACTION</i>	100µL	Substrate solution (PNPP)		15 min at 37°C
<i>STOP REACTION</i>	50µL	NaOH 1N		

## SYMBOLS USED

	EC Declaration of conformity		Number of test
	ELISA Test		Consult Instructions
	<i>In Vitro</i> Diagnostic Device		Manufactured by
	Catalogue number		Temperature limitation
	Lot Number		Biological risk
	Expiry Date		Reconstitute with

***BioMédical Diagnostics SA***

**Office**  
Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)

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Rheumatoid Factors. *In : Autoantibodies. JB Peter and Y. Schoenfeld Eds. Elsevier Sciences B.V. 1996 ; 706-715.*