**REF** HM 037 2 x 0,5ml

**REF** HM 137 100 x 0,5ml

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**REF**    **HM 037**    2 x 0,5mL    **REF**    **HM 137**    100 x 0,5mL

#### INTENDED USE

This reagent, prepared from human serum, is intended for use as an unassayed precision control with **DOT, ELISA and FIDIS procedures** for determination of the following autoantibodies:

- anti-thyropoxidase
- anti-thyroglobulin
- anti-Cardiolipin IgM isotype
- anti-Phospholipides IgM Isotype
- anti-Beta 2 - glycoprotein 1 IgM and IgG isotypes
- anti-gliadin IgA and IgG isotypes
- anti-transglutaminase IgA isotype

Codes
TPO
TG
aCL IgM
aPL IgM
$\beta_2$ GP1 IgG and IgM
Glia IgA and IgG
tTG IgA

#### DESCRIPTION

**IMMUNO-TROL** reagents have been designed for use with **DOT, ELISA and FIDIS** methods for monitoring assay performance and maintaining quality of assays.

Each **IMMUNO-TROL** reagent lot is designed to produce a positive reaction for each autoantibodies above listed.

**IMMUNO-TROL** reagents do not have assigned values.

Each laboratory should establish its own target value for each marker and for each new lot on a minimum of 30 assays.

#### PACKAGING

**HM 037 : 2 vials de 0,5mL**  
**HM 137 : 100 vials de 0,5mL**

**CONTROL +**

The antigenic specificities are confirmed on the vial label. The antibody codes are specified above.

#### WARNING

- Human sources for the preparation of this control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus or hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.
- The reagent contains 0.02% ProClin® or less of 0.1% (w/v) NaN<sub>3</sub>. Do not eat and avoid contact with skin and eyes.

#### PROCEDURE

**IMMUNO-TROL** reagents should be analyzed in the same manner as patient samples, according to the instructions supplied by manufacturers of the test kits to be used.

#### STABILITY AND STORAGE

- Store the reagents at +2°C to +8°C when not in use.
- Do not use reagents beyond the expiration date.
- Open vial stability : 60 days at +2°C to +8°C.
- Do not freeze.

#### LIMITS

- Performances and levels of reactivity of the **IMMUNO-TROL** reagents may vary with different manufacturers' test and assay procedures.
- **IMMUNO-TROL** reagents must not be substituted for the positive and the negative control reagents provided with manufacturers' test kits.
- Because **IMMUNO-TROL** reagents are formulated with human sera and non human derived constituents, incompatibilities between this reagent and test kits may exist.
- **IMMUNO-TROL** reagents are provide for quality assurance and must not be used for calibration.
- It is the responsibility of each laboratory to implement its own quality assurance program, to determine the suitability of **IMMUNO-TROL** reagents for its particular use and to establish guidelines for interpretation of results obtained with **IMMUNO-TROL** reagents.
- Adverse storage conditions or use of outdated reagents mau produce erroneous results.
- **IMMUNO-TROL** reagents have been tested and showed a reactivity in test procedure of the following bmd kits: DOT, ELISA and FIDIS.

#### SYMBOLS USED

	Biological risk		Temperature limitation
	Read instructions for use	<b>REF</b>	Catalog Number
<b>IVD</b>	In Vitro Diagnostic Use	<b>LOT</b>	Lot Number
<b>CE</b>	EC Declaration of Conformity		Use by
<b>CONTROL +</b>		Positive control	

**BioMédical Diagnostics SA**

**IVD**

**CE**

**Office**  
 Actipole 25 - 4 bld de Beaubourg  
 77435 Marne la Vallée Cx2  
 France

Tél : 33 1 64 62 10 12  
 Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
 Internet : [www.bmd-net.com](http://www.bmd-net.com)