ENA-DOT 7

REF HM 042

DÉFINITION

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BIBLIOGRAPHIE

CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance.

Clin. And Exp. Rheumatol. 1997 ; 15, 457-465

ERMENS A. et al. Simple Dot-Blot method evaluated for detection of autoantibodies against extractable antigens. *Clin. Chemistry* 1997 ; 43, 12, 2420-2422.

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MEWIS A. et al.

Evaluation of a DOT-Blot method for identification of antibodies against extractable nuclear antigens and anti-cytoplasmic antibodies. *Clin Chem., 1999, 45 (2), 311-312.*

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SPENCER GREEN G. et al.

Test performance in systemic sclerosis : anti-centromère and anti-Scl70 antibodies. *Am. J. of Med., 1997, 18, 291-300.*

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REF **HM 042**

DEFINITION

The **ENA-DOT 7** kit is a qualitative enzyme immunoassay (EIA) for the simultaneous detection of 7 antibody specificities:

- ↪ against Extractable Nuclear Antigens : SS-A, SS-B, Sm, Sm/RNP and Scl-70,
- ↪ and against 2 additional antigens : centromere and Jo-1.

DIAGNOSTIC VALUE

The detection and characterisation of anti-ENA (Extractable Nuclear Antigens) antibodies are an important element for the clinician in the diagnosis and classification of sub-groups of connective tissue diseases.

Many elaborate serological studies have demonstrated significant associations between these autoantibodies and the various connective diseases:

Antibodies associated with different autoimmune pathologies

- **Anti-SS-A and anti-SS-B antibodies**

Either one or both of these are observed, in the same pathological conditions Sjögren's syndrome and Systemic Lupus Erythematosus (SLE) (6). Anti-SSA are also found in mothers who are carrying a baby with neonatal lupus syndrome including heart block.

- **Anti-RNP antibody**

Found primarily in Mixed Connective Tissue Disease or Sharp's Syndrome, for which they are the biological marker. They are also observed in SLE.

Specific antibodies for a single connective tissue disease

- **Anti-Sm antibodies** characterize severe forms of SLE.
- **Anti-Scl70 antibodies** are specific for diffuse Scleroderma.
- **Anti-Jo1 antibodies** are found exclusively in sera from patients suffering from Polymyositis.
- **Anti-Centromere antibodies** are most often found in sera from patients having Scleroderma with limited skin involvement **CREST**, called **CREST syndrome**.

SAMPLES COLLECTION AND HANDLING


- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific bindings, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

ASSAY PRINCIPLE

The antigens are coated on distinct membranes on the assay strip.

- First, the assay strip is incubated in a reaction tube containing the diluted sample of the patient. If this sample contains at least one of the antibodies, these antibodies will recognize and bind to the corresponding antigen. After incubation, a first wash step removes all of the unbound proteins.
- An alkaline phosphatase labelled protein A conjugate binds to the captured antibodies. The excess of unbound conjugate is removed by a second wash step.
- The bound conjugate is visualized with an insoluble substrate for the phosphatase enzyme (BCIP/NBT). This will give rise to a blue colored circle, confirming the presence of autoantibodies in the sample tested.

REAGENTS

		HM042
Assay strips	 <p>Coated antigens: purified Ro/SS-A protein (60kD), purified La/SS-B protein (47kD), purified Sm protein (mix of several polypeptides of differing molecular weights), purified Sm/RNP protein (mix of several polypeptides of differing molecular weights), purified Scl-70 protein (70kD), purified Jo-1 protein (58kD) and recombinant CENP-B (80kD).</p> <p>STRIP</p>	25
1 vial of PAL conjugate (alkaline phosphatase labelled protein A) To be diluted	CONJ IgG	650 µL
1 vial of Substrate (BCIP/NBT) Ready to use	SUBS BCIP-NBT	35 mL
1 vial of Phosphate-Tween buffer (concentrate 10x) To be diluted in distilled water	BUF WASH 10x	30 mL
Reaction tubes (1,5mL)		100
Self-stick labels for returned results		25

MATERIAL REQUIRED BUT NOT PROVIDED

- Pipettes
- Timer
- Distilled water (to dilute the concentrated Phosphate-Tween buffer)
- Adsorbent towels (to blot dry the assay strip)
- Support for 96 reaction tubes (bmd, cat. HM501)

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

PRECAUTIONS

Be sure that the strips are dry after each wash.

If all specificities are found with the same patient, it is advisable to control this sample with another method.

Reagents in solution (except for substrate buffer) contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

ENA-DOT 7 has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labelling of dangerous preparations.

ENA-DOT 7 has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

SETUP

All reagents must be at room temperature before use.

1. Preparation of the wash and dilution buffer

- Dilute 10 times the concentrated Phosphate-Tween buffer in distilled water.
- Storage period: 3 months at +2°C/+8°C

2. Preparation of samples

- Mark one strip per sample.
- Mark 3 reaction tubes per sample and number from 1 to 3:
 - ≈ tube 1: sample incubation.
 - ≈ tube 2: conjugate incubation.
 - ≈ tube 3: substrate incubation.
- Dispense 1.2mL of wash and dilution buffer in tubes 1 and 2.
- Mark another microtube reserved for the two first washing steps and dispenses 1.2mL of wash and dilution buffer.

Remarks and precautions:

- ❖ The substrate is light sensitive. It requires to be dispense immediately at the time of his use in the tubes 3.
- ❖ Certain environmental conditions as the temperature may influence the final result. For temperatures > +30°C, it is recommended to be vigilant on the coherence of the results obtained.
- ❖ Certain operational conditions may influence the intensity of the blue color and should therefore be avoided:
 - the agitation time in each wash step is less than a minute
 - the membrane sticking to the wall of the vessel during incubation steps.
 - incomplete mixing of sample and conjugate

INTERPRETATION OF RESULTS

1. Validation criterion for the manipulation:

- The positive control should give a blue colored circle with a well defined outline.
- If 4 or more specificities are found with the same patient (example: SS-A – SS-B - Sm/RNP – Scl-70), it is recommended to control this sample with another method.

2. Reading of the strips:

In order to interpret correctly make sure that the strips are completely dry.

For each specificity, interpret according to this scheme:

APPEARANCE OF THE WELL	RESULTS
no circle	absence of auto-antibodies NEGATIVE
presence of a circle of which the staining is variable from bright to dark	presence of auto-antibodies POSITIVE

Each membrane showing a difficult to distinguish circle should be controlled with a second blood sample, drawn a few weeks later. The interpretation of the results should be done in relation to other examinations and in the clinical context.

QUALITY CONTROL

bmd offers a multiparametric quality controls (Immunotrol IV, Cat. No HM051, for determination of SS-A, SS-B, Sm and Sm-RNP and Immunotrol V, Cat. No HM052, for determination of ribosome, centromere, Jo-1, Scl-70 and mitochondrial antigen). They contains antibodies directed against one of the 7 specificities; and must be handled in the same way as samples.

CHARACTERISTICS AND PERFORMANCE OF THE TEST

1. INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH AN OTHER METHOD

ENA-DOT 7 has been compared:

- with ENA-LISA (cat. N° HM008) commercialized by bmd for quantitative determination by ELISA method of: SS-A, SS-B, Sm, Sm/RNP, Scl-70 and Jo-1 specificities.
- with HEp-2000™ (cat. N° IC 2040 RO), commercialized by bmd for centromere, nuclear antigen specificity by indirect immunofluorescence.

The study was performed on 211 samples:

- ⇒ 120 samples with at least one of the evaluated specificities
- ⇒ 91 samples with potential interferences factors: aspect (turbidity, globular, plasma) or biological markers (cryoglobulinemia, hypergammaglobulinemia, monoclonal immunoglobulins, circulating immune complexes)

	Relative sensitivity (%) *	Relative specificity (%) **	Concordance (%) ***
SS-A	97,9 (1/47)	99,4 (1/164)	99,0 (2/211)
SS-B	100 (0/14)	100 (0/197)	100 (0/211)
Sm	90,0 (1/9)	99,5 (1/202)	99,0 (2/211)
Sm/RNP	100 (0/27)	100 (0/184)	100 (0/211)
Scl-70	100 (0/13)	100 (0/198)	100 (0/211)
Jo-1	100 (0/13)	100 (0/198)	100 (0/211)
Centromere	97,3 (1/36)	99,4 (1/175)	99,0 (2/211)

*: discordant samples / total number of positive samples detected

** : discordant samples / total number of negative samples detected

***: discordant samples / total number of samples tested

2. ESTIMATION OF PREDICTIVE VALUE

Calculated on 1477 tests:

- Negative predictive value : 99,8% (1318/1321)
- Positive predictive value : 98,1% (159/162)

3. VERIFICATION OF POSITIVE THRESHOLD

Calculated on the 91 samples described above:

- These samples have been found by the three methods: ENA-DOT 7, ENA-LISA (threshold: 40 arbitrary units) or HEp-2000™ (threshold: 1/80).

4. REPRODUCIBILITY

Calculated for 3 specificities of ENA-DOT 7 kit, with different signals of positivity: positive, weakly positive or borderline and negative.

Specificity	Titre ENA-LISA Or HEp-2000™	intra-assay		between-assays	
		n	ENA-DOT 7 Appearance	n	ENA-DOT 7 Appearance
Anti-SS-A	74UA	5	dark circle	5	dark circle
	50UA		bright circle		bright circle
	4UA		no circle		no circle
Anti-Sm/RNP	90UA	5	dark circle	5	dark circle
	47UA		bright circle		bright circle
	3UA		no circle		no circle
Anti-centromere	>1/160	5	dark circle	5	dark circle
	1/40		bright circle		bright circle
	negative		no circle		no circle

Performance comparison of ENA-DOT 7, new version (short procedure) to ENA-DOT 7, old version (long procedure)

The study was performed on 183 samples:

- ⇒ 91 samples with at least one of the evaluated specificities,
- ⇒ 34 samples with potential interferences factors: aspect (turbidity, globular, plasma) or biological markers (cryoglobulinemia, hypergammaglobulinemia, monoclonal immunoglobulins, circulating immune complexes),
- ⇒ 58 samples of blood donors.



	Relative sensitivity (%) *	Relative specificity (%) **	Concordance (%) ***
SS-A	100 (0/30)	100 (0/153)	100 (0/183)
SS-B	100 (0/4)	100 (0/179)	100 (0/183)
Sm	100 (0/9)	99,4 (1/174)	99,5 (1/183)
Sm/RNP	100 (0/31)	100 (0/152)	100 (0/183)
Scl-70	100 (0/7)	99,4 (1/176)	99,5 (1/183)
Jo-1	100 (0/10)	100 (0/173)	100 (0/183)
Centromere	100 (0/18)	100 (0/165)	100 (0/183)

*: discordant samples / total number of positive samples detected

** : discordant samples / total number of negative samples detected

***: discordant samples / total number of samples tested

SYMBOLS USED

	Catalog number
	Batch code
	Number of tests
	Rapid test
	In vitro diagnostic device
	Use by
	Storage temperature limitation
	Read instructions for use
	EC Declaration of Conformity
	Reconstitute with water
	Contains sodium azide
	Contains Proclin

BIBLIOGRAPHY

CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance. *Clin. and Exp. Rheumatol.* 1997 ;15, 457-465

ERMENS A. et al. Simple Dot-Blot method evaluated for detection of autoantibodies against extractable antigens. *Clin. Chemistry* 1997 ; 43, 12, 2420-2422.

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MEWIS A. et al.
Evaluation of a DOT-Blot method for identification of antibodies against extractable nuclear antigens and anti-cytoplasmic antibodies. *Clin Chem., 1999, 45 (2), 311-312.*

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SPENCER GREEN G. et al.
Test performance in systemic sclerosis: anti-centromère and anti-Scl70 antibodies. *Am. J. of Med., 1997, 18, 291-300.*

SUMMARY OF METHOD

	OPERATIONS STEPS	INCUBATION TIME
INCUBATION OF SAMPLES	<ul style="list-style-type: none"> - Add 20µL sample to tube 1 containing wash and dilution buffer. - Put in the identified strip. - Start the alarm. - Move* the strip manually 20 times up and down (during 15s) - Then incubate <u>without shaking</u>. 	10 minutes at R.T
During sample incubation, add 20µL conjugate to tube 2 containing wash and dilution buffer		
WASH	<ul style="list-style-type: none"> - Put the strip in the corresponding wash tube with wash and dilution buffer. - Move* manually. - Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel. 	1 minute at R.T
INCUBATION WITH CONJUGATE	<ul style="list-style-type: none"> - Put the strip in tube 2. - Start the alarm. - Move* the strip manually 20 times up and down (during 15s) - Then incubate <u>without shaking</u>. 	10 minutes at R.T
At the end of the conjugate incubation, dispense 1.2mL of substrate in tube 3		
WASH	<ul style="list-style-type: none"> - Put the strip in the corresponding wash tube with wash and dilution buffer. - Move* manually. - Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel. 	1 minute at R.T
INCUBATION OF SUBSTRATE	<ul style="list-style-type: none"> - Put the strip in tube 3. - Start the alarm. - Move* the strip manually 20 times up and down (during 15s) - Then incubate <u>without shaking</u>. 	5 minutes at R.T
WASH	<ul style="list-style-type: none"> - Wash the strip under tap water. 	
DRY	<ul style="list-style-type: none"> - Remove water by blotting the strip with an adsorbent towel and let dry completely before reading the result. 	about 30 minutes at R.T

R.T : Room Temperature between +18°C and + 25°C. * : Up and down agitations should be slow.

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com