# CYTO-DOT

**REF** HM 043

**Français**

## DÉFINITION

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**REF** HM 043

## DEFINITION

The **CYTO-DOT 4** kit is a qualitative enzyme immunoassay (EIA) for the simultaneous detection of 4 cytoplasmic specificities: ribosome, mitochondria (M2), Jo-1 and LKM-1.

## DIAGNOSTIC VALUE

- Mitochondrial antibodies type M2**  
 Among the 9 types of antibodies directed against different components of mitochondria (M1 to M9), the antibody of type M2 is the most frequent. It is a sensitive and specific marker for primary biliary cirrhosis (present in 90% of the cases). It is found rarely in certain chronic hepatic diseases.
- Anti-ribosomal antibodies**  
 These antibodies are found in association with anti-nuclear antibodies in 10 to 20% of lupus. *Their interest lays in their strong association with central neurological abnormalities in lupus.*
- Anti-Jo-1 antibodies**  
 These antibodies recognize the histidyl-tRNA synthetase, an enzyme responsible for complexing the amino acid histidine to its cognate transfer RNA (50kD). It is present, with a variable frequency, in serum of patients with polymyositis, dermatomyositis or overlap syndromes. These antibodies are characteristic particularly in the forms with diffuse interstitial lung fibrosis.
- Anti-LKM-1 antibodies**  
 Anti-LKM-1 antibodies react with cytochrome P450 II D6. They are markers of autoimmune hepatitis type II, for which they recognize a major epitope of 9 amino acids (262-270) on linear form. Other etiologic factors can induce the production of anti-LKM-1 antibodies: drug (halothane), virus (hepatitis C), and multifactorial autoimmune primitive disorder.

## MATERIAL REQUIRED BUT NOT PROVIDED


- Pipettes
- Timer
- Distilled water (to dilute the concentrated Phosphate-Tween buffer)
- Adsorbent towels (to blot dry the assay strip)
- Support for 96 reaction tubes (bmd, cat. HM501)

## ASSAY PRINCIPLE

The antigens are coated on distinct membranes on the assay strip.

- First, the assay strip is incubated in a reaction tube containing the diluted sample of the patient. If this sample contains the antibodies, these antibodies will recognize and bind to the corresponding antigen.
- After incubation, a first wash step removes all of the unbound proteins.
- An alkaline phosphatase labeled protein A conjugate binds to the captured antibodies. The excess of unbound conjugate is removed by a second wash step.
- The bound conjugate is visualised with an insoluble substrate for the phosphatase enzyme (BCIP/NBT). This will give rise to a blue colored circle, confirming the presence of auto-antibodies in the sample tested.

## REAGENTS

		HM043
Assay strips		10
	Antigens coated: ribosomal proteins (mix of several purified phosphoproteins of different molecular weights), mitochondrial M2 proteins (mix of several purified proteins of differing molecular weights), purified Jo1 protein (58kD) and recombinant LKM-1 protein (60kD). <b>STRIP</b>	
1 vial of PAL Conjugate (Phosphatase alkaline labeled protein A) To be diluted	CONJ   IgG	300µL
1 vial of Substrate (BCIP/NBT) Ready to use	SUBS   BCIP-NBT	12mL
1 vial of Phosphate-Tween buffer (concentrate 10x) To be diluted in distilled water	BUF   WASH   10x	30mL
Reaction tubes (1,5mL)		40
Self-stick labels for returned results		10

## SAMPLES COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

## PRECAUTIONS

Be sure that the strips are dry after each wash.

If all specificities are found with the same patient, it is advisable to control this sample with an other method.

Reagents in solution (except for substrate buffer) contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**CYTO-DOT 4** has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labelling of dangerous preparations.

**CYTO-DOT 4** has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

## STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## SETUP

All reagents must be at room temperature before use.

### 1. Preparation of the wash and dilution buffer

- Dilute 10 times the concentrated Phosphate-Tween buffer in distilled water.
- Storage period: 3 months at +2°C/+8°C

### 2. Preparation of samples

- Mark one strip per sample.
- Mark 3 reaction tubes per sample and number from 1 to 3:
  - ≠ tube 1: sample incubation.
  - ≠ tube 2: conjugate incubation.
  - ≠ tube 3: substrate incubation.
- Dispense 1mL of wash and dilution buffer in tubes 1 and 2.
- Mark another microtube reserved for the two first washing steps and dispenses 1mL of wash and dilution buffer.

## Remarks and precautions:

- ❖ The substrate is light sensitive. It requires to be dispense immediately at the time of his use in the tubes 3.
- ❖ Certain environmental conditions as the temperature may influence the final result. For temperatures > +30°C, it is recommended to be vigilant on the coherence of the results obtained.
- ❖ Certain operational conditions may influence the intensity of the blue color and should therefore be avoided:
  - the agitation time in each wash step is less than a minute
  - the membrane sticking to the wall of the vessel during incubation step.
  - incomplete mixing of sample and conjugate

## INTERPRETATION OF RESULTS

### 1. Validation criterion for the manipulation:

- The positive control should give a blue colored circle with a well defined outline.
- If all specificities are found with the same patient, it is recommended to control this sample with an other method.

### 2. Reading of the strips:

In order to interpret correctly make sure that the strips are completely dry.

For each specificity interpret this scheme:

APPEARANCE OF THE WELL	RESULTS
no circle	absence of auto-antibodies NEGATIVE
presence of a circle of which the staining is variable from bright to dark	presence of auto-antibodies POSITIVE

*Each membrane showing a difficult to distinguish circle should be controlled with a second blood sample, drawn a few weeks later. The interpretation of the results should be done in relation to other examinations and in the clinical context.*

## QUALITY CONTROL

bmd offers a multiparametric quality control (Immunotrol V, Cat. No: HM052) which contains antibodies directed against one of the 4 specificities; it must be handled in the same way as samples.

## CHARACTERISTICS AND PERFORMANCE OF THE TEST

### 1. INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH AN OTHER METHOD

**CYTO-DOT 4** has been compared :

- With RAT TRIPLE SUBSTRATE (Cat. ME1248) commercialized by bmd for screening and titering of anti-ribosomal, anti-mitochondrial subtype M2 and anti-LKM-1 autoantibodies by indirect immunofluorescence on rat kidney, stomach and liver substrate.
- With ENA-LISA (Cat. HM008) commercialized by bmd for semi-quantitative determination by ELISA method of anti-ENA antibodies and especially for Jo-1 specificity.

Study was performed on 187 samples :

- ⇒ 91 samples with at least one of the evaluated specificities,
- ⇒ 49 samples of blood donor,
- ⇒ 47 samples with potential interference factors : aspect (turbidity, globular lysis, plasma) or biological markers (cryoglobulinemia, hypergammaglobulinemia, monoclonal immunoglobulins, circulating immune complexes).

	Relative sensibility		Relative specificity		Concordance	
	(%)	*	(%)	**	(%)	***
<b>Ribosome</b>	100	(0/10)	100	(0/177)	100	(0/187)
<b>M2</b>	96.8	(2/61)	100	(0/124)	98.9	(2/187)
<b>Jo-1</b>	100	(0/11)	100	(0/176)	100	(0/187)
<b>LKM-1</b>	100	(0/9)	100	(0/178)	100	(0/187)

- \* : Discordant samples / total number of positive samples detected
- \*\* : Discordant samples / total number of negative samples detected
- \*\*\* : Discordant samples / total number of samples tested

## 2. ESTIMATION OF PREDICTIVE VALUE

Calculated on 748 tests :

- Negative predictive value : 99.7% (655/657)
- Positive predictive value : 100% (93/93)

## 3. VERIFICATION OF POSITIVE THRESHOLD

Calculated on the 145 samples :

- ⇒ 49 samples of blood donors,
- ⇒ 47 samples with potential interference factors : aspect (turbidity, globular lysis, plasma) or biological markers (cryoglobulinemia, hypergammaglobulinemia, monoclonal immunoglobulins, circulating immune complexes).

These samples have been found negative.

## 4. REPRODUCIBILITY

Calculated for 2 specificities of **CYTO-DOT 4** kit, with different positivity degrees : positive, weakly positive or borderline and negative.

Specificity	Triple substrate Titer	INTRA-ASSAY		INTER-ASSAYS	
		n	CYTO-DOT 4 appearance	n	CYTO-DOT 4 appearance
<b>Anti-ribosomes</b>	>1/320 1/40 négatif	5	Dark circle Bright circle No circle	5	Dark circle Bright circle No circle
<b>Anti-LKM-1</b>	>1/160 1/40 négatif	5	Dark circle Bright circle No circle	5	Dark circle Bright circle No circle

Precision of the test couldn't be evaluated by a coefficient of variation, because **CYTO-DOT 4** is not a quantitative method.

But all tests have been compared together; they have the same appearance. The reproducibility is excellent.

## 5. PERFORMANCE COMPARISON OF CYTO-DOT 4, NEW VERSION (SHORT PROCEDURE) TO CYTO-DOT 4, PREVIOUS VERSION (LONG PROCEDURE)



Study was performed on 193 samples :

- ⇒ 69 samples with at least one of the evaluated specificities,
- ⇒ 62 samples of blood donor,
- ⇒ 62 samples with potential interference factors : aspect (turbidity, globular lysis, plasma) or biological markers (cryoglobulinemia, hypergammaglobulinemia, monoclonal immunoglobulins, circulating immune complexes).

	Relative sensibility		Relative specificity		Concordance	
	(%)	*	(%)	**	(%)	***
<b>Ribosome</b>	100	(0/9)	100	(0/184)	100	(0/193)
<b>M2</b>	100	(0/42)	99.3	(1/151)	99.5	(1/193)
<b>Jo-1</b>	100	(0/10)	100	(0/59)	100	(0/193)
<b>LKM-1</b>	100	(0/8)	99.5	(1/185)	99.5	(1/193)

- \* : Discordant samples / total number of positive samples detected
- \*\* : Discordant samples / total number of negative samples detected
- \*\*\* : Discordant samples / total number of samples tested

## SYMBOLS USED

	Catalog number
	Batch code
	Number of tests
	Rapid test
	In vitro diagnostic device
	Use by
	Storage temperature limitation
	Read instructions for use
	EC Declaration of Conformity
	Reconstitute with water
	Contains sodium azide
	Contains Proclin

## BIBLIOGRAPHY

- ABUAF N. et al.  
Non-organ-specific autoantibodies associated with chronic C virus hepatitis. J. Hepatol., 1993, 18, 359-364.
- ARNET F.C. et al.  
Ribosomal P autoantibodies in systemic lupus erythematosus. Arthritis Rheum., 1996, 39, 1833-183.
- BASSENDINE M.F.  
Antimitochondrial and other autoantibodies in primary biliary cirrhosis. In: Krawitt EL, Wiesner RH and Nishioka M. Autoimmune liver disease, second edition. Elsevier Science B.V., 1998, Amsterdam. 287-304, 1998
- ERMENS A. et al.  
Simple dot-Blot method evaluated for detection of antibodies against extractable nuclear antigens. Clin Chem.; 43 (12), 1997, 2420-2422.

HOMBERG J.C. et al.  
Antibodies to endoplasmic reticulum and cytosol in auto-immune hepatitis. In: Krawitt EL, Wiesner RH and Nishioka M. Autoimmune liver disease, second edition. Elsevier Science B.V., Amsterdam. 1998, 257-271.

LEUNG P.S.C. et al.  
Mitochondrial autoantibodies. In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam. 494-500.

MADDISON P.D.  
Aminoacyl-tRNA histidyl (Jo-1) synthetase autoantibodies. In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam. 31-35.

MEWIS A. et al.  
Evaluation of a DOT-Blot method for identification of antibodies against extractable nuclear antigens and anti-cytoplasmic antibodies. Clin Chem., 1999, 45 (2), 311-312.

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YOSHIO T. et al.  
A close temporal relationship of liver disease to antiribosomal P0 protein antibodies and central nervous system disease in patients with systemic lupus erythematosus. J. Rheumatol., 1998, 25, 681-688.

## SUMMARY OF METHOD

	OPERATIONS STEPS	INCUBATION TIME
<b>INCUBATION OF SAMPLES</b>	<ul style="list-style-type: none"> <li>- Add 20µL sample to tube 1 containing wash and dilution buffer.</li> <li>- Put in the identified strip.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	10 minutes at R.T
During sample incubation, add 20µL conjugate to tube 2 containing wash and dilution buffer		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>INCUBATION OF CONJUGATE</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 2.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	10 minutes at R.T
At the end of the conjugate incubation, dispense 1mL of substrate in tube 3		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>INCUBATION OF SUBSTRATE</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 3.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	5 minutes at R.T
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Wash the strip under tap water.</li> </ul>	
<b>DRY</b>	<ul style="list-style-type: none"> <li>- Remove water by blotting the strip with an adsorbent towel and <b>let dry completely</b> before reading the result.</li> </ul>	about 30 minutes at R.T

R.T : Room Temperature between +18°C and + 25°C.

\* : Up and down agitations should be slow.

**BioMédical Diagnostics SA**

**Office**  
Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)