# CELIAC-DOT

**REF** HM 047 IgG 
**REF** HM 048 IgA 
**Français**

## DÉFINITION

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### BIBLIOGRAPHIE

**DIETERICH W et al**, Identification of tissue transglutaminase as the autoantigen of celiac disease. Nat Med 1997; 3: 797-801

**FOTOUAKI M et al**, Clinical application of immunological markers as monitoring tests in celiac disease. Dig Dis Sci. 1999 Oct;44(10):2133-8.

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English

## DEFINITION

The **CELIAC-DOT** kit is a qualitative enzyme immunoassay (EIA) based on membrane assay strip for autoantibodies detection.

- ↳ **CELIAC-DOT IgA** is designed for the simultaneous detection of human IgA isotype autoantibodies directed against Gliadin and Tissue Transglutaminase Enzyme.
- ↳ **CELIAC-DOT IgG** is designed for the detection of human IgG isotype autoantibodies directed against Gliadin and Tissue Transglutaminase Enzyme.

## DIAGNOSTIC VALUE

Coeliac disease is an absorption deficiency syndrome characterized by weight loss, abdominal pain, and dietary gluten intolerance. Many affected patients have antibodies directed against gliadin (anti-gliadin), endomysium/tissue transglutaminase (anti-tTG) and reticulin (anti-reticulin). Anti-reticulin antibodies may be considered to be less sensitive than anti-gliadin in the assessment of patients with celiac disease. The detection of these antibodies is valuable in identifying patients with celiac disease, and to assure appropriate treatment and follow-up

Anti-gliadin antibodies (IgA and IgG) appear to be less specific than anti-endomysium or anti-tTG antibodies. Nevertheless, they may offer a benefit when assessing children under two years of age since affected persons in this age group may develop IgA anti-gliadin before developing IgA anti-endomysium. Anti-gliadin IgG may also be important for use in patients with celiac disease who also suffer from a variety of IgA deficiencies. After several months of treatment with a gluten free diet, antibody titers frequently decrease; though antibody levels may increase again on discontinuation of the gluten free diet.

The combined anti-gliadin and anti-tTG assay offers adult and children patients an easy, rapid, sensitive and specific technique for the detection and the treatment follow-up of Coeliac disease.


- \* The simultaneous presence of IgA anti-gliadin and anti-tTG antibodies, with or without IgG antibodies, will enable to confirm the quasi evidence of a celiac disease diagnosis. It is also an indication that a duodenojejunal biopsy should be performed.
- \* With reference made to the patient clinical data, the presence of IgA and IgG anti-gliadin antibodies can in some cases confirm the diagnosis of celiac disease.
- \* IgG anti-tTG antibodies are highly reliable marker of Coeliac Disease in IgA deficient subjects.
- \* The presence of IgG anti-gliadin antibodies alone is not indicative of a celiac disease, nevertheless, it is recommended to take into account an eventual IgA deficiency.

## ASSAY PRINCIPLE

The antigens are coated on distinct membranes on the assay strip.

- First, the assay strip is incubated in a reaction tube containing the diluted sample of the patient. If this sample contains at least one of the antibodies, these antibodies will recognize and bind to the corresponding antigen. After incubation, a first wash step removes all of the unbound proteins.
- Second step, the assay strip is incubated in a reaction tube containing the enhancer. After incubation, a wash step removes all of the unbound enhancer.
- An alkaline phosphatase labeled to the conjugate binds to the captured antibodies. The excess of unbound conjugate is removed by a third wash step.
- The bound conjugate is visualized with an insoluble substrate for the phosphatase enzyme (BCIP/NBT). This will result into a blue colored spot, confirming the presence of autoantibodies in the sample tested.

## REAGENTS

	HM047	HM048
Assay strips		
	Antigens coated : Purified Gliadin and recombinant human tissue transglutaminase( tTG)	25
1 vial of IgA-PAL Conjugate (Phosphatase alkanin labeled IgA) To be diluted	CONJ   IgA	650 µL
1 vial of IgG-PAL Conjugate (Phosphatase alkanin labeled IgG) To be diluted	CONJ   IgG	650 µL
1 vial of Substrate (BCIP/NBT) Ready to use	SUBS   BCIP-NBT	18 mL
1 vial of Enhancer corresponding for each isotype detection Ready to use	SOLN   ENH	650 µL
1 vial of Phosphate-Tween buffer (concentrate 10x) To be diluted in distilled water	BUF   WASH   10x	30 mL
Reaction tubes (1,5mL)	125	125
Self-stick labels for returned results	25	25

## MATERIAL REQUIRED BUT NOT PROVIDED

- Pipettes
- Timer
- Distilled water (to dilute the concentrated Phosphate-Tween buffer)
- Adsorbent towels (to blot dry the assay strip)
- Support for 96 microtubes (bmd, cat. HM501)

## SAMPLES COLLECTION AND HANDLING

- The tests are performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

## STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## PRECAUTIONS

Be sure that the strips are dry after each wash.

Reagents in solution (except for substrate buffer) contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**CELIAC-DOT** kits have been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

**CELIAC-DOT IgG and IgA** have been optimized for the use as describe in this procedure. Do not substitute other manufacture's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

## QUALITY CONTROL

bmd offers a multiparametric quality control (Immunotrol II, Cat. No: HM037) which contains antibodies directed against one of both specificities; it must be handled in the same way as samples.

## SETUP

All reagents must be at room temperature before use.

### 1. Preparation of the wash and dilution buffer

- Dilute 10 times the concentrated Phosphate-Tween buffer in distilled water.
- Storage period: 3 months at +2°C/+8°C

### 2. Preparation of samples

- Mark one strip per sample.
- Mark 4 reaction tubes per sample and number from 1 to 4:
  - ≠ tube 1: sample incubation.
  - ≠ tube 2: enhancer incubation.
  - ≠ tube 3: conjugate incubation.
  - ≠ tube 4: substrate incubation
- Dispense 600µL of wash and dilution buffer in tubes 1, 2 and 3.
- Mark an other microtube reserved for the three first washing steps and dispenses 600µL of wash and dilution buffer.

### Remarks and precautions:

- ❖ *The substrate is light sensitive. It requires to be dispense immediately at the time of his use in the tubes 4.*
- ❖ *Certain environmental conditions as the temperature may influence the final result. For temperatures > +30°C, it is recommended to be vigilant on the coherence of the results obtained.*
- ❖ *Certain operational conditions may influence the intensity of the blue color and should therefore be avoided:*
  - *the agitation time in each wash step is less than a minute*
  - *the membrane sticking to the wall of the vessel during incubation steps.*
  - *incomplete mixing of sample and conjugate*

## INTERPRETATION OF RESULTS

### 1. Validation criterion for the manipulation:

- The positive control should give a blue colored circle with a well defined outline.

### 2. Reading of the strips:

In order to interpret correctly make sure that the strips are completely dry.

For each specificity interpret this scheme:

APPEARANCE OF THE WELL	RESULTS
no circle	absence of auto-antibodies NEGATIVE
presence of a circle of which the staining is variable from bright to dark	presence of auto-antibodies POSITIVE

According to the CELIAC DOT evaluation of performances, samples giving a negative dot pattern correspond to titles lower than 20AU/mL with bmd GLIA-LISA and FIDIS™ CELIAC.

*Each membrane showing a difficult to distinguish circle should be controlled with a second blood sample, drawn a few weeks later. The interpretation of the results should be done in relation to other examinations and in the clinical context.*

**CHARACTERISTICS AND PERFORMANCE OF THE TEST**

Specificities	Glia	tTG	Glia	tTG
	IgG	IgG	IgA	IgA
Predicate test	ELISA tests			
Population evaluated	* 72 positive samples for one or more parameters related to Coeliac Disease. * 104 samples of healthy individuals and patients suffering of other diseases.		* 71 positive samples for one or more parameters related to Coeliac Disease. * 104 samples of healthy individuals and patients suffering of other diseases.	
	All the discordant results between the two methods are been controlled by FIDIS™ CELIAC (ref. MX009 / MX010) commercialised by bmd.			

**IgG Isotype**

Anti-Gliadin IgG			
176		ELISA	
		Positive	Negative
CELIAC-DOT	Positive	57	15
	Negative	3	101

Anti-tTG IgG			
176		ELISA	
		Positive	Negative
CELIAC-DOT	Positive	29	9
	Negative	1	137

**Anti-Gliadin IgG**

**Anti-tTG IgG**

*Indication of the diagnostic value*

Relative sensitivity: 95.0% (57/60)      Relative sensitivity: 96.7% (29/30)  
 Relative specific affinity: 87.1% (101/116)      Relative specific affinity: 93.8% (137/146)  
 Concordance: 89.8% (158/176)      Concordance: 94.3% (166/176)

*Controlled by FIDIS™ CELIAC IgG*

<p>↪ 15 samples positive on the DOT and negative by ELISA, 12 are positive with FIDIS™ CELIAC.</p> <p>↪ 3 samples negative on the DOT and positive by ELISA, 2 are negative with FIDIS™ CELIAC.</p>	<p>↪ 9 samples positive on the DOT and negative by ELISA, 3 are positive with FIDIS™ CELIAC.</p> <p>↪ 1 sample negative on the DOT and positive by ELISA, and positive with FIDIS™ CELIAC.</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Isotype IgA**

Anti-Gliadin IgA			
175		ELISA	
		Positive	Negative
CELIAC-DOT	Positive	58	16
	Negative	1	100

Anti-tTG IgA			
175		ELISA	
		Positive	Negative
CELIAC-DOT	Positive	57	3
	Negative	0	115

**Anti- Gliadin IgA**

**Anti-tTG IgA**

*Indication of the diagnostic value*

Relative sensitivity: 98.3% (58/59)      Relative sensitivity: 100% (57/57)  
 Relative specific affinity: 86.2% (100/116)      Relative specific affinity: 97.5% (115/118)  
 Concordance: 90.3% (158/175)      Concordance: 98.3% (172/175)

*Controlled by FIDIS™ CELIAC IgG:*

<p>↪ 16 samples positive on the DOT and negative by ELISA, 5 are positive with FIDIS™ CELIAC.</p> <p>↪ 1 sample negative on the DOT and positive by ELISA, and negative with FIDIS™ CELIAC.</p>	<p>↪ 3 samples positive on the DOT and negative by ELISA, all are positive with FIDIS™ CELIAC.</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------

**SYMBOLS USED**

	Catalog number
	Batch code
	Number of tests
	Rapid test
	In vitro diagnostic device
	Use by
	Storage temperature limitation
	Read instructions for use
	EC Declaration of Conformity
	Reconstitute with water
	Contains sodium azide
	Contains Proclin

**BIBLIOGRAPHY**

**DIETERICH W et al**, Identification of tissue transglutaminase as the autoantigen of celiac disease. Nat Med 1997; 3: 797-801

**FOTOULAKI M et al**, Clinical application of immunological markers as monitoring tests in celiac disease. Dig Dis Sci. 1999 Oct;44(10):2133-8.

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**SARDY M et al**, Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy. Clin Chem. 1999 Dec;45 (12):2142-9.

**VITORIA JC. et al**, Antibodies to gliadin, endomysium, and tissue transglutaminase for the diagnosis of celiac disease. J Pediatr Gastroenterol Nutr. 1999 Nov;29(5):571-4.

**WALTER-SMITH JA et al**, Working of European society of pediatric gastroenterology and nutrition (ESPGAN). Revised criteria for diagnosis of celiac disease Arch. Dis of Child 1990; 65: 909

## SUMMARY OF METHOD

	OPERATIONS STEPS	INCUBATION TIME
<b>INCUBATION OF SAMPLES</b>	<ul style="list-style-type: none"> <li>- Add 20µL sample to tube 1 containing wash and dilution buffer.</li> <li>- Put in the identified strip.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	20 minutes at R.T
During sample incubation, add 20µL enhancer to tube 2 containing wash and dilution buffer		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>ENHANCER</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 2.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	15 minutes at R.T
During sample incubation, add 20µL conjugate to tube 3 containing wash and dilution buffer		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>INCUBATION OF CONJUGATE</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 3.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	15 minutes at R.T
At the end of the conjugate incubation, dispense 600µL of substrate in tube 4		
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Put the strip in the corresponding wash tube with wash and dilution buffer.</li> <li>- Move* manually.</li> <li>- Take the strip out of the tube and remove residual buffer by blotting the strip with an adsorbent towel.</li> </ul>	1 minute at R.T
<b>INCUBATION OF SUBSTRATE</b>	<ul style="list-style-type: none"> <li>- Put the strip in tube 4.</li> <li>- Start the alarm.</li> <li>- Move* the strip manually 20 times up and down (during 15s)</li> <li>- Then incubate <u>without shaking</u>.</li> </ul>	5 minutes at R.T
<b>WASH</b>	<ul style="list-style-type: none"> <li>- Wash the strip under tap water.</li> </ul>	
<b>DRY</b>	<ul style="list-style-type: none"> <li>- Remove water by blotting the strip with an adsorbent towel and <b>let dry completely</b> before reading the result.</li> </ul>	about 30 minutes at R.T

R.T : Room Temperature between +18°C and + 25°C.

\* : Up and down agitations should be slow.

**BioMédical Diagnostics SA**

**Office**  
Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)