REF	HM 051	2 x 0,5ml
REF	HM 151	100 x 0,5ml

IMMUNO-TROL IV

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« dsDNA, SS-A 60, SS-A 52, SS-B, Sm, Sm/RNP »

INTENDED USE

This reagent, prepared from human serum, is intended for use as an unassayed precision control with **DOT, ELISA and FIDIS** procedures for determination of the following autoantibodies:

	Codes
• Anti-double stranded deoxyribonucleic acid	DNA
• Anti-Sjögren's Syndrom antigen A - 60 kDa.	SS-A 60 kDa
• Anti-Sjögren's Syndrom antigen A - 52 kDa	SS-A 52 kDa
• Anti-Sjögren's Syndrom antigen B	SS-B
• Anti-Smith antigen	Sm
• Anti-smith and ribonucleoproteins antigens	Sm/RNP

DESCRIPTION

IMMUNO-TROL reagents have been designed for use with **DOT, ELISA and FIDIS** methods for monitoring assay performance and maintaining quality of assays.

Each **IMMUNO-TROL** reagent lot is designed to produce a positive reaction for each autoantibodies above listed.


IMMUNO-TROL reagents do not have assigned values.

Each laboratory should establish its own target value for each marker and for each new lot on a minimum of 30 assays.

PACKAGING

HM 051: 2 vials of 0,5mL

HM 151: 100 vials of 0,5mL

The antigenic specificities are confirmed on the vial label. The antibody codes are specified above.

WARNING

- Human sources for the preparation of this control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus or hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.
- The reagent contains 0.02% ProClin® or less of 0.1% (w/v) NaN₃. Do not eat and avoid contact with skin and eyes.

PROCEDURE

IMMUNO-TROL reagents should be analyzed in the same manner as patient samples, according to the instructions supplied by manufacturers of the test kits to be used.


STABILITY AND STORAGE

- Store the reagents at +2°C to +8°C when not in use.
- Do not use reagents beyond the expiration date.
- Open vial stability : 60 days at +2°C to +8°C.
- Do not freeze.

LIMITS

- Performances and levels of reactivity of the **IMMUNO-TROL** reagents may vary with different manufacturers' test and assay procedures.
- **IMMUNO-TROL** reagents must not be substituted for the positive and the negative control reagents provided with manufacturers' test kits.
- Because **IMMUNO-TROL** reagents are formulated with human sera and non human derived constituents, incompatibilities between this reagent and test kits may exist.
- **IMMUNO-TROL** reagents are provide for quality assurance and must not be used for calibration.
- It is the responsibility of each laboratory to implement its own quality assurance program, to determine the suitability of **IMMUNO-TROL** reagents for its particular use and to establish guidelines for interpretation of results obtained with **IMMUNO-TROL** reagents.
- Adverse storage conditions or use of outdated reagents may produce erroneous results.
- **IMMUNO-TROL** reagents have been tested and showed a reactivity in test procedure of the following bmd kits: DOT, ELISA and FIDIS.

SYMBOLS USED

	Biological risk		Temperature limitation
	Read instructions for use		Catalog Number
	In Vitro Diagnostic Use		Lot Number
	EC Declaration of Conformity		Use by
		Positive control	