REF	HM 052	<i>2 x 0,5ml</i>
REF	HM 152	<i>100 x 0,5ml</i>

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Multiparameters quality control

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INTENDED USE

This reagent, prepared from human serum, is intended for use as an unassayed precision control with **DOT, ELISA, FIDIS** procedures for determination of the following autoantibodies:

- Anti-ribosomal antigen
- Anti-centromere antigen
- Anti-histidyl tRNA synthetase antigen
- Anti-scleroderma antigen - 70 kDA
- Anti-mitochondrial antigen

Codes
Ribo
CENP-B
Jo-1
Scl-70
Mito

DESCRIPTION

IMMUNO-TROL reagents have been designed for use with **DOT, ELISA, FIDIS** methods for monitoring assay performance and maintaining quality of assays.


Each **IMMUNO-TROL** reagent lot is designed to produce a positive reaction for each autoantibodies above listed.

IMMUNO-TROL reagents do not have assigned values.

Each laboratory should establish its own target value for each marker and for each nex lot on a minimum of 30 assays.

PACKAGING

HM 052: 2 vials of 0,5mL
HM 152: 100 vials of 0,5mL

The antigenic specificities are confirmed on the vial label. The antibody codes are specified above.

WARNING

- Human sources for the preparation of this control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus or hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepate B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.
- The reagent contains 0.02% ProClin® or less of 0.1% (w/v) NaN₃. Do not eat and avoid contact with skin and eyes.

PROCEDURE

IMMUNO-TROL reagents should be analyzed in the same manner as patient samples, according to the instructions supplied by manufacturers of the test kits to be used.

STABILITY AND STORAGE

- Store the reagents at +2°C to +8°C when not in use.
- Do not use reagents beyond the expiration date.
- Open vial stability : 60 days at +2°C to +8°C.
- Do not freeze.

LIMITS

- Performaces and levels of reactivity of the **IMMUNO-TROL** reagents may vary with different manufacturers' test and assay procedures.
- **IMMUNO-TROL** reagents must not be substituted for the positive and the negative control reagents provided with manufacturers' test kits.
- Because **IMMUNO-TROL** reagents are formulated with human sera and non human derived constituents, incompatibilities between this reagent and test kits may exist.
- **IMMUNO-TROL** reagents are provide for quality assurance and must not be used for calibration.
- It is the responsibility of each laboratory to implement its own quality assurance program, to determine the suitability of **IMMUNO-TROL** reagents for its particular use and to establish guidelines for interpretation of results obtained with **IMMUNO-TROL** reagents.
- Adverse storage conditions or use of outdated reagents may produce erroneous results.
- **IMMUNO-TROL** reagents have been tested and showed a reactivity in test procedure of the following bmd kits: DOT, ELISA and FIDIS.

SYMBOLS USED

	Biological risk		Temperature limitation
	Read instructions for use	REF	Catalog Number
IVD	In Vitro Diagnostic Use	LOT	Lot Number
CE	EC Declaration of Conformity		Use by
CONTROL +	Positive control		