ANA HEp-2

REF
HME 0060 Coffret
60
HME 0120 Coffret
120
HME 1600 Coffret
1600
HME 1006 10 Lames HEp-2 (6 puits)
HME 1012
10 Lames HEp-2 (12 puits)
HME 1016
10 Lames HEp-2 (16 puits)

DEFINITION

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BIBLIOGRAPHIE

- Tan, E., Antinuclear Antibodies: Diagnostic Markers for Autoimmune Diseases and Probes for Cell Biology, *Adv. Immuno.*, 44:93-151, 1989.
- Nakamura, R., C. Peebles, R. Rubin, et al, Autoantibodies to Nuclear Antigens, 2nd éd., ASCP, Chicago, 1985.
- Tan, E., E. Chan, K. Sullivan, et al, Antinuclear Antibodies (ANA's): Diagnostically Specific Immune Markers and Clues Toward the Understanding of Systemic Autoimmunity, *Clin. Immunol. Immunopathol.*, 47:121-141, 1988.
- Barnett, E.V., Antinuclear Antibodies and Nuclear Antigens, *California Medicine*, 104:463-469, 1966.
- Casais, S.P., G.J. Friou, L.L. Meyers, Significance of Antibody to DNA in Systemic Lupus Erythematosus, *Arth. Rheum.*, 7:379-390, 1964.
- Douvas, A.S., M. Achten, E.M. Tan, Identification of a Nuclear Protein (Scl-70) as a Unique Target of Human Antinuclear Antibodies in Scleroderma, *J. Biol. Chem.*, 244:10514-10522, 1979.






- Moroi, Y., C. Peebles, M.J. Fritzler, et al, Autoantibody to Centromere (Kinetochore) in Scleroderma Sera, *Proc. Natl. Acad. Sci.*, 77:1627-1631, 1980.
- Sharp, G.C., W.S. Irwin, E.M. Tan, et al, Mixed Connective Tissue Disease – An Apparently Distinct Rheumatic Disease Syndrome Associated with a Specific Antibody to an Extractable Nuclear Antigen (ENA), *Am. J. Med.*, 52:148-159, 1972.
- Alsaugh, M.A., W. Talal, E.M. Tan, Differentiation and Characterization of Autoantibodies and their Antigens in Sjogrens Syndrome, *Arth. Rheum.*, 19:216-222, 1976.
- Wolfe, J.F., E. Adelstein, G.C. Sharp, Antinuclear Antibody with Distinct Specificity for Polymyositis, *J. Clin. Invest.*, 59:176-178, 1977.
- Nakamura, R.M., CL. Peebles, D.P. Molden, E.M. Tan, Advances in Laboratory Tests for Autoantibodies to Nuclear Antigens in Systemic Rheumatic Diseases, *Lab. Med.*, 15:190-198, 1984.
- Moore, A.E., L. Sabachewsky, H.W. Toolan, Culture Characteristics of Four Permanent Lines of Human Cells, *Cancer Res.*, 15:598-602, 1955.
- McCarty, G.A., J.R. Rice, Characterization and Comparison of Commercially Available Antinuclear Antibody Kits Using Single Pattern Index Sera, *J. Rheum.*, 7:339-347, 1980.
- McCarty, G.A., et al, A Unique Antinuclear Antibody Staining Only the Mitotic Spindle Apparatus, *New Engl. J. Med.*, 305:703, 1981.
- Miyachi, K., M.J. Fritzler, E.M. Tan, Autoantibody to a Nuclear Antigen in Proliferating Cells, *J. Immunol.*, 121:2228-2234, 1978.
- Schur, P.H., M. Moore, N. Rothfield, The Subclass of Antinuclear and Antinucleic Acid Antibodies, *Arth. Rheum.*, 15:174-180, 1972.
- Weller, T.H., A.H. Coons, Fluorescent Antibody Studies with Agents of Varicella and Herpès Zoster Propagated In Vitro, *Proc. Soc. Exp. Biol. Med.*, 86:789-794, 1954.
- Friou, F.J., Clinical Application of Lupus Serum-Nucleoprotein Reaction Using the Fluorescent Antibody Technique, *J. Clin. Invest.*, 36:890, 1957.
- Holman, H.R., H.G. Kunkel, Affinity between the Lupus Erythematosus Serum Factor and Cell Nuclei and Nucleoprotein, *Science*, 126:162-163, 1957.
- Lyerla, H.C., F.T. Forrester, The Immunofluorescence (IF) Test, in *Immunofluorescence Methods in Virology*, USDHHS, Georgia, 71-81, 1979.
- Peter, J.B., R.L. Dawkins, Evaluating Autoimmune Diseases, *Diag. Med.*, 1-10, Sept/Oct 1979.
- Notman, D.D., N. Kurata, E.M. Tan, Profiles of Antinuclear Antibodies in Systemic Rheumatic Diseases, *Annals Intern. Med.*, 83:464-469, 1975.
- Nakamura, R.M., CL. Peebles, R.L. Rubin, et al, Autoantibodies to Nuclear Antigens (ANA), ASCP, Chicago, 1985.
- Molden, D.P., ANA Profiles in Systemic Rheumatic Disease, *Diag. Med.*, 12-18, June 1985.
- Fritzler, M.J., D.W. Valencia, G.A. McCarty, Speckled Pattern Antinuclear Antibodies Resembling Anticentromere Antibodies, *Arth. Rheum.*, 27(1):92-96, 1984.
- Fritzler, M.J., T.D. Kinsella, E. Garbutt, The CREST Syndrome: A Distinct Sérologie Entity with Anticentromere Antibodies, *Am. J. Med.*, 69:520-526, 1980.
- Tan, E.M., G.P. Rodnan, I. Garcia, et al, Diversity of Antinuclear Antibodies in Progressive Systemic Sclerosis, *Arth. Rheum.*, 23:617-625, 1980.
- Penning, C, V. Steen, G. Reimer, et al, An Analysis of Systemic Sclerosis Patients with Autoantibodies to the Nucleolar Proteins PM-Scl, RNA Polymerase I and Fibrillarin, *A.R.A. Mtg.*, D29, 1987.
- Pinnas, J.L., J.D. Northway, E.M. Tan, Antinucleolar Antibodies in Human Sera, *J. Immunol.*, 111:996-1004, 1973.
- Ritchie, R.F., Antinucleolar Antibodies: Their Frequency and Diagnostic Association, *New Engl. J. Med.*, 282:1174-1178, 1970.
- Bianchi, F.B., M. Rizzetto, P. Penfold, et al, Ultrastructural Localization and Characterization of a Ribosomal Antibody Detected by Immunofluorescence in Systemic Lupus Erythematosus, *Clin. Exp. Immunol.*, 17:629-636, 1974.
- Rothfield, N.F., Detection of Antibodies to Nuclear Antigens by Immunofluorescence, *Manual of Clinical Immunology*, ASM, Washington D.C., 85:647-651, 1976.
- Pollack, V.E., Antinuclear Antibodies in Families of Patients with Systemic Lupus Erythematosus, *New Engl. J. Med.*, 271:165-171, 1964.
- Tan, E.M., M.J. Fritzler, P.A. Reyes-Lopez, Autoantibodies as Biological Markers of Neoplasia: Basic and Applied Aspects, Elsevier-North Holland, New York, 1978.
- Cleymaet, J.E., R.M. Nakamura, Indirect Immunofluorescent Antinuclear Antibody Tests: Comparison of Sensitivity and Specificity of Different Substrates, *Am. J. Clin.Path.*, 58:388-393, 19

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DEFINITION

The **ANA HEp-2 (Anti-Nuclear Antibody)** () test is an indirect fluorescent antibody assay utilizing HEp-2 tissue culture cells as a substrate for the qualitative and/or semi-quantitative determination of antinuclear antibodies in human serum. The **ANA HEp-2** () test is intended for use as an aid in the diagnosis of certain autoimmune diseases.

SUMMARY AND EXPLANATION

Antinuclear antibodies (ANA) are a group of auto-antibodies characterized by specificity for numerous antigenic determinants of cell nuclei. While the role of ANA's in the pathogenesis of autoimmune disease is controversial, they are quite useful as disease markers, primarily for diagnostic screening and also to monitor the course of connective tissue diseases. ^(1,2,3)

Because of the high correlation of positive antinuclear antibodies with SLE a negative ANA essentially rules out this disease. ⁽⁴⁾ Although antibodies specific to DNA have a high correlation with SLE, ⁽⁵⁾ antibodies to a number of other nuclear antigens appear to be of diagnostic and/or prognostic significance in diseases such as Progressive Systemic Sclerosis, ^(6,7) Mixed Connective Tissue Disease, ⁽⁸⁾ Sjogren's Syndrome, ⁽⁹⁾ and Polymyositis; ⁽¹⁰⁾ making ANA testing useful not only for SLE, but as a general screening tool for connective tissue diseases. ⁽¹¹⁾

Among the methodologies available to detect ANA's are EIA, ELISA, Dot Blot and the Indirect Fluorescent Antibody (IFA) technique. The antigen for the first three methods can either be a spectrum of clinically significant, specific auto-antigens, a single mixture of auto-antigens from a cell lysate, or a combination of the two. These methods are not as sensitive as IFA, nor can they detect the variety of autoantibodies. They also do not have the pattern recognition quality of the IFA. The IFA test is sensitive, screens for a wide variety of known and unknown autoantibodies and, through pattern recognition, offers insights into the probable identity of the antigen and associated autoimmune disorder. It is the dominant methodology in clinical laboratories at this time ⁽²⁾ and the method of choice for ANA screening and semi-quantitation.

The antigen of the **ANA HEp-2** () substrate is a human epithelial cell (HEp-2) line established by Moore, Sabachensky and Toolan. ⁽¹²⁾ HEp-2 cells have been shown to have greater sensitivity than tissue sections and yield sharper pattern recognition. ⁽¹³⁾ The presence of mitotic figures aids in differential pattern recognition as well as detecting previously unreported nuclear antibodies. ^(14,15) Antinuclear antibodies can be found in all major immunoglobulin classes (IgG, IgA or IgM), therefore, antihuman gammaglobulin conjugate that detects all classes is recommended for use in routine ANA testing. ⁽¹⁶⁾

REAGENTS

		REF HME 0060	REF HME 0120	REF HME 1600
SORB SLD n Stable in sealed foil pouch at 8° C, or lower, until labeled expiration date.	Individually foil wrapped slides with HEp-2 tissue culture cells fixed onto each well.	10 slides x 6 wells	10 slides x 12 wells	100 slides x 16 wells + coverslip
CONJ FITC IgG Stable at 2-8° C away from direct light until labeled expiration date.	Dropper vials containing fluorescein isothiocyanate labeled goat antihuman immunoglobulins with 0.001% Evans Blue counterstain, protein stabilizer, less than 0.1% sodium azide and 0.001% thimerosal added.	1 x 3,5mL	2 x 3,5mL	4 x 20mL
CONTROL + Stable at 2-8° C until labeled expiration date.	Vial of ANA positive human control serum (homogeneous pattern) with protein stabilizer and 0.005% thimerosal.	1 x 0,5mL	1 x 0,5mL	1 x 0,5mL
When used undiluted, as provided, specific fluorescent intensity of 3+ or greater should be seen. Optionnaly, the positive control can be tittered to endpoint. If tittered, the control should be serially diluted in PBS; When the control has been tested for the endpoint titer by bmd. Enterptises, Ltd. An endpoint titer is printed on the positive control vial. Due to variations within each laboratory (fluorescent microscope, etc.) each laboratory should establish its own mean titer for each lot of positive control (generally ± one dilution from stated endpoint).				
CONTROL - Stable at 2-8°C until labeled expiration date. The control is intended to be used undiluted as provided. The staining reaction should exhibit less than 1+ fluorescence.	Vial of ANA negative human control serum with protein stabilizer and 0.005% thimerosal.	1 x 0,5mL	1 x 0,5mL	
MM Stable at 2-8° C until labeled expiration date.	Vial of Mounting Medium (phosphate buffered glycerol of pH 7.4 ± 0.2)	1 x 3,5mL	1 x 3,5mL	
DIL SPE Stable at 2-8° C until the labeled expiration date provided no gross contamination is seen.	Bottle containing sample diluent, with less than 0.1% thimerosal, formulated to reduce nonspecific staining.	1 x 60mL	1 x 60mL	3 x 250mL
Do not use if the solution turns cloudy, or if a precipitate forms.				
BUF WASH RCNS 1L H₂O Stable in sealed packet at +25° C, or lower, until labeled expiration date.	One-liter packet of dry PBS for the reconstitution of 1 liter of wash buffer	2	2	

ASSAY PRINCIPLE

The ANA HEp-2 (bmd) test utilizes the indirect fluorescent antibody assay method first described by Weller and Coons,⁽¹⁷⁾ and applied to ANA detection by Friou,⁽¹⁸⁾ and Holman and Kunkel.⁽¹⁹⁾

The procedure is carried out in two basic reaction steps:

Step 1 - Human serum is reacted with the antigen substrate. Antibodies, if present, will bind to the antigen forming stable antigen-antibody complexes. If no antibodies are present, the complex will not be formed and serum components will be washed away.

Step 2 - Fluorescein labeled antihuman antibody is added to the reaction site which binds with the complexes formed in step one. This results in a positive reaction of bright apple-green fluorescence when viewed with a properly equipped fluorescence microscope. If no complexes are formed in step one, the fluorescein labeled antibody will be washed away, exhibiting a negative result.

MATERIAL REQUIRED BUT NOT PROVIDED

- One-liter volumetric flask or one liter graduated cylinder
- Distilled water - CAP Type one or equivalent
- One-liter screw capped container
- Disposable test tubes and rack
- Disposable serological pipettes
- Calibrated pipettes to deliver 10 µl and 100 µl, with disposable pipette tips
- Pasteur pipettes and bulbs
- Moist chambers
- Plastic squeeze wash bottle
- Coplin jars or staining dishes with slide racks
- Coverslips, No. 1, 24 x 60 mm
- Permanent black felt tip marking pen
- Fluorescence microscope equipped with a mercury or tungsten-halogen light source, a 390-490nm excitation filter and 515-520nm barrier filter, and optics to give a total magnification of 400X. The excitation wavelength of FITC is 490nm and the emission wavelength is 520nm.

REAGENTS AVAILABLE INDIVIDUALLY

HME 1006	HEp-2 slides (10 x 6 wells)
HME 1012	HEp-2 slides (10 x 12 wells)
HME 1016	HEp-2 slides (10 x 16 wells)
HME 9920	FITC conjugate (20 mL)
HME DIL250	IFA Diluent (250 mL)
HME MM	Mounting medium (3,5 mL)
HME PBS10	Dry PBS packet (10x1L)
HME NEG	ANA Negative control (0,5mL)
HME HOM	Homogeneous positive control (0,5mL)
HME NUC	Nucleolar positive control (0,5mL)
HME RNP	Speckled/RNP positive control (0,5mL)
HME ACA	Centromere positive control (0,5mL)
HME SSA	SS-A positive control (0,5mL)
HME SSB	SS-B positive control (0,5mL)
HME SCL	Scl-70 positive control (0,5mL)
HME 12CS	Coverslips 24x60 (1x100 pieces)
HME 16CS	Coverslips 24x70 (1x100 pieces)

STABILITY AND STORAGE

- All reagents and substrate slides are stable if they are stored between +2°C and +8°C in their own packaging.
- Do not use any kits or reagents beyond the stated expiration date.

PREPARATION OF REAGENTS

BUFFER PREPARATION

Place contents of a one-liter PBS packet into a one-liter volumetric flask, add *distilled water to the one-liter mark, mix and leave several hours or overnight to dissolve. Reconstituted buffer should have a pH of 7.4 ± 0.2.

Adjust with 1N NaOH or 1N HCl if pH value is outside the stated range. Store in a clean screw capped bottle at 25°C or lower. Stable during 1 month if no gross contamination is seen. Do not use if pH changes, if the solution turns cloudy, or if a precipitate forms.

*Use deionized water with caution, as pH of this type of water may vary causing the pH of PBS to become unstable upon prolonged storage.

PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Do not remove slides from pouch until ready for testing. Do not use if pouch has been punctured, as indicated by a flat pouch.
3. All reagents should be brought to room temperature (+20°C-+25°C) prior to use.
4. Abnormal test results may be seen if the antigen substrate slides are allowed to dry during the staining procedure.
5. Refrigeration (+2°C-+8°C) of kit immediately upon arrival will insure stability until labeled expiration date.
6. Reagents should not be used beyond stated expiration date.
7. Substitution of components other than those provided may yield inconsistent results.
8. Do not expose conjugate to strong light during storage or use.
9. Avoid microbial contamination of all reagents involved in the testing procedure or incorrect results may occur.
10. Deviation from the defined test procedure, such as incubation times or temperatures, may give erroneous results.
11. Lipemic, hemolyzed or contaminated sera may yield erroneous results.
12. Previously frozen specimens after thawing should be thoroughly mixed prior to testing.
13. Reusable glassware must be washed and thoroughly rinsed free of detergents.
14. Care should be taken to avoid splashing and generation of aerosols.
15. Patient samples, as well as ail materials coming into contact with them, should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the CDC/NIH manual "Biosafety in Microbiological and Biomedical Laboratories", 1984 Edition. Never pipette by mouth. Avoid contact with skin and mucous membranes.

16. Sera used to prepare positive and negative controls have been tested by an FDA approved method and found to be negative (were not repeatedly reactive) for the presence of Hepatitis B surface Antigen (HBsAg) and antibodies to Hepatitis C (HepCAb) and HIV 1 & 2. However, because no test method can offer complete assurance of the absence of these or other infectious agents, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the CDC/NIH manual "Biosafety in Microbiological and Biomedical Laboratories", 1984 Edition.
17. The preservatives used in conjugates and controls are toxic if ingested. Azides may react with copper or lead plumbing to form explosive metal azides. When disposing, flush drains with water to minimize buildup of azide and metal compounds.
18. **IFA diluent** () should be used ONLY as a diluent for patient specimens. Do NOT prepare serial dilutions for endpoint titer with the IFA diluent. Do NOT use in any of the wash steps.

TEST PROCEDURE

NOTE: Bring slides and reagents to room temperature (+20°C-+25° C) before use.

1. SPECIMEN PREPARATION

Screening:

Prepare a 1:40 dilution of each patient's serum by adding 0.01 ml (10 µl) of patient's serum to 0.39 ml of IFA Diluent.

NOTE: Remove only the amount of IFA Diluent needed to perform each test run to reduce the possibility of product contamination.

Semi-quantitation:

The significance of a positive screening result should be confirmed by repeating the test with dilutions of the serum. Each laboratory should establish its own titrating protocol; however, the following fourfold serial titration is suggested:

- a. Prepare a 1:40 dilution of each patient's serum by adding 0.01 ml (10µl) of patient's serum to 0.39 ml of IFA Diluent.
- b. Add 0.3 ml PBS to tubes #2, #3, #4 and #5. (Do NOT use IFA Diluent for serial dilutions.)
- c. Using a 100 µl pipette, transfer 0.1 ml (100 µl) from tube #1 to tube #2. Mix. Using a new tip for each dilution, transfer 0.1 ml (100 µl) from the second tube to the third, from the third tube to the fourth, and from the fourth tube to the fifth, mixing after each transfer. This will give a fourfold titration with the following dilutions:

-  Tube #1 = 1:40
-  Tube #2 = 1:160
-  Tube #3 = 1:640
-  Tube #4 = 1:2560
-  Tube #5 = 1:10240

NOTE: The Positive and Negative Controls are intended to be used undiluted. However, if performing a semi-quantitative test, the Positive Control should be diluted as suggested above. However, DO NOT use IFA Diluent to make the 1:40 or any other dilution of the Positive Control.

2. SLIDE PREPARATION

Remove as many slides as are required from the refrigerator and allow to equilibrate to room temperature (+20°C-+25°C) for at least five minutes. Remove slides from sealed foil pouches, being careful not to touch antigen surface. Identify each slide using a permanent black marking pen.

3. SPECIMEN APPLICATION

Using separate Pasteur pipettes, apply one drop (20-30µl) of each control, and one drop (20-30µl) of each patient serum dilution to individual wells of the slide. Do not touch the antigen surface with the pipette while dropping. Do not allow drops to mix, as cross contamination of samples between wells could cause erroneous results.

4. INCUBATION 1

Incubate in a moist chamber at room temperature (+20°C-+25°C) for 30 minutes.

NOTE: THE ANTIGEN MUST NOT BE ALLOWED TO DRY DURING ANY OF THE FOLLOWING STEPS. Nonspecific binding may occur if the reagent is allowed to dry on the slide.

5. RINSE 1

Remove slides from moist chamber one at a time and rinse GENTLY with PBS using a squeeze wash bottle. Do not focus the PBS stream directly onto the wells. To prevent cross contamination tilt slide first toward wells (1-6 or 1-8) and, running PBS stream along the midline of the slide, allow the PBS to run off the top edge of the slide. Then tilt the slide toward wells (7-12 or 9-16), and repeat this procedure, allowing the PBS to run off the bottom edge of the slide. For six well slides, tilt slide down and run the PBS stream across the slide above the wells, allowing the PBS to run off the bottom edge of the slide.

6. WASH 1

Place slides in Coplin jars or staining dishes and wash in two changes of PBS for 5 to 6 minutes each wash, agitating gently at entry and prior to removal.

7. CONJUGATE APPLICATION

Remove slides from the wash one at a time, shake off excess PBS, dry around outside edges if necessary and return each slide to the moist chamber. Using dropper bottle provided, apply one drop of conjugate (30-50µl) to each well of each slide, making sure each well is completely covered.

8. INCUBATION 2

Incubate in a moist chamber at room temperature (+20°C-+25°C) for 30 minutes. Protect slides from excessive light.

9. RINSE 2

Remove slides from moist chamber and rinse GENTLY with PBS using a squeeze wash bottle. As suggested in step 5, do not focus PBS stream directly onto the wells.

10. WASH 2

Place slides in Coplin jars or staining dishes and wash in two changes of PBS for 5 to 6 minutes each wash, agitating gently at entry and prior to removal.

11. COVERSLIP

Remove slides one at a time from last PBS wash, shake off excess PBS and immediately add two to four drops of mounting medium across the slide. Tilt slide and rest the edge of the coverslip against the bottom of the slide allowing the mounting medium to form a continuous bead between the coverslip and slide. Gently lower the coverslip from the bottom of the slide to the top, being careful to avoid air bubbles. Drain excess mounting medium by holding the edge of the slide against absorbent paper. Wipe off back of slide.

12. READ

Examine stained slides as soon as possible using a properly equipped fluorescence microscope. It is recommended that slides be examined on the same day they are stained. If any delay is anticipated, store slides in the refrigerator (+2°C-+8° C) away from direct light and read the following day. Do not allow mounting medium to dry between slide and coverslip. If drying should occur, add additional mounting medium or recoverslip slide.

FLUORESCENT INTENSITY GRADING

Fluorescent intensity may be semi-quantitated by following the guidelines established by the Centers for Disease Control, Atlanta, Georgia. ⁽²⁰⁾

4+ = Maximal fluorescence; brilliant yellow-green.

3+ = Less brilliant yellow-green fluorescence.

2+ = Definite but dull yellow-green fluorescence.

1+ = Very dim subdued fluorescence.

The degree of fluorescent intensity is not clinically relevant and has only limited value as an indicator of titer. Differences in fluorescence microscope optics, filters and light sources may result in differences of 1+ or more fluorescent intensity when observing the same slide using different microscopes.

QUALITY CONTROL

Specificity control

Both a positive and negative antibody control must be included with each run. These controls must be examined prior to reading test samples and should demonstrate the following results:

Negative Control

Using the ANA NEGATIVE CONTROL SERUM ([bmd](#)) as provided with the ANA HEp-2 ([bmd](#)) test, the cells should exhibit less than 1+ fluorescence and appear reddish-orange due to the counterstain.

Positive Control

Using the ANA POSITIVE CONTROL SERUM ([bmd](#)) as provided with the ANA HEp-2 ([bmd](#)) test, the cells should exhibit a homogeneous staining pattern with a fluorescent intensity of 3+ to 4+.

Each control must demonstrate the expected reaction in order to validate the test. If the controls fail to appear as described above, the test results should not be reported and the test should be repeated.

If upon repeat testing the controls still fail to show the proper reaction, do not report test results. The specificity of the antigen substrate can further be tested by running a panel of various types and patterns of antinuclear antibodies. (These are available separately from [bmd](#) enterprises, Ltd.)

Sensitivity control

A titrated control included with each run tests substrate sensitivity, as well as, checks technique, conjugate quality and the microscope optical System. The endpoint titer of each lot of ANA POSITIVE CONTROL SERUM ([bmd](#)) must be determined. There must not be more than a twofold difference (+/-) in titer from the stated endpoint. Each run should include the endpoint dilution, one fourfold dilution above and one fourfold dilution below the endpoint dilution. The more concentrated dilution should be positive and the less concentrated negative. If the control does not behave as described, the test results are invalid and the test should be repeated. If the control fails to show the proper reaction upon repeat testing, do not report the test results.

[bmd](#) can provide you ready to use controls (see § Reagent available individually) which can be tested at the same time.

RESULTS

NEGATIVE

A serum dilution is considered negative for antinuclear antibodies if the cells exhibit less than 1+ fluorescence and lack a clearly discernible pattern. Cells will appear reddish-orange due to the Evans Blue counterstain.

A sample is considered negative for antinuclear antibodies if it exhibits less than 1+ fluorescence at a serum dilution of 1:40 and all greater dilutions, or if the fluorescence observed is not a discernible ANA pattern.

... Negative samples may exhibit fluorescent staining slightly greater than the Negative Control, but less than 1+.

... Some sera may show a low degree of nuclear or cytoplasmic fluorescence with no clearly discernible staining pattern. This phenomenon is generally due to heterophile antibodies and should be reported as negative. ⁽²¹⁾

... Intense non-nuclear staining may be observed in some sera containing Anti-Mitochondrial, Anti-Smooth Muscle or other cytoplasmic antibodies.

POSITIVE

A serum dilution is considered positive for antinuclear antibodies if the fluorescent staining is at an intensity of 1+ or greater with a clearly discernible pattern of fluorescence. A sample is considered positive for antinuclear antibodies if it exhibits a characteristic ANA staining pattern with a fluorescent intensity of 1+ or greater at a serum dilution of 1:40 or greater. Multiple antinuclear antibodies may be present in a given specimen; one masking the other.

Serially diluting the specimen will aid in distinguishing these patterns.

Report all titers and patterns seen.

TITRATION

If a semi-quantitative titration is performed, the result should be reported as the reciprocal of the last dilution in which 1+ apple-green fluorescent intensity with a clearly discernible staining pattern is detected. When reading fourfold serial dilutions, endpoints can be extrapolated where necessary.

EXAMPLE OF ENDPOINT EXTRAPOLATION:

1:40 = 3+
1:160 = 2+
1:640 = +/-
1:2560 = Neg

The extrapolated endpoint is reported as 320.

Report all titers and patterns seen, extrapolating the titer where necessary.

EXAMPLE:

1:40 = 4+ Peripheral and 3+ Homogeneous
1:160 = 3+ Homogeneous
1:640 = 2+ Homogeneous and 3+ Speckled
1:2560 = 1+ Homogeneous and 2+ Speckled
1:10,240 = +/- Speckled

Report: 40 Peripheral, 2560 Homogeneous and 5120 Speckled.

INTERPRETATION OF RESULTS

Four major staining patterns which may occur singly or in combinations have been described:

1. Peripheral (shaggy, rim, membranous)
2. Homogeneous (diffuse, solid)
3. Speckled (including ACA)
4. Nucleolar

Other patterns less frequently seen include Spindle and Ribosomal RNP.

1. PERIPHERAL - smooth staining primarily around the outer region of the nucleus with weaker staining in the center. Not all cells within a well may appear peripheral; some may appear homogeneous. The chromosome region of mitotic cells will exhibit a bright positive staining pattern.

Note: A very thin distinct line around the nucleus is not a peripheral ANA, but a nuclear membrane antibody. It differs in that the chromosome region within the mitotic cells is negative.

NUCLEAR ANTIGENS: nDNA or Histones.

DISEASE ASSOCIATION: High titers to nDNA are suggestive of active SLE. Lower titers are suggestive of SLE or other connective tissue diseases. ⁽²²⁾

2. HOMOGENEOUS - diffuse staining of the entire nucleus, with or without apparent masking of the nucleoli. Some specimens may exhibit moderately large irregular areas of more intense staining; or the pattern may appear granular, especially as the antibody reaches its endpoint. The chromosome region of mitotic cells will exhibit a bright positive staining pattern.

NUCLEAR ANTIGENS: nDNA, DNP or Histones.

DISEASE ASSOCIATION: High titers are suggestive of SLE, while low titers may be found in SLE, and Rheumatoid Arthritis. ⁽⁵⁾

- a. Antibodies to DNP have the same specificity as the LE Cell Factor. ⁽¹¹⁾
- b. Antibodies to Histone alone have a high association with drug-induced lupus. ⁽¹¹⁾

3. SPECKLED - fluorescent aggregates throughout the nucleus which can be very fine to very coarse depending on the type of antibody present. More than one type of speckle may be seen in any one specimen. The chromosome region of mitotic cells is usually negative.

- a. Sm and nRNP antibodies usually present as a coarse speckle with the chromosome region of mitotic cells negative.
- b. SS-A/Ro and SS-B/La antibodies present as small uniform speckles in a uniform distribution with the chromosome region of mitotic cells negative.
- c. Fine dense speckles with positive staining of the nucleoli and chromosome region of mitotic cells may indicate an Scl-70 antibody.
- d. Varying fluorescence of fine to coarse speckles, in approximately 30-60% of cells with positive or negative staining of the mitotic cells, may represent a PCNA (Proliferating Cell Nuclear Antigen) antibody.

NUCLEAR ANTIGENS: Sm, nRNP, SS-A/Ro, SS-B/La, Scl-70 (DNA-Topoisomerase 1) or PCNA (Cyclin).

DISEASE ASSOCIATION:

- a. Sm antibodies are highly specific for SLE and appear to be a "marker" for this disease. ⁽²³⁾
- b. nRNP along with other types of ANA's have been found in SLE, Rheumatoid Arthritis (RA), and Progressive Systemic Sclerosis (Scleroderma). High levels of nRNP antibodies alone are characteristic of Mixed Connective Tissue Disease (MCTD). ⁽²³⁾

- c. SS-A/Ro and SS-B/La antibodies are frequently present in patients with Sjogren's Syndrome without associated RA. Both may be found less frequently in patients with SLE. SS-A/Ro antibodies are found in a high percentage of infants with congenital heart block, neonatal lupus or both. ⁽²⁴⁾
- d. Scl-70 antibodies appear to be a "marker" for Progressive Systemic Sclerosis (Scleroderma). ⁽²³⁾
- e. Anti-PCNA has been found in a small percentage of patients with SLE. ⁽¹⁵⁾

4. ANTI-CENTROMERE ANTIBODY (ACA) - discrete uniform speckles throughout the nucleus, the number of which corresponds to a multiple of the normal chromosome number. The staining pattern of the mitotic cells will follow that of the chromosomes, with pairs of dots arranging themselves in an equatorial plane during metaphase and then moving towards their respective centrosomes during anaphase.

Note: An antibody closely resembling the ACA is the Pseudo Centromere (or NSp-1) Antibody. This antibody can be differentiated from the ACA in that the chromosome region of mitotic cells does not stain. ⁽²⁵⁾

NUCLEAR ANTIGEN: Centromere/kinetochore portion of the chromosome.

DISEASE ASSOCIATION: This antibody is considered a "marker" for the CREST variant of Progressive Systemic Sclerosis. ⁽²⁶⁾ It is infrequently found in diffuse Scleroderma and Raynaud's disease. ⁽²⁷⁾

5. NUCLEOLAR - fluorescent staining of the nucleoli within the nucleus, sharply separated from the unstained nucleoplasm. The nucleolar fluorescence may be homogeneous, speckled or clumpy. Frequently accompanied by a speckled pattern.

NUCLEAR ANTIGENS: PM/Scl (PM-1), RNA Polymerase 1 or Fibrillarin (U3). ⁽²⁸⁾

DISEASE ASSOCIATION: High titers are highly specific for Progressive Systemic Sclerosis (PSS), with lower titers found in PSS, SLE, Sjogren's Syndrome and Raynaud's disease. ^(29,30)

6. ANTI-SPINDLE ANTIBODY - a network of threads connecting the centrosomes to each other in the mitotic cells.

NUCLEAR ANTIGEN: Spindle apparatus in cells undergoing mitosis.

DISEASE ASSOCIATION. There may be some association with Carpal Tunnel Syndrome. ⁽¹⁴⁾

7. CYTOPLASMIC FLUORESCENCE:

- a. Ribosomal RNP (rRNP) - diffuse granular fluorescence throughout the cytoplasm. May be confirmed by placing the specimen on a stomach section, where it will stain the chief cells. Frequently accompanied by a nucleolar pattern.

NUCLEAR ANTIGEN: Ribosomal or cytoplasmic RNP.

DISEASE ASSOCIATION: This antibody has been found in a small percentage of cases of SLE. ⁽³¹⁾

- b. Non-ANA Cytoplasmic Fluorescence: The two most common autoantibodies are:

1) *Anti-Mitochondrial Antibody (AMA)* - discrete speckles throughout the cytoplasm in a fibrous network, with denser speckling in the perinuclear region. May be confirmed by testing on appropriate tissue substrate.

2) *Anti-Smooth Muscle Antibody (ASMA)* - fluorescent strands in the cytoplasm in a spidery network, with fibrils extending from the cell membrane. May be confirmed by testing on appropriate tissue substrate.

**CORRELATION OF ANA'S
WITH IFA STAINING PATTERNS**

Antibody Against:	Staining Pattern:	Mitotic Cells:
nDNA	Peripheral & Homogeneous	Positive
Histones	Peripheral & Homogeneous	Positive
DNP	Homogeneous	Positive
Sm	Coarse Speckle	Negative
nRNP	Coarse Speckle	Negative
SS-A	Small Uniform Speckle	Negative
SS-B	Small Uniform Speckle	Negative
Scl-70	Fine Dense Speckle & Nucleolar	Positive
PCNA	Variable Speckle	Negative or Positive
ACA	Discrete Uniform Speckle	Positive Centromeres
PM/Scl	Nucleolar (homogeneous)	Negative
RNA Polymerase 1	Nucleolar (speckled)	Few Discrete Speckles
Fibrillar	Nucleolar (clumpy)	Positive fibers
Spindle	Spindle Apparatus	
rRNP	Cytoplasmic	Negative

LIMITATIONS OF PROCEDURE

- ANA serological test results should be used in conjunction with information available from the clinical evaluation and other diagnostic information.
- Two to ten percent of a normal adult population have antinuclear antibodies.⁽³²⁾
- Antinuclear antibodies are known to be age and sex related. With increasing age there is an increased incidence of ANA's; therefore, a positive low titer result may be normal for certain individuals in the absence of other clinical signs and symptoms. Antinuclear antibodies are not usually found, however, in normal young individuals.
- Some positive reactions have been reported in relatives of patients suffering from a connective tissue disease who may develop such a disease at a later time.⁽³³⁾
- Positive test results from cord blood or neonates should be interpreted with caution. The presence of antinuclear antibodies in cord blood is usually the result of passive transfer from mother to the fetus. A negative test, however, may be useful in excluding a possible autoimmune process.
- Positive ANA results may be seen in a small percentage of patients with Infectious and/or Neoplastic diseases, and also in diseases of drug etiology.^(22,34)
- SLE patients undergoing steroid therapy or in remission may have a negative ANA.⁽¹¹⁾
- Positive test results may not be valid in persons who have received blood transfusions or various blood products within the past several months.
- Test results on specimens from immunosuppressed patients and pregnant women may be difficult to interpret.
- Endpoint reactions may vary between laboratories due to differences in type or condition of fluorescence microscope employed or assay procedure used.⁽³⁵⁾
- If both the positive and negative control substrate cells are not visible when viewed using the fluorescence microscope, it may be necessary to replace or realign the light source and check the specific filters.
- Cell culture substrate slides may exhibit nonspecific fluorescence due to contamination of antibodies or PBS rinse-wash solutions with bacteria or fungi. It is very important that personnel reading the staining results have experience in fluorescence microscopy.
- In general titers of 1:40 and 1:80 are considered low titers, 1:160 and 1:320 are considered medium titers, and 1:640 and greater are considered high titers. It is recommended that each laboratory establish its own reference ranges.

EXPECTED VALUES

The following chart presents the incidence of antinuclear antibodies utilizing a HEp-2 cell ANA substrate in patient population studies performed at the Duke University Medical Center Division of Rheumatic and Genetic Disease laboratories over a two-year period. This represents a study of over 9,000 control sera and over 4,500 abnormal sera.

Clinical Diagnosis	% Positive	Diagnostic clinique	% Positive
Controls :		Vasculitides	20.0%
20-60 ans	2.0%	Childhood SLE	64.0%
70-80 ans	3.5%	JRA	
SLE	95.0%	Systemic	14.0%
RA	40.0%	Polyarthricular	6.0%
MCTD	99.0%	Pauciarticular	
PSS (Diffuse)	85.0%	HLA B27 pos.	0.0%
PSS (CREST variant)	93.0%	HLA B27 neg.	26.0%
PM/DM	25.0%		

SPECIFIC PERFORMANCE CHARACTERISTICS

To investigate the relative specificity and sensitivity of the ANA HEp-2 test, one hundred twenty serum specimens were compared qualitatively and semi-quantitatively with another commercially available indirect fluorescent ANA HEp-2 kit. All tables represent averaged results from two independent readers.

The relative sensitivity and specificity are summarized in TABLE 1. The one specimen in which the two kits did not agree was a CDC Reference Serum having high levels of antibodies to SS-A.

TABLE 1 – SUMMARY OF RELATIVE COMPARISON TESTING

OTHER KIT	ANA HEp-2		Relative Sensitivity	Relative Specificity
	Positive	Negative		
	71	0	100%	100%
	1	48		

In the following tables H = Homogeneous pattern, S = Speckled pattern, C = Centromere antibody, N = Nucleolar pattern, cR = cRNP antibody (cytoplasmic pattern). TABLE 2 represents the titers and patterns obtained with the 72 positive specimens in TABLE 1.

TABLE 2 – SUMMARY OF RELATIVE SENSITIVITY TESTING

Spec #	BMD	Other	Spec #	BMD	Other	Spec #	BMD	Other
1	80 H	80 H	25	40 S	40 S	49	320 S	160 S
2	40 H	40 H	26	80 S	80 S	50	640 S	640 S
3	80 H	40 H	27	80 S	40 S	51	640 S	160 S
4	40 H	40 H	28	40 S	40 S	52	160 S	40 S
5	160 H	80 H	29	40 S	40 S	53	2560 S	320 S
6	320 H	160 H	30	80 S	40 S	54	1280 S	640 S
7	1280 H	640 H	31	2560 S	2560 S	55	640 S	640 S
8	5120 H	2560 H	32	10240 S	5120 S	56	2560 S	2560 S
9	5120 H	5120 H	33	2560 S	40 S	57	160 S	< 40
10	640 H	320 H	34	640 S	640 S	58	640 S	320 S
11	10240 H	10240 H	35	2560 S	640 S	59	80 S	40 S
12	2560 H	1280 H	36	320 S	320 S	60	40 S	40 S
13	640 H	320 H	37	1280 S	1280 S	61	640 S	320 S
14	1280 H	640 H	38	10240 S	5120 S	62	640 C	640 C
15	640 H	320 H	39	5120 S	40 S	63	10240 C	5120 C
16	640 H	320 H	40	5120 S	640 S	64	320 C	320 C
17	1280 H	1280 H	41	10240 S	5120 S	65	640 C	320 C
18	1280 H	640 H	42	640 S	1280 S	66	1280 N	1280 N
19	640 H	320 H	43	10240 S	5120 S	67	640 N	320 N
20	80 H	40 H	44	160 S	160 S	68	320 N	160 N
21	1280 H	640 H	45	10240 S	10240 S	69	160 N	160 N
22	80 S	80 S	46	5120 S	5120 S	70	80 N	80 N
23	40 S	40 S	47	160 S	160 S	71	320cR	640cR
24	320 S	320 S	48	160 S	80 S	72	640S/160H	160 H

All but nine of the positive specimens agreed within a twofold dilution. The nine serum specimens were characterized by additional testing as follows:

- . #’s 33, 39, 57 (CDC Reference Serum) and 72 were characterized as having antibodies against SS-A.
- . #35 was characterized as having antibodies against Scl-70.
- . #40 was characterized as having antibodies to both SS-A and SS-B.
- . #53 (CDC Reference Serum) was characterized as having antibodies against SS-B.
- . #’s 51 and 52 were unable to be characterized by additional testing.

Interlot precision of the **ANA HEp-2 SYSTEM**, represented in TABLE 3, was evaluated by testing eleven serum specimens (2 negative and 9 positive over a range of titers) on three successive days using three different lot numbers. There was no more than a twofold difference (+/-) in titer between any of the comparison testings, which is within the confidence limits of this methodology.³⁵

Intralot precision of the **ANA HEp-2 SYSTEM**, represented in TABLE 4, was evaluated by running nine different types of antinuclear antibodies three times within one run using three different slides from the same lot. Again, there was no more than a twofold difference in titer between any of the comparison testings.³⁵

TABLE 3 – SUMMARY OF INTERLOT PRECISION

Spec #	Lot 1	Lot 2	Lot 3
1	160 N	160 N	160 N
2	160 S	160 S	80 S
3	320 S	640 S	320 S
4	2560 S	2560 S	1280 S
5	640 H	1280 H	1280 H
6	<40	<40	<40
7	320 C	640 C	640 C
8	40 S	80 S	40 S
9	<40	<40	<40
10	2560 S	5120 S	2560 S
11	2560 S	5120 S	2560 S

TABLE 4 – Summary of Intralot Precision

Spec #	Test #1	Test #2	Test #3
1	1280 N	1280 N	1280 N
2	320 H	320 H	320 H
3	1280 S	1280 S	1280 S
4	1280 C	1280 C	1280 C
5	2560 S	5120 S	5120 S
6	320 S	320 S	320 S
7	320 S	320 S	320 S
8	80 cR	80 cR	80 cR
9	320 S	320 S	320 S

BIBLIOGRAPHY

1. Tan, E., Antinuclear Antibodies: Diagnostic Markers for Autoimmune Diseases and Probes for Cell Biology, *Adv. Immuno.*, 44:93-151, 1989.
2. Nakamura, R., C. Peebles, R. Rubin, et al, Autoantibodies to Nuclear Antigens, 2nd éd., ASCP, Chicago, 1985.
3. Tan, E., E. Chan, K. Sullivan, et al, Antinuclear Antibodies (ANA's): Diagnostically Specific Immune Markers and Clues Toward the Understanding of Systemic Autoimmunity, *Clin. Immunol. Immunopathol.*, 47:121-141, 1988.
4. Barnett, E.V., Antinuclear Antibodies and Nuclear Antigens, *California Medicine*, 104:463-469, 1966.
5. Casais, S.P., G.J. Friou, L.L. Meyers, Significance of Antibody to DNA in Systemic Lupus Erythematosus, *Arth. Rheum.*, 7:379-390, 1964.
6. Douvas, A.S., M. Achten, E.M. Tan, Identification of a Nuclear Protein (Scl-70) as a Unique Target of Human Antinuclear Antibodies in Scleroderma, *J. Biol. Chem.*, 244:10514-10522, 1979.





7. Moroi, Y., C. Peebles, M.J. Fritzler, et al, Autoantibody to Centromere (Kinetochores) in Scleroderma Sera, *Proc. Natl. Acad. Sci.*, 77:1627-1631, 1980.
8. Sharp, G.C., W.S. Irwin, E.M. Tan, et al, Mixed Connective Tissue Disease – An Apparently Distinct Rheumatic Disease Syndrome Associated with a Specific Antibody to an Extractable Nuclear Antigen (ENA), *Am. J. Med.*, 52:148-159, 1972.
9. Alsough, M.A., W. Talal, E.M. Tan, Differentiation and Characterization of Autoantibodies and their Antigens in Sjogrens Syndrome, *Arth. Rheum.*, 19:216-222, 1976.
10. Wolfe, J.F., E. Adelstein, G.C. Sharp, Antinuclear Antibody with Distinct Specificity for Polymyositis, *J. Clin. Invest.*, 59:176-178, 1977.
11. Nakamura, R.M., CL. Peebles, D.P. Molden, E.M. Tan, Advances in Laboratory Tests for Autoantibodies to Nuclear Antigens in Systemic Rheumatic Diseases, *Lab. Med.*, 15:190-198, 1984.
12. Moore, A.E., L. Sabachewsky, H.W. Toolan, Culture Characteristics of Four Permanent Lines of Human Cells, *Cancer Res.*, 15:598-602, 1955.
13. McCarty, G.A., J.R. Rice, Characterization and Comparison of Commercially Available Antinuclear Antibody Kits Using Single Pattern Index Sera, *J. Rheum.*, 7:339-347, 1980.
14. McCarty, G.A., et al, A Unique Antinuclear Antibody Staining Only the Mitotic Spindle Apparatus, *New Engl. J. Med.*, 305:703, 1981.
15. Miyachi, K., M.J. Fritzler, E.M. Tan, Autoantibody to a Nuclear Antigen in Proliferating Cells, *J. Immunol.*, 121:2228-2234, 1978.
16. Schur, P.H., M. Moore, N. Rothfield, The Subclass of Antinuclear and Antinucleic Acid Antibodies, *Arth. Rheum.*, 15:174-180, 1972.
17. Weller, T.H., A.H. Coons, Fluorescent Antibody Studies with Agents of Varicella and Herpes Zoster Propagated In Vitro, *Proc. Soc. Exp. Biol. Med.*, 86:789-794, 1954.
18. Friou, F.J., Clinical Application of Lupus Serum-Nucleoprotein Reaction Using the Fluorescent Antibody Technique, *J. Clin. Invest.*, 36:890, 1957.
19. Holman, H.R., H.G. Kunkel, Affinity between the Lupus Erythematosus Serum Factor and Cell Nuclei and Nucleoprotein, *Science*, 126:162-163, 1957.
20. Lyster, H.C., F.T. Forrester, The Immunofluorescence (IF) Test, in *Immunofluorescence Methods in Virology*, USDHHS, Georgia, 71-81, 1979.
21. Peter, J.B., R.L. Dawkins, Evaluating Autoimmune Diseases, *Diag. Med.*, 1-10, Sept/Oct/1979.
22. Notman, D.D., N. Kurata, E.M. Tan, Profiles of Antinuclear Antibodies in Systemic Rheumatic Diseases, *Annals Intern. Med.*, 83:464-469, 1975.
23. Nakamura, R.M., CL. Peebles, R.L. Rubin, et al, Autoantibodies to Nuclear Antigens (ANA), ASCP, Chicago, 1985.
24. Molden, D.P., ANA Profiles in Systemic Rheumatic Disease, *Diag. Med.*, 12-18, June 1985.
25. Fritzler, M.J., D.W. Valencia, G.A. McCarty, Speckled Pattern Antinuclear Antibodies Resembling Anticentromere Antibodies, *Arth. Rheum.*, 27(1):92-96, 1984.
26. Fritzler, M.J., T.D. Kinsella, E. Garbutt, The CREST Syndrome: A Distinct Sérologie Entity with Anticentromere Antibodies, *Am. J. Med.*, 69:520-526, 1980.
27. Tan, E.M., G.P. Rodnan, I. Garcia, et al, Diversity of Antinuclear Antibodies in Progressive Systemic Sclerosis, *Arth. Rheum.*, 23:617-625, 1980.
28. Penning, C, V. Steen, G. Reimer, et al, An Analysis of Systemic Sclerosis Patients with Autoantibodies to the Nucleolar Proteins PM-Scl, RNA Polymerase I and Fibrillar, *A.R.A. Mtg.*, D29, 1987.
29. Pinnas, J.L., J.D. Northway, E.M. Tan, Antinucleolar Antibodies in Human Sera, *J. Immunol.*, 111:996-1004, 1973.
30. Ritchie, R.F., Antinucleolar Antibodies: Their Frequency and Diagnostic Association, *New Engl. J. Med.*, 282:1174-1178, 1970.
31. Bianchi, F.B., M. Rizzetto, P. Penfold, et al, Ultrastructural Localization and Characterization of a Ribosomal Antibody Detected by Immunofluorescence in Systemic Lupus Erythematosus, *Clin. Exp. Immunol.*, 17:629-636, 1974.
32. Rothfield, N.F., Detection of Antibodies to Nuclear Antigens by Immunofluorescence, *Manual of Clinical Immunology*, ASM, Washington D.C., 85:647-651, 1976.
33. Pollack, V.E., Antinuclear Antibodies in Families of Patients with Systemic Lupus Erythematosus, *New Engl. J. Med.*, 271:165-171, 1964.
34. Tan, E.M., M.J. Fritzler, P.A. Reyes-Lopez, Autoantibodies as Biological Markers of Neoplasia: Basic and Applied Aspects, Elsevier-North Holland, New York, 1978.
35. Cleymaet, J.E., R.M. Nakamura, Indirect Immunofluorescent Antinuclear Antibody Tests: Comparison of Sensitivity and Specificity of Different Substrates, *Am. J. Clin. Path.*, 58:388-393, 19.

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
SPECIMEN APPLICATION	1 drop (~20-30µl) 1 drop (~20-30µl) 1 drop (~20-30µl)	Negative control Positive Control Diluted sample	30 min at R.T <u>in a moist chamber</u>
RINSE 1	Using a squeeze wash bottle, rinse GENTLY with PBS		
WASH 1	Wash slides in Coplin jars full of PBS ☞ PBS must be change 2 times		5 to 6 minutes
CONJUGATE APPLICATION	1 drop (~20-30µl)	FITC conjugate	30 min at R.T <u>in a moist chamber</u>
RINSE 2	Using a squeeze wash bottle, rinse GENTLY with PBS		
WASH 2	Wash slides in Coplin jars full of PBS ☞ PBS must be change 2 times		5 to 6 minutes
READ	2 to 4 drops	Mounting Medium	

R.T: Room Temperature (+20°C/+25°C)

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Slide with n wells		Negative control		Mounting Medium
	Fluorescent antibody conjugate		Positive control		IFA Diluent
	Wash Buffer		Reconstitute with 1L distilled water		Immunofluorescent assay
	Contains sodium azide				Contains thimerosal

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com