Fr. Mars 2009

#### DÉFINITION

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Ag. March 2009

**INTENDED USE**

The ANA HEp-2 Optional Reagents (bmd) may be utilized in the indirect fluorescent antibody assay method for the qualitative and/or semi-quantitative determination of specific IgG antibodies in human serum.

The ANA HEp-2 Optional Reagents (bmd) are intended for use in the ANA HEp-2 (bmd) Test Kits.

**DESCRIPTION**

**FITC Conjugate** (REF HME 9920) (CONJ) (FITC) (IgG)

Dropper vial containing 20 mL fluorescein isothiocyanate labeled goat antihuman immunoglobulins with Evans Blue counterstain, protein stabilizer, less than 0.1% sodium azide and 0.001% thimerosal. Stable at +2°C/+8°C away from direct light until labeled expiration date.

**Mounting medium** (REF HME MM) (MM)

Dropper vial containing 3.5 mL phosphate buffered glycerol of pH 7.4 ± 0.2.

**IFA diluent** (REF HME DIL 250) (DIL) (SPE)

Vial containing 250 mL of sample diluent, with less than 0.1% thimerosal, formulated to reduce nonspecific staining.

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One-liter packet of dry PBS (10 packets)

**Coverslips 24x60** (REF HME 12CS)

Box of 100 cover glass for slides of 4, 6, 8 or 12 wells.

**Coverslips 24x70** (REF HME 16CS)

Box of 100 coverslips for slides of 16 wells.

**WARNINGS AND PRECAUTIONS**

- For *in vitro* diagnostic use only.
- Store all reagents at appropriate temperature. Do not use if stored improperly. All reagents should be brought to room temperature (+20°C/+25°C) prior to use.
- Substitution of components other than those provided may yield inconsistent results.
- Do not expose conjugate to strong light during storage or use. The preservatives used in conjugates are toxic if ingested. Azides may react with copper or lead plumbing to form explosive metal azides. When disposing, flush drains with water to minimize buildup of azide and metal compounds.
- Avoid microbial contamination of all reagents involved in the testing procedure or incorrect results may occur.
- Incubation times or temperatures other than those specified may give erroneous results.
- Reusable glassware must be washed and thoroughly rinsed free of detergents.
- IFA diluent** (bmd) should be used ONLY as a diluent for patient specimens. Do NOT prepare serial dilutions for endpoint titer with the IFA diluent. Do NOT use in any of the wash steps.

**PACKAGING**

HME 9920:	FITC Conjugate (1x20mL)
HME MM:	Mounting Medium (1 x 3.5mL)
HME DIL250:	IFA Diluent (1 x 250mL)
HME PBS10:	PBS packet (10 packets, each for 1L)
HME 12CS or 16CS :	Coverslips (box of 100 pieces)

**STABILITY AND STORAGE**

- FITC Conjugate** is stable at +2°C/+8°C away from direct light until labeled expiration date.
- Mounting Medium** is stable at +2°C/+8°C C until labeled expiration date.
- IFA diluent** is stable at +2°C/+8°C until the labeled expiration date provided no gross contamination is seen. Do not use if the solution turns cloudy, or if a precipitate forms.
- PBS** is stable in sealed packet at +25°C or lower until labeled expiration date.
- Do not freeze.
- Do not used beyond stated expiration date.

**PREPARATION OF REAGENTS**

The ANA HEp-2 Optional Reagents (bmd) are ready to use except the PBS.

**BUFFER PREPARATION**

Place contents of a one-liter PBS packet into a one-liter volumetric flask, add\* distilled water to the one-liter mark, mix and leave several hours or overnight to dissolve. Reconstituted buffer should have a pH of 7.4 ± 0.2.

Adjust with 1N NaOH or 1N HCl if pH value is outside the stated range. Store in a clean screw capped bottle at +25°C or lower. Stable during 1 month if no gross contamination is seen. Do not use if pH changes, if the solution turns cloudy, or if a precipitate forms.

\*Use deionized water with caution, as pH of this type of water may vary causing the pH of PBS to become unstable upon prolonged storage.


**PERFORMANCE CHARACTERISTICS**

These reagents have been tested with ANA HEp-2 (bmd) Test kits and found to meet performance requirements.

**LIMITS**

They have not been tested for performance in other manufacturer's test systems.

**SYMBOLS USED**

	Temperature limitation	(REF)	Catalogue Number
	Upper limit of temperature	(LOT)	Batch Code
(IVD)	In Vitro Diagnostic Medical Device		Used by
	Read instructions for use	CE	EC Declaration of Conformity
(RCNS) 1L H <sub>2</sub> O	Reconstitute with 1L distilled water		
(CONT) NaN <sub>3</sub>	Contains sodium azide		
(CONT) C9H9HgNaO2S	Contains thimerosal		