ELISIS ENA Profile

LIS 1101

Français

DÉFINITION

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REF LIS 1101

DEFINITION

ELISIS ENA Profile kit ([bmd](#)) is an enzyme linked immunoassay (ELISA) for the qualitative determination of 8 antibodies in human samples. The wells are separately coated with U1-snRNP, snRNP/Sm, Sm, SS-A (52 and 60 kDa), SS-B, Scl-70, CENP-B and Jo-1.

The results of the **ELISIS ENA Profile kit** ([bmd](#)) are to be used in conjunction with the clinical findings and the other laboratory tests to aid in the diagnosis of connective diseases (systemic lupus erythematosus (SLE), Sjogren's syndrome, mixed connective tissue disease (MCTD), scleroderma, dermatomyositis and CREST syndrome).

ELISIS ENA Profile kit ([bmd](#)) may be used with the **CARIS™ system** (diluting and dispensing device).

DIAGNOSTIC VALUE

Connective diseases are systemic auto-immune diseases and each one is, most often, selectively directed against a specific organ. Their classification is based upon clinical and biological data. The detection of anti-nuclear antibodies, which often appears in the patient serum before any definite symptomatic sign, allows the identification of the specific connective disease.

Many elaborate serological studies have demonstrated significant associations between these autoantibodies and the various connective diseases:

Antibodies associated with different autoimmune pathologies

• **Anti-SSA and anti-SSB antibodies**

Either one or both of these are observed in two pathologies: Sjögren's syndrome and Systemic Lupus Erythematosus (SLE). Anti-SSA antibodies are also found in mothers who are carrying a baby with neonatal lupus syndrome including heart block.

• **Anti-RNP antibodies**

- ✓ Antibody against U1 snRNP is directed to the 70 kDa protein of U1 snRNP. They are pathognomic for MCTD but do also occur in SLE. A high titer of antibodies against this antigen is typical for the Sharp-Syndrome.
- ✓ Antibody against snRNP/Sm are directed against Sm and U1 -snRNP proteins (70 kDa, A and C). They occur in SLE, Sjögren's syndrome, scleroderma and polymyositis.

Antibodies specific to a single connective tissue disease

- **Anti-Sm antibodies** are directed against core proteins (B,B\ D1-D3, E, F, G) of small nuclear ribonucleoproteins (snRNPs). Anti-Sm as well as antibodies against double stranded DNA (dsDNA) are highly specific for SLE and thus are included in diagnostic and classification criteria for SLE.
- **Anti-Scl-70 antibodies** are directed against DNA-topoisomerase I. They are highly specific for systemic scleroderma and give a hint for a severe course.
- **Anti-Jo1 antibodies** are directed against histidyl-tRNA synthetase (cytoplasmic protein involved in protein biosynthesis) and are found in 20-40 % of patients with polymyositis and dermatomyositis.
- **Anti-Centromere antibodies (CENP-B)** are most often found in sera from patients having Scleroderma with limited skin involvement CREST, called CREST syndrome.

ASSAY PRINCIPLE

Each different extractable nuclear antigen is coated onto specific wells of a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well for each specific antigen, which allows to bind. After incubation, unbound proteins are removed by washing.
- Then, a monoclonal anti-human IgG conjugated to horseradish peroxidase is added to each of the wells. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of TMB substrate (3,3',5,5' - tetramethylbenzidine). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution HCL (1 M) serves to stop the enzymatic reaction.
- After stopping the reaction by HCL (1 M) the optical density is read by a spectrophotometer at 450nm
- A control calibrator allows to express in arbitrary units the titre of the patient sample.
- The rate of color formation from the chromogen is a function of the amount of conjugate bound to the antigen-antibody complex and this is proportional to the initial concentration of the respective antibodies in the patient sample.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum. Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

REAGENTS

96 wells microplate coated with recombinant 70 kDa U1-snRNP, SS-B, SS-A 52 kDa, Scl-70, CENP-B, Jo-1 and purified native human snRNP/Sm, Sm and SS-A 60 kDa (12x8 individual breakaway microwells)	12 strips
2 vials of Cut-off Calibrator (capped blue: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use	2 x 1.8mL
2 vials of Negative control (capped green: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use	2 x 1.8mL
1 vial of anti-human IgG-HRP Conjugate (capped blue: blue solution) Ready to use	1 x 15mL
1 vial of TMB Substrate (capped black) containing: stabilized TMB/H ₂ O ₂ Ready to use	1 x 15mL
1 vial of Sample Diluent - 5x concentrated (capped white: yellow solution) containing: Tris, NaCl, BSA, sodium azide < 0.1% (preservative) To dilute	1 x 20mL
1 vial of Wash Buffer - 50x concentrated (capped white: green solution) Containing: Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative) To dilute	1 x 20mL
1 vial of Stop Solution (capped white: colorless solution) containing: 1M Hydrochloric Acid Ready to use	1 x 15mL

ELISIS ENA Profile kit ([bmd](#)) is calibrated against reference sera from the CDC (Centres for Disease Control and Prevention) Atlanta. The results are expressed in AU/mL.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

- Purified water (according to the definition of the United States Pharmacopeia and the European Pharmacopeia)
- Microtiter plate reader 450nm reading filter and optional 620nm reference filter (600-690nm)
- Glass ware (cylinder 100-1000ml)
- Test tubes for serum dilutions
- Vortex mixer
- Precision pipettes (10, 100, 200, 500, 1000µl) or adjustable multipipette (100-1000ml)
- Microplate washing device (300µl repeating or multichannel pipette or automated System)
- Adsorbent paper

SETUP

All the reagents should be prepared as required:

1. Preparation of Sample diluent solution

- Dilute concentrated Sample diluent 1/5 in distilled water.
Ex: 20mL plus 80mL
- Storage period: 1 month at +2°C/+8°C.

2. Preparation of wash buffer

- Dilute concentrated Wash buffer 1/50 in distilled water.
Ex: 20mL plus 980mL
- Prepare 20mL of diluted wash buffer for 8 wells or 200mL for 96 wells
- Storage period: 1 month at +2°C/+8°C.

3. Preparation of samples

- Dilute to 1/101 in Sample diluent (1x)
Ex: 1000µl sample diluent (1x) + 10 µl serum
- Shake vigorously at the vortex.

4. Automated washing:

- Consider excess volumes required for setting up the instrument and dead volume of robot pipette.

5. Manual washing:

- Discard liquid from wells by inverting the plate.
- Knock the microwell frame with wells downside vigorously on clean adsorbent paper.
- Pipette 300µl of diluted wash buffer into each well, wait for 20 seconds.
- Repeat the whole procedure twice again.

6. Microplates

- Calculate the number of wells required for the test.
- Remove unused wells from the frame, replace and store in the provided plastic bag, together with desiccant, seal tightly (+2°C/+8°C).

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C / +25°C) before handling.

Check that all plates are well drained after each wash.

Do not smoke, eat or drink when manipulating the kit.

Do not pipette by mouth.

Human sources for the preparation of calibrator and control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

This kit contains potentially hazardous components. Though kit reagents are not classified being irritant to eyes and skin we recommend to avoid contact with eyes and skin and wear disposable gloves.

Calibrator, Control and Buffer contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

ELISIS ENA Profile kit ([bmd](#)) has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

ELISIS ENA Profile kit ([bmd](#)) has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.

Incubation: We recommend test performance at +30°C for automated Systems.

Never expose components to higher temperature than +37°C.

Always pipette substrate solution with brand new tips only. Protect this reagent from light. Never pipette conjugate with tips used with other reagents prior.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

Set out:

- Cut-off calibrator (Ctrl Ref)
- Negative control (Ctrl -)
- Patient samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
U1-RNP	A	Ctrl Ref	Ctrl -	S1	S2	S3						
SnRNP/Sm	B	Ctrl Ref	Ctrl -	S1	S2	S3						
Sm	C	Ctrl Ref	Ctrl -	S1	S2	S3						
SS-A 52/60	D	Ctrl Ref	Ctrl -	S1	S2	S3						
SS-B	E	Ctrl Ref	Ctrl -	S1	S2	S3						
Scl-70	F	Ctrl Ref	Ctrl -	S1	S2	S3						
CenpB	G	Ctrl Ref	Ctrl -	S1	S2	S3						
Jo-1	H	Ctrl Ref	Ctrl -	S1	S2	S3						

2. Sample incubation

Add 100µL of Cut-off calibrator and Negative control into the designated microwells

Add 100µL of Samples.

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing buffer

3. Incubation of conjugate

Add 100µL of Conjugate

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 300µL of washing buffer

4. Incubation of substrate

Add 100µL TMB Substrate into each well

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Protect the wells from intense light

5. Stop solution

Add 100µL of Stop solution into each well

Incubate for 5 minutes minimum

Agitate plate carefully for 5 sec.

6. Reading

Read the optical density of each well at 450nm (optionally 450/620 nm) within 30 minutes.

INTERPRETATION OF RESULTS

Read the optical density of the cut-off control and the patient samples. Multiply the OD of the cut-off calibrator by the parameter specific factor, provided with the lot specific QC certificate. Compare patient ODs with the calculated parameter OD cut-off value. All samples which are higher than cut-off are considered positive.

OD Cut-off calibrator x antigen specific factor = OD Cut-off parameter.

Negative:	OD	Patient	<	OD cut-off
Positive:	OD	Patient	>	OD cut-off

Calibrators	OD 450/620 nm	CV % (Variation)
Negative Control	0.047	2.6
Cut-off Calibrator	0.350	1.8

Example of interpretation

We recommend pipetting cut-off calibrator in parallel for each run.

U1-snRNP70kD Factor= 1.0

OD Cut-off calibrator 0.550 x 1.0 = 0.550

OD patient 0.999 > 0.550 = positive

Do not use this example for interpreting patients results!

bmd recommend to retest samples that are equivocal.

For lot specific data, see enclosed quality control leaflet.

Each laboratory should establish its own normal range based upon its own techniques, controls, equipment and patient population according to their own established procedures.

For **semi-quantification of the results**, each patient-OD value can be expressed by the Index-Value. The Index-Value is calculated by dividing the patient-OD by the cut-off OD:

$$\text{Index Value} = \frac{\text{OD (patient sample)}}{\text{OD (cut-off calibrator)}}$$

Negative => Index Value < 1.0

Positive => Index Value > 1.0

A definite clinical diagnosis should not be based on the results of the performed test only, but should be made by the physician after all clinical and laboratory findings have been evaluated.

The diagnosis is to be verified using different diagnostic methods.

LIMITATION

Hemolytic, lipemic, icteric samples or contaminated samples with may confound the results of this assay.

QUALITY CONTROL

The kit contains a negative control. The negative control should be included with each run of the test.

bmd offers an external multiparameter quality control (IMMUNO-TROL IV, Cat: HM051; IMMUNO-TROL V, Cat: HM052). These controls should be handled as samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

Sensitivity, specific affinity and agreement to other methods

ELISIS ENA Profile kit has been compared to one predicate ELISA kit.

The study was performed on 340 samples from patients with associated systemic auto-immune diseases.

Results

Antigenic specificity	Relative sensitivity (%)	Relative specificity (%)	Concordance (%)
U1-snRNP70kD	87.0	100	95.8
snRNP/Sm	96.2	97.9	97.0
Sm	96.0	94.2	95.1
SS-A 52 + 60 kDa	80.3	95.8	84.7
SS-B	84.7	100	91.1
Scl-70	100	94.3	97.1
Centromère	96.1	98.1	97.1
Jo-1	100	100	100

Test Precision

Estimated for each parameter on three positive samples : strongly, moderate and weakly positive.

	ELISIS ENA Profile			
	Intra-assay (24 tests in a same assay)		Inter-assay (18 tests in 3 differents assays)	
	Titer mean AU/mL	CV (%)	Titer mean AU/mL	CV (%)
U1-snRNP70kD	151	2.3	148	2.6
	72	3.6	70	4.1
	26	4.5	24	3.9
snRNP/Sm	85.1	6.2	83.6	6.4
	50.5	3.3	51.3	3.7
	19.1	1.0	19.1	1.6
Sm	64.2	4.6	50.4	3.3
	41.5	3.3	29.6	2.3
	20.4	1.9	12.7	1.4
SS-A 52 + 60 kDa	78.2	8.1	102.4	3.8
	44.3	2.8	62.5	5.5
	22.9	2.0	26.3	1.9
SS-B	87.0	6.0	88.5	5.2
	55.5	3.1	58.9	4.2
	20.7	1.8	21.6	1.8
Scl-70	66.5	7.7	57.5	4.4
	37.7	4.3	52.2	4.2
	19.4	1.9	17.7	1.3
Centromere	265	1.5	261	1.4
	122	2.3	119	2.5
	44	4.1	46	3.8
Jo-1	112.9	9.3	112.9	9.3
	63.3	7.5	63.3	7.5
	15.3	1.0	15.3	1.0

Threshold value

Calculated on 40 samples from « blood donors ».

	ELISIS ENA Profile
Cut-off : Mean value +3SD	15 AU/mL

REFERENCES

CONRAD K. et al.

Autoantibodies – Diagnostic, pathogenic and prognostic relevance.

Clin. and Exp. Rheumatol. 1997 ;15, 457-465

GAL I et al.

Comparison of anti-Ro/SSA antibody profile between patients with primary and secondary syndrome.

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HOMBURGER H.A.

Cascade testing for autoantibodies in connective tissue disease.

Mayo Clin. Proc., 1995, 70, 183-184.

KEECH CL et al.

SS-B (La) autoantibodies. Autoantibodies.

Peter JB and Shoenfeld Y, editors. Ed Elsevier 1996, 789-797.

MEWIS A. et al.

Evaluation of a DOT-Blot method for identification of antibodies against extractable nuclear antigens and anti-cytoplasmic antibodies.

Clin Chem., 1999, 45 (2), 311-312.

RAHMAN MA et al.

Autoantibodies in systemic lupus erythematosus.

Curr Opin Rheumatol 1994, (6), 468-473.

SMEENK RJ et al.

dsDNA autoantibodies.

In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234.

SPENCER GREEN G. et al.

Test performance in systemic sclerosis: anti-centromere and anti-Scl70 antibodies.

Am. J. of Med., 1997, 18, 291-300.

THE LS et al.

Antiribosomal P protein antibodies in systemic lupus erythematosus.


Arthritis Rheum. 1994, (37), 307-315.

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	100µL 100µL 100µL	Cut-off calibrator (in duplicate) Negative Control (in duplicate) Diluted samples (in duplicate)	30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)		
INCUBATION OF CONJUGATE	100µL	Anti-human IgG-HRP conjugate	30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)		
ENZYME REACTION	100µL	TMB Substrate	30 min at RT protected from intense light
STOP REACTION	100µL	Stop solution (1M HCl)	5 min

RT: room temperature (+20°C/+32°C)

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Sample Diluent		Stop Solution		Cut-off calibrator
	Wash Buffer		Calibrator		Microplate
	Control		Substrate		Conjugate

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com