ELISIS RF IgG/IgM/IgA

REF LIS 1161 

DÉFINITION

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SHMERLING RH, DELBANCO TL.

The rheumatoid factor : an analysis of clinical utility.
Am J Med 1991 ;91 :528-534.

BORRETZEN M et al.

Rheumatoid Factors.
In : Autoantibodies. JB Peter and Y. Schoenfeld Eds. Elsevier Sciences B.V. 1996 ;706-715.

ARNETT FC et al.

The American Rheumatism Association 1987 revised criteria for classification of rheumatoid arthritis.
Arthritis Rheum 1988 ;31 :315-324.

MOORE TT et al.

Rheumatoid Factors.
Clin Biochem 1993 ;26 :75-84.

CORPER AL et al.

Structure of human IgM rheumatoid factor Fab bound to its autoantigen IgG Fc reveals a novel topology of antibody-antigen interaction.
Nat Struct Biol. 1997 May;4(5):374-81.

WITTE T, HARTUNG K, SACHSE C, MATTHIAS T, FRICKE M, KALDEN JR, LAKOMEK HJ, PETER HH, SCHMIDT RE.

Rheumatoid factors in systemic lupus erythematosus: Association with clinical and laboratory parameters. SLE study group. *Rheumatol Int 19: 107-111 ; 2000.*

SMITH JB, HAYNES MK.

Rheumatoid arthritis - a molecular understanding.
Ann Intern Med 136: 908 - 922 ; 2002.

SCHELLEKENS GA, VISSER H, de JONG BAW and al.

The diagnostic properties of rheumatoid arthritis antibodies recognizing a cyclic citrullinated peptide.
Arthritis Rheum 43: 155-163 ; 2000.

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

Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies.
J Clin Invest 101: 273-281.

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English

DEFINITION

ELISIS RF IgG/IgM/IgA kit ([bmd](#)) is an enzyme linked immunoassay (ELISA) for the qualitative or quantitative determination of IgG, IgM and IgA class rheumatoid factors directed against Fc determinants from human G Immunoglobulins in human serum samples.

The results of the **ELISIS RF IgG/IgM/IgA kit** ([bmd](#)) are to be used in conjunction with the clinical findings and the other laboratory tests to aid in the diagnosis of rheumatoid arthritis (RA).

ELISIS RF IgG/IgM/IgA kit ([bmd](#)) may be used with the **CARIS™ system** (diluting and dispensing device).

DIAGNOSTIC VALUE

Rheumatoid arthritis (RA) is one of the most common autoimmune diseases, affecting approximately 1% of the population worldwide. It is a chronic, systemic, inflammatory disease that predominantly affects the synovial membranes of multiple joints in the body.

Up to now, rheumatoid factors (RF) have been the main serological markers for this disease. Rheumatoid factors, first described in 1940 as antibodies reacting with gamma globulins, are autoantibodies directed against the C-terminal part of the constant region of the IgG heavy chain, the IgG Fc. Although named after the disease they were initially associated with, RFs are found both in the healthy population and several diseases.

The diseases commonly associated with high RF concentrations are rheumatoid arthritis (50-90%) and Sjögren's syndrome (75-95%). They are also found in systemic lupus erythematosus (SLE; 15-35%), systemic sclerosis (20-30%) polymyositis / dermatomyositis (5-10%), cryoglobulinemia (40-100%) and mixed connective tissue diseases (MCTD; 50-60%). Therefore, despite being a sensitive marker for RA, rheumatoid factor lack specificity for the disease. It has long been known that anti-perinuclear autoantibodies, also called anti-keratin antibodies, are found in patients with RA. Recently it has been discovered that these autoantibodies target an epitope containing the deimidated form of arginine called citrulline. Antibodies against synthetic cyclic citrullinated peptides (CCP) have been described as being highly specific markers for RA with sensitivities comparable to RF.

ASSAY PRINCIPLE

Fc fragment of Human G Immunoglobulins antigen(s) antibodies are coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Then, a monoclonal anti-human IgG, IgM or IgA conjugated to horseradish peroxidase is added to each well. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of TMB substrate (3,3',5,5' – tetramethylbenzidine). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution HCl (1 M) serves to stop the enzymatic reaction.
- After stopping the reaction by HCl (1 M) the optical density is read by a spectrophotometer at 450nm
- The rate of color formation from the chromogen is a function of the amount of conjugate bound to the antigen-antibody complex and this is proportional to the initial concentration of the respective antibodies in the patient sample.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum. Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

REAGENTS

96 wells microplate coated with Fc fragment of Human G Immunoglobulins antigen(s). (12x8 individual breakaway microwells)	MP	12 strips
1 vial of Cut-off Calibrator (capped blue: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use	CONTROL REF	1 x 1.5mL
6 vials of Calibrators (color increasing with concentration: yellow solution) 0, 3, 10, 30, 100, 300 IU/mL containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use	CAL	6 x 1.5mL
1 vial of Positive control (capped red: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use	CONTROL +	1 x 1.5mL
1 vial of Negative control (capped green: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use	CONTROL -	1 x 1.5mL
1 vial of anti-human IgA-HRP Conjugate (capped red: red solution) Ready to use	CONJ IgA	1 x 15mL
1 vial of anti-human IgG-HRP Conjugate (capped blue: blue solution) Ready to use	CONJ IgG	1 x 15mL
1 vial of anti-human IgM-HRP Conjugate (capped green: green solution) Ready to use	CONJ IgM	1 x 15mL
1 vial of TMB Substrate (capped black) Containing: stabilized TMB/H ₂ O ₂ Ready to use	SUBS TMB	1 x 15mL
1 vial of Sample Diluent - 5x concentrated (capped white: yellow solution) Containing: Tris, NaCl, BSA, sodium azide < 0.1% (preservative) To dilute	DIL SPE 5X	1 x 20mL
1 vial of Wash Buffer - 50x concentrated (capped white: green solution) Containing: Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative) To dilute	BUF WASH 50x	1 x 20mL
1 vial of Stop Solution (capped white: colorless solution) containing: 1M Hydrochloric Acid Ready to use	SOLN STOP	1 x 15mL

ELISIS RF IgG/IgM/IgA kit ([bmd](#)) is calibrated in arbitrary units (AU/mL) for IgG and IgA rheumatoid factors and for IgM rheumatoid factors the results are expressed in IU/ml (International WHO Standard).

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

- Purified water (according to the definition of the United States Pharmacopeia and the European Pharmacopeia)
- Microtiter plate reader 450nm reading filter and optional 620 nm reference filter (600-690nm)
- Glass ware (cylinder 100-1000mL)
- Test tubes for serum dilutions
- Vortex mixer
- Precision pipettes (10, 100, 200, 500, 1000µL) or adjustable multipipette (100-1000µL)
- Microplate washing device (300µL repeating or multichannel pipette or automated System)
- Adsorbent paper

SETUP

All the reagents should be prepared as required:

1. Preparation of Sample diluent solution

- Dilute concentrated Sample diluent 1/5 in distilled water.
Ex: 20mL plus 80mL
- Storage period: 1 month at +2°C/+8°C.

2. Preparation of wash buffer

- Dilute concentrated Wash buffer 1/50 in distilled water.
Ex: 20mL plus 980mL
- Prepare 20mL of diluted wash buffer for 8 wells or 200mL for 96 wells
- Storage period: 1 month at +2°C/+8°C.

3. Preparation of samples

- Dilute to 1/101 in Sample diluent (1x)
Ex: 1000µL sample buffer (1x) + 10µL serum
- Shake vigorously at the vortex.

4. Washing Step

- ❖ Automated washing:
 - Consider excess volumes required for setting up the instrument and dead volume of robot pipette.
- ❖ Manual washing:
 - Discard liquid from wells by inverting the plate.
 - Knock the microwell frame with wells downside vigorously on clean adsorbent paper.
 - Pipette 300µL of diluted wash buffer into each well, wait for 20 seconds.
 - Repeat the whole procedure twice again.

5. Microplates

- Calculate the number of wells required for the test.
- Remove unused wells from the frame, replace and store in the provided plastic bag, together with desiccant, seal tightly (+2°C/+8°C).

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C / +25°C) before handling.

Check that all plates are well drained after each wash.

Do not smoke, eat or drink when manipulating the kit.

Do not pipette by mouth.

Human sources for the preparation of calibrators and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

This kit contains potentially hazardous components. Though kit reagents are not classified being irritant to eyes and skin we recommend to avoid contact with eyes and skin and wear disposable gloves.

Calibrators, Controls and Buffers contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

ELISIS RF IgG/IgM/IgA kit ([bmd](#)) has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

ELISIS RF IgG/IgM/IgA kit ([bmd](#)) has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.

Incubation: We recommend test performance at +30°C for automated Systems.

Never expose components to higher temperature than +37°C.

Always pipette Substrate and Conjugate with brand new tips only.

Protect Substrate from light.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

For quantitative interpretation use calibrators to establish a standard curve

Set out in duplicate:

- Calibrators (Cal A, Cal B, Cal C, Cal D, Cal E and Cal F)
- Positive control (Ctrl +)
- Negative control (Ctrl -)
- Diluted samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
A	Cal A	Cal E	S1	S9								
B	Cal A	Cal E	S2	S10								
C	Cal B	Cal F	S3	S11								
D	Cal B	Cal F	S4	S12								
E	Cal C	Ctrl +	S5	S13								
F	Cal C	Ctrl +	S6	S14								
G	Cal D	Ctrl -	S7									
H	Cal D	Ctrl -	S8									

For qualitative interpretation use cut-off calibrators

Set out in duplicate:

- Negative control (Ctrl -)
- Cut-off calibrator (Ctrl Ref)
- Positive control (Ctrl +)
- Diluted samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
A	Ctrl -	S3	S11	S19								
B	Ctrl -	S4	S12	S20								
C	Ctrl Ref	S5	S13	S21								
D	Ctrl Ref	S6	S14	S22								
E	Ctrl +	S7	S15	S23								
F	Ctrl +	S8	S16	S24								
G	S1	S9	S17									
H	S2	S10	S18									

Note: If IgA, IgG and IgM are determined in parallel calibrators, controls and samples have to be done for each subclass separately.

2. Sample incubation

Add 100µL of calibrators OR cut-off calibrator and controls into the designated microwells

Add 100µL of diluted samples.

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing buffer.

3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 300µL of washing buffer

4. Incubation of substrate

Add 100µL TMB substrate into each well

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Protect the wells from intense light

5. Stop solution

Add 100µL of Stop solution to each well

Incubate for 5 minutes minimum

Agitate plate carefully for 5 sec.

6. Reading

Read the optical density of each well at 450 nm (optionally 450/620nm) within 30 minutes.

INTERPRETATION OF RESULTS

For **quantitative interpretation** establish the standard curve by plotting the *optical density (OD) of each calibrator (y-axis)* with respect to the corresponding concentration values in *IU/mL or AU/mL*. For best results [bmd](#) recommend log (x) / lin (y) coordinates and 4-Parameter Fit. From the OD of each sample, read the corresponding antibody concentrations expressed in *IU/mL or AU/mL*.

Negative results	Positive Results
≤ 15 IU/mL or AU/mL	> 15 IU/mL or AU/mL

Example of a standard curve

We recommend pipetting calibrators in parallel for each run.

Calibrators IgA/IgG/IgM	OD 450/620 nm	CV % (Variation)
0 IU/mL or AU/mL	0.035	2.3
3 IU/mL or AU/mL	0.138	2.6
10 IU/mL or AU/mL	0.342	3.2
30 IU/mL or AU/mL	0.632	3.2
100 IU/mL or AU/mL	1.216	0.5
300 IU/mL or AU/mL	2.178	0.1

Example of calculation

Patient	Replicate (OD)	Mean (OD)	Result (IU/mL or AU/mL)
P01	0.872/0.922	0.897	54.7
P02	1.159/1.188	1.174	86.3

Do not use this example for interpreting patient's results!

For lot specific data, see enclosed quality control certificate.

Each laboratory should establish its own normal range based upon its own techniques, controls, equipment and patient population according to their own established procedures.

For **qualitative interpretation** read the optical density of the cut-off calibrator and the patient samples. Compare patient's OD with the OD of the cut-off calibrator. All samples which are higher than cut-off are considered positive.

Negative: OD Patient < OD cut-off

Positive: OD Patient > OD cut-off

A definite clinical diagnosis should not be based on the results of the performed test only, but should be made by the physician after all clinical and laboratory findings have been evaluated. The diagnosis is to be verified using different diagnostic methods.

LIMITATION

Hemolytic, lipemic, icteric samples or contaminated samples with may confound the results of this assay.

QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive and negative controls concentration should fall within the range printed on the quality control certificate. If either reagent control is invalid, the test should be repeated.

[bmd](#) offers an external multiparameter quality control (IMMUNO-TROL III, Cat: HM050). This control should be handled as samples.

REFERENCES

WILLIAMS DG.

Rheumatoid arthritis : autoimmunity in rheumatoid arthritis.
In : Klippel JH, Dieppe PA, eds. Rheumatology, 1st edition. St louis 1994 ;91 :9-11.

SHMERLING RH, DELBANCO TL.

The rheumatoid factor : an analysis of clinical utility.
Am J Med 1991 ;91 :528-534.

BORRETZEN M et al.

Rheumatoid Factors.
In : Autoantibodies. JB Peter and Y. Schoenfeld Eds. Elsevier Sciences B.V. 1996 ;706-715.

ARNETT FC et al.

The American Rheumatism Association 1987 revised criteria for classification of rheumatoid arthritis.
Arthritis Rheum 1988 ;31 :315-324.

MOORE TT et al.

Rheumatoid Factors.
Clin Biochem 1993 ;26 :75-84.

CORPER AL et al.

Structure of human IgM rheumatoid factor Fab bound to its autoantigen IgG Fc reveals a novel topology of antibody-antigen interaction. *Nat Struct Biol. 1997 May;4(5):374-81.*

WITTE T, HARTUNG K, SACHSE C, MATTHIAS T, FRICKE M, KALDEN JR, LAKOMEK HJ, PETER HH, SCHMIDT RE.

Rheumatoid factors in systemic lupus erythematosus: Association with clinical and laboratory parameters. SLE study group. *Rheumatol Int 19: 107-111 ; 2000.*

SMITH JB, HAYNES MK.

Rheumatoid arthritis - a molecular understanding.
Ann Intern Med 136: 908 - 922 ; 2002.

SCHELLEKENS GA, VISSER H, de JONG BAW and al.

The diagnostic properties of rheumatoid arthritis antibodies recognizing a cyclic citrullinated peptide.
Arthritis Rheum 43: 155-163 ; 2000.

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Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies.
J Clin Invest 101: 273-281 ; 1998

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS		INCUBATION CONDITIONS
		QUANTITATIVE METHOD	QUALITATIVE METHOD	
INCUBATION OF SAMPLES	100µL	Calibrators A, B, C, D, E and F (in duplicate)	Negative Control (in duplicate)	30 min at RT
	100µL	Positive Control (in duplicate)	Cut-off calibrator (in duplicate)	
	100µL	Negative Control (in duplicate)	Positive Control (in duplicate)	
	100µL	Diluted samples (in duplicate)	Diluted samples (in duplicate)	
WASHING	Wash 3 times with washing buffer (300µL/ well)			
INCUBATION OF CONJUGATE	100µL	Anti-human IgG-HRP conjugate OR Anti-human IgM-HRP conjugate OR Anti-human IgA-HRP conjugate		30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)			
ENZYME REACTION	100µL	TMB Substrate		30 min at RT protected from intense light
STOP REACTION	100µL	Stop solution (1M HCl)		5 min

RT: room temperature (+20°C/+32°C)

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Sample Diluent		Stop Solution		Cut-off calibrator
	Wash Buffer		Calibrator		Microplate
	Control		Substrate		Conjugate

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com