# ELISIS $\beta$ 2-GPI IgG / IgM

REF LIS 1206 

## DÉFINITION

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REF LIS 1206

## DEFINITION

**ELISIS  $\beta$ 2-GPI IgG/IgM kit** () is an enzyme linked immunoassay (ELISA) for the qualitative or quantitative determination of IgG and/or IgM anti- $\beta$ 2-GPI antibodies in human serum samples. Anti- $\beta$ 2-GPI antibodies recognize specific epitopes on human  $\beta$ 2-glycoprotein I which are only expressed when  $\beta$ 2-glycoprotein I interacts with lipid membranes.

The results of the **ELISIS  $\beta$ 2-GPI IgG/IgM kit** () are to be used in conjunction with the clinical findings and the other laboratory tests to aid in the diagnosis and risk estimation of primary and secondary AntiPhospholipid Syndrome (APS).

**ELISIS  $\beta$ 2-GPI IgG/IgM kit** () may be used with the **CARIS™ system** (diluting and dispensing device).

## DIAGNOSTIC VALUE

The antiphospholipid antibodies (aPL) are an extremely homogenous family of antibodies, detected by the *in vitro* elongation of phospholipid-dependent coagulation tests and by ELISA techniques which employ cardiolipin or other anionic or neutral phospholipids. The illnesses most often associated with aPL are systemic lupus erythematosus, drug-induced lupus, thrombo-embolic disease and thromboses in general, repeated abortion, autoimmune peripheral cytopenia, primitive biliary cirrhosis, a number of acute infections, retroviral infections, AIDS and HTLV. The expression, antiphospholipid syndrome denotes the association of one or several clinical criteria, venous and/or arterial thrombosis, consecutive abortions, premature births, fetal deaths, to the presence of a circulating anticoagulant, an aPL and/or an anti-  $\beta$ 2-GPI detected on two or more occasions.

The detection of some anionic aPLs found in autoimmune diseases, by ELISA technique calls for the presence of a serum protein,  $\beta$ 2-glycoprotein I ( $\beta$ 2-GPI) or apolipoprotein H. Conversely, infections which are not generally accompanied by thrombotic elements, do not require its presence.

$\beta$ 2-GPI is a 50 kDa plasma protein, occurring at a concentration of 200 $\mu$ g/ml, in the form of a polypeptide chain of 326 amino-acids. Its great affinity for the phospholipids is the basis of its *in vitro* anticoagulant effect on coagulation and platelet aggregation.

Today most authors consider that the  $\beta$ 2-GPI-dependent aPLs are probably directed against the protein which has been structurally modified by linkage with the PLs. Their *in vitro* recognition of  $\beta$ 2-GPI requires the use of an irradiated solid phase in order to reproduce the change in structure.

They confirm the value of detecting the aPL directed against the anionic phospholipids, by testing their antigen specific affinity to the  $\beta$ 2-GPI, in order to increase their ability to predict the risk of thrombosis.

## ASSAY PRINCIPLE

The purified  $\beta$ 2-GPI antigen (human origin) is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Then, a monoclonal anti-human IgG or IgM conjugated to horseradish peroxidase is added to each well. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of TMB substrate (3,3',5,5' – tetramethylbenzidine). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution HCl (1 M) serves to stop the enzymatic reaction.
- After stopping the reaction by HCl (1 M) the optical density is read by a spectrophotometer at 450nm
- The rate of color formation from the chromogen is a function of the amount of conjugate bound to the antigen-antibody complex and this is proportional to the initial concentration of the respective antibodies in the patient sample.

## SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum. Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

## STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

## REAGENTS

96 wells microplate coated with highly purified native human $\beta$ 2-glycoprotein I (12x8 individual breakaway microwells) <span style="float: right;">MP</span>	12 strips
1 vial of Cut-off Calibrator (capped blue: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use <span style="float: right;">CONTROL REF</span>	1 x 1.5mL
6 vials of Calibrators (color increasing with concentration: yellow solution) 0, 3, 10, 30, 100, 300 U/mL containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use <span style="float: right;">CAL</span>	6 x 1.5mL
1 vial of Positive control (capped red: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use <span style="float: right;">CONTROL +</span>	1 x 1.5mL
1 vial of Negative control (capped green: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use <span style="float: right;">CONTROL -</span>	1 x 1.5mL
1 vial of anti-human IgG-HRP Conjugate (capped blue: blue solution) Ready to use <span style="float: right;">CONJ IgG</span>	1 x 15mL
1 vial of anti-human IgM-HRP Conjugate (capped green: green solution) Ready to use <span style="float: right;">CONJ IgM</span>	1 x 15mL
1 vial of TMB Substrate (capped black) containing: stabilized TMB/ H <sub>2</sub> O <sub>2</sub> Ready to use <span style="float: right;">SUBS TMB</span>	1 x 15mL
1 vial of Sample Diluent - 5x concentrated (capped white: yellow solution) containing: Tris, NaCl, BSA, sodium azide < 0.1% (preservative) To dilute <span style="float: right;">DIL SPE 5x</span>	1 x 20mL
1 vial of Wash Buffer - 50x concentrated (capped white: green solution) containing Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative) To dilute <span style="float: right;">BUF WASH 50x</span>	1 x 20mL
1 vial of Stop Solution (capped white: colorless solution) containing: 1M Hydrochloric Acid Ready to use <span style="float: right;">SOLN STOP</span>	1 x 15mL

**ELISIS  $\beta$ 2-GPI IgG/IgM kit** ([bmd](#)) is calibrated in arbitrary units (AU/mL).

**ELISIS  $\beta$ 2-GPI IgG/IgM kit** ([bmd](#)) is standardized using the Sapporo-Standards: HCAL for IgG and EY2C9 for IgM.

## MATERIAL REQUIRED BUT NOT PROVIDED

- Purified water (according to the definition of the United States Pharmacopeia and the European Pharmacopeia)
- Microtiter plate reader 450nm reading filter and optional 620nm reference filter (600-690nm)
- Glass ware (cylinder 100-1000mL)
- Test tubes for serum dilutions
- Vortex mixer
- Precision pipettes (10, 100, 200, 500, 1000 $\mu$ L) or adjustable multipipette (100-1000 $\mu$ L)
- Microplate washing device (300 $\mu$ L repeating or multichannel pipette or automated System)
- Adsorbent paper

## SETUP

All the reagents should be prepared as required:

### 1. Preparation of Sample diluent solution

- Dilute concentrated Sample diluent 1/5 in distilled water. Ex: 20mL plus 80mL
- Storage period: 1 month at +2°C/+8°C.

### 2. Preparation of wash buffer

- Dilute concentrated Wash buffer 1/50 in distilled water. Ex: 20mL plus 980mL
- Prepare 20mL of diluted wash buffer for 8 wells or 200mL for 96 wells.
- Storage period: 1 month at +2°C/+8°C.

### 3. Preparation of samples

- Dilute to 1/101 in Sample diluent (1x) Ex: 1000 $\mu$ L sample buffer (1x) + 10 $\mu$ L serum
- Shake vigorously at the vortex.

### 4. Washing Step

- ❖ Automated washing:
  - Consider excess volumes required for setting up the instrument and dead volume of robot pipette.
- ❖ Manual washing:
  - Discard liquid from wells by inverting the plate.
  - Knock the microwell frame with wells downside vigorously on clean adsorbent paper.
  - Pipette 300 $\mu$ L of diluted wash buffer into each well, wait for 20 seconds.
  - Repeat the whole procedure twice again.

### 5. Microplates

- Calculate the number of wells required for the test.
- Remove unused wells from the frame, replace and store in the provided plastic bag, together with desiccant, seal tightly (+2°C/+8°C).

## PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C / +25°C) before handling.

Check that all plates are well drained after each wash.

Do not smoke, eat or drink when manipulating the kit.

Do not pipette by mouth.

*Human sources for the preparation of calibrators and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.*

This kit contains potentially hazardous components. Though kit reagents are not classified being irritant to eyes and skin we recommend to avoid contact with eyes and skin and wear disposable gloves.

Calibrators, Controls and Buffers contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

**ELISIS β2-GPI IgG/IgM kit** ([bmd](#)) has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

**ELISIS β2-GPI IgG/IgM kit** ([bmd](#)) has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.

**Incubation: We recommend test performance at +30°C for automated Systems.**

Never expose components to higher temperature than +37°C. Always pipette Substrate and Conjugate with brand new tips only.

Protect Substrate from light.

## METHOD

### 1. Preparing the test

Use the work sheet to note the sample locations.

Detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

**For quantitative interpretation use calibrators to establish a standard curve**

Set out in duplicate:

- Calibrators (Cal A, Cal B, Cal C, Cal D, Cal E and Cal F)
- Positive control (Ctrl +)
- Negative control (Ctrl -)
- Diluted samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
A	Cal A	Cal E	S1	S9								
B	Cal A	Cal E	S2	S10								
C	Cal B	Cal F	S3	S11								
D	Cal B	Cal F	S4	S12								
E	Cal C	Ctrl +	S5	S13								
F	Cal C	Ctrl +	S6	S14								
G	Cal D	Ctrl -	S7									
H	Cal D	Ctrl -	S8									

**For qualitative interpretation use cut-off calibrators**

Set out in duplicate:

- Negative control (Ctrl -)
- Cut-off calibrator (Ctrl Ref)
- Positive control (Ctrl +)
- Diluted samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
A	Ctrl -	S3	S11	S19								
B	Ctrl -	S4	S12	S20								
C	Ctrl Ref	S5	S13	S21								
D	Ctrl Ref	S6	S14	S22								
E	Ctrl +	S7	S15	S23								
F	Ctrl +	S8	S16	S24								
G	S1	S9	S17									
H	S2	S10	S18									

### 2. Sample incubation

Add 100µL of calibrators OR cut-off calibrator and controls into the designated microwells

Add 100µL of diluted samples

Incubate for 30 minutes at room temperature (+20°C/+32°C)

*Wash step:*

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing buffer.

### 3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature (+20°C/+32°C)

*Wash step:*

Remove the contents of the wells by rapid inversion

Wash 3 times with 300µL of washing buffer

### 4. Incubation of substrate

Add 100µL TMB substrate into each well

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Protect the wells from intense light

### 5. Stop solution

Add 100µL of Stop solution to each well

Incubate for 5 minutes minimum

Agitate plate carefully for 5 sec.

### 6. Reading

Read the optical density of each well at 450nm (optionally 450/620nm) within 30 minutes.

## INTERPRETATION OF RESULTS

For **quantitative interpretation** establish the standard curve by plotting the *optical density (OD) of each calibrator (y-axis)* with respect to the corresponding concentration values in *AU/mL (x-axis)*. For best results [bmd](#) recommend log (x) / lin (y) coordinates and 4-Parameter Fit. From the OD of each sample, read the corresponding antibody concentrations expressed in *AU/mL*.

Negative Results	Positive Results
≤ 15 AU/mL	> 15 AU/mL

### Example of a standard curve

We recommend pipetting calibrators in parallel for each run.

Calibrators IgG	OD 450/620 nm	CV % (Variation)
0 AU/ml	0.041	1.7
3 AU/ml	0.132	0.0
10 AU/ml	0.280	2.6
30 AU/ml	0.584	2.1
100 AU/ml	1.211	0.0
300 AU/ml	2.042	0.6

### Example of calculation

Patient	Replicate (OD)	Mean (OD)	Result (AU/mL)
P01	0.756/0.757	0.757	45.9
P02	1.344/1.352	1.348	123.3

### Do not use this example for interpreting patients results!

For lot specific data, see enclosed quality control certificate.

Each laboratory should establish its own normal range based upon its own techniques, controls, equipment and patient population according to their own established procedures.

For **qualitative interpretation** read the optical density of the cut-off calibrator and the patient samples. Compare patient's OD with the OD of the cut-off calibrator. All samples which are higher than cut-off are considered positive.

<b>Negative:</b>	OD Patient < OD cut-off
<b>Positive:</b>	OD Patient > OD cut-off

**A definite clinical diagnosis should not be based on the results of the performed test only, but should be made by the physician after all clinical and laboratory findings have been evaluated. The diagnosis is to be verified using different diagnostic methods.**

## LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration may confound the results of this assay.

## QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive and negative controls concentration should fall within the range printed on the quality control certificate. If either reagent control is invalid, the test should be repeated.

[bmd](#) offers an external multiparameter quality control (IMMUNO-TROL II, Cat: HM037). This control should be handled as samples.

## CHARACTERISTICS AND PERFORMANCES OF THE TEST

### Sensitivity, specific affinity and agreement to other methods

**ELISIS β<sub>2</sub>-GPI IgG/IgM kit** ([bmd](#)) has been compared to one predicate ELISA kit.

The study was performed on 115 samples of:

⇒ 104 samples related to a pathological population (APS, SLE and other autoimmune diseases)

⇒ 11 samples selected from healthy donors

### Results

		β <sub>2</sub> -GPI IgG predicate kit		
		+	-	
ELISIS β <sub>2</sub> -GPI IgG	+	30	18	Relative sensibility : 96.8% Relative specificity : 78.6% Accuracy : 83.5%
	-	1	66	

		β <sub>2</sub> -GPI IgM predicate kit		
		+	-	
ELISIS β <sub>2</sub> -GPI IgM	+	14	11	Relative sensibility : 93.3% Relative specificity : 89.0% Accuracy : 89.6%
	-	1	89	

### Test Precision

Estimated for IgG and IgM isotype on three positive samples : strongly, moderate and weakly positive.

ELISIS β <sub>2</sub> -GPI IgG/IgM		ELISIS β <sub>2</sub> -GPI IgG/IgM	
Intra-assay (24 tests in a same assay)		Inter-assay (18 tests in 3 differents assays)	
Titer mean AU/mL	CV (%)	Titer mean AU/mL	CV (%)
228.5	6.6	172.1	6.5
54.3	8.8	70.3	6.9
17.2	5.1	19.6	4.4

### Threshold value

Calculated on 70 samples from « blood donors ».

ELISIS β <sub>2</sub> -LISA IgG/IgM	
<b>Cut-off :</b> Mean value +3SD	15 AU/mL

## REFERENCES

---

### **Abuaf N et al**

Prevalence of anti-  $\beta_2$ -glycoprotein I and anti-prothrombin antibodies in 2174 patients and 107 healthy controls.

*Nouv Rev Fr Hematol* 1995, 37: S53-S56.

### **Alarcon-Segovia D, Cabral A**

The concept and classification of antiphospholipid/cofactor syndromes.

*Lupus* 1996, 5: 364-367.

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Measurement of antiphospholipid antibodies by ELISA using  $\beta_2$ -glycoprotein I as an antigen.

*J Immunol Methods* 1991, 143: 223-229.

### **Boffa MC**

Detection of antiphospholipid-protein antibodies and management of the antiphospholipid syndrome.

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*Sang Thrombose Vaisseaux* 1994, 6: 687-692.

### **Galli M et al**

Anticardiolipin antibodies directed not to cardiolipin but to a plasma protein cofactor.

*Lancet* 1990, 335: 1544-1551.

### **Hunt JE, Krilis SA**

The fifth domain of  $\beta_2$ -glycoprotein I contains a phospholipid binding site (Cys 281-Cys 288) and a region recognized by anticardiolipin antibodies.

*J Immunol* 1994, 152: 653-659.

### **Martinuzzo ME et al**

Anti- $\beta_2$ -glycoprotein I antibodies: detection and association with thrombosis.

*BJ Haematol* 1995, 89: 397-402.

### **Matsuura E et al**

Anticardiolipin cofactor (s) and differential diagnosis of autoimmune disease.

*Lancet* 1990, 336: 177-178.

### **Matsuura E et al**

Anticardiolipin antibodies recognize  $\beta_2$ -glycoprotein I structure altered by interacting with an oxygen modified solid phase surface.

*J Exp Med* 1994, 179: 457-462.

### **McNeil HP et al**

Antiphospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation;  $\beta_2$ -glycoprotein I (apolipoprotein H).

*Proc Natl Acad Sci USA* 1990; 87: 4120-4124.

### **Miyakis S. et al**

International consensus statement on an update of the classification criteria for definite anti-phospholipid syndrome (APS).

*J. Thromb. Haemost* 2005, 3: 1-12.

### **Sorice M et al**

Anticardiolipin and anti  $\beta_2$ -GPI are two distinct populations of autoantibodies.

*Thromb Haemost* 1996, 75: 393-398.

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Anticardiolipin autoantibodies recognize  $\beta_2$ -glycoprotein I in the absence of phospholipid. Importance of antigen density and bivalent binding.



*J Immunol* 1995, 154: 954-960.

## SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS		INCUBATION CONDITIONS
		QUANTITATIVE METHOD	QUALITATIVE METHOD	
INCUBATION OF SAMPLES	100µL	Calibrators A, B, C, D, E and F (in duplicate)	Negative Control (in duplicate)	30 min at RT
	100µL	Positive Control (in duplicate)	Cut-off calibrator (in duplicate)	
	100µL	Negative Control (in duplicate)	Positive Control (in duplicate)	
	100µL	Diluted samples (in duplicate)	Diluted samples (in duplicate)	
WASHING	Wash 3 times with washing buffer (300µL/ well)			
INCUBATION OF CONJUGATE	100µL	Anti-human IgG-HRP conjugate OR Anti-human IgM-HRP conjugate		30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)			
ENZYME REACTION	100µL	TMB Substrate		30 min at RT protected from intense light
STOP REACTION	100µL	Stop solution (1M HCl)		5 min

RT: room temperature (+20°C/+32°C)

## SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Sample Diluent		Stop Solution		Cut-off calibrator
	Wash Buffer		Calibrator		Microplate
	Control		Substrate		Conjugate

***BioMédical Diagnostics SA***

**Office**  
Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)