ELISIS P-Ethanolamine IgG / IgM

REF LIS 1209 

DÉFINITION

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BIBLIOGRAPHIE

ALARCON-SEGOVIA D, CABRAL A.R.

Functional and immunochemical heterogeneity of antiphospholipid antibodies : A classification. *J. Rheum.*, 19, 1166-1169, 1992.

ASHERSON R.A., KHAMASTA M.A., ROS J.

The "primary" antiphospholipid syndrome : major clinical and serological features. *Med.*, 68, 366-373, 1989.

GALLI M., BARBUI T., CONFURIUS P. ZWALL R.F.A., BEVERS E.M.

Anticoagulant mechanisms of antiphospholipids antibodies. Lupus Anticoagulants : Recent advances in diagnosis and understanding *ISTH Congress, Satellite Symposium- NEW-YORK 1993*

HARRIS E.N.

Annotation - Antiphospholipid antibodies. *Br. J. Haematol.*, 74, 1-9, 1990

LAROCHE P, BERARD M, ROUQUETTE AM, DESGRUELLE, C. BOFFA MC.

Advantage of using both anionic and switterionic phospholipid antigens for the detection of antiphospholipid antibodies. *Am J Clin Pathol* 106, 549-554, 1996

LOIZOU S., MACKWORTH-YOUNG C.G., COFINER C., WALPORT M.J.

Heterogeneity of binding reactivity to different phospholipids of antibodies from patients with systemic lupus erythematosus. *Clin. Exp. Immunol.*, 80, 171-176, 1990

Mc INTYRE J.A., WAGENKNECHT D.R., SUGI T.

Role of lipoprotein cofactors in the detection of antiphospholipid antibodies. Lupus Anticoagulants : Recent Advances in diagnosis and understanding – *ISTH Congress, Satellite Symposium- NEW-YORK 1993*

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TRIPLETT D.A., BRANDT J.T., MUSGRAVE K.A., ORR C.A.

The relationship between lupus anticoagulants and antibodies to phospholipids. *J. Am. Med. Ass.*, 259, 550-554, 1988.

BOEY, M.L., COLACO, C.B., GHARAVI, A.E., and al.

Thrombosis in systemic lupus erythematosus: striking association with the presence of circulating lupus anticoagulant. *Br. Med. J.* 287: 1021-1023, 1983.

GASTINEAU, D.A., KAZMIER, F.J., NICHOLS, W.L., BOWIE, E.J.

Lupus anticoagulant: an analysis of the clinical and laboratory features of 219 cases. *Am. J. Hematol.* 19:265-267, 1985.

McNEIL HP, SIMPSON RJ, CHESTERMAN CN, KIRILIS SA

Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β 2-Glycoprotein I (apolipoprotein H). *Proc Natl Acad Sci USA* 87: 4120-4124, 1990.

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Clinical relevance of antibodies against different phospholipids. *Journal of Autoimmunity* 15: A60, 2000.

E. BALADA, J. ORDI-ROS, F. PAREDES, J. VILLARREAL, M. MAURI, M. VILARDELL-TARRES

Anti-phosphatidylethanolamine antibodies contribute to the diagnosis of antiphospholipid syndrome in patients with systemic lupus erythematosus. *Scand J Rheumatol* 30: 235-241; 2001.

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


REF LIS 1209

DEFINITION

ELISIS P-Ethanolamine IgG/IgM kit ([bmd](#)) is an enzyme linked immunoassay (ELISA) for the quantitative and qualitative detection of IgG and/or IgM antibodies directed against phosphatidyl-ethanolamin in human serum samples.

Anti-phosphatidyl-ethanolamin antibodies recognize specific epitopes on a complex composed out of phosphatidyl-ethanolamin and β_2 -glycoprotein I which are only present when these two molecules interact.

The results of **ELISIS P-Ethanolamine IgG/IgM** ([bmd](#)) are to be used in conjunction with the clinical findings and the other laboratory tests to aid in the diagnosis and risk estimation of thrombosis in patients with systemic lupus erythematosus (SLE) and the antiphospholipid syndrome (APS).

ELISIS P-Ethanolamine IgG/IgM kit ([bmd](#)) may be used with the **CARIS™ system** (diluting and dispensing device).

DIAGNOSTIC VALUE

Phospholipids are components of biological membranes. Antibodies against phosphatidyl-ethanolamin, an phospholipid derived from glycerol, belong to the group of anti-phospholipid antibodies specific for phospholipids such as Cardiolipin, Phosphatidyl-serine, Phosphatidyl-inositol, Phosphatidyl-choline, Sphingomyelin and phosphatidic acid.

Anti-phospholipid antibodies are frequently found in sera of patients with systemic lupus erythematosus (SLE) and related diseases. The occurrence of anti-phospholipid antibodies in patients with SLE and related diseases is typical for a secondary anti-phospholipid syndrome (APS). In contrast, anti-phospholipid antibodies in patients with no other autoimmune diseases characterize the primary APS.

Many studies have shown a correlation between these autoantibodies and an enhanced incidence of thrombosis, thrombocytopenia and habitual abortions (as a consequence of placental infarct). The exact mechanisms by which pathogenic anti-phospholipid antibodies induce thrombosis is not yet revealed fully.

Reports of anti-phosphatidyl-ethanolamin antibodies with similar or identical pathogenic associations as those described for anti-cardiolipin and anti-phosphatidylserine have repeatedly been described in the literature. PE is a zwitterionic phospholipid present in both internal and external sides of cell membranes, thought to be involved in the development of thrombotic events via its role in the protein C pathway and the inactivation of factor Va by activated protein C. Antibodies directed to PE could play an important role in the pathogenesis of the APS by inhibiting activated protein C. In some instances, anti-PE antibodies are the sole antiphospholipid antibodies observed thus being an important diagnostic tool in addition to cardiolipin antibodies.

ASSAY PRINCIPLE

Highly purified phosphatidyl-ethanolamin bound to the native human β_2 -glycoprotein I is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Then, a monoclonal anti-human IgG or IgM conjugated to horseradish peroxidase is added to each well. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of TMB substrate (3,3',5,5' – tetramethylbenzidine). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution HCl (1 M) serves to stop the enzymatic reaction.
- After stopping the reaction by HCl (1 M) the optical density is read by a spectrophotometer at 450nm
- The rate of color formation from the chromogen is a function of the amount of conjugate bound to the antigen-antibody complex and this is proportional to the initial concentration of the respective antibodies in the patient sample.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum. Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

REAGENTS

96 wells microplate individually coated with highly purified phosphatidyl-ethanolamin (bovine origin) and native human β_2 -glycoprotein I (12x8 individual breakaway microwells) MP	12 strips
1 vial of Cut-off Calibrator (capped blue: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CONTROL REF	1 x 1.5mL
6 vials of Calibrators (color increasing with concentration: yellow solution) 0, 3, 10, 30, 100, 300 U/mL containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CAL	6 x 1.5mL
1 vial of Negative control (capped green: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CONTROL -	1 x 1.5mL
1 vial of Positive control (capped red: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CONTROL +	1 x 1.5mL
1 vial of anti-human IgG-HRP Conjugate (capped blue: blue solution) <u>Ready to use</u> CONJ IgG	1 x 15mL
1 vial of anti-human IgM-HRP Conjugate (capped green: green solution) <u>Ready to use</u> CONJ IgM	1 x 15mL
1 vial of TMB Substrate (capped black) containing: stabilized TMB/ H ₂ O ₂ <u>Ready to use</u> SUBS TMB	1 x 15mL
1 vial of Sample Diluent - 5x concentrated (capped white: yellow solution) containing: Tris, NaCl, BSA, sodium azide < 0.1% (preservative) <u>To dilute</u> DIL SPE 5x	1 x 20mL
1 vial of Wash Buffer - 50x concentrated (capped white: green solution) containing Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative) <u>To dilute</u> BUF WASH 50x	1 x 20mL
1 vial of Stop Solution (capped white: colorless solution) containing: 1M Hydrochloric Acid <u>Ready to use</u> SOLN STOP	1 x 15mL

ELISIS P-Ethanolamine IgG/IgM kit ([bmd](#)) is calibrated in arbitrary units (AU/mL).

MATERIAL REQUIRED BUT NOT PROVIDED

- Purified water (according to the definition of the United States Pharmacopeia and the European Pharmacopeia)
- Microtiter plate reader 450nm reading filter and optional 620nm reference filter (600-690nm)
- Glass ware (cylinder 100-1000mL)
- Test tubes for serum dilutions
- Vortex mixer
- Precision pipettes (10, 100, 200, 500, 1000 μ L) or adjustable multipipette (100-1000 μ L)
- Microplate washing device (300 μ L repeating or multichannel pipette or automated System)
- Adsorbent paper

SETUP

All the reagents should be prepared as required:

1. Preparation of Sample diluent solution

- Dilute concentrated Sample diluent 1/5 in distilled water.
Ex: 20mL plus 80mL
- Storage period: 1 month at +2°C/+8°C.

2. Preparation of wash buffer

- Dilute concentrated Wash buffer 1/50 in distilled water.
Ex: 20mL plus 980mL
- Prepare 20mL of diluted wash buffer for 8 wells or 200mL for 96 wells.
- Storage period: 1 month at +2°C/+8°C.

3. Preparation of samples

- Dilute to 1/101 in Sample diluent (1x)
Ex: 1000 μ L sample buffer (1x) + 10 μ L serum
- Shake vigorously at the vortex.

4. Washing Step

- ❖ Automated washing:
 - Consider excess volumes required for setting up the instrument and dead volume of robot pipette.
- ❖ Manual washing:
 - Discard liquid from wells by inverting the plate.
 - Knock the microwell frame with wells downside vigorously on clean adsorbent paper.
 - Pipette 300 μ l of diluted wash buffer into each well, wait for 20 seconds.
 - Repeat the whole procedure twice again.

5. Microplates

- Calculate the number of wells required for the test.
- Remove unused wells from the frame, replace and store in the provided plastic bag, together with desiccant, seal tightly (+2°C/+8°C).

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C / +25°C) before handling.

Check that all plates are well drained after each wash.

Do not smoke, eat or drink when manipulating the kit.

Do not pipette by mouth.

Human sources for the preparation of calibrators and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

This kit contains potentially hazardous components. Though kit reagents are not classified being irritant to eyes and skin we recommend to avoid contact with eyes and skin and wear disposable gloves.

Calibrators, Controls and Buffers contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

ELISIS P-Ethanolamine IgG/IgM kit ([bmd](#)) has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

ELISIS P-Ethanolamine IgG/IgM ([bmd](#)) has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.

Incubation: We recommend test performance at +30°C for automated Systems.

Never expose components to higher temperature than +37°C. Always pipette Substrate and Conjugate with brand new tips only.

Protect Substrate from light.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

For quantitative interpretation use calibrators to establish a standard curve

Set out in duplicate:

- Calibrators (Cal A, Cal B, Cal C, Cal D, Cal E and Cal F)
- Positive control (Ctrl +)
- Negative control (Ctrl -)
- Diluted samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
A	Cal A	Cal E	S1	S9								
B	Cal A	Cal E	S2	S10								
C	Cal B	Cal F	S3	S11								
D	Cal B	Cal F	S4	S12								
E	Cal	Ctrl+	S5	S13								
F	Cal C	Ctrl+	S6	S14								
G	Cal D	Ctrl-	S7									
H	Cal D	Ctrl-	S8									

For qualitative interpretation use cut-off calibrators

Set out in duplicate:

- Negative control (Ctrl -)
- Cut-off calibrator (Ctrl Ref)
- Positive control (Ctrl +)
- Diluted samples (S)

	1	2	3	4	5	6	7	8	9	10	11	12
A	Ctrl-	S3	S11	S19								
B	Ctrl-	S4	S12	S20								
C	Ctrl Ref	S5	S13	S21								
D	Ctrl Ref	S6	S14	S22								
E	Ctrl+	S7	S15	S23								
F	Ctrl+	S8	S16	S24								
G	S1	S9	S17									
H	S2	S10	S18									

Note: If IgG and IgM are determined in parallel calibrators, controls and samples have to be done for each subclass separately.

2. Sample incubation

Add 100µL of calibrators OR cut-off calibrator and controls into the designated microwells

Add 100µL of diluted samples

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing buffer.

3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 300µL of washing buffer

4. Incubation of substrate

Add 100µL TMB substrate into each well

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Protect the wells from intense light

5. Stop solution

Add 100µL of Stop solution to each well

Incubate for 5 minutes minimum

Agitate plate carefully for 5 sec.

6. Reading

Read the optical density of each well at 450nm (optionally 450/620nm) within 30 minutes.

INTERPRETATION OF RESULTS

For **quantitative interpretation** establish the standard curve by plotting the **optical density (OD) of each calibrator (y-axis)** with respect to the corresponding concentration values in **AU/mL (x-axis)**. For best results [bmd](#) recommend log (x) / lin (y) coordinates and 4-Parameter Fit. From the OD of each sample, read the corresponding antibody concentrations expressed in **AU/mL**.

Negative Results	Positive Results
≤ 15 AU/mL	> 15 AU/mL

Example of a standard curve

We recommend pipetting calibrators in parallel for each run.

Calibrators IgG/IgM	OD 450/620 nm	CV % (Variation)
0 AU/ml	0.039	3.2
3 AU/ml	0.172	5.4
10 AU/ml	0.330	3.2
30 AU/ml	0.663	0.6
100 AU/ml	1.302	2.0
300 AU/ml	2.115	0.6

Example of calculation

Patient	Replicate (OD)	Mean (OD)	Result (AU/mL)
P01	0.917/0.910	0.914	52.0
P02	0.443/0.454	0.449	15.8

Do not use this example for interpreting patients results!.

For lot specific data, see enclosed quality control certificate.

Each laboratory should establish its own normal range based upon its own techniques, controls, equipment and patient population according to their own established procedures.

For **qualitative interpretation** read the optical density of the cut-off calibrator and the patient samples. Compare patient's OD with the OD of the cut-off calibrator. All samples which are higher than cut-off are considered positive.

Negative:	OD Patient < OD cut-off
Positive:	OD Patient > OD cut-off

A definite clinical diagnosis should not be based on the results of the performed test only, but should be made by the physician after all clinical and laboratory findings have been evaluated. The diagnosis is to be verified using different diagnostic methods.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration may confound the results of this assay.

QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive and negative controls concentration should fall within the range printed on the quality control certificate. If either reagent control is invalid, the test should be repeated.

REFERENCES

ALARCON-SEGOVIA D, CABRAL A.R.

Functional and immunochemical heterogeneity of antiphospholipid antibodies : A classification.
J. Rheum., 19, 1166-1169, 1992.

ASHERSON R.A., KHAMASTA M.A., ROS J.

The "primary" antiphospholipid syndrome : major clinical and serological features.
Med., 68, 366-373, 1989.

GALLI M., BARBUI T., CONFURIUS P. ZWALL R.F.A., BEVERS E.M.

Anticoagulant mechanisms of antiphospholipids antibodies. Lupus Anticoagulants : Recent advances in diagnosis and understanding
ISTH Congress, Satellite Symposium- NEW-YORK 1993

HARRIS E.N.

Annotation - Antiphospholipid antibodies.
Br. J. Haematol., 74, 1-9, 1990

LAROCHE P, BERARD M, ROUQUETTE AM, DESGRUELLE, C. BOFFA MC.

Advantage of using both anionic andwitterionic phospholipid antigens for the detection of antiphospholipid antibodies.
Am J Clin Pathol 106, 549-554, 1996

LOIZOU S., MACKWORTH-YOUNG C.G., COFINER C., WALPORT M.J.

Heterogeneity of binding reactivity to different phospholipids of antibodies from patients with systemic lupus erythematosus.
Clin. Exp. Immunol., 80, 171-176, 1990

Mc INTYRE J.A., WAGENKNECHT D.R., SUGI T.

Role of lipoprotein cofactors in the detection of antiphospholipid antibodies. Lupus Anticoagulants : Recent Advances in diagnosis and understanding –
ISTH Congress, Satellite Symposium- NEW-YORK 1993

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TRIPLETT D.A., BRANDT J.T., MUSGRAVE K.A., ORR C.A.

The relationship between lupus anticoagulants and antibodies to phospholipids.
J. Am. Med. Ass., 259, 550-554, 1988.

BOEY, M.L., COLACO, C.B., GHARAVI, A.E., and al.

Thrombosis in systemic lupus erythematosus: striking association with the presence of circulating lupus anticoagulant.
Br. Med. J. 287: 1021-1023, 1983.

GASTINEAU, D.A., KAZMIER, F.J., NICHOLS, W.L., BOWIE, E.J.

Lupus anticoagulant: an analysis of the clinical and laboratory features of 219 cases.
Am. J. Hematol. 19:265-267, 1985.

McNEIL HP, SIMPSON RJ, CHESTERMAN CN, KIRILIS SA

Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β 2-Glycoprotein I (apolipoprotein H).
Proc Natl Acad Sci USA 87: 4120-4124, 1990.

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Clinical relevance of antibodies against different phospholipids.
Journal of Autoimmunity 15: A60, 2000.

E. BALADA, J. ORDI-ROS, F. PAREDES, J. VILLARREAL, M. MAURI, M. VILARDELL-TARRES

Anti-phosphatidylethanolamine antibodies contribute to the diagnosis of antiphospholipid syndrome in patients with systemic lupus erythematosus.
Scand J Rheumatol 30: 235-241; 2001

Mc INTYRE J.A., WAGENKNECHT D.R.




Anti-phosphatidylethanolamine (aPE) antibodies: a survey.
JAutoimmun 15(2): 185-93; 2000.

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS		INCUBATION CONDITIONS
		QUANTITATIVE METHOD	QUALITATIVE METHOD	
INCUBATION OF SAMPLES	100µL	Calibrators A, B, C, D, E and F (in duplicate)	Negative Control (in duplicate)	30 min at RT
	100µL	Positive Control (in duplicate)	Cut-off calibrator (in duplicate)	
	100µL	Negative Control (in duplicate)	Positive Control (in duplicate)	
	100µL	Diluted samples (in duplicate)	Diluted samples (in duplicate)	
WASHING	Wash 3 times with washing buffer (300µL/ well)			
INCUBATION OF CONJUGATE	100µL	Anti-human IgG-HRP conjugate OR Anti-human IgM-HRP conjugate		30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)			
ENZYME REACTION	100µL	TMB Substrate		30 min at RT protected from intense light
STOP REACTION	100µL	Stop solution (1M HCl)		5 min

RT: room temperature (+20°C/+32°C)

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Sample Diluent		Stop Solution		Cut-off calibrator
	Wash Buffer		Calibrator		Microplate
	Control		Substrate		Conjugate

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com