ELISIS APS Profile IgG / IgM

REF LIS 1234 

DÉFINITION

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A new anti-phospholipid-antibody ELISA - sensitivity and specificity for cerebrovascular insults in autoimmune diseases.

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



REF LIS 1234

DEFINITION

ELISIS APS Profile IgG/IgM kit ([bmd](#)) is an enzyme linked immunoassay (ELISA) for the individual qualitative IgG and/or IgM screening of 8 antibodies against Prothrombin, Thrombin, Cardiolipin, Phosphatidyl-serine, Phosphatidyl-inositol, Phosphatidyl-ethanolamin, Phosphatidyl-choline and Sphingomyelin in human serum samples.

The results of **ELISIS APS Profile IgG/IgM** ([bmd](#)) are to be used in conjunction with the clinical findings and the other laboratory tests to aid in the diagnosis and risk estimation of thrombosis in patients with systemic lupus erythematosus (SLE).

ELISIS APS Profile IgG/IgM kit ([bmd](#)) may be used with the **CARIS™ system** (diluting and dispensing device).

DIAGNOSTIC VALUE

Antibodies against phospholipids, components of the biological membranes, are specific for phospholipids such as Cardiolipin, Phosphatidyl -inositol, -ethanolamine, -choline and Sphingomyelin.

Anti-phospholipid antibodies are frequently found in sera of patients with systemic lupus erythematosus (SLE) and related diseases. The occurrence of anti-phospholipid antibodies in patients with SLE and related diseases is typical for a secondary anti-phospholipid syndrome (APS). In contrast, anti-phospholipid antibodies in patients with no other autoimmune diseases characterize the primary APS.

Many studies have shown a correlation between these autoantibodies and an enhanced incidence of thrombosis, thrombocytopenia and habitual abortions (as a consequence of placental infarct). The exact mechanisms by which pathogenic anti-phospholipid antibodies induce thrombosis is not yet revealed fully.

Antibodies targeting prothrombin alone have recently been described being highly associated with fetal loss in APS patients. Prothrombin antibodies are the the first marker for this serious complication as all other known antiphospholipid antibodies failed to correlate with fetal loss.

The presence of anti-thrombin antibodies has been associated with the clinical features of the antiphospholipid syndrome (APS). Anti-Thrombin antibodies seem to correlate with deep venous thrombosis.

ASSAY PRINCIPLE

8 highly purified phospholipids (Prothrombin, Thrombin, Cardiolipin, Phosphatidyl-serine, Phosphatidyl-inositol, Phosphatidyl-ethanolamin, Phosphatidyl-choline and Sphingomyelin) are individually coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Then, a monoclonal anti-human IgG or IgM conjugated to horseradish peroxidase is added to each well. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of TMB substrate (3,3',5,5' – tetramethylbenzidine). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution HCl (1 M) serves to stop the enzymatic reaction.
- After stopping the reaction by HCl (1 M) the optical density is read by a spectrophotometer at 450nm
- The rate of color formation from the chromogen is a function of the amount of conjugate bound to the antigen-antibody complex and this is proportional to the initial concentration of the respective antibodies in the patient sample.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum. Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

REAGENTS

96 wells microplate individually coated with 8 highly purified phospholipids (Prothrombin, Thrombin, Cardiolipin, Phosphatidyl-serine, Phosphatidyl-inositol, Phosphatidyl-ethanolamin, Phosphatidyl-choline and Sphingomyelin) <input type="checkbox"/> MP (12x8 individual breakaway microwells)	12 strips
8 vials of Cut-off Calibrator (capped white: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <input type="checkbox"/> CONTROL <input type="checkbox"/> REF Ready to use	8 x 1.5mL
1 vial of Negative control (capped green: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) Ready to use <input type="checkbox"/> CONTROL <input type="checkbox"/> -	2 x 1.5mL
1 vial of anti-human IgG-HRP Conjugate (capped blue: blue solution) Ready to use <input type="checkbox"/> CONJ <input type="checkbox"/> IgG	1 x 15mL
1 vial of anti-human IgM-HRP Conjugate (capped green: green solution) Ready to use <input type="checkbox"/> CONJ <input type="checkbox"/> IgM	1 x 15mL
1 vial of TMB Substrate (capped black) containing: stabilized TMB/ H ₂ O ₂ Ready to use <input type="checkbox"/> SUBS <input type="checkbox"/> TMB	1 x 15mL
1 vial of Sample Diluent - 5x concentrated (capped white: yellow solution) containing: Tris, NaCl, BSA, sodium azide < 0.1% (preservative) To dilute <input type="checkbox"/> DIL <input type="checkbox"/> SPE <input type="checkbox"/> 5x	1 x 20mL
1 vial of Wash Buffer - 50x concentrated (capped white: green solution) containing Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative) To dilute <input type="checkbox"/> BUF <input type="checkbox"/> WASH <input type="checkbox"/> 50x	1 x 20mL
1 vial of Stop Solution (capped white: colorless solution) containing: 1M Hydrochloric Acid Ready to use <input type="checkbox"/> SOLN <input type="checkbox"/> STOP	1 x 15mL

MATERIAL REQUIRED BUT NOT PROVIDED

- Purified water (according to the definition of the United States Pharmacopeia and the European Pharmacopeia)
- Microtiter plate reader 450nm reading filter and optional 620nm reference filter (600-690nm)
- Glass ware (cylinder 100-1000mL)
- Test tubes for serum dilutions
- Vortex mixer
- Precision pipettes (10, 100, 200, 500, 1000µL) or adjustable multipipette (100-1000µL)
- Microplate washing device (300µL repeating or multichannel pipette or automated System)
- Adsorbent paper

SETUP

All the reagents should be prepared as required:

1. Preparation of Sample diluent solution

- Dilute concentrated Sample diluent 1/5 in distilled water.
Ex: 20mL plus 80mL
- Storage period: 1 month at +2°C/+8°C.

2. Preparation of wash buffer

- Dilute concentrated Wash buffer 1/50 in distilled water.
Ex: 20mL plus 980mL
- Prepare 20mL of diluted wash buffer for 8 wells or 200mL for 96 wells.
- Storage period: 1 month at +2°C/+8°C.

3. Preparation of samples

- Dilute to 1/101 in Sample diluent (1x)
Ex: 1000µL sample buffer (1x) + 10µL serum
- Shake vigorously at the vortex.

4. Washing Step

- ❖ Automated washing:
 - Consider excess volumes required for setting up the instrument and dead volume of robot pipette.
- ❖ Manual washing:
 - Discard liquid from wells by inverting the plate.
 - Knock the microwell frame with wells downside vigorously on clean adsorbent paper.
 - Pipette 300µl of diluted wash buffer into each well, wait for 20 seconds.
 - Repeat the whole procedure twice again.

5. Microplates

- Calculate the number of wells required for the test.
- Remove unused wells from the frame, replace and store in the provided plastic bag, together with desiccant, seal tightly (+2°C/+8°C).

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C / +25°C) before handling.

Check that all plates are well drained after each wash.

Do not smoke, eat or drink when manipulating the kit.

Do not pipette by mouth.

Human sources for the preparation of calibrator and control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

This kit contains potentially hazardous components. Though kit reagents are not classified being irritant to eyes and skin we recommend to avoid contact with eyes and skin and wear disposable gloves.

Calibrator, Control and Buffers contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

ELISIS APS Profile IgG/IgM kit ([bmd](#)) has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

ELISIS APS Profile IgG/IgM kit ([bmd](#)) has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.

Incubation: We recommend test performance at +30°C for automated Systems.

Never expose components to higher temperature than +37°C. Always pipette Substrate and Conjugate with brand new tips only.

Protect Substrate from light.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

For qualitative interpretation use cut-off calibrators

Set out in duplicate:

- Cut-off calibrator (Ctrl Ref)
- Negative control (Ctrl -)
- Diluted samples (S)

2. Sample incubation

Add 100µL of cut-off calibrator and control into the designated microwells

Add 100µL of diluted samples

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing buffer.

3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 300µL of washing buffer

4. Incubation of substrate

Add 100µL TMB substrate into each well

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Protect the wells from intense light

5. Stop solution

Add 100µL of Stop solution to each well

Incubate for 5 minutes minimum

Agitate plate carefully for 5 sec.

6. Reading

Read the optical density of each well at 450nm (optionally 450/620nm) within 30 minutes.

INTERPRETATION OF RESULTS

Read the optical density of the cut-off control and the patient samples. Multiply the OD of the cut-off calibrator by the parameter specific factor, provided with the lot specific QC certificate. Compare patient ODs with the calculated parameter OD cut-off value. All samples which are higher than cut-off are considered positive.

OD Cut-off calibrator x antigen specific factor = OD Cut-off parameter.

Negative: OD Patient < OD cut-off

Positive: OD Patient > OD cut-off

Calibrators	OD 450/620 nm
Negative Control	0.034
Cut-off Calibrator	0.450

Example of interpretation

We recommend pipetting cut-off calibrator in parallel for each run.

Cardiolipin Factor= 1.0

OD Cut-off calibrator 0.450 x 1.0 = 0.550

OD patient 0.999 > 0.450 = positive

Do not use this example for interpreting patients results!

[bmd](#) recommend to retest samples that are equivocal.

For lot specific data, see enclosed quality control leaflet.

Each laboratory should establish its own normal range based upon its own techniques, controls, equipment and patient population according to their own established procedures.

For **semi-quantification of the results**, each patient-OD value can be expressed by the Index-Value. The Index-Value is calculated by dividing the patient-OD by the cut-off OD:

$$\text{Index Value} = \frac{\text{OD (patient sample)}}{\text{OD (cut-off calibrator)}}$$

Negative => Index Value < 1.0

Positive => Index Value > 1.0

A definite clinical diagnosis should not be based on the results of the performed test only, but should be made by the physician after all clinical and laboratory findings have been evaluated.

The diagnosis is to be verified using different diagnostic methods.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration may confound the results of this assay.

QUALITY CONTROL

The kit contains a negative control. The negative control should be included with each run of the test.

[bmd](#) offers external multiparameter quality controls (IMMUNO-TROL I, Cat: HM036; IMMUNO-TROL II, Cat: HM037). These controls should be handled as samples.

REFERENCES

Boey, M.L., Colaco, C.B., Gharavi, A.E., et al. (1983):

Thrombosis in systemic lupus erythematosus: striking association with the presence of circulating lupus anticoagulant. Br. Med. J. 287: 1021-1023.

Gastineau, D.A., Kazmier, F.J., Nichols, W.L., Bowie, E.J. (1985):

Lupus anticoagulant: an analysis of the clinical and laboratory features of 219 cases. Am. J. Hematol. 19:265-267.

P. von Landenberg, T. Matthias, J. Zaech, M. Schultz, M. Lorber, M. Blank, Y. Shoenfeld (2003).

Antiprothrombin antibodies are associated with pregnancy loss in patients with the antiphospholipid syndrome.

Am J Reprod Immunol 49: 51-56.

Wöhrle R, Matthias T, von Landenberg P, Oppermann M, Helmke K, Förger F (2000):

Clinical relevance of antibodies against different phospholipids.

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



Journal of Autoimmunity Vol 15: A 60.

SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	100µL 100µL 100µL	Cut-off calibrator (in duplicate) Negative Control (in duplicate) Diluted samples (in duplicate)	30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)		
INCUBATION OF CONJUGATE	100µL	Anti-human IgG-HRP conjugate OR Anti-human IgM-HRP conjugate	30 min at RT
WASHING	Wash 3 times with washing buffer (300µL/ well)		
ENZYME REACTION	100µL	TMB Substrate	30 min at RT protected from intense light
STOP REACTION	100µL	Stop solution (1M HCl)	5 min

RT: room temperature (+20°C/+32°C)

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Sample Diluent		Stop Solution		Cut-off calibrator
	Wash Buffer		Calibrator		Microplate
	Control		Substrate		Conjugate

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com