ELISIS CELIAC GliA IgA

REF LIS 1501

DÉFINITION

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REF LIS 1501

DEFINITION

ELISIS CELIAC Glia IgA kit ([bmd](#)) is an enzyme linked immunoassay (ELISA) for the qualitative or quantitative determination of IgA anti-gliadin (Glia) antibodies in human serum samples.

The results of the **ELISIS CELIAC Glia IgA kit** ([bmd](#)) are to be used in conjunction with the clinical findings and the other laboratory tests to aid in the diagnosis and monitoring of Celiac disease (gluten-sensitive enteropathy).

ELISIS CELIAC Glia IgA kit ([bmd](#)) may be used with the **CARIST™ system** (diluting and dispensing device).

DIAGNOSTIC VALUE

Celiac disease is an absorption deficiency syndrome characterized by weight loss, abdominal pain, and dietary gluten intolerance. Many affected patients have antibodies directed against gliadin (anti-gliadin), endomysium/tissue transglutaminase (anti-tTG) and reticulin (anti-reticulin). Anti-reticulin antibodies may be considered to be less sensitive than anti-gliadin in the assessment of patients with celiac disease.

Anti-tTG antibodies appear to be more sensitive and specific than anti-gliadin antibodies. Nevertheless, anti-gliadin antibodies may offer a benefit when assessing children under two years of age since affected persons in this age group may develop IgA anti-gliadin before developing IgA anti-endomysium. Anti-gliadin IgG may also be important for use in patients with celiac disease who also suffer from a variety of IgA deficiencies.

After several months of treatment with a gluten free diet, antibody titers frequently decrease; though antibody levels may increase again on discontinuation of the gluten free diet.

The combined anti-gliadin and anti-tTG assay offers adult and children patients an easy, rapid, sensitive and specific technique for the detection and the treatment monitoring of celiac disease.

* The simultaneous presence of IgA anti-gliadin and anti-tTG antibodies, with or without IgG antibodies, will enable to confirm the quasi evidence of a celiac disease diagnosis. It is also an indication that a duodenojejunal biopsy should be performed.

* IgG anti-tTG antibodies are highly reliable marker of Celiac Disease in IgA deficient subjects.

* With reference made to the patient clinical data, the presence of IgA and IgG anti-gliadin antibodies can in some cases confirm the diagnosis of celiac disease.

* The presence of IgG anti-gliadin antibodies alone is not indicative of a celiac disease, nevertheless, it is recommended to take into account an eventual IgA deficiency.

ASSAY PRINCIPLE

Gliadin is coated onto a polystyrene microtiter plate (12 strips of 8 wells).

- First, the diluted sample is added to the coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Then, a monoclonal anti-human IgA conjugated to horseradish peroxidase is added to each well. The conjugate binds to the antigen-antibody complex. After incubation, the excess conjugate is removed by a second wash.
- The chromogenic step is performed by the addition of TMB substrate (3,3',5,5' – tetramethylbenzidine). During this step, a colour will develop in proportion to the amount of antibodies in the sample.
- Addition of Stop Solution HCl (1 M) serves to stop the enzymatic reaction.
- After stopping the reaction by HCl (1 M) the optical density is read by a spectrophotometer at 450nm
- The rate of color formation from the chromogen is a function of the amount of conjugate bound to the antigen-antibody complex and this is proportional to the initial concentration of the respective antibodies in the patient sample.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on serum. Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If determination is not performed immediately, samples should be stored at +2°C/+8°C for no longer than a week or frozen.
- To avoid any non-specific fixation, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

STABILITY AND STORAGE

- Store reagents and assay strips at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

REAGENTS

| | |
|---|-----------|
| 96 wells microplate coated with gliadin (12x8 individual breakaway microwells) MP | 12 strips |
| 1 vial of Cut-off Calibrator (capped blue: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CONTROL REF | 1 x 1.5mL |
| 6 vials of Calibrators (color increasing with concentration: yellow solution) 0, 3, 10, 30, 100, 300 U/mL containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CAL | 6 x 1.5mL |
| 1 vial of Positive control (capped red: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CONTROL + | 1 x 1.5mL |
| 1 vial of Negative control (capped green: yellow solution) containing: human serum (diluted), sodium azide < 0.1% (preservative) <u>Ready to use</u> CONTROL - | 1 x 1.5mL |
| 1 vial of anti-human IgA-HRP Conjugate (capped blue: blue solution) <u>Ready to use</u> CONJ IgA | 1 x 15mL |
| 1 vial of TMB Substrate (capped black) containing: stabilized TMB/ H ₂ O ₂ <u>Ready to use</u> SUBS TMB | 1 x 15mL |
| 1 vial of Sample Diluent - 5x concentrated (capped white: yellow solution) containing: Tris, NaCl, BSA, sodium azide < 0.1% (preservative) <u>To dilute</u> DIL SPE 5x | 1 x 20mL |
| 1 vial of Wash Buffer - 50x concentrated (capped white: green solution) Containing: Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative) <u>To dilute</u> BUF WASH 50x | 1 x 20mL |
| 1 vial of Stop Solution (capped white: colorless solution) containing: 1M Hydrochloric Acid) <u>Ready to use</u> SOLN STOP | 1 x 15mL |

ELISIS CELIAC Glia IgA kit () is calibrated in arbitrary units (AU/mL).

MATERIAL REQUIRED BUT NOT PROVIDED

- Purified water (according to the definition of the United States Pharmacopeia and the European Pharmacopeia)
- Microtiter plate reader 450nm reading filter and optional 620nm reference filter (600-690nm)
- Glass ware (cylinder 100-1000mL)
- Test tubes for serum dilutions
- Vortex mixer
- Precision pipettes (10, 100, 200, 500, 1000µL) or adjustable multipipette (100-1000µL)
- Microplate washing device (300µL repeating or multichannel pipette or automated System)
- Adsorbent paper

SETUP

All the reagents should be prepared as required:

1. Preparation of Sample Diluent Solution

- Dilute concentrated Sample Diluent 1/5 in distilled water.
Ex: 20mL plus 80mL
- Storage period: 1 month at +2°C/+8°C.

2. Preparation of wash buffer

- Dilute concentrated Wash buffer 1/50 in distilled water.
Ex: 20mL plus 980mL
- Prepare 20mL of diluted wash buffer for 8 wells or 200mL for 96 wells.
- Storage period: 1 month at +2°C/+8°C.

3. Preparation of samples

- Dilute to 1/101 in Sample Diluent (1x)
Ex: 1000µL Sample Diluent (1x) + 10µL serum
- Shake vigorously at the vortex.

4. Washing Step

- ❖ Automated washing:
 - Consider excess volumes required for setting up the instrument and dead volume of robot pipette.
- ❖ Manual washing:
 - Discard liquid from wells by inverting the plate.
 - Knock the microwell frame with wells downside vigorously on clean adsorbent paper.
 - Pipette 300µl of diluted wash buffer into each well, wait for 20 seconds.
 - Repeat the whole procedure twice again.

5. Microplates

- Calculate the number of wells required for the test.
- Remove unused wells from the frame, replace and store in the provided plastic bag, together with desiccant, seal tightly (+2°C/+8°C).

PRECAUTIONS

Allow all reagents and samples to come to room temperature (+18°C / +25°C) before handling.

Check that all plates are well drained after each wash.

Do not smoke, eat or drink when manipulating the kit.

Do not pipette by mouth.

Human sources for the preparation of calibrators and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

This kit contains potentially hazardous components. Though kit reagents are not classified being irritant to eyes and skin we recommend to avoid contact with eyes and skin and wear disposable gloves.

Calibrators, Controls and Buffers contain less of 0.1% (w/v) sodium azide. Do not eat and avoid contact with skin and eyes. This substance can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

ELISIS CELIAC GliA IgA kit ([bmd](#)) has been developed according CE Directives 67/548/EEC and 1999/45/EC relating to the classification, packaging and labeling of dangerous preparations.

ELISIS CELIAC GliA IgA kit ([bmd](#)) has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.

Incubation: We recommend test performance at +30°C for automated Systems.

Never expose components to higher temperature than +37°C.

Always pipette Substrate and Conjugate swith brand new tips only.

Protect Substrate from light.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Detach the exact number of wells needed and return unused wells to plastic pouch provided in the kit, with the desiccant bag.

For quantitative interpretation use calibrators to establish a standard curve

Set out in duplicate:

- Calibrators (Cal A, Cal B, Cal C, Cal D, Cal E and Cal F)
- Positive control (Ctrl +)
- Negative control (Ctrl -)
- Diluted samples (S)

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------|--------|----|-----|---|---|---|---|---|----|----|----|
| A | Cal A | Cal E | S1 | S9 | | | | | | | | |
| B | Cal A | Cal E | S2 | S10 | | | | | | | | |
| C | Cal B | Cal F | S3 | S11 | | | | | | | | |
| D | Cal B | Cal F | S4 | S12 | | | | | | | | |
| E | Cal C | Ctrl + | S5 | S13 | | | | | | | | |
| F | Cal C | Ctrl + | S6 | S14 | | | | | | | | |
| G | Cal D | Ctrl - | S7 | | | | | | | | | |
| H | Cal D | Ctrl - | S8 | | | | | | | | | |

For qualitative interpretation use cut-off calibrators

Set out in duplicate:

- Negative control (Ctrl -)
- Cut-off calibrator (Ctrl Ref)
- Positive control (Ctrl +)
- Diluted samples (S)

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|----------|-----|-----|-----|---|---|---|---|---|----|----|----|
| A | Ctrl - | S3 | S11 | S19 | | | | | | | | |
| B | Ctrl - | S4 | S12 | S20 | | | | | | | | |
| C | Ctrl Ref | S5 | S13 | S21 | | | | | | | | |
| D | Ctrl Ref | S6 | S14 | S22 | | | | | | | | |
| E | Ctrl + | S7 | S15 | S23 | | | | | | | | |
| F | Ctrl + | S8 | S16 | S24 | | | | | | | | |
| G | S1 | S9 | S17 | | | | | | | | | |
| H | S2 | S10 | S18 | | | | | | | | | |

2. Sample incubation

Add 100µL of Calibrators OR Cut-off Calibrator and controls into the designated microwells

Add 100µL of diluted samples.

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion.

Wash 3 times with 300µL of washing buffer.

3. Incubation of conjugate

Add 100µL of conjugate

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Wash step:

Remove the contents of the wells by rapid inversion

Wash 3 times with 300µL of washing buffer

4. Incubation of substrate

Add 100µL TMB substrate into each well

Incubate for 30 minutes at room temperature (+20°C/+32°C)

Protect the wells from intense light

5. Stop solution

Add 100µL of Stop solution to each well

Incubate for 5 minutes minimum

Agitate plate carefully for 5 sec.

6. Reading

Read the optical density of each well at 450nm (optionally 450/620nm) within 30 minutes.

INTERPRETATION OF RESULTS

For **quantitative interpretation** establish the standard curve by plotting the *optical density (OD) of each calibrator (y-axis)* with respect to the corresponding concentration values in *AU/mL (x-axis)*. For best results [bmd](#) recommend log (x) / lin (y) coordinates and 4-Parameter Fit. From the OD of each sample, read the corresponding antibody concentrations expressed in *AU/mL*.

| Negative Results | Positive Results |
|------------------|------------------|
| ≤ 15 AU/mL | > 15 AU/mL |

Example of a standard curve

We recommend pipetting calibrators in parallel for each run.

| Calibrators IgA | OD 450/620nm | CV % (Variation) |
|-----------------|--------------|------------------|
| 0 AU/ml | 0.059 | 1.4 |
| 3 AU/ml | 0.182 | 1.2 |
| 10 AU/ml | 0.323 | 2.2 |
| 30 AU/ml | 0.667 | 0.7 |
| 100 AU/ml | 1.316 | 0.9 |
| 300 AU/ml | 2.203 | 0.1 |

Example of calculation

| Patient | Replicate (OD) | Mean (OD) | Result (AU/mL) |
|---------|----------------|-----------|----------------|
| P01 | 0.654/0.633 | 0.644 | 27.6 |
| P02 | 1.284/1.263 | 1.274 | 89.9 |

Do not use this example for interpreting patients results!

For lot specific data, see enclosed quality control certificate.

Each laboratory should establish its own normal range based upon its own techniques, controls, equipment and patient population according to their own established procedures.

For **qualitative interpretation** read the optical density of the cut-off calibrator and the patient samples. Compare patient's OD with the OD of the cut-off calibrator. All samples which are higher than cut-off are considered positive.

| | |
|------------------|-------------------------|
| Negative: | OD Patient < OD cut-off |
| Positive: | OD Patient > OD cut-off |

A definite clinical diagnosis should not be based on the results of the performed test only, but should be made by the physician after all clinical and laboratory findings have been evaluated. The diagnosis is to be verified using different diagnostic methods.

LIMITATION

Hemolytic, lipemic, icteric samples or contaminated samples with may confound the results of this assay.

QUALITY CONTROL

The kit contains positive and negative controls. Both positive and negative control should be included with each run of the test. Positive and negative controls concentration should fall within the range printed on the quality control certificate. If either reagent control is invalid, the test should be repeated.

[bmd](#) offers an external multiparameter quality control (IMMUNO-TROL II, Cat: HM037). This control should be handled as samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

Sensitivity, specific affinity and agreement to other methods

ELISIS CELIAC Glia IgA kit ([bmd](#)) has been evaluated clinically.

The study was performed on 204 samples clinically identified:

⇒ 195 samples from patients with confirmed Celiac Disease, Crohn's disease, Ulcerative Colitis or other related autoimmune Diseases

⇒ 9 samples from healthy patients

Results

| | | Gliadin IgA predicate kit | | |
|------------------------|---|---------------------------|-----|---|
| | | + | - | |
| ELISIS CELIAC (TG IgA) | + | 14 | 20 | Relative sensibility : 87.5% Relative specificity : 89.4% Accuracy: 89.2% |
| | - | 2 | 168 | |

Test precision

Estimated for IgA isotype on three positive samples : strongly, moderate and weakly positive.

| ELISIS CELIAC Glia IgA | | | |
|---|--------|---|--------|
| Intra-essai (18 tests in a same assay) | | Inter-essais (18 tests in 3 differents assays) | |
| Titer mean AU/mL | CV (%) | Titer mean AU/mL | CV (%) |
| 50.4 | 2.8 | 46.2 | 3.0 |
| 32.2 | 5.9 | 29.0 | 4.6 |
| 13.4 | 4.7 | 14.6 | 8.8 |

Threshold value

Calculated on 76 samples from « blood donors ».

| ELISIS CELIAC Glia IgA | |
|------------------------|------------------|
| Cut-off : 15AU/mL | 97.4th percentil |

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DIETERICH W et al

Identification of tissue transglutaminase as the autoantigen of celiac disease.

Nat Med 1997; 3: 797-801

FOTOULAKI M et al.

Clinical application of immunological markers as monitoring tests in celiac disease.

Dig Dis Sci. 1999 Oct;44(10):2133-8.

KUMAR V et al.

Celiac Disease and Immunoglobulin A Deficiency: How effective are the serological method and diagnosis.

Clin. Diagn. Lab. Immunol. 2002 Nov, 9: 1295-1300.

SARDY M et al.

Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy.

Clin Chem. 1999 Dec;45 (12):2142-9.

VITORIA JC. et al.

Antibodies to gliadin, endomysium, and tissue transglutaminase for the diagnosis of celiac disease.

J Pediatr Gastroenterol Nutr. 1999 Nov;29(5):571-4.

WALTER-SMITH JA et al.

Working of European society of pediatric gastroenterology and nutrition (ESPGAN). Revised criteria for diagnosis of celiac disease





Arch. Dis of Child 1990; 65: 909

SUMMARY OF METHOD

| | AMOUNT TO BE DISTRIBUTED | REAGENTS | | INCUBATION CONDITIONS |
|-------------------------|--|--|-----------------------------------|---|
| | | QUANTITATIVE METHOD | QUALITATIVE METHOD | |
| INCUBATION OF SAMPLES | 100µL | Calibrators A, B, C, D, E and F (in duplicate) | Negative Control (in duplicate) | 30 min at RT |
| | 100µL | Positive Control (in duplicate) | Cut-off calibrator (in duplicate) | |
| | 100µL | Negative Control (in duplicate) | Positive Control (in duplicate) | |
| | 100µL | Diluted samples (in duplicate) | Diluted samples (in duplicate) | |
| WASHING | Wash 3 times with washing buffer (300µL/ well) | | | |
| INCUBATION OF CONJUGATE | 100µL | Anti-human IgA-HRP conjugate | | 30 min at RT |
| WASHING | Wash 3 times with washing buffer (300µL/ well) | | | |
| ENZYME REACTION | 100µL | TMB Substrate | | 30 min at RT protected from intense light |
| STOP REACTION | 100µL | Stop solution (1M HCl) | | 5 min |

RT: room temperature (+20°C/+32°C)

SYMBOLS USED

| | | | | | |
|---|---------------------------|---|-------------------------|---|------------------------------|
|  | Biological risk |  | Temperature limitation |  | Catalog Number |
|  | Read instructions for use |  | In Vitro Diagnostic Use |  | Lot Number |
|  | Number of tests |  | Use by |  | EC Declaration of Conformity |
|  | Sample Diluent |  | Stop Solution |  | Cut-off calibrator |
|  | Wash Buffer |  | Calibrator |  | Microplate |
|  | Control |  | Substrate |  | Conjugate |

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com