### DÉFINITION

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### DEFINITION

This reagent, prepared from human serum, is intended for use as an unassayed precision control with Immunofluorescence Assay (IFA).

REF	Parameter
ME 0022 CONTROL DNA	Positive control dsDNA
ME 0032 CONTROL AMA	Positive control anti-mitochondrial
ME 0032 KLM CONTROL LKM	Positive control anti-LKM1
ME 0042 CONTROL ASMA	Positive control anti-smooth muscle
ME 0052 CONTROL APCA	Positive control anti-parietals cells
ME 0062 ES CONTROL BM	Positive control anti-basal membrane (ASA-Pemphigoid)
ME 0062 ESD CONTROL AEmA	Positive control anti-endomysial
ME 0063 ES CONTROL ICS	Positive control anti-intercellular substance (ASA-Pemphigus)
ME 0082 CONTROL MIX	Multiple controls (Positive ANA-AMA-ASMA-APCA-Negative)

**bmd CAN PROVIDE YOU OTHERS SPECIFIC CONTROLS.  
PLEASE CONTACT US.**

### DESCRIPTION

**IFA control serums** have been designed for use with I.F.A method for monitoring assay performance and maintaining quality of assays.

Each **IFA control serum** lot is designed to produce a positive reaction for each autoantibodies above listed.

**IFA control Serums** do not have assigned values.

### PACKAGING

1 vial of 500µL. The antigenic specificity is confirmed on the vial label.

### WARNING

- Human sources for the preparation of this control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus or hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.
- The reagent contains less than 0.1% sodium Azide. Do not eat and avoid contact with skin and eyes.
- Reagents should be brought to room temperature (+20°C-+25°C) prior to use.

### PROCEDURE

**IFA control serums** are ready to use.

Due to variations within each laboratory (fluorescent microscopy, etc.) each laboratory should establish its own mean titer for each lot of positive control (generally ± one dilution from stated endpoint).

#### Negative Control:

Using the **NEGATIVE CONTROL SERUM** (**bmd**) as provided should exhibit less than 1+ fluorescence and appears reddish-orange due to the counterstain.

#### Positive Control:

Using the **POSITIVE CONTROL SERUM** (**bmd**) as provided should exhibit well defined specific fluorescent staining at an intensity of 3+ or greater.

#### Fluorescent intensity grading:

Fluorescent intensity may be semi-quantitated by following the guidelines established by the Centers for Disease Control, Atlanta, Georgia:

- 4+ = Maximal fluorescence; brilliant yellow-green.
- 3+ = Less brilliant yellow-green fluorescence.
- 2+ = Definite but dull yellow-green fluorescence.
- 1+ = Very dim subdued fluorescence.

The degree of fluorescent intensity is not clinically relevant and has only limited value as an indicator of titer. Differences in fluorescence microscope optics, filters and light sources may result in differences of 1+ or more fluorescent intensity when observing the same slide using different microscopes.

### STABILITY AND STORAGE

- Store the reagents at +2°C to +8°C when not in use.
- Do not use reagents beyond the expiration date.
- Open vial stability: 90 days at +2°C to +8°C.
- Do not freeze.

### LIMITS

- Performances and levels of reactivity of the **IFA control serums** may vary with different manufacturers' test and assay procedures.
- **IFA control serums** must not be substituted for the positive and the negative control reagents provided with manufacturers' test kits.
- Because **IFA control serums** are formulated with human sera and non human derived constituents, incompatibilities between this reagent and test kits may exist.
- **IFA control serums** are provide for quality assurance and must not be used for calibration.
- It is the responsibility of each laboratory to implement its own quality assurance program, to determine the suitability of **IFA control serums** for its particular use and to establish guidelines for interpretation of results obtained with **IFA control serums**.
- Adverse storage conditions or use of outdated reagents may produce erroneous results.
- These controls should be only use on the **bmd** slide ME series.

**BioMédical Diagnostics SA**

**Office**  
Actipole 25  
4 bld de Beaubourg  
77435 Marne la Vallée Cx2  
France

Tel : 33 1 64 62 10 12  
Fax : 33 1 64 62 09 66

E-mail : [support@bmd-net.com](mailto:support@bmd-net.com)  
Internet : [www.bmd-net.com](http://www.bmd-net.com)