Crithidia luciliae (IF)

REF

ME 0248

Coffret

ME 0296

Coffret

ME 0201

Lame de substrat (4 puits)

ME 0202

Lame de substrat (8 puits)

Français

DEFINITION

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BIBLIOGRAPHIE

1. Aarden LA, DeGroot ER, Feltkamp TEW. Immunology of DNA III
2. Ceppellini R, Polli E, Celada E. A-DNA reacting factor in serum of patients with lupus erythematosus diffusus. Proc Soc Exp Biol Med 96:572 (1957)
3. Cohen SA, et al. Character of anti-DAN antibodies in systemic lupus erythematosus. Clin Exp Immunol 8:551 (1971)
4. Crowe W, Kushner I. An immunofluorescent method using *Crithidia luciliae* to detect antibodies to double-stranded DNA. Arthritis Rheum 20:811 (1977)
5. Edmonds JP, et al. The value of tests for antibodies to DNA in monitoring the clinical course of SLE. Clin Exp Immunol 22:9 (1975)
6. Gershwin ME, Steinberg AD. Qualitative characteristics of anti-DNA antibodies in lupus nephritis. Arthritis Rheum 17:947 (1974)
7. Hughes GRV, Cohen SA. Anti-DNA activity in systemic lupus erythematosus. A diagnostic and therapeutic guide. Ann Rheum Dis 30:259 (1971)
8. Slater NGP, Cameron JS, Lessof MH. The *crithidia luciliae* kinetoplast immunoplast immunofluorescence test in systemic lupus erythematosus. Clin Exp Immunol 25:480
9. Stingl G, et al. An immunofluorescence procedure for the demonstration of antibodies to native, double-stranded DNA and f circulating DNA-anti -DNA complexes. Clin Immunol Immunopathol 6:131 (1976)
10. Sennekamp J, et al. Determination of antibodies against native double-stranded DNA with an indirect immunofluorescence technique using flagellate *crithidia luciliae*. A Rheumatol 37(Suppl 5):113 (1978)

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

Anglais

INTENDED USE

Crithidia luciliae (IF) kit is a screening test for the semi-quantitative detection of auto-antibodies against nDNA in human serum by immunofluorescence.

DIAGNOSTIC VALUE

The initial finding that auto-antibodies against native, double-stranded deoxyribonucleic acid (nDNA) are present in sera of patients with systemic lupus erythematosus (SLE) was made more than 40 years ago⁽²⁾.

Through the years much evidence has accumulated that auto-antibodies against nDNA are almost specific for sera of SLE patients especially with renal involvement. The specificity of anti-nDNA antibodies is much greater than that of antinuclear antibodies (ANA)⁽⁵⁾⁽⁶⁾⁽⁹⁾.

In recent years, with the availability of a well defined substrate, a simple immunofluorescence procedure has been developed for the detection of anti-nDNA antibodies⁽¹⁾⁽⁴⁾⁽⁹⁾. It is presently being used as an aid to differentiate disease group.

PRINCIPLE OF THE TEST

Antibodies against nDNA do not exhibit species or organ specificity. The **Crithidia luciliae (IF)** kit utilizes as substrate the kinetoplast (a large mitochondrion) of *Crithidia luciliae*.

This hemoflagellate has high amounts of DNA concentrated in the kinetoplast, an organelle detectable between the nucleus and the base of the flagellum.

- First, C Luciliae substrate on microscope slides is overlaid with the diluted patient's serum and incubated to allow the anti-nDNA antibodies (A-nDNA-A) present in positive sera, to bind of the DNA.
- After a wash step, sites where A-nDNA-A (human immunoglobulins) reacted with DNA (in the kinetoplast) are made visible by the binding of fluorescently labelled anti-human antibodies (the conjugate).
- A second wash eliminates the free conjugate.
- After mounting, the slide is observed with a fluorescence microscope. A positive reaction results in a fluorescence of the kinetoplast.

REAGENTS AVAILABLE INDIVIDUALLY

ME 0201	<i>Crithidia luciliae</i> slide of 4 wells
ME 0202	<i>Crithidia luciliae</i> slide of 8 wells
ME 902.2	Fluorescent anti-IgG conjugate (2 mL)
ME 902.8	Fluorescent anti-IgG conjugate (8 mL)
ME 0022	Positive Control dsDNA (500 µL)
HME NEG	Negative control (500 µL)
HME MM	Mounting Medium buffered glycerol (3,5 mL)
HME PBS10	One-liter packet of dry PBS 10g (10 x 1L)
HME 12CS	24x60mm glass coverslips (100 units)

COMPONENTS

	ME 0248	ME 0296
Crithidia luciliae slides <div style="border: 1px solid black; display: inline-block; padding: 2px;">SORB</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">SLD</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">n</div>	12 x 4 wells	12 x 8 wells
Positive control dsDNA <i>Dropper vial containing positive patient sample producing fluorescent staining in the kinetoplast of C.luciliae only or kinetoplast + nuclear.</i> Ready to use <div style="border: 1px solid black; display: inline-block; padding: 2px;">CONTROL</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">DNA</div>	1 x 0,5 mL	1 x 0,5 mL
Negative Control <i>Dropper vial containing negative patient sample.</i> Ready to use <div style="border: 1px solid black; display: inline-block; padding: 2px;">CONTROL</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">-</div>	1 x 0,5 mL	1 x 0,5 mL
Anti-IgG conjugate <i>Fluorescently labelled anti-human IgG globulin in phosphate buffered saline (PBS) with 1% bovine serum albumin and with evans blue.</i> Ready to use <div style="border: 1px solid black; display: inline-block; padding: 2px;">CONJ</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">DTAF</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">IgG</div>	1 x 2 mL	2 x 2 mL
Phosphate buffer saline <i>Each PBS Powder pouch contains sufficient buffer powder to make 1 liter.(10g/pouch)</i> <div style="border: 1px solid black; display: inline-block; padding: 2px;">BUF</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">WASH</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">RCNS</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">1L H₂O</div> To dilute	2 pouch	2 pouch
Mounting Medium <i>Dropper vial containing Buffer glycerol.</i> Ready to use <div style="border: 1px solid black; display: inline-block; padding: 2px;">MM</div>	3 mL	3 mL

MATERIAL REQUIRED BUT NOT PROVIDED I

- Volumetric pipettes to deliver 10-15 µL
- Small test tubes and rack or U-bottom microtitration plates
- Humid chambers – Petri dishes with moist filter paper
- Coplin staining jars with slide holders
- One liter volumetric flask or graduated cylinder
- Squeeze bottle for PBS
- Cover glasses #1, 24 x 60 mm.
- Distilled water
- Fluorescence microscope
- Timer
- Forceps
- Positive and negative controls

SPECIMEN COLLECTION

- Collect blood specimens in plain untreated tubes.
- After clot formation, immediately harvest the serum.
- The serum specimens may be stored in a refrigerator (2°C-8°C) for up to 48 hours. If longer storage is desired, store at -20°C.
- Turbidity, haemolysis, bacterial growth or drugs capable of fluorescing in the serum may adversely affect the performance of the test. Do not use hemolyzed or lipemic sera.

STABILITY AND STORAGE

- Expiration dates for substrates (slides), conjugate and controls are given for storage at -20°C or below.
- Once thawed the conjugate and controls are stable for 90 days when stored at +2° to +8°C.
- Repeated freezing and thawing should be avoided since it is detrimental to the products.
- The mounting medium is stable until the expiration date when stored in the refrigerator.
- The PBS powder can be stored at room temperature until its expiration date.

PRECAUTIONS

All bmd Test Kit components are for *in vitro* diagnostic use ONLY.

All human serum components have been tested by FDA approved method and found negative for antibody to HIV, HCV and non-reactive for HBsAg. However, no test method can offer complete assurance of the absence of these antigens. Care should be exercised when handling ALL human sera containing specimens.
Avoid mouth pipetting.

All controls, conjugate and mounting medium contain sodium azide as preservative in concentration 1:1,000. This preservative is toxic if ingested. In addition, sodium azide may react with lead or copper plumbing to form explosive azides. On disposal of sodium azide containing solutions, flush with a large volume of water to prevent azide build-up.

Do not substituted others reagent's manufacturers.

Do not use components or reagents beyond their expiration dates. Extended incubation during the wash step or improper washing may result in poor morphology of substrate and increased background fluorescence.

Avoid using pen to identify the slide

Remove slides from freezer and allow them to reach room temperature – **Don't remove from protective envelopes before use.**

Use each slide immediately after opening the envelope to avoid damaging the substrate by drying.

All reagents must be at room temperature before use (at least 15 mn)

SET-UP

All reagents must be at room temperature before use.

1. Preparation of the wash and dilution buffer

- Reconstitute buffer pack with one liter distilled water.
Wait 10 mn before use
- Storage period : 4 weeks at +2°C/+8°C

2. Sera dilution

- Dilute freshly collected and properly stored patient serum 1:10 in PBS.
Ex : 20µL serum + 180µL PBS
For titrations, set up doubling dilutions (i.e 1/20, 1/40 ...)
- Mix with vortex

TEST PROCEDURE

1. Incubation of samples

- Identify wells which are going to contain positive and negative control sera and patient's serum dilutions.
- Remove slide from protective envelopes and put on the humid chamber
- In the identified well of the slide, add
 - **20µL sample by avoiding to touch the substrate with the tip of the pipette.**
 - 1 drop of positive control ready to use
 - 1 drop of negative control ready to use
- Do the same with each slide.
- Incubate for 30 minutes at room temperature.

Do not let the well dry out

2. Wash with PBS

- Remove slides from humid chamber.
- Tap each slide quickly on its side to avoid contaminations.
- Dip the slides several times in a beaker containing fresh PBS. **Do not use wash bottle to rinse the slides.**
- Immediately place slides in a slide holder or a coplin jar and wash with PBS and incubate for 10 minutes at room temperature. Agitation is advised but not useful.
Note: If more than 10 slides are tested, wash two time in PBS bath (5 minutes each).

⚠ Extended incubation and improper washing may result in poor morphology of substrate and increased background fluorescence.

3. Incubation of Conjugate

- Remove the slide from PBS
- Eliminate excess of PBS on Teflon with the blotter
- Immediately apply 1 drop of conjugate to each well to cover substrates
- Incubate for 30 minutes at room temperature and **far from light.**

Do not let the wells dry out

4. PBS rinse

- Repeat Step 2 "Wash with PBS".

5. Mount Cover glasses

- Apply 3 or 4 drops mounting medium to each slide.
- Cover with cover glasses – Lateral movement and pressure will damage the substrate

6. Reading

It is recommended to examine the day of the test with an appropriately equipped fluorescence microscope (x10) and use microscopes (x25 or x40) to confirm a positive fluorescence.

Note: If delay in examining slides is anticipated, seal coverglass with clear nail polish and store in refrigerator.

QUALITY CONTROL

Inside the kit of 48 and 96 tests, [bmd](#) provide:

- 1 positive control: Positive patient sample producing fluorescent staining in the kinetoplast of *C. luciliae* only or kinetoplast + nuclear.
- 1 negative control: patient sample producing no fluorescent staining pattern in the kinetoplast of *C. luciliae*.

ME 0022 Positive control can be used as external quality control.

INTERPRETATION OF RESULTS

The end point is the highest dilution that produces a positive result.

All positive sera should be titrated to the highest titre using a series of doubling dilutions.

Read only single well defined organism in each field. Some of the organisms will look deformed, due to their stage of growth and orientation on the viewing interference.

• Positive result

Patient's sera giving bright apple green appearance of the kinetoplast alone or the kinetoplast and the nucleus of dilutions 1:10 and higher are considered positive. Being highly specific for SLE, the immunofluorescence test is most useful in primary diagnosis of this disease⁽⁸⁾. Antibodies to native DNA are thought to correlate with clinical activity of the disease⁽⁶⁾. Mean titre of SLE patients is 1:160 but titres as high as 5,120 have been reported⁽⁹⁾.

• Negative result

No fluorescence on the kinetoplast whatever the fluorescence of the nuclear or the flagella.

EXPECTED VALUES

Normal population is negative with the screening dilution 1/10e

LIMITATIONS OF THE TEST

The *Crithidia luciliae* kit is a valuable diagnostic aid. However, it should not be considered diagnostic by itself only. Positive nDNA test results may be obtained with sera of healthy people as a result of many factors (e.g. age, drug treatment, etc). Therefore, all nDNA test results should be interpreted by a medical authority in light of all known patient's conditions.

PERFORMANCE CHARACTERISTICS

Assay units:

Indirect immunofluorescence (IIF) was assayed using accepted arbitrary units of intensity under fluorescence microscopy. The units from greatest to least intensity are; 4+, 3+, 2+, 1+, +/-, -.

Titer reproducibility:

Four serum samples were tested on three different lot numbers of [bmd](#) kits on three different occasions. All positive values titered to negative at the same dilution for each of the three replicates.

Normal samples:

Sera from 743 healthy blood donors with no history of autoimmune disease were tested on [bmd](#) Crithidia luciliae kits. Sera were diluted 1:20. A negative reading at 1:20 was considered negative. No positive readings of 1+ or higher at dilutions of 1:20 or greater were scored as positive. No sera of health subjects is reported as positive with [bmd](#) Crithidia luciliae kit.

Sera from patients with diagnosed disorders:

Sera from 1019 patients with variety of diagnosed disorders were tested with [bmd](#) Crithidia luciliae kits. Sera were diluted from 1:10 to 1:320. No positive readings of 1+ or higher at dilutions of 1:40 were scored as positive for A-nDNA-A.

REFERENCES

1. Aarden LA, DeGroot ER, Feltkamp TEW. Immunology of DNA III
2. Ceppellini R, Polli E, Celada E. A-DNA reacting factor in serum of patients with lupus erythematosus diffusus. Proc Soc Exp Biol Med 96:572 (1957)
3. Cohen SA, et al. Character of anti-DNA antibodies in systemic lupus erythematosus. Clin Exp Immunol 8:551 (1971)
4. Crowe W, Kushner I. An immunofluorescent method using Crithidia luciliae to detect antibodies to double-stranded DNA. Arthritis Rheum 20:811 (1977)
5. Edmonds JP, et al. The value of tests for antibodies to DNA in monitoring the clinical course of SLE. Clin Exp Immunol 22:9 (1975)
6. Gershwin ME, Steinberg AD. Qualitative characteristics of anti-DNA antibodies in lupus nephritis. Arthritis Rheum 17:947 (1974)
7. Hughes GRV, Cohen SA. Anti-DNA activity in systemic lupus erythematosus. A diagnostic and therapeutic guide. Ann Rheum Dis 30:259 (1971)
8. Slater NGP, Cameron JS, Lessof MH. The crithidia luciliae kinetoplast immunoplast immunofluorescence test in systemic lupus erythematosus. Clin Exp Immunol 25:480
9. Stingl G, et al. An immunofluorescence procedure for the demonstration of antibodies to native, double-stranded DNA and of circulating DNA-anti-DNA complexes. Clin Immunol Immunopathol 6:131 (1976)
10. Sennekamp J, et al. Determination of antibodies against native double-stranded DNA with an indirect immunofluorescence technique using flagellate crithidia luciliae. A Rheumatol 37(Suppl 5):113 (1978)





SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	20 µL 1 drop 1 drop	Diluted samples Positive control (ready to use) Negative control (ready to use)	30 mn at R.T
<i>Do the same with each slide</i>			
PBS WASH	⇒ Tap each slide quickly on its side to avoid contaminations ⇒ Dip the slide several times in a beaker containing fresh PBS ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. ⇒ N.B : Agitation is not useful.		10 mn at R.T
<i>Eliminate excess of PBS with the blotters provided (*) – Dry uniformly and without friction using a roller</i>			
INCUBATION OF CONJUGATE	1 drop to each well	Fluorescein conjugated Affinipure Goat anti-human IgG (ready to use)	30 mn at R.T in the dark
<i>Verify that all wells are perfectly covered</i>			
PBS RINSE	⇒ Tap quickly each slide on its side to eliminate the conjugate. ⇒ Dip the slide several times in a beaker containing fresh PBS. ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>DO NOT ELIMINATE EXCESS OF PBS (Teflon) (*)</i>			
MOUNT COVER GLASSES	3-4 drops to each slide	Mounting Medium (ready to use)	-
READING	Examine with an appropriately equipped fluorescence microscope (x 10 – x 25 - x 40).		

(*) **Do not let wells dry out.** Every deposit must be made at once.

R.T : Room Temperature

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Slide with "n" wells		Negative control		Mounting Medium
	Fluorescent antibody conjugate		Positive control		Wash Buffer
	Immunofluorescent assay				Reconstitute with 1L distilled water

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com