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REF | ME 0600 | **Lame de substrat (4 puits)**

DÉFINITION

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BIBLIOGRAPHIE

- 1 Nicholson, G.A, McCleod.JG, Griffiths, LR, comparison of diagnostic tests in myasthenia Gravis, Clin Exp Neurol. 18:45, 1982.
- 2 Gotti C , Mantegazza R. And Clementi F : new Antigen for Antibody Detection in Myastheniae Gravis, Neurol 34:374-77, 1984.
- 3 Gihlus N.E, Aarli JA, Matre R, Myasthenia Gravis, The specificities of skeletal Muscle and thymus antibodies.
- 4 Irvine, WJ; Kalden, JR: Muscle in allergic disease. In Clinical Aspects of immunology.Ch 51. Gell, PGH, coombs, RRA and Lackmann, PJ (Eds) Blackwell Scientific Publication , Oxford Engand 1975
- 5 Yamamoto AM et Bach JF (1999) – Autoanticorps anti-récepteurs de l'Acétylcholine (RACH) et myasthénie - Autoanticorps marqueurs de maladies auto-immune. bmd Editions p.393
- 6 Humbel RL (1994) – Endocrinopathie auto-immune – Auto-anticorps et maladies auto-immunes Editions Option Bio Elsevier version N°1 – p.207.

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
Substrat slide (4 wells)

INTENDED USE

The **sections of specific-organs** of primate monkey () are used for the qualitative and semi-quantitative detection of human autoantibodies by indirect immunofluorescence method. These substrate slides are useful in the diagnosis of certain autoimmune disorders listed below.

Catalog No	Substrate	Autoimmune disorders
ME 0600 AG	Adrenal gland	Idiopathic Addison's disease
ME 0600 BM	Bone marrow	Hemolytic anemia
ME 0600 BR	Brain	Schizophrenics
ME 0600 CO	Colon	Chronic ulcerative colitis
ME 0600 HE	Heart	Myocardial infarction
ME 0600 LI	Lip	Pemphigus, pemphigoid
ME 0600 LU	Lung	Goodpasture's syndrome Hamman-Rich Syndrome
ME 0600 OV	Ovary	Gonadal failure, Premature ovarian
ME 0600 PG	Pituitary Gland	Endocrine diseases
ME 0600 SG	Salivary Gland	Sjögren's syndrome
ME 0600 SN	Sciatic Nerve	Sclerosis Guillian-Barre syndrome
ME 0600 TE	Testes	Testicular failure, Infertility
ME 0600 SM	Striated Muscle	Multiple Sclerosis Guillian-Barre syndrome

DIAGNOSTIC VALUE

Example:

Anti-Skeletal muscle antibodies (ASkMA): the myasthenia gravis is a neuromuscular disease leading to fluctuating muscle weakness and fatigability.

Clinical studies on the age and sex distribution clearly show a prevalence of females in patients below the age of 40 years in contrast a higher incidence of males in patients over 60 years.⁽¹⁾

Circulating anti-acetylcholine receptor (AChR) antibodies have been detected in 75% to 90% of myasthenic patients.^(2,5)

Antibodies to the acetylcholine receptor are important in the pathogenesis of myasthenia gravis (MG) but recent studies indicate they are not solely responsible for the clinical features of the disease, since thymic epithelial factors may also be implicated.⁽³⁾

Mc Farlin, Johnson & Seymor (1968) found that the skeletal muscle antibodies (ASkMA) in patients with MG belong to the IgG class⁽⁴⁾.

Anti-adrenal antibodies (AAdA): Antibodies anti-adrenal cortex are found in patients with Addison's disease⁽⁶⁾.

Anti-salivary gland antibodies: the development of anti-salivary gland antibodies which occur in Sjogren's syndrome. These antibodies bind to excretion canal of epithelial cells.

However non-organ-specific auto-antibodies can bind to the salivary glands and especially the antimyochondrial antibodies of type 2 (AMA-M2).

PRINCIPLE OF THE TEST

Indirect immunofluorescence is a method usually used in the detection of autoantibodies.

First:

The substrate sections of monkey on microscope slides are overlaid with diluted patient's serum and incubated to allow the antibodies present in positive sera.

Second:

Sites that bind to antibodies are made visible with a fluorescently labelled antihuman antibodies (the conjugate).

In the fluorescence microscope, a positive reaction results in an apple green fluorescence pattern of the characteristic structural features of the substrate sections.

COMPONENTS

	ME 0600
Monkey substrate slides (specific sections of monkey in multiple well-slides) <div style="border: 1px solid black; display: inline-block; padding: 2px;">SORB SLD 4</div>	1 slide of 4 wells
Blotter for drying the substrate slides around the wells.	1

REAGENTS AVAILABLE INDIVIDUALLY

ME 902.2	Fluorescent anti-IgG conjugate (2 mL)
ME 902.8	Fluorescent anti-IgG conjugate (8 mL)
ME 0904	Fluorescent anti-IgA conjugate (2 mL)
HME NEG	Negative control (500 µL)
HME MM	Mounting Medium buffered glycerol (3,5 mL)
HME PBS10	One-liter packet of dry PBS 10g (10x1L)
HME 12CS	24x60mm glass coverslips (100 units)

 can provide you **positive controls**, please contact us.

MATERIAL REQUIRED BUT NOT PROVIDED

- Volumetric pipettes to deliver 10-15 µL
- Small test tubes and rack or U-bottom microtiter plates
- Humid chambers – Petri dishes with moist filter paper
- Coplin staining jars with slide holders
- One liter volumetric flask or graduated cylinder
- Beaker for PBS
- 24 x 60 mm #1 glass coverslips
- Distilled water Cover glasses
- Fluorescence microscope
- Timer
- Forceps
- Absorbent paper

SPECIMEN COLLECTION

- Collect blood specimens in plain untreated tubes.
- After clot formation, immediately harvest the serum.
- The serum specimens may be stored in a refrigerator (+2°C-+8°C) for up to 48 hours. If longer storage is desired, store at -20°C.
- Turbidity, haemolysis, bacterial growth or drugs capable of fluorescing in the serum may adversely affect the performance of the test. Do not use hemolyzed or lipemic sera.
- Specimen should be handled as though they contained agents infectious to humans

STABILITY AND STORAGE

- Expiration date for slides (tissue substrates) is given for storage at -20°C or below.

We recommend to use with this substrate, the reagents listed above § " REAGENT AVAILABLE INDIVIDUALLY ".

- Once thawed the conjugate and controls are stable for 90 days when stored at +2° to +8°C. Do not alicote the conjugate at -20°C.
- Repeated freezing and thawing should be avoided since it is detrimental to the products.
- The mounting medium is stable until the expiration date when stored in the refrigerator.
- The PBS powder can be stored at room temperature until its expiration date.

PRECAUTIONS

All bmd substrate slides are for *in vitro* diagnostic use ONLY. Do not use components or reagents beyond their expiration dates. All reagents must be at room temperature before use (at least 15 mn).

After freezing, use each slide immediately after opening the envelope to avoid the drying of the wells.

In order to avoid damaging the substrate:

- Apply each drop on the surround of wells in order to touch the substrate with the tip of the pipette.
- During the wash step, do not rince the slides with a pipette or squeeze (wash) bottle.
- Do not apply any pressure on the slide during to mount the cover glasses.
- Do not move the substrate slide after adding the mounting medium buffered.
- Avoid direct pressure.
- Avoid using pen to identify the slide.
- **Always avoid mouth pipetting.**

SET-UP

All reagents must be at room temperature before use.

1. Preparation of the wash and dilution buffer

- Reconstitute buffer pack with one liter distilled water. Wait 10 mn before use
- Storage period : 4 weeks at +2°C/+8°C

2. Sera dilution

- Dilute freshly collected and properly stored patient serum 1:10 in PBS.
Ex : 20µL serum + 180µL PBS
- Mix with vortex

TEST PROCEDURE

1. Incubation of samples

- Identify wells which are going to contain positive and negative control sera and patient's serum dilutions.
- Remove slide from protective envelopes and put into the humid chamber
- In the identified well of the slide, add
 - 40µL of diluted sample, **by avoiding to touch the substrate with the tip of the pipette.**
 - 1 drop of ready to use positive control
 - 1 drop of ready to use negative control
- Do the same with each slide.
- Incubate for 30 minutes at room temperature.

Do not let the wells dry out

2. Wash with PBS

- Remove slides from humid chamber.
 - Tap each slide quickly on its side to avoid contaminations and let the sera run off onto absorbent paper.
 - Dip the slides several times in a beaker containing fresh PBS. **Do not use wash bottle to rinse the slides.**
 - Immediately place slides in a slide holder or a coplin jar and wash with PBS and incubate for 10 minutes at room temperature. Agitation is not useful.
- Note: If more than 10 slides are tested, wash two time in PBS bath (5 minutes each).

⚠ Extended incubation and improper washing may result in poor morphology of substrate and increased background fluorescence.

3. Incubation of Conjugate

- Remove the slide from PBS
- Eliminate excess of PBS onto some absorbent paper.
- Blot the slide with blotter provided and put into the humid chamber
- Immediately apply 1 drop of conjugate to each well to cover substrates.
- Incubate for 30 minutes at room temperature and **far from light.**

Do not let the wells dry out

4. PBS rinse

- Repeat Step 2 "Wash with PBS".

5. Mount Cover glasses

- Apply 3 or 4 drops of mounting medium to each slide
- Cover with a #1 glass coverslip – Pressing or laterally moving the coverslip will damage the substrate.

6. Reading

It is recommended to examine the day of the test with an appropriately equipped fluorescence microscope (x 10) and to use (x25 or x40) microscopes to confirm a positive fluorescence.

Note: If delay in examining slides is anticipated, seal coverglass with clear nail polish and store in refrigerator.

INTERPRETATION OF RESULTS

All positive sera should be titrated to the highest titer using a series of doubling dilutions. If more than one antibody is present in the patient's serum, interference may occur. Titration of such serum may reveal first one and then another autoantibody.

No definitive association between the fluorescence pattern and any specific disease should be made. As a reference ONLY, certain correlations that have been cited are shown in the table below.

Examples:

Substrate	Ig class	Cut Off	Major disease association(s)
Adrenal gland	IgG	1/5	Addison's disease Adrenatitis Idiopathic disease
Heart	IgG	1/10	Cardiomyopathy Dilated myocarditis Infarctions Acute rheumatic fever
Skeletal muscle	IgG	1/60	Myasthenia gravis disease
Ovary	IgG	1/10	Polyendocrinopathy syndrome type 1
Pituitary Gland	IgG	1/10	Hypophyseal fealure primary Cushing disease Diabete mellitus

QUALITY CONTROL

On each run, both positive and negative controls, should be included.

When results obtained are distant from expected values, test run should be repeated. Should atypical observations continue, we recommend to investigate if :

- adequate amount of serum is applied on the wells
- wells were allowed to dry out during incubation of serum or conjugate
- there is any sign of microbial contamination of the liquid reagents due to incorrect storage and handling. If turbidity appears on any of the liquid reagents, discard and use another batch.
- Fluorescence microscope is properly aligned, using correct bulb, correct filters, etc...

LIMITATIONS OF THE TEST

This test is a valuable diagnostic aid. However, it should not be considered diagnostic by itself alone. Positive results may be obtained with sera of apparently healthy people as a result of many factors. Therefore, all test results should be interpreted by a medical authority in light of all known patient's conditions.

BIBLIOGRAPHY

- 1 Nicholson, G.A, McCleod.JG, Griffiths, LR, comparison of diagnostic tests in myastenia Gravis, Clin Exp Neurol. 18:45, 1982.
- 2 Gotti C, Mantegazza R. And Clementi F : new Antigen for Antibody Detection in Myasteniae Gravis, Neurol 34:374-77, 1984.
- 3 Gihlus N.E, Aarli JA, Matre R, Myastenia Gravis, The specificities of skeletal Muscle and thymus antibodies.
- 4 Irvine, WJ; Kalden, JR: Muscle in allergic disease. In Clinical Aspects of immunology.Ch 51. Gell, PGH, coombs, RRA and Lackmann, PJ (Eds) Blackwell Scientific Publication , Oxford Engand 1975
- 5 Yamamoto AM et Bach JF (1999) – Autoanticorps anti-récepteurs de l'Acétylcholine (RACH) et myasthénie - Autoanticorps marqueurs de maladies auto-immune. bmd Editions p.393
- 6 Humbel RL (1994) – Endocrinopathie auto-immune – Auto-anticorps et maladies auto-immunes Editions Option Bio Elsevier version N°1 – p.207.



SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	40 µL	Diluted samples	30 mn at R.T
	1 drop	Positive control (ready to use)	
	1 drop	Negative control (ready to use)	
<i>Do the same with each slide</i>			
PBS WASH	⇒ Tap each slide quickly on its side to avoid contaminations ⇒ Dip the slide several times in a beaker containing fresh PBS ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>Eliminate excess of PBS on Teflon with the blotter(*) - Dry uniformly and without friction using a roller</i>			
INCUBATION OF CONJUGATE	1 drop to each well	Fluorescein conjugated Affinipure (ready to use)	30 mn at R.T in the dark
<i>Verify that all wells are perfectly covered</i>			
PBS RINSE	⇒ Tap quickly each slide on its side to eliminate the conjugate. ⇒ Dip the slide several times in a beaker containing fresh PBS. ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>DO NOT ELIMINATE EXCESS OF PBS (Teflon) (*)</i>			
MOUNT COVER GLASSES	3-4 drops to each slide	Mounting Medium (ready to use)	-
READING	Examine with an appropriately equipped fluorescence microscope (x 10 – x25 - x 40).		

(*) Do not let wells dry out. Every deposite must be made at once.

R.T : Room Temperature

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Slide of 4 wells		Immunofluorescent assay		

BioMédical Diagnostics SA

Office
 Actipole 25
 4 bld de Beaubourg
 77435 Marne la Vallée Cx2
 France

Tel : 33 1 64 62 10 12
 Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
 Internet : www.bmd-net.com