Endomysium (AEmA) IF

REF

ME 0648 ESD

Coffret

48

ME 0604 ESD

Lame de substrat (4 puits)

ME 0696 ESD

Coffret

96

ME 0608 ESD

Lame de substrat (8 puits)

DEFINITION

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BIBLIOGRAPHIE

1. Beutner EH, Chorzelski TP, Kumar V, Leonard J, Krasny S : Sensitivity and specificity of IgA-class anti-endomysial antibodies for dermatitis herpetiformis and findings relevant to their pathogenic significance. Journal of the American Academy of Dermatology, Vol 15, N°3, p 464 (1986).
2. Kapuscinska A, Zalewski T, Chorzelski T, Sulej J, Beutner EH, Kumar V, Rossi T: Disease specificity, sensitivity and dynamics of changes in IgA class anti-endomysial antibodies in celiac disease. J Pediatric Gastroenterol Nutr, Vol 6, N°4, p529 (1987).
3. Blazer S, Naveh Y, Berant M et al. : Serum IgG antibodies to gliadin in children with celiac disease as measured by an immunofluorescence method. J Pediatric Gastroenterol Nutr, Vol 3, p205 (1984).
4. Chorzelski TP, Beutner EH, Sulej J, Tchorzewska H, Jablonski S, Kumar V, Kapuscinska A : IgA anti-endomysium antibody. A new immunological marker for herpetiformis and coeliac disease. British Journal of Dermatology Vol 111, p 395 (1984).
5. Leonard JN, Haffenden GP, Fry L : Dermatitis herpetiformis in “ Immunopathology of the skin ”, Beutner EH, Chorzelski TP and Kumar V. Eds, John Wiley and Sons , New York, 3rd Ed. p 433 (1987).
6. Weinstein WM : Latent celiac sprue. Gastroenterology , Vol 66 N° 4, p 489 (1974).
7. Unsworth DJ, Leonard JN, Fry L : Anti-reticulin and anti-gliadin antibodies in dermatitis herpetiformis and celiac disease. “ Immunopathology of the skin ”, Beutner EH, Chorzelski TP and Kumar V. Eds, John Wiley and Sons , New York, 3rd Ed. p 455 (1987).
8. Volta U, Molinario N, Fusconi M et al. IgA anti-endomysial antibody test a step forward in celiac disease screening, Diag. Disease and Sci. Vol 36, N°6, 1991.
9. Kumar V, Jarzabek-Chorzelska M, Sulej J, Rajadhyaksha M., and Jablonska S : Tissue transglutaminase and endomysial antibodies - Diagnostic markers of gluten-sensitive enteropathy in dermatitis herpetiformis, Clinical Immunology, Vol 98, N°3, p378-382 (2001).
10. Johanet C. Marqueurs sérologiques de la maladie coeliaque. Spectra Biologie, Vol 17, N°99 (1998).
11. Ghedira I; et al. : Intérêt des anticorps anti-endomysium, anti-réticuline et anti-gliadine dans le diagnostic de la maladie coeliaque chez l'adulte. Ann Biol Clin, 57 : 717-9 (1999).

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

ME 0608 ESD Substrat slide (8 puits)

INTENDED USE

The **Endomysium (AEmA)** kit is an indirect Immunofluorescence test used for the qualitative and semi-quantitative detection of IgA class endomysial autoantibodies in human serum as an aid in the diagnosis of Celiac Disease and other related disorders such as Dermatitis Herpetiformis.

DIAGNOSTIC VALUE

Celiac Disease (CD), considered to be an autoimmune disease, is a common gastrointestinal disorder, primarily of children, which is induced by gluten. Virtually all active cases of CD and 70% of a related disorder the dermatitis herpetiformis ⁽¹⁾ exhibit anti-endomysial antibodies of the IgA class. In most CD patients, the symptoms and anti-endomysial antibody titers decline once patients are put on a gluten-free diet ⁽³⁾. Thus, the detection of anti-endomysial antibodies can be used for the diagnosis of CD and for monitoring of the disease process ^(10,11).

Dermatitis Herpetiformis is an intensely pruritic vesicular skin disease characterized clinically by urticarial lesions and blisters. Any part may be involved but the lesions are usually distributed symmetrically over the elbows, knees, upper back and buttocks ^(5,6,9).

Celiac Disease entails intolerance to wheat, barley, rye and oat protein (gluten and its fractions), resulting in severe structural damages to the small intestine mucosa. The incidence of celiac disease is considered to range between 1/300 and 1/6500. A definitive diagnosis usually demands repeated intestinal biopsies, a controlled diet and follow up over a prolonged period ⁽²⁾.

PRINCIPLE OF THE TEST

Anti-endomysial antibodies, exhibiting some organ and species specificity ⁽³⁾, are detected by indirect immunofluorescence on sections of primate oesophagus (distal), where they react with the endomysium or the reticulin components of the smooth muscle filaments.

First:


Sections of monkey oesophagus distal on microscope slides are overlaid with diluted patient's serum and incubated to allow the anti-endomysial antibodies – present in positive sera – to bind to the antigen sites in the oesophagus tissue.

Second:

Sites where the anti-endomysial antibodies reacted with the substrate are made visible by the binding of the fluorescently labelled anti-human globulin (the conjugate).

This microphotography shows in a stronger magnification the typical pattern of esophagus in longitudinal section.

- E :** Stratified squamous epithelium
- C :** Chorion
- MM :** Muscularis mucosae
- S :** Submucosa
- CM :** Muscularis : circular layer
- ML :** Longitudinal muscularis
- A :** Adventitia

COMPONENTS

	ME 1248	ME 1296
Endomysium Substrate slides (monkey esophagus distal sections in multiple well-slides) SORB SLD n	12 x 4 wells	12 x 8 wells
Anti-IgA conjugate Fluorescently labelled anti-human IgA globulin in phosphate buffered saline (PBS) with 1% bovine serum albumin. Ready to use CONJ DTAF IgA	1 x 3 mL	1 x 2 mL + 1 x 3 mL
Positive control anti-endomysial (AEmA) Dropper vial containing positive patient sample producing a characteristic fluorescent pattern in the endomysium Ready to use CONTROL AEmA	1 x 0,5 mL	1 x 0,5 mL
Negative Control Dropper vial containing negative patient sample producing no fluorescent staining of the endomysium of oesophagus substrate. Ready to use CONTROL -	1 x 0,5 mL	1 x 0,5 mL
Phosphate buffered saline (PBS) Each PBS Powder pouch contains sufficient buffer powder to make 1 liter (10g/pouch). BUF WASH RCNS 1L H ₂ O To dilute	2 pouches	2 pouches
Mounting Medium Dropper vial containing buffered glycerol. Ready to use MM	3 mL	3 mL
Blotters for drying the substrate slides around the wells.	12	12

MATERIAL REQUIRED BUT NOT PROVIDED

- Volumetric pipettes to deliver 10-15 µL
- Small test tubes and rack or U-bottom microtitration plates
- Humid chambers – Petri dishes with moist filter paper
- Coplin staining jars with slide holders
- One liter volumetric flask or graduated cylinder
- Squeeze bottle for PBS
- Cover glasses #1, 24 x 60 mm.
- Distilled water
- Fluorescence microscope
- Timer
- Forceps

REAGENTS AVAILABLE INDIVIDUALLY

ME 0604 ESD	Endomysium Slide of 4 wells
ME 0608 ESD	Endomysium Slide of 8 wells
ME 0904	Anti-IgA conjugate (2 mL)
ME 0062 ESD	Positive Control anti-endomysial AEmA (0,5mL)
HME NEG	Negative control (500 µL)
HME MM	Mounting Medium buffered glycerol (3,5 mL)
HME PBS10	One-liter packet of dry PBS 10g (10x1L)
HME 12CS	24x60mm glass coverslips (100 units)

SPECIMEN COLLECTION

- Collect blood specimens in plain untreated tubes.
- After clot formation, immediately harvest the serum.
- The serum specimens may be stored in a refrigerator (+2°C-+8°C) for up to 48 hours. If longer storage is desired, store at -20°C.
- Turbidity, haemolysis, bacterial growth or drugs capable of fluorescing in the serum may adversely affect the performance of the test. Do not use hemolyzed or lipemic sera.

STABILITY AND STORAGE

- Expiration dates for substrates (slides), conjugate and controls are given for storage at -20°C or below.
- Once thawed the conjugate and controls are stable for 90 days when stored at +2° to +8°C.
- Repeated freezing and thawing should be avoided since it is detrimental to the products.
- The mounting medium is stable until the expiration date when stored in the refrigerator.
- The PBS powder can be stored at room temperature until its expiration date.

PRECAUTIONS

All  Test Kit components are for in vitro diagnostic use ONLY.

All human serum components have been tested by FDA approved method and found negative for antibody to HIV, HCV and non-reactive for HBsAg. However, no test method can offer complete assurance of the absence of these antigens. Care should be exercised when handling ALL human sera containing specimens.

Avoid mouth pipetting.

All controls, conjugate and mounting medium contain sodium azide as preservative in concentration 1:1,000. This preservative is toxic if ingested. In addition, sodium azide may react with lead or copper plumbing to form explosive azides.

On disposal of sodium azide containing solutions, flush with a large volume of water to prevent azide build-up.

Do not substituted others reagent's manufacturers.

Do not use components or reagents beyond their expiration dates. Extended incubation during the wash step or improper washing may result in poor morphology of substrate and increased background fluorescence.

Avoid using pen to identify the slide.

Remove slides from freezer and allow them to reach room temperature – **Don't remove from protective envelopes before use.**

Use each slide immediately after opening the envelope to avoid damaging the substrate by drying.

All reagents must be at room temperature before use (at least 15 mn).

SET-UP

All reagents must be at room temperature before use.

1. Preparation of the wash and dilution buffer

- Reconstitute buffer pack with one liter distilled water. Wait 10 mn before use
- Storage period: 4 weeks at +2°C/+8°C

2. Sera dilution

- Dilute freshly collected and properly stored patient serum 1:10 in PBS. For titrations, set up doubling dilutions, i.e. 1:20, 1:40, etc.
- Mix with vortex

TEST PROCEDURE

1. Incubation of samples

- Identify wells which are going to contain positive and negative control sera and patient's serum dilutions.
- Remove slide from protective envelopes and put into the humid chamber
- In the identified well of the slide, add
 - **40µL of diluted sample, by avoiding to touch the substrate with the tip of the pipette.**
 - 1 drop of ready to use positive control
 - 1 drop of ready to use negative control
- Do the same with each slide.
- Incubate for 30 minutes at room temperature.

Do not let the wells dry out

2. Wash with PBS

- Remove slides from humid chamber.
- Tap each slide quickly on its side to avoid contaminations.
- Dip the slides several times in a beaker containing fresh PBS. **Do not use wash bottle to rinse the slides.**
- Immediately place slides in a slide holder or a coplin jar and wash with PBS and incubate for 10 minutes at room temperature. Agitation is not useful.

Note: If more than 10 slides are tested, wash two time in PBS bath (5 minutes each).

⚠ Extended incubation and improper washing may result in poor morphology of substrate and increased background fluorescence.

3. Incubation of Conjugate

- Remove the slide from PBS
- Eliminate excess of PBS onto some absorbent paper.
- Blot the slide with blotter provided and put into the humid chamber
- Immediately apply 1 drop of conjugate to each well to cover substrates.
- Incubate for 30 minutes at room temperature and **far from light**.

Do not let the wells dry out

4. PBS rinse

- Repeat Step 2 "Wash with PBS".

5. Mount Cover glasses

- Apply 3 or 4 drops mounting medium to each slide
- Cover with cover glasses – Lateral movement and pressure will damage the substrate

6. Reading

It is recommended to examine the day of the test with an appropriately equipped fluorescence microscope (x 10) and to use (x25 or x40) microscopes to confirm a positive fluorescence.

Note: If delay in examining slides is anticipated, seal coverglass with clear nail polish and store in refrigerator.

INTERPRETATION OF RESULTS

Negative reaction

Produces no fluorescence staining of the endomysium of the oesophagus distal substrate. Samples may produce background staining or other reaction due to anti-nuclear or anti-smooth muscle antibodies interference.

Positive reaction

Produces an apple green fluorescence staining of the endomysium or the reticulin components of the smooth muscle filaments of the distal oesophagus substrate.

QUALITY CONTROL

On each run, both positive and negative controls, provided in each kits, should be included.

The negative control should produce no fluorescence staining on the endomysium of oesophagus distal substrate.

Anti-Endomysial antibody (AEmA) positive control should result in an apple green fluorescence staining of the endomysium or the reticulin components of the smooth muscle filaments of distal oesophagus substrate at a fluorescence value of 2+ or higher.

When results obtained are distant from expected values, test run should be repeated. Should atypical observations continue, we recommend to investigate if :

- adequate amount of serum is applied on the wells
- wells were allowed to dry out during incubation of serum or conjugate
- there is any sign of microbial contamination of the liquid reagents due to incorrect storage and handling. If turbidity appears on any of the liquid reagents, discard and use another batch.
- Fluorescence microscope is properly aligned, using correct bulb, correct filters, etc...

LIMITATIONS OF THE TEST

The **Endomysium (AEmA)** test kit is a valuable diagnostic aid. However, it should not be considered diagnostic by itself alone. All tests should be interpreted by a medical authority in light of all patients conditions.

Interference or other auto-antibodies such as anti-nuclear antibodies (ANA) and anti-smooth muscle antibodies (ASMA) may result in the reaction of the muscle layers of distal oesophagus substrate. [bmd](#) offers Multiple Antibody Test Kit to further evaluate the presence and titre of these antibodies.

IgA endomysial antibody titer decreases with time, when patients are put on gluten free diet.

Weak or negative reactions to IgA class endomysial antibodies may be exhibited by a small percentage of celiac patients who are IgA deficient. Testing for IgG class endomysial and or IgG and IgA gliadin antibodies may be beneficial.

Sensitivity of the test kit may be influenced by several factors aside from the components. We recommend user to take note of factors such as optical system, improper storage and handling, technician / observer, etc.

Samples that were screened positive at 1:10 should be further diluted. All positive sera should be titrated to the highest titre using a series of doubling dilutions, i.e. 1:20, 1:40, etc. The end point titre is the highest dilution that produces a positive result.

EXPECTED VALUES

Endomysial antibodies are present in almost 100% patients with celiac disease, 70% of patients with dermatitis herpetiformis, and in about 20% of patients with other enteropathies^(3,4).

There is a small percentage of IgG AEmA that will be negative IgA when screened. A negative result exhibited by a patient with overt clinical symptoms may to be considered for IgG testing.

PERFORMANCE CHARACTERISTICS

Serums of patients tested with the Endomysium (AEmA) [bmd](#) kit:

- patients whose diagnosed disorders are compatible with high rate of anti-endomysium antibodies,
- patients with no symptom of coeliac disease and dermatitis herpetiformis,
- healthy adults

The **sensitivity** of the **Endomysium (AEmA)** [bmd](#) kit was estimated at **95%** in the population of patients with coeliac disease (interval of values reported in the literature: 90%-100%).

The **sensitivity** of the **Endomysium (AEmA)** [bmd](#) kit was estimated at **64%** in the population of patients with dermatitis herpetiformis (interval of values reported in the literature: 60%-70%).

The **specificity** of the **Endomysium (AEmA)** [bmd](#) kit was estimated between **97% and 99%** in the population of patients with other autoimmune diseases (mean of values reported in the literature: 98%).

The **specificity** of the **Endomysium (AEmA)** [bmd](#) kit was estimated at **98%** in the control group (mean of values reported in the literature: 98%).

The comparative studies between Endomysium (AEmA)  kit and the kit of reference showed almost all agreement between the two techniques with comparable performances for sensitivity, specificity and reproducibility on the population tested.

The results of this evaluation on human serums of different specificities show that the **Endomysium (AEmA) ** kit has very high sensitivity and specificity for detection of the anti-endomysium antibodies, and thus it is useful to the diagnosis of the coeliac disease and associated pathologies. The evaluation confirm that the kit is in agreement frequencies of positive and negative results from the scientific literature.

BIBLIOGRAPHY

1. Beutner EH, Chorzelski TP, Kumar V, Leonard J, Krasny S : Sensitivity and specificity of IgA-class anti-endomysial antibodies for dermatitis herpetiformis and findings relevant to their pathogenic significance. *Journal of the American Academy of Dermatology*, Vol 15, N°3, p 464 (1986).
2. Kapuscinska A, Zalewski T, Chorzelski T, Sulej J, Beutner EH, Kumar V, Rossi T: Disease specificity, sensitivity and dynamics of changes in IgA class anti-endomysial antibodies in celiac disease. *J Pediatric Gastroenterol Nutr*, Vol 6, N°4, p529 (1987).
3. Blazer S, Naveh Y, Berant M et al. : Serum IgG antibodies to gliadin in children with celiac disease as measured by an immunofluorescence method. *J Pediatric Gastroenterol Nutr*, Vol 3, p205 (1984).
4. Chorzelski TP, Beutner EH, Sulej J, Tchorzewska H, Jablonski S, Kumar V, Kapuscinska A : IgA anti-endomysium antibody. A new immunological marker for herpetiformis and coeliac disease. *British Journal of Dermatology* Vol 111, p 395 (1984).
5. Leonard JN, Haffenden GP, Fry L : Dermatitis herpetiformis in “ Immunopathology of the skin ”, Beutner EH, Chorzelski TP and Kumar V. Eds, John Wiley and Sons , New York, 3rd Ed. p 433 (1987).
6. Weinstein WM : Latent celiac sprue. *Gastroenterology* , Vol 66 N° 4, p 489 (1974).
7. Unsworth DJ, Leonard JN, Fry L : Anti-reticulín and anti-gliadin antibodies in dermatitis herpetiformis and celiac disease. “ Immunopathology of the skin ”, Beutner EH, Chorzelski TP and Kumar V. Eds, John Wiley and Sons , New York, 3rd Ed. p 455 (1987).
8. Volta U, Molinaro N, Fusconi M et al. IgA anti-endomysial antibody test a step forward in celiac disease screening. *Diag. Disease and Sci.* Vol 36, N°6, 1991.
9. Kumar V, Jarzabek-Chorzelska M, Sulej J, Rajadhyaksha M., and Jablonska S : Tissue transglutaminase and endomysial antibodies - Diagnostic markers of gluten-sensitive enteropathy in dermatitis herpetiformis, *Clinical Immunology*, Vol 98, N°3, p378-382 (2001).
10. Johanet C. Marqueurs sérologiques de la maladie coeliaque. *Spectra Biologie*, Vol 17, N°99 (1998).
11. Ghedira I; et al. : Intérêt des anticorps anti-endomysium, anti-réticuline et anti-gliadine dans le diagnostic de la maladie coeliaque chez l'adulte. *Ann Biol Clin*, 57 : 717-9 (1999).



SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	40 µL	Diluted samples	30 mn at R.T
	1 drop	Positive control (ready to use)	
	1 drop	Negative control (ready to use)	
<i>Do the same with each slide</i>			
PBS WASH	⇒ Tap each slide quickly on its side to avoid contaminations ⇒ Dip the slide several times in a beaker containing fresh PBS ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>Eliminate excess of PBS on Teflon with the blotter(*) - Dry uniformly and without friction using a roller</i>			
INCUBATION OF CONJUGATE	1 drop to each well	Fluorescein conjugated Affinipure Goat anti-human IgA (ready to use)	30 mn at R.T in the dark
<i>Verify that all wells are perfectly covered</i>			
PBS RINSE	⇒ Tap quickly each slide on its side to eliminate the conjugate. ⇒ Dip the slide several times in a beaker containing fresh PBS. ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>DO NOT ELIMINATE EXCESS OF PBS (Teflon) (*)</i>			
MOUNT COVER GLASSES	3-4 drops to each slide	Mounting Medium (ready to use)	-
READING	Examine with an appropriately equipped fluorescence microscope (x 10 – x25 - x 40).		

(*)Do not let wells dry out. Every deposit must be made at once.

R.T : Room Temperature

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Slide with "n" wells		Negative control		Mounting Medium
	Fluorescent antibody conjugate		Positive control		Wash Buffer
	Immunofluorescent assay				Reconstitute with 1L distilled water

BioMedical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com