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BIBLIOGRAPHIE

1. Bean SF, et al
Cicatrical pemphigoid, immunofluorescent studies. Arch Dermatol 106,195 (1972)
2. Beutner EH, Jordon RE
Immunofluorescence staining. Proc Soc Exp Biol Med 117, 505 (1964)
3. Beutner EH, et al
Immunopathology of skin. Churchill
4. Chorzelski TP, et al
Clinical significance of autoantibodies in pemphigus. Arch Dermatol 93,570 (1966)
5. Chorzelski T, et al
Autoantibodies in pemphigoid. Dermatologica 136,325 (1968)
6. Cooperative Study
Uses of immunofluorescence tests of skin and sera. Arch Dermatol 111, 371 (1975)
7. Fellner MJ
Immunology of skin diseases. Elsevier, New York, 1980
8. Fry L, Seah PP
Immunological aspects of skin diseases. John Wiley & Sons, New York, 1974
9. Hodge L, et al
Bullous pemphigoid : The frequency of mucosal involvement and concurrent malignancy related to indirect immunofluorescence findings. Brit J Dermatol 105,65 (1981)
10. Nishikawa T, et al
Binding of bullous pemphigoid antibodies to basal cells. J Invest Dermatol 74,389 (1975)
11. Tufanelli DL
Cutaneous immunopathology : Recent Observations. J Invest Dermatol 65,143 (1975)
12. Valenzuela R, Bergfeld WF, Deodhar SD
The interpretation of immunofluorescence patterns in skin disease. Amer Soc Clin Pathol Press, Chicago, 1984
13. Richard C ; Joly P (1999)
Anticorps anti-peau et maladies bulleuses auto-immune - Autoanticorps marqueurs de maladies auto-immune. [bnd](#) Editions p.349.
14. Humbel RL (1994)
Maladies auto-immunes de la peau – Auto-anticorps et maladies auto-immunes Editions Option Bio Elsevier version N°1 – p.211.

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




Coffret

INTENDED USE

The **Anti-Skin (IF) – Monkey esophagus** kit is an indirect Immunofluorescence test used for the qualitative and semi-quantitative detection of skin disease's auto-antibodies using microtome sections of the monkey oesophagus as substrate in human serum ^(2,6,10,13). This test is a particular value as an aid in the diagnosis and management of patients with one of the above-mentioned skin diseases, i.e. pemphigus vulgaris or bullous pemphigoid.

DIAGNOSTIC VALUE

Skin lesions may reflect a wide variety of systemic diseases ^(3, 7, 8, 9, 12). The skin diseases with a primary immune disorders of the pathogenesis are referred to as autoimmune disorders of the skin and include pemphigus vulgaris and bullous pemphigoid.

Pemphigus vulgaris is a blistering disorder of the skin and mucous membranes affecting mostly individuals between 40 and 60 years of age. Circulating pemphigus antibodies, reacting with the intercellular membranes of the epithelial substrate, can be detected in 80% to 90% of the patients ⁽⁵⁾.

Bullous pemphigoid is a chronic bullous disorder of the elderly. Circulating pemphigoid antibodies, reacting with the basement membrane of the squamous epithelial cell layer, can be found in approximately 70% of the patients ^(1,5,10). Circulating serum antibodies reacting with epithelial tissues, are not seen in other vesiculo-bullous dermatoses. Detection of circulating antibodies is therefore valuable in the differential diagnosis of these disorders ⁽⁴⁾.

The circulating antibodies in both of the above mentioned autoimmune disorders have been shown to belong, principally, to the IgG class. There is a positive relationship between the antibody titer and the disease activity.

This microphotography shows in a stronger magnification the typical pattern of esophagus in longitudinal section.

- E :** Stratified squamous epithelium
- C :** Chorion
- MM :** Muscularis mucosae
- S :** Submucosa
- CM :** Muscularis : circular layer
- ML :** Longitudinal muscularis
- A :** Adventitia

PRINCIPLE OF THE TEST

Both biopsy and serum studies are of critical value in the diagnosis. The **Anti-Skin (IF) – Monkey esophagus** kit is an immunofluorescence test.

First:

The substrate sections of monkey oesophagus on microscope slides are overlaid with diluted patient's serum and incubated to allow the anti-skin (anti-epidermal) antibodies present in positive sera – to bind to the antigen sites in the oesophageal tissue.

Second:

Sites where the anti-skin interacted with the antigenic constituents of the substrate are made visible by binding to them fluorescently labelled antihuman antibodies (the conjugate).

In the fluorescence microscope, a positive reaction results in an apple green fluorescence pattern of the characteristic structural features of the substrate sections.

COMPONENTS

	ME 0748 ES
Monkey esophagus substrate slides (sections of monkey oesophagus in multiple well-slides) <div style="text-align: right;">SORB SLD 4</div>	12 x 4 wells
Anti-IgG conjugate Fluorescently labelled anti-human IgG globulin in phosphate buffered saline (PBS) with 1% bovine serum albumin and with evans blue. Ready to use <div style="text-align: right;">CONJ DTAF IgG</div>	1 x 2 mL
Positive control anti-intercellular substance Dropper vial containing positive patient sample exhibiting fluorescence pattern, characterizing intercellular substance reaction (ASA-Pemphigus), using immunofluorescence method. Ready to use <div style="text-align: right;">CONTROL ICS</div>	1 x 0,5 mL
Negative Control Dropper vial containing negative patient sample producing no fluorescence staining with the oesophagus substrate. Ready to use <div style="text-align: right;">CONTROL -</div>	1 x 0,5 mL
Phosphate buffered saline (PBS) Each PBS Powder pouch contains sufficient buffer powder to make 1 liter (10g/pouch). To dilute <div style="text-align: right;">BUF WASH RCNS 1L H₂O</div>	2 pouches
Mounting Medium Dropper vial containing buffered glycerol. Ready to use <div style="text-align: right;">MM</div>	3 mL
Blotters for drying the substrate slides around the wells.	12

MATERIAL REQUIRED BUT NOT PROVIDED

- Volumetric pipettes to deliver 10-15 µL
- Small test tubes and rack or U-bottom microtitration plates
- Humid chambers – Petri dishes with moist filter paper
- Coplin staining jars with slide holders
- One liter volumetric flask or graduated cylinder
- Squeeze bottle for PBS
- Cover glasses #1, 24 x 60 mm.
- Distilled water
- Fluorescence microscope
- Timer
- Forceps

REAGENTS AVAILABLE INDIVIDUALLY

ME 0063 ES	Positive control anti-intercellular substance ASA-Pemphigus (0,5 mL)
ME 0062 ES	Positive control anti-basal membrane ASA-Pemphigoide (0,5 mL)
ME 902.2	Anti-IgG conjugate (2 mL)
HME NEG	Negative control (500 µL)
HME MM	Mounting Medium buffered glycerol (3,5 mL)
HME PBS10	One-liter packet of dry PBS 10g (10x1L)
HME 12CS	24x60mm glass coverslips (100 units)

SPECIMEN COLLECTION

- Collect blood specimens in plain untreated tubes.
- After clot formation, immediately harvest the serum.
- The serum specimens may be stored in a refrigerator (+2°C-+8°C) for up to 48 hours. If longer storage is desired, store at -20°C.
- Turbidity, haemolysis, bacterial growth or drugs capable of fluorescing in the serum may adversely affect the performance of the test. Do not use hemolyzed or lipemic sera.

STABILITY AND STORAGE

- Expiration dates for substrates (slides), conjugate and controls are given for storage at -20°C or below.
- Once thawed the conjugate and controls are stable for 90 days when stored at +2° to +8°C.
- Repeated freezing and thawing should be avoided since it is detrimental to the products.
- The mounting medium is stable until the expiration date when stored in the refrigerator.
- The PBS powder can be stored at room temperature until its expiration date.

PRECAUTIONS

All [Lmd](#) Test Kit components are for in vitro diagnostic use ONLY.

All human serum components have been tested by FDA approved method and found negative for antibody to HIV, HCV and non-reactive for HBsAg. However, no test method can offer complete assurance of the absence of these antigens. Care should be exercised when handling ALL human sera containing specimens. Avoid mouth pipetting.

All controls, conjugate and mounting medium contain sodium azide as preservative in concentration 1:1,000. This preservative is toxic if ingested. In addition, sodium azide may react with lead or copper plumbing to form explosive azides.

On disposal of sodium azide containing solutions, flush with a large volume of water to prevent azide build-up.

Do not substituted others reagent's manufacturers.

Do not use components or reagents beyond their expiration dates. Extended incubation during the wash step or improper washing may result in poor morphology of substrate and increased background fluorescence.

Avoid using pen to identify the slide.

Remove slides from freezer and allow them to reach room temperature – **Don't remove from protective envelopes before use.**

Use each slide immediately after opening the envelope to avoid damaging the substrate by drying.

All reagents must be at room temperature before use (at least 15 mn).

SET-UP

All reagents must be at room temperature before use.

1. Preparation of the wash and dilution buffer

- Reconstitute buffer pack with one liter distilled water. Wait 10 mn before use
- Storage period: 4 weeks at +2°C/+8°C

2. Sera dilution

- Dilute freshly collected and properly stored patient serum 1:10 in PBS.
Ex : 20µL serum + 180 µL PBS
For titrations, set up doubling dilutions, i.e. 1:20, 1:40, etc.
- Mix with vortex

TEST PROCEDURE

1. Incubation of samples

- Identify wells which are going to contain positive and negative control sera and patient's serum dilutions.
- Remove slide from protective envelopes and put into the humid chamber
- In the identified well of the slide, add
 - **40µL** of diluted sample, **by avoiding to touch the substrate with the tip of the pipette.**
 - 1 drop of ready to use positive control
 - 1 drop of ready to use negative control
- Do the same with each slide.
- Incubate for 30 minutes at room temperature.

Do not let the wells dry out

2. Wash with PBS

- Remove slides from humid chamber.
- Tap each slide quickly on its side to avoid contaminations.
- Dip the slides several times in a beaker containing fresh PBS. **Do not use wash bottle to rinse the slides.**
- Immediately place slides in a slide holder or a coplin jar and wash with PBS and incubate for 10 minutes at room temperature. Agitation is not useful.

Note: If more than 10 slides are tested, wash two time in PBS bath (5 minutes each).

⚠ Extended incubation and improper washing may result in poor morphology of substrate and increased background fluorescence.

3. **Incubation of Conjugate**

- Remove the slide from PBS
- Eliminate excess of PBS onto some absorbent paper.
- Blot the slide with blotter provided and put into the humid chamber
- Immediately apply 1 drop of conjugate to each well to cover substrates.
- Incubate for 30 minutes at room temperature and **far from light.**

Do not let the wells dry out

4. **PBS rinse**

- Repeat Step 2 “Wash with PBS”.

5. **Mount Cover glasses**

- Apply 3 or 4 drops mounting medium to each slide
- Cover with cover glasses – Lateral movement and pressure will damage the substrate

6. **Reading**

It is recommended to examine the day of the test with an appropriately equipped fluorescence microscope (x 10) and to use (x25 or x40) microscopes to confirm a positive fluorescence.

Note: If delay in examining slides is anticipated, seal coverglass with clear nail polish and store in refrigerator.

INTERPRETATION OF RESULTS

All positive sera should be titrated to the highest titer using a series of doubling dilutions. If more than one antibody is present in the patient's serum, interference may occur. Titration of such serum may reveal first one and then another autoantibody.

Specific staining of the intercellular substance between the squamous epithelial cells is considered a positive test for pemphigus antibodies.

Specific staining of the basement membrane is considered a positive test for bullous pemphigoid.

Auto-antibodies detected	Cut off
Intercellular membranes	> 1/10
Basement membrane	> 1/10

QUALITY CONTROL

On each run, both positive and negative controls, provided in each kits, should be included.

When results obtained are distant from expected values, test run should be repeated. Should atypical observations continue, we recommend to investigate if :

- adequate amount of serum is applied on the wells
- wells were allowed to dry out during incubation of serum or conjugate
- there is any sign of microbial contamination of the liquid reagents due to incorrect storage and handling. If turbidity appears on any of the liquid reagents, discard and use another batch.
- Fluorescence microscope is properly aligned, using correct bulb, correct filters, etc...

LIMITATIONS OF THE TEST

This test is a valuable diagnostic aid. However, it should not be considered diagnostic by itself alone. Positive results may be obtained with sera of apparently healthy people as a result of many factors. Therefore, all test results should be interpreted by a medical authority in light of all known patient’s conditions.

PERFORMANCE CHARACTERISTICS

Assay units:

Indirect immunofluorescence (IIF) was assayed using accepted arbitrary units of intensity under fluorescence microscopy. The units from greatest to least intensity are; 4+, 3+, 2+, 1+, +/-, -.

Titer reproducibility:

Four serum samples were tested on three different lot numbers of bmd kits on three different occasions. All positive values titred to negative at the same dilution for each of the three replicates.

Normal samples:

Sera from 743 healthy blood donors with no history of autoimmune disease were tested on bmd **Anti-Skin (IF) – Monkey esophagus** kits. Sera were diluted 1:20. A negative reading at 1:20 was considered negative. Positive readings of 1+ or higher at dilutions of 1:20 or greater were scored as positive. These data are reported in the table below as the observed percent of sera tested that scored as positive for healthy subjects.

Sera from patients with diagnosed disorders:

Sera from 1019 patients with variety of diagnosed disorders were tested with bmd **Anti-Skin (IF) – Monkey esophagus** kits. Sera were diluted from 1:10 to 1:320. The observed and expected values of sera are notices in the hereafter tables.

<i>Observed percent (%) of sera scored as positive</i>	
<i>Pemphigoid, bullous</i>	68
<i>Pemphigus vulgaris</i>	87
Healthy subjects	4

<i>Expected percent (%) of sera scored as positive</i>	
<i>Pemphigoid, bullous</i>	70
<i>Pemphigus vulgaris</i>	85
Healthy subjects	2

BIBLIOGRAPHY

1. Bean SF, et al
Cicatrical pemphigoid, immunofluorescent studies. Arch Dermatol 106,195 (1972)
2. Beutner EH, Jordon RE
Immunofluorescence staining. Proc Soc Exp Biol Med 117, 505 (1964)
3. Beutner EH, et al
Immunopathology of skin. Churchill
4. Chorzelski TP, et al
Clinical significance of autoantibodies in pemphigus. Arch Dermatol 93,570 (1966)
5. Chorzelski T, et al
Autoantibodies in pemphigoid. Dermatologica 136,325 (1968)
6. Cooperative Study
Uses of immunofluorescence tests of skin and sera. Arch Dermatol 111, 371 (1975)
7. Fellner MJ
Immunology of skin diseases. Elsevier, New York, 1980
8. Fry L, Seah PP
Immunological aspects of skin diseases. John Wiley & Sons, New York, 1974
9. Hodge L, et al
Bullous pemphigoid : The frequency of mucosal involvement and concurrent malignancy related to indirect immunofluorescence findings. Brit J Dermatol 105,65 (1981)
10. Nishikawa T, et al
Binding of bullous pemphigoid antibodies to basal cells. J Invest Dermatol 74,389 (1975)
11. Tufanelli DL
Cutaneous immunopathology : Recent Observations. J Invest Dermatol 65,143 (1975)
12. Valenzuela R, Bergfeld WF, Deodhar SD
The interpretation of immunofluorescence patterns in skin disease. Amer Soc Clin Pathol Press, Chicago, 1984
13. Richard C ; Joly P (1999)
Anticorps anti-peau et maladies bulleuses auto-immune - Autoanticorps marqueurs de maladies auto-immune. bmd Editions p.349.
14. Humbel RL (1994)
Maladies auto-immunes de la peau – Auto-anticorps et maladies auto-immunes Editions Option Bio Elsevier version N°1 – p.211.





SUMMARY OF METHOD

	AMOUNT TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
INCUBATION OF SAMPLES	40 µL	Diluted samples	30 mn at R.T
	1 drop	Positive control (ready to use)	
	1 drop	Negative control (ready to use)	
<i>Do the same with each slide</i>			
PBS WASH	⇒ Tap each slide quickly on its side to avoid contaminations ⇒ Dip the slide several times in a beaker containing fresh PBS ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>Eliminate excess of PBS on Teflon with the blotter(*) - Dry uniformly and without friction using a roller</i>			
INCUBATION OF CONJUGATE	1 drop to each well	Fluorescein conjugated Affinipure Goat anti-human IgG (ready to use)	30 mn at R.T in the dark
<i>Verify that all wells are perfectly covered</i>			
PBS RINSE	⇒ Tap quickly each slide on its side to eliminate the conjugate. ⇒ Dip the slide several times in a beaker containing fresh PBS. ⇒ Place slides in a slide holder or a coplin jar and wash with PBS. N.B : Agitation is not useful.		10 mn at R.T
<i>DO NOT ELIMINATE EXCESS OF PBS (Teflon) (*)</i>			
MOUNT COVER GLASSES	3-4 drops to each slide	Mounting Medium (ready to use)	-
READING	Examine with an appropriately equipped fluorescence microscope (x 10 – x25 - x 40).		

(*) Do not let wells dry out. Every deposit must be made at once.

R.T : Room Temperature

SYMBOLS USED

	Biological risks		Temperature limitation		Catalogue Number
	Consult instructions for use		In Vitro Diagnostic Medical Device		Batch Code
	Number of tests		Use by		EC Declaration of Conformity
	Slide of 4 wells		Negative control		Mounting Medium
	Fluorescent antibody conjugate		Positive control		Wash Buffer
	Immunofluorescent assay				Reconstitute with 1L distilled water

BioMédical Diagnostics SA

Office
 Actipole 25
 4 bld de Beaubourg
 77435 Marne la Vallée Cx2
 France

Tel : 33 1 64 62 10 12
 Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
 Internet : www.bmd-net.com