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REF

ME 902.2 2 ml
ME 902.8 8 ml

Français

DÉFINITION

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Anglais

INTENDED USE

Fluorescein conjugated Affinipure Goat anti-human IgG.

SPECIFICITY

The antibody reacts with the heavy chains on human IgG and with light chains common to most human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins, but antibodies may cross-react with immunoglobulins from other species.

COMPONENTS

Appearance: Nonsterile liquid, ready to use
Fluorophore: Dichlorotriazinyl Amino Fluorescein (DTAF)
Antibody concentration: 3,9 mg/ml
Buffer solution: 0,02M Phosphate Buffer + 0,25M NaCl pH 7.6
Stabilizer: 15mg Bovine Serum Albumin (BSA) added per ml
Preservative: < 0,1 % sodium azide
Fluorescein/Protein ratio: 14.2 ug per mg antibody

PROCEDURE

This reagent must be at room temperature before use. Apply this one in the appropriate wells corresponding to the chosen substrate (cells or organ sections), according to the procedure recommended for every substrate. Generally, the appropriate volume is 1 drop/well

PRECAUTIONS

- This reagent contains sodium azide; upon disposal of this reagent, flush with large volumes of water to avoid metal azides formation with lead and copper plumbing.
- Do not use this product beyond their expiration date.
- Reagent to protect from light.
- Do not use this reagent if signs of contamination are present.

STABILITY AND STORAGE

Expiration dates for conjugate is given for storage at -20°C or below. Once thawed, it is stable for 90 days when stored between +2°C and +8°C.

LIMITATION OF USE

It should be used for in vitro diagnostic. All tests should be interpreted by a medical authority in light of all patients' conditions.

BioMédical Diagnostics SA

Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tél : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com