FIDIS™ dsDNA

REF

MX 005

96

DEFINITION

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ISENBERG et al. Clinical laboratory assays for measuring anti-dsDNA antibodies. Were are we now? Lupus 2002, 11, 797-800.

RAHMAN A et al. Anti-DNA antibodies – structure and function. Lupus 2002, 11, 776-779.

RAHMAN A et al. Molecular expression systems for anti-DNA antibodies-1–2. Lupus 2002, 11, 824-832, 833-842.







RAVIRAJAN CT et al. An analysis of clinical disease activity and nephritis associated serum autoantibody profiles in patients with systemic lupus erythematosus: a cross-sectional study. Rheumatology 2001, 40, 1405-1412.

SMEENK RJ et al. dsDNA autoantibodies. In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234.

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REF **MX 005**

INTENDED USE

The **FIDIS™ dsDNA** ([bmd](#)) is a homogeneous fluorescent-based microparticles immunoassay using flow cytometry readings. It is designed for the quantitative detection of autoantibodies directed against **double stranded DNA (dsDNA)**.

The **FIDIS™ dsDNA** ([bmd](#)) kits can be used with a dispensing/diluting device, the **CARIS™ system**.

SUMMARY AND EXPLANATION

Connective diseases are systemic auto-immune diseases which are, most often, selectively directed against a specific organ. Their classification is based upon clinical and biological data. The detection of anti-nuclear antibodies, which often appears in the patient serum before any definite symptomatic signs, will allow the identification of the specific connective disease.

Among them, Anti-dsDNA antibodies are recognized as the major serologic marker of SLE. Their specificity and their sensitivity give them a high diagnostic value. Therefore they are a part of the clinical and biological criteria established in 1982 by the American Rheumatism Association (ARA) for the diagnosis of the SLE.

The detection of anti-dsDNA antibodies is particularly useful in two different ways : as an aid to the diagnosis of SLE and as a tool to monitor the course of the disease.

For the second purpose repeated serum sampling of individual patients can be very informative about the clinical course of the disease because a clear-cut relationship exists between anti-DNA and diseases activity : flares of SLE are generally preceded by a rise in anti-dsDNA levels, followed by a steep drop during the exacerbation (particularly in nephritis).

Furthermore different treatments of patients have varying influences on anti-dsDNA levels and can be adapted by a regular follow-up of these antibodies.

ASSAY PRINCIPLE

FIDIS™ dsDNA ([bmd](#)) kits are based on the use of distinct uniform size color-coded microspheres and a benchtop flow cytometer interfaced to digital signal processing hardware and software. A red diode laser beam in the flow cytometer classifies each set of microspheres on the basis of its unique fluorescence intensity (red to orange) thus identifying which analyte is being tested. At the same time, a green laser beam illuminates the external second molecule

fluorescence to quantify the reaction related to the specific analyte.

FIDIS™ dsDNA ([bmd](#)) kits allow the detection of **anti-dsDNA** antibodies.

The antigen, recombinant dsDNA, is covalently coupled to an individual set of microspheres through its surface functional groups.

The test is performed in a 96 well microplate with a filtering membrane at the bottom of the wells.

- In the first step, the sample is distributed in each well containing the microspheres mixture. If this sample contains the suspected antibodies, this antibody bind to the dsDNA antigen on the microspheres.
- After incubation, a wash step using a filtration process removes the unbound antibodies.
- A phycoerythrin labeled **anti-human IgG conjugate** is then added that binds to the previously bound antibodies.
- A final wash step allows to stop the reaction.
- The reaction is then directly measured by the flow cytometer, which differentiates each set of microspheres according to its fluorescence color while simultaneously measuring the average fluorescence emitted by the conjugate.
- A calibration system allows the determination of the titer (IU/mL) of each sample by interpolation.

ADDITIONAL MATERIAL – NOT SUPPLIED

- Precision pipettes capable of accurately delivering 5µL to 1000µL
- Multichannel Pipettes or dispensers capable of accurately delivering 40µL to 300µL or 5µL to 2 mL
- Vortex mixer
- Laboratory timer to monitor incubation steps
- Ultrasonic bath
- Absorbent towel
- Serological pipettes
- Microplate sealing films

SPECIMEN COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

REAGENTS SUPPLIED

96 wells microplate with filtering membrane and lid. MP	1 plate
Vial (A) of color-coded microsphere set sensitized by dsDNA. <u>Lyophilized</u> (to be reconstituted with the buffer named D) MICROSPHERES	sq 6mL
Vial (B1) of sample dilution buffer (white vial) <u>Ready to use</u> DIL SPE	2 x 115mL
Vial of calibrator* <u>Ready to use</u> Cal <i>The titer is printed on the vial label</i>	1 x 1,5mL
Vial of positive control concentrated. This control has a standard reactivity, which provides evidence of the proper reagents activity and proper assay performance. <u>To be diluted</u> CONTROL + <i>Expected values are printed on the vial label.</i>	1 x 250 µL
Vial of negative control concentrate <u>To be diluted</u> CONTROL -	1 x 250µL
Vial of anti-human IgG coupled to phycoerythrin <u>Ready to use</u> CONJ IgG	1 x 12mL
Vial (C1) of washing buffer (black vial) <u>Ready to use</u> BUF WASH	1 x 100mL
Vial (D) of reconstitution buffer for the microsphere set <u>Ready to use</u> BUF MICROSPHERES	1 x 6mL

Calibrator titer for dsDNA (IU/ml) is standardized against the WHO international reference: Wo/80.

STABILITY AND STORAGE

- Store reagents in their original packaging at +2°C to +8°C.
- Do not freeze reagents.
- Do not use kits beyond the expiration date.
- After use, store all components immediately back at +2°C to +8°C.

CAUTION

- Reagents in solution contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
- Avoid to use reagents if signs of contamination or other visible changes occur
- The **FIDIS™ dsDNA** (cmd) kit has been developed according to CE Directives 67/548/EEC and 1999/45/EC regarding classification, packaging and labelling of dangerous preparations.
- The **FIDIS™ dsDNA** kit (cmd) has been optimized for use as described in this procedure. Do not substitute other manufacturer reagents. Dilution or alteration of these reagents may also alter the performance of the test. Follow thoroughly the test procedure to insure optimal performance.
- *Calibrators and controls are from human origin. The human sera used in the preparation of these products were tested and found non-reactive for antibodies to HIV-1, HIV-2, anti-HCV and Hepatitis B antigen. Because no test method can offer complete assurance that infectious agents are absent, handle as if capable of transmitting infectious diseases.*

TEST SET-UP

FIDIS™ Washer

☞ Check the tubing connecting the pump to the stand and the manometer setting (wheel completely closed).


Using FIDIS™ Analyzer and MLX-BOOSTER Software

☞ See the User's Manual provided with the **FIDIS™ Instrument** for detailed instructions on running the equipment and for calculation. For additional information and/or trouble shooting problems, please contact cmd subsidiaries or distributors.

1. Start the run as described in the User's Manual.
NOTE: The FIDIS™ takes 30 minutes to warm up after being turned on. A new warm up is necessary after 4 hours of system inactivity.
2. Calibration and controls are described in the User's Manual. Calibration and controls should be routinely performed once per month and each time a new lot of Sheath fluid reagent is used in order to insure optimal instrument performance.
Calibration should also be performed when the temperature is out of the determined range shown on the run batch screen.
3. Programming a batch or a multibatch of test protocol on the FIDIS™.


Example 1:

- The microplate is totally blank (no well has yet been used): Use A1 as the 1st well.

Example 2:

- 16 wells have already been used: Use A3 as the 1st well.

4. Load the microplate into the plate holder of the FIDIS™ as described in the User’s Manual.
5. Analyze the results according to the User’s Manual.
6. When finished for the day, perform the sanitizing and soak operations, prior to turning the Analyzer off according to the shutting down procedure described in the User’s Manual.

Reagents Preparation

☞ Remove the individual components from storage, allow them to warm up to **room temperature** (+18°C to +25°C), and mix them well.

1. Microsphere preparation

- Reconstitute the microsphere solution by adding the vial **D** of reconstitution buffer to the vial **A** of lyophilized microspheres. Wait 5 minutes and vortex.
- Stable 2 months at +2°C to +8°C, after reconstitution.

2. Preparation of samples and controls.

- Dilute the samples and controls at **1:201** in dilution buffer (**B1**).
Ex.: 10µL sample in 2000µL dilution buffer (**B1**).
- **Vortex vigorously.**

3. Assay configuration

Use the work-sheet included in the kit to identify the location of the samples.

a. When setting up the test, systematically take into account the following well requirements: ⇒ See examples below

- ✗ 1 “reagent-blank” well
- ✗ 1 well for negative control
- ✗ 1 well for positive control
- ✗ 2 “calibrator” wells

b. Calculation of the correct number of wells necessary and their location.

In the following examples, 2 different configurations are described, according to the number of wells needed, and the availability of unused wells. The sample dispensing must be systematically carried out in a column. Do not leave empty wells.

The same microplate can be used for more than one serie of tests if unused wells are protected by microplate sealing film. Any single well can only be used once.

To avoid any error in using a well more than once, identify the used wells with a marker (see example 2 below).

c. Programming of the test protocol on the FIDIS™

☞ Refer to the FIDIS™ user’s manual: “Creating a batch or a multi-batch”.

ASSAY PROCEDURE

1. Microspheres dispensing

- Before to start the assay and if it is necessary, protect the unused wells with a microplate sealing film.
- **Vortex the microspheres reagent vigorously for 20s** and dispense 50µL in each well to be used.

PLEASE NOTE:

- If CARIS™ is used, the microspheres should be sonicated for 5 minutes.
- If a dispensing/diluting device is used (like CARIS™), the procedure steps 1 and 2 should be reversed:
⇒ dispense diluted sera in first and then, dispense microspheres.

2. Addition and Incubation of the samples

- Dispense 100µL of sample dilution buffer (**B1**) for reagent blank.
- Dispense 100µL of prediluted controls, of ready to use calibrator and of prediluted samples.
- Cover the plate and incubate 30 minutes at room temperature, away from direct sunlight and without shaking.

3. Wash step 1

Wash the plate 2 times with washing buffer (**C1**) using filtration unit:

- a. Remove the microplate lid (verify if the vacuum break valve is closed) and place it on the washer.
- b. Start the pump. Stop it when the whole incubation liquid has been removed in totality and open the vacuum break valve. Close the vacuum break valve.

- c. Dispense 300µL of washing buffer **(C1)** in all the used wells.
- d. Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- e. Dispense 300µL of washing buffer **(C1)** in all the used wells.
- f. Start the pump. Stop it after 5 additional seconds when all the buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- g. Remove the microplate from the washer and remove the residual buffer by tapping vigorously the plate 10 times on an absorbent towel.
- h. Place again the microplate on the washer and start the pump for 5 seconds. Stop it and quickly open the vacuum break valve.
- i. Remove the microplate from the washer and remove the residual buffer by blotting vigorously the plate on an absorbent towel.
- j. Place the plate on a totally dry surface before starting the conjugate incubation.

4. Incubation of the conjugate

- Dispense 100µL of ready to use conjugate in each well used.
- Cover the plate and incubate 30 minutes at room temperature away from direct sunlight and without shaking. The incubation time starts after the conjugate has been added to all wells. If this timing is not followed, the results might be erroneous.

Note : the conjugate is photosensitive.
 ⇒ After using, close the vial carefully.

5. Wash step 2

Wash the plate using the filtration unit on 1 cycle with dilution buffer **(B1)**.

Use two different reagent reservoirs for the 2 washing steps. Either clean and dry correctly the reservoir between both wash steps (avoid to mix both B1 and C1 buffers).

- a. Remove the microplate lid (verify that the vacuum break valve is closed) and place it on the washer.
- b. Start the pump. Stop it when the whole incubation liquid has been removed and open the vacuum break valve. Close the vacuum break valve.
- c. Dispense 100µL of dilution buffer **(B1)** in all the used wells.
- d. Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- e. Remove the microplate from the washer and remove the residual buffer by blotting the plate on an absorbent towel.

- f. Dispense 100µL of dilution buffer **(B1)** in all the used wells. Place the plate on a totally dry surface at room temperature.

The reading must be doing in the hour following adding the solution (B1). During this time, keep the plate at room temperature away from direct sunlight.

6. Reading of the test

☞ Follow the **MLX-BOOSTER™** and **FIDIS™** user's manual: Processing batches.

VALIDATION CRITERIA OF THE RESULTS

Calibrator, negative and positive controls have to be run with every batch of samples to ensure that all reagents and procedures performed properly.

In order to validate the results, all the criteria listed below must be met. Otherwise, the test is invalid and must be repeated.

- a. The positive controls should show a value within the limits printed on the corresponding vial labels.
- b. The negative controls should be less than 30 IU/mL.

The microsphere set contains Internal standard beads allowing to verify the presence of serum in wells, as well as the good respect of the test procedure. Wells presenting a not corresponding fluorescent signal will be no valid with MLX BOOSTER™ software.

INTERPRETATION OF RESULTS

International Units (IU/ml) for the dsDNA	< 30 IU/mL	30 – 40 IU/mL	> 40 IU/mL
Interpretation	Negative	Borderline(*)	Positive

(*) *Borderline results should be controlled on a second sample and the interpretation of the results should be done in the frame of additional testing and taking into account the clinical status of the patient.*

Each laboratory should establish and maintain its own references (normal) range values, based on the patient population and other local factors.

CALCUL OF RESULTS

☞ The results are automatically calculated by MLX-BOOSTER™ software and can be printed for each assay.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

The serums from Immunodeficient patients will give a no valid result.

QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities. Immunotrol I and IV multiparametric controls ([bmd](#), Cat. n° HM036 and HM051) contain human auto-antibodies directed against dsDNA specificity. To be assayed in the same manner as the unknown samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

FIDIS™ dsDNA ([bmd](#)) is a derived product from FIDIS™ Connective ([bmd](#)). These two kits are manufactured under the same process, particularly for coupling procedure of microspheres to dsDNA.

FIDIS™ dsDNA ([bmd](#)) is designed for the specific follow-up of anti-dsDNA autoantibodies levels after detection of a positive sample by FIDIS™ Connective. Interpretation values are equal for the two kits.

INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH OTHER METHODS

FIDIS™ dsDNA kit ([bmd](#)) has been compared to FIDIS Connective kit ([bmd](#)):

The study was performed on 182 samples :

- 54 samples related to systemic autoimmune diseases
- 128 samples from blood donors


Results

dsDNA		FIDIS™ Connective	
		Positive	Negative
FIDIS™ dsDNA	Positive	54	0
	Negative	0	128

Relative sensitivity : 100%
 Relative specificity : 100%
 Concordance : 100%

CORRELATION

The linear regression analysis of the two products showed that both are equivalent.

Reproducibility

Estimated on three samples: weakly positive, positive and strongly positive.

Within-run (10 tests in the same run)		Between-run (5 tests in 5 different runs)	
Mean value	CV (%)	Mean value	CV (%)
44	3.4	41	5.5
104	4.0	102	5.9
233	3.9	221	2.5

BIBLIOGRAPHY

CONRAD K. et al. Autoantibodies – Diagnostic, pathogenic and prognostic relevance. Clin. and Exp. Rheumatol. 1997 ;15, 457-465

HOLLINGSWORTH P.N. et al. Antinuclear antibodies In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.74-90.

HOMBURGER H.A. Cascade testing for autoantibodies in connective tissue disease. Mayo Clin. Proc., 1995, 70, 183-184.

ISENBERG et al. Clinical laboratory assays for measuring anti-dsDNA antibodies. Were are we now? Lupus 2002, 11, 797-800.

RAHMAN A et al. Anti-DNA antibodies – structure and function. Lupus 2002, 11, 776-779.

AHMAN A et al. Molecular expression systems for anti-DNA antibodies - 1 – 2. Lupus 2002, 11, 824-832, 833-842.







RAVIRAJAN CT et al. An analysis of clinical disease activity and nephritis associated serum autoantibody profiles in patients with systemic lupus erythematosus: a cross-sectional study. Rheumatology 2001, 40, 1405-1412.

SMEENK RJ et al. dsDNA autoantibodies. In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234.

SUMMARY OF THE TEST PROCEDURE

	VOLUME TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
Dispense 50µL of microspheres reagent in each well to be used			
Incubation of the samples	100µL	Dilution buffer (B1) : Reagent blank	30 min. Room temperature
	100µL	Diluted negative control (1:201)	
	100µL	Diluted positive control (1:201)	
	100µL	Calibrator (in duplicate)	
	100µL	Diluted samples (1:201)	
Washing 1	Wash twice in washing buffer (C1) (300µL/well)		
Incubation of the conjugate	100µL	Ready to use Anti-IgG conjugate	30 min Room temperature
Washing 2	Wash 1 time in Dilution buffer (B1) (100µL/well)		
Reading	<p>Add 100µL/well of Dilution buffer (B1) and proceed immediately with the reading in FIDIS™ Instrument</p> <p>The reading must be doing in the hour following adding the solution (B1). During this time, keep the plate at room temperature away from direct sunlight.</p>		

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Microplate		Negative control		Microspheres
	Specimen diluent		Positive control		IgG Conjugate
	Calibrator		Wash Buffer		Microsphere reconstitution buffer
	FIDIS test		Contains sodium azide		Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com