FIDIS™ Connective 10

REF

MX 006

96

DÉFINITION

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BIBLIOGRAPHIE

SUSAN S. COPLLE, THOMAS B. MARTINS, C. MASTERSON, E. JOLY, HARRY R. HILL

Comparison of Three Multiplex Immunoassays for Detection of Antibodies to Extractable Nuclear Antigen Using Clinically Defined Sera.

Ann. N.Y Acad. Sci. 1109: 464-472 (2007).

BULIARD A, FORTENFANT F, GHILLANI-DALBIN P, MUSSET L, OKSMAN F, OLSSON N.O.

Related Articles, Links [Analysis of nine autoantibodies associated with systemic autoimmune diseases using the Luminex technology. Results of a multicenter study].

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Evaluation of the new multiplexed immunoassay, FIDIS, for simultaneous quantitative determination of antinuclear antibodies and comparison with conventional methods. Hematology Laboratory, Tenon Hospital, 75020 Paris, France

Am J Clin Pathol., 2003; 120 (5): 676-81.

GAL I et al.

Comparison of anti-Ro/SSA antibody profile between patients with primary and secondary syndrome.

Autoimmunity 2000; Sep 32, 2, 89-92.

MEWIS A. et al.

Evaluation of a DOT-Blot method for identification of antibodies against extractable nuclear antigens and anti-cytoplasmic antibodies.

Clin Chem., 1999, 45 (2), 311-312.

CONRAD K. et al.

Autoantibodies – Diagnostic, pathogenic and prognostic relevance.

Clin. and Exp. Rheumatol. 1997 ;15, 457-465

ERMENS A. et al.

Simple Dot-Blot method evaluated for detection of autoantibodies against extractable antigens.

Clin. Chemistry 1997; 43, 12, 2420-2422.

SPENCER GREEN G. et al.

Test performance in systemic sclerosis: anti-centromere and anti-Sci70 antibodies.

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KEECH CL et al.

SS-B (La) autoantibodies. Autoantibodies. Peter JB and Shoenfeld Y, editors.

Ed Elsevier 1996, 789-797.

MEENK RJ et al.

dsDNA autoantibodies.

In: Peter JB and Shoenfeld Y. Autoantibodies. Elsevier Science, 1996, Amsterdam.227-234.

HOMBURGER H.A.

Cascade testing for autoantibodies in connective tissue disease. Mayo Clin. Proc., 1995, 70, 183-184.

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Ann. Biol. Clin., 1995, 53, 487-490.

THE LS et al.

Antiribosomal P protein antibodies in systemic lupus erythematosus.

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KEEBLE AH, KHAN Z, FORSTER A, JAMES LC.

TRIM21 is an IgG receptor that is structurally, thermodynamically, and kinetically conserved.

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





REF **MX 006**

INTENDED USE

The **FIDIS™ Connective 10** kit is a semi-quantitative homogeneous fluorescent-based microparticles immunoassay using flow cytometry readings. It is designed for the simultaneous detection of 10 autoantibodies specificities: **double stranded DNA (dsDNA), SS-A (60 kDa), TRIM21** (*previously named **SS-A 52 kDa**), **SS-B, Sm, Sm/RNP, Scl-70, Jo-1, ribosomes and centromere**.

SUMMARY AND EXPLANATION

Connective diseases are systemic auto-immune diseases which are, most often, selectively directed against a specific organ. Their classification is based upon clinical and biological data. The detection of anti-nuclear antibodies, which often appears in the patient serum before any definite symptomatic signs, will allow the identification of the specific connective disease. The titer of some autoantibodies, in particular for anti-dsDNA, provides vital information on progression of the disease and for monitoring treatment.

Many elaborate serological studies have demonstrated significant associations between these autoantibodies and the various connective diseases:

Antibodies associated with different autoimmune pathologies

Anti-SS-A (SS-A 60kDa, TRIM21) and anti-SS-B antibodies
Either one or both of these are observed, in the same pathological conditions: Sjögren's syndrome and Systemic Lupus Erythematosus (SLE). Anti-SS-A are also found in mothers who are carrying a baby with neonatal lupus syndrome including heart block.

Anti-RNP antibodies

Found primarily in Mixed Connective Tissue Disease or Sharp's Syndrome, for which they are the biological marker. They are also observed in SLE.

Antibodies specific to a single connective tissue disease

Anti-dsDNA antibodies are the most frequently detected and the most specific serological marker for SLE.

Anti-Sm antibodies characterize severe forms of SLE.

Anti-ribosomal antibodies are strongly associated with central neurological abnormalities in Lupus.

Anti-Scl-70 antibodies are specific for diffuse Scleroderma.

Anti-Jo-1 antibodies are found exclusively in sera from patients suffering from Polymyositis.

Anti-Centromere antibodies are most often found in sera from patients having Scleroderma with limited skin involvement CREST, called CREST syndrome.

ASSAY PRINCIPLE

FIDIS™ Connective 10 kits are based on the use of distinct uniform size color-coded microspheres and a benchtop flow cytometer interfaced to digital signal processing hardware and software. A red diode laser beam in the flow cytometer classifies each set of microspheres on the basis of its unique fluorescence intensity (red to orange) thus identifying which analyte is being tested. At the same time, a green laser beam illuminates the external second molecule fluorescence to quantify the reaction related to the specific analyte.

Each antigen, recombinant dsDNA, SS-B, TRIM21, CENP-B, Jo-1 and native purified Scl-70, SmRNP, Sm, ribosomal and SS-A 60 kDa, required for the assay is covalently coupled to an individual set of microspheres through its surface functional groups. The different antigen coupled microspheres are mixed together, constituting the final microsphere reagent.

FIDIS™ Connective 10 kits allow the detection of 10 autoantibodies specificities: **anti-double stranded DNA (dsDNA), anti-SS-A (60 kDa), anti-TRIM21** (*previously named **SS-A 52 kDa**), **anti-SS-B, anti-Sm, anti-Sm/RNP, anti-Scl-70, anti-Jo-1, anti-ribosomes and anti-centromere**.

The test is performed in a 96 well microplate with a filtering membrane at the bottom of the wells.

- In the first step, the sample is distributed in each well containing the microspheres mixture. If this sample contains one or more of the suspected antibodies, this(ese) antibody(ies) bind to the corresponding antigen(s) on the various sets of microspheres.
- After incubation, a wash step using a filtration process removes the unbound antibodies.
- A phycoerythrin labeled **anti-human IgG conjugate** is then added that binds to the previously bound antibodies.
- A final wash step allows to stop the reaction.
- The reaction is then directly measured by the flow cytometer, which differentiates each set of microspheres according to its fluorescence color while simultaneously measuring the average fluorescence emitted by the conjugate.
- A calibration system allows the determination of the titer (AU/mL) of each sample by interpolation for each following antigenic specificity : SS-A 60 kDa, TRIM21, SS-B, Sm, Sm/RNP, Scl-70, Jo-1, ribosomes and centromère and the titer (IU/mL) of each sample by interpolation for the following antigenic specificity: dsDNA. The **FIDIS™ Connective 10** kits can be used with a dispensing/diluting device, the **CARIS™ system**.

SPECIMEN COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

REAGENTS SUPPLIED

| | |
|--|------------|
| 96 wells microplate with filtering membrane and lid. MP | 1 plate |
| Vial (A) of color-coded microsphere set of 10 sensitized by dsDNA, SS-A 60 kDa, TRIM21, SS-B, Sm, Sm/RNP, Scl-70, Jo-1, ribosomes and centromere antigen(s). MICROSPHERES <u>Lyophilized</u> (to be reconstituted with the buffer named D) | sq 6mL |
| Vial (B1) of sample dilution buffer (white vial) <u>Ready to use</u> DIL SPE | 2 x 115mL |
| Vial of calibrator* <u>Ready to use</u> Each titer is printed on the vial label CAL | 1 x 1,5mL |
| Vial of positive control concentrated. This control has a standard reactivity, which provides evidence of the proper reagents activity and proper assay performance. <u>To be diluted</u> CONTROL + Expected values are printed on the vial label. | 1 x 250 µL |
| Vial of negative control concentrate <u>To be diluted</u> CONTROL - | 1 x 250µL |
| Vial of anti-human IgG coupled to phycoerythrin <u>Ready to use</u> CONJ IgG | 1 x 12mL |
| Vial (C1) of washing buffer (black vial) <u>Ready to use</u> BUF WASH | 1 x 100mL |
| Vial (D) of reconstitution buffer for the microsphere set <u>Ready to use</u> BUF MICROSPHERES | 1 x 6mL |

Calibrator titers* are expressed in arbitrary units per ml (AU/ml) except for dsDNA.

Calibrator titer for dsDNA (IU/ml) is standardized against the international reference: Wo/80 WHO.

ADDITIONAL MATERIAL – NOT SUPPLIED

- Precision pipettes capable of accurately delivering 5µL to 1000µL
- Multichannel Pipettes or dispensers capable of accurately delivering 40µL to 300µL or 5µL to 2 mL
- Vortex mixer

- Laboratory timer to monitor incubation steps
- Ultrasonic bath
- Absorbent towel
- Serological pipettes
- Microplate sealing films

STABILITY AND STORAGE

- Store reagents in their original packaging at +2°C to +8°C.
- Do not freeze reagents.
- Do not use kits beyond the expiration date.
- After use, store all components immediately back at +2°C to +8°C.

CAUTION

- Reagents in solution contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
- Avoid to use reagents if signs of contamination or other visible changes occur
- The **FIDIS™ Connective 10** kit has been developed according to CE Directives 67/548/EEC and 1999/45/EC regarding classification, packaging and labelling of dangerous preparations.
- The **FIDIS™ Connective 10** kit has been optimized for use as described in this procedure. Do not substitute other manufacturer reagents. Dilution or alteration of these reagents may also alter the performance of the test. Follow thoroughly the test procedure to insure optimal performance.
- *Calibrators and controls are from human origin. The human sera used in the preparation of these products were tested and found non-reactive for antibodies to HIV-1, HIV-2, anti-HCV and Hepatitis B antigen. Because no test method can offer complete assurance that infectious agents are absent, handle as if capable of transmitting infectious diseases.*
- If all specificities are found in the same patient sample, it is advisable to control this sample by another method.

TEST SET-UP

FIDIS™ Washer

- ☞ Check the tubing connecting the pump to the stand and the manometer setting (wheel completely closed).

Using FIDIS™ Analyzer and MLX-BOOSTER Software

- ☞ See the User's Manual provided with the **FIDIS™ Instrument** for detailed instructions on running the equipment and for calculation. For additional information and/or trouble shooting problems, please contact bmd subsidiaries or distributors.

1. Start the run as described in the User's Manual.
NOTE: The FIDIS™ takes 30 minutes to warm up after being turned on. A new warm up is necessary after 4 hours of system inactivity.
2. Calibration and controls are described in the User's Manual. Calibration and controls should be routinely performed once per month and each time a new lot of Sheath fluid reagent is used in order to insure optimal instrument performance. Calibration should also be performed when the temperature is out of the determined range shown on the run batch screen.
3. Programming a batch or a multibatch of test protocol on the FIDIS™.
4. Load the microplate into the plate holder of the FIDIS™ as described in the User's Manual.
5. Analyze the results according to the User's Manual.
6. When finished for the day, perform the sanitizing and soak operations, prior to turning the Analyzer off according to the shutting down procedure described in the User's Manual.

Reagents Preparation

☞ Remove the individual components from storage, allow them to warm up to **room temperature** (+18°C to +25°C), and mix them well.

1. Microsphere preparation

- Reconstitute the microsphere solution by adding the vial **D** of reconstitution buffer to the vial **A** of lyophilized microspheres. Wait 5 minutes and vortex.
- Stable 2 months at +2°C to +8°C, after reconstitution.

2. Preparation of samples and controls.

- Dilute the samples and controls at **1:201** in dilution buffer (**B1**).
Ex.: 10µL sample in 2000µL dilution buffer (B1).
- **Vortex vigorously.**

3. Assay configuration

Use the work-sheet included in the kit to identify the location of the samples.

a. When setting up the test, systematically take into account the following well requirements: ⇒ See examples below

- ✗ 1 "reagent-blank" well
- ✗ 1 well for negative control
- ✗ 1 well for positive control
- ✗ 2 "calibrator" wells

b. Calculation of the correct number of wells necessary and their location.

In the following examples, 2 different configurations are described, according to the number of wells needed, and the availability of unused wells. The sample dispensing must be systematically carried out in a column. Do not leave empty wells.

The same microplate can be used for more than one serie of tests if unused wells are protected by

microplate sealing film. Any single well can only be used once.

To avoid any error in using a well more than once, identify the used wells with a marker (see example 2 below).

c. Programming of the test protocol on the FIDIS™

☞ Refer to the FIDIS™ user's manual: "Creating a batch or a multi-batch".

Example 1:

- The microplate is totally blank (no well has yet been used): Use A1 as the 1st well.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------|-----|-----|-----|---|---|---|---|---|----|----|----|
| A | Blank | S4 | S12 | S20 | | | | | | | | |
| B | Ctrl- | S5 | S13 | S21 | | | | | | | | |
| C | Ctrl+ | S6 | S14 | S22 | | | | | | | | |
| D | Cal | S7 | S15 | S23 | | | | | | | | |
| E | Cal | S8 | S16 | S24 | | | | | | | | |
| F | S1 | S9 | S17 | S25 | | | | | | | | |
| G | S2 | S10 | S18 | | | | | | | | | |
| H | S3 | S11 | S19 | | | | | | | | | |

Example 2:

- 16 wells have already been used: Use A3 as the 1st well.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|-------|-----|-----|-----|-----|-----|-----|-----|----|----|
| A | | | Blank | S4 | S12 | S20 | S28 | S36 | S44 | S52 | | |
| B | | | Ctrl- | S5 | S13 | S21 | S29 | S37 | S45 | S53 | | |
| C | | | Ctrl+ | S6 | S14 | S22 | S30 | S38 | S46 | | | |
| D | | | Cal | S7 | S15 | S23 | S31 | S39 | S47 | | | |
| E | | | Cal | S8 | S16 | S24 | S32 | S40 | S48 | | | |
| F | | | S1 | S9 | S17 | S25 | S33 | S41 | S49 | | | |
| G | | | S2 | S10 | S18 | S26 | S34 | S42 | S50 | | | |
| H | | | S3 | S11 | S19 | S27 | S35 | S43 | S51 | | | |

ASSAY PROCEDURE

1. Microspheres dispensing

- Before to start the assay and if it is necessary, protect the unused wells with a microplate sealing film.
- **Vortex the microspheres reagent vigorously for 20s** and dispense 50µL in each well to be used.

PLEASE NOTE:

- ☛ If CARIS™ is used, the microspheres should be sonicated for 5 minutes.
- ☛ If a dispensing/diluting device is used (like CARIS™), the procedure steps 1 and 2 should be reversed:
⇒ dispense diluted sera in first and then, dispense microspheres

2. Addition and Incubation of the samples

- Dispense 100µL of sample dilution buffer (**B1**) for reagent blank.
- Dispense 100µL of prediluted controls, of ready to use calibrator and of prediluted samples.
- Cover the plate and incubate 30 minutes at room temperature, away from direct sunlight and without shaking.

3. Wash step 1

Wash the plate 2 times with washing buffer (C1) using filtration unit:

- Remove the microplate lid (verify if the vacuum break valve is closed) and place it on the washer.
- Start the pump. Stop it when the whole incubation liquid has been removed in totality and open the vacuum break valve. Close the vacuum break valve.
- Dispense 300µL of washing buffer (C1) in all the used wells.
- Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- Dispense 300µL of washing buffer (C1) in all the used wells.
- Start the pump. Stop it after 5 additional seconds when all the buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by tapping vigorously the plate 10 times on an absorbent towel.
- Place again the microplate on the washer and start the pump for 5 seconds. Stop it and quickly open the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by blotting vigorously the plate on an absorbent towel.
- Place the plate on a totally dry surface before starting the conjugate incubation.

4. Incubation of the conjugate

- Dispense 100µL of ready to use conjugate in each well used.
- Cover the plate and incubate 30 minutes at room temperature away from direct sunlight and without shaking. The incubation time starts after the conjugate has been added to all wells. If this timing is not followed, the results might be erroneous.

Note : the conjugate is photosensitive.
 ⇒ **After using, close the vial carefully.**

5. Wash step 2

Wash the plate using the filtration unit on 1 cycle with dilution buffer (B1).

Use two different reagent reservoirs for the 2 washing steps. Either clean and dry correctly the reservoir between both wash steps (avoid to mix both B1 and C1 buffers).

- Remove the microplate lid (verify that the vacuum break valve is closed) and place it on the washer.
- Start the pump. Stop it when the whole incubation liquid has been removed and open the vacuum break valve.
Close the vacuum break valve.

- Dispense 100µL of dilution buffer (B1) in all the used wells.
- Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by blotting the plate on an absorbent towel.
- Dispense 100µL of dilution buffer (B1) in all the used wells. Place the plate on a totally dry surface at room temperature.
The reading must be doing in the hour following adding the solution (B1). During this time, keep the plate at room temperature away from direct sunlight.

6. Reading of the test

☞ Follow the **MLX-BOOSTER™** and **FIDIS™** user's manual: Processing batches.

VALIDATION CRITERIA OF THE RESULTS

Calibrator, negative and positive controls have to be run with every batch of samples to ensure that all reagents and procedures performed properly.

In order to validate the results, all the criteria listed below must be met. Otherwise, the test is invalid and must be repeated.

- The positive controls should show a value within the limits printed on the corresponding vial labels.
- The negative controls should be less than 30 AU/mL or 30 IU/mL.

The microsphere set contains Internal standard beads allowing to verify the presence of serum in wells, as well as the good respect of the test procedure. Wells presenting a not corresponding fluorescent signal will be no valid with MLX BOOSTER™ software.

INTERPRETATION OF RESULTS

| | | | |
|--|------------|---------------|------------|
| Arbitrary Units (AU/ml) except for the dsDNA | < 30 AU/ml | 30 – 40 AU/ml | > 40 AU/ml |
| Interpretation | Negative | Borderline(*) | Positive |

| | | | |
|---|------------|---------------|------------|
| International Units (IU/ml) for the dsDNA | < 30 IU/ml | 30 – 40 IU/ml | > 40 IU/ml |
| Interpretation | Negative | Borderline(*) | Positive |

(*) *Borderline results should be controlled on a second sample and the interpretation of the results should be done in the frame of additional testing and taking into account the clinical status of the patient.*

Each laboratory should establish and maintain its own references (normal) range values, based on the patient population and other local factors.

CALCUL OF RESULTS

The results are automatically calculated by MLX-BOOSTER™ software and can be printed for each assay.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

The serums from Immunodeficient patients will give a no valid result.

QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities. Immunotrol I, IV and V multiparametric controls (bmd, Cat. n°: HM036, HM051 and HM052) contain human auto-antibodies directed against researched specificities. To be assayed in the same manner as the unknown samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH OTHER METHODS.

FIDIS™ Connective 10 kit has been compared to:

- ⇒ ELISA method for the detection of antibodies anti-SS-A 60 kDa, SS-B, Sm, Sm/RNP, Scl-70, Jo-1 and dsDNA and anti-TRIM21.
- ⇒ Dot-Blot method for the detection of antibodies anti-centromere and anti-ribosome.

Confirmation methods for discrepant results:

- ⇒ Dot-Blot method for antibodies anti-SS-A 60 kDa, SS-B, Sm, Sm/RNP, Scl-70, Jo-1, dsDNA and TRIM21.
- ⇒ Immunofluorescence on HEp2000™ for antibodies anti-centromere.
- ⇒ Immunofluorescence on triple substrate for antibodies anti-ribosome.

The study was performed on 209 samples:

- ⇒ 112 samples related to systemic autoimmune diseases
- ⇒ 48 samples from blood donors
- ⇒ 49 samples selected for their potential biological interferences (cryoglobulinemia, hypergammaglobulinemia, IgG and IgM monoclonal immunoglobulins, complement, rheumatoid factor, plasma, haemolysed and turbid samples).

Results

| Antigenic Specificity | Nber of positive results FIDIS™ / ELISA or DOT* | Relative sensitivity (%) | Relative specificity (%) | Accuracy (%) |
|-----------------------|---|--------------------------|--------------------------|--------------|
| dsDNA | 40/42 | 95.5 | 100 | 99 |
| SS-A 60kDa | 58/58 | 100 | 100 | 100 |
| TRIM21 | 50/51 | 98.1 | 100 | 99.5 |
| SS-B | 21/20 | 95.2 | 99 | 98.6 |
| Sm | 10/11 | 91.7 | 100 | 99.5 |
| Sm/RNP | 22/21 | 100 | 99.5 | 99.5 |
| Scl-70 | 9/9 | 100 | 100 | 100 |
| Jo-1 | 12/12 | 100 | 100 | 100 |
| Centromere | 9/9* | 100 | 100 | 100 |
| Ribosome | 5/5* | 100 | 100 | 100 |

Discrepant results analysis (8/209)

| Antigenic specificity | FIDIS™ Results | ELISA Results | Dot-Blot result | Associated specificities |
|-----------------------|----------------|-----------------|-----------------|--------------------------|
| dsDNA | negative | positive | negative | SSA 60 kDa |
| dsDNA | negative | Weakly positive | negative | SSA 60 kDa TRIM21 |
| TRIM21 | negative | Weakly positive | positive | none |
| SS-B | negative | Weakly positive | positive | SSA 52 kDa |

| Antigenic specificity | FIDIS™ Results | ELISA Results | Dot-Blot result | Associated specificities |
|-----------------------|-----------------|-----------------|-----------------|---------------------------|
| SS-B | Weakly positive | negative | positive | SS-A 60kDa |
| SS-B | positive | borderline | negative | SS-A 60kDa TRIM21 dsDNA |
| Sm | negative | Weakly positive | positive | Sm/RNP, SS-A 60kDa, dsDNA |
| Sm/RNP | positive | negative | positive | SS-A 60kDa, dsDNA |

All these 8 results are in limit of positivity, which can lead to no significant discrepancies.

⇒ 4/8 FIDIS results were confirmed by Dot-Blot.

Method for determining threshold values

Estimated on the 97 samples from « blood donors » and « biological interferences » populations and according standard Wo/80 for dsDNA specificity.

The positive thresholds (40AU/ml or 40IU/ml) correspond to the 100th percentile of the values distribution for all the 10 specificities.

The negative thresholds (30AU/ml or 30IU/ml) correspond to the 96th percentile for all the 10 specificities.

Between these two thresholds, samples are considered borderline.

REPRODUCIBILITY / PRECISION

| Antigenic Specificity | Within-run (10 tests in the same run) | | Between-run (5 tests in 5 different runs) | |
|-----------------------|--|--------|--|--------|
| | Mean value | CV (%) | Mean value | CV (%) |
| dsDNA | 150 | 3.8 | 61 | 5.0 |
| SS-A 60kDa | 136 | 6.2 | 36 | 3.6 |
| TRIM21 | 171 | 6.1 | 55 | 4.2 |
| SS-B | 89 | 5.2 | 160 | 3.1 |
| Sm | 78 | 6.0 | 67 | 6.7 |
| Sm/RNP | 124 | 4.3 | 83 | 6.3 |
| Scl-70 | 185 | 3.5 | 123 | 3.9 |
| Jo-1 | 135 | 6.6 | 80 | 4.4 |
| Centromere | 283 | 4.5 | 167 | 2.8 |
| Ribosome | 68 | 2.2 | 120 | 3.0 |

BIBLIOGRAPHY

SUSAN S. COPLLE, THOMAS B. MARTINS, C. MASTERSON, E. JOLY, HARRY R. HILL

Comparison of Three Multiplex Immunoassays for Detection of Antibodies to Extractable Nuclear Antigen Using Clinically Defined Sera.
Ann. N.Y Acad. Sci. 1109: 464-472 (2007).

BULIARD A, FORTENFANT F, GHILLANI-DALBIN P, MUSSET L, OKSMAN F, OLSSON N.O.

Related Articles, Links [Analysis of nine autoantibodies associated with systemic autoimmune diseases using the Luminex technology. Results of a multicenter study].
Ann Biol Clin (Paris), 2005; Jan-Feb, 63 (1): 51-8.

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Evaluation of the new multiplexed immunoassay, FIDIS, for simultaneous quantitative determination of antinuclear antibodies and comparison with conventional methods. Hematology Laboratory, Tenon Hospital, 75020 Paris, France
Am J Clin Pathol., 2003; 120 (5): 676-81.

GAL I et al.

Comparison of anti-Ro/SSA antibody profile between patients with primary and secondary syndrome. Autoimmunity 2000; Sep 32, 2, 89-92.

MEWIS A. et al.

Evaluation of a DOT-Blot method for identification of antibodies against extractable nuclear antigens and anti-cytoplasmic antibodies.
Clin Chem., 1999, 45 (2), 311-312.

CONRAD K. et al.

Autoantibodies – Diagnostic, pathogenic and prognostic relevance.
Clin. and Exp. Rheumatol. 1997 ;15, 457-465

ERMENS A. et al.

Simple Dot-Blot method evaluated for detection of autoantibodies against extractable antigens.
Clin. Chemistry 1997; 43, 12, 2420-2422.

SPENCER GREEN G. et al.

Test performance in systemic sclerosis: anti-centromere and anti-Scl70 antibodies.
Am. J. of Med., 1997, 18, 291-300.

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KEEBLE AH, KHAN Z, FORSTER A, JAMES LC.

TRIM21 is an IgG receptor that is structurally, thermodynamically, and kinetically conserved.
Proc Natl Acad Sci U S A. 2008 Apr 22;105(16):6045-50. Epub 2008 Apr 17.

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





Latest update on the RO/SS-A autoantibody system.
Autoimmun Rev. 2009 Feb 11.

SUMMARY OF THE TEST PROCEDURE

| | VOLUME TO BE DISTRIBUTED | REAGENTS | INCUBATION CONDITIONS |
|---|--|---|-----------------------------|
| Dispense 50µL of microspheres reagent in each well to be used | | | |
| Incubation of the samples | 100µL | Dilution buffer (B1) : Reagent blank | 30 min. Room temperature |
| | 100µL | Diluted negative control | |
| | 100µL | Diluted positive control | |
| | 100µL | Calibrator (in duplicate) | |
| | 100µL | Diluted samples | |
| Washing 1 | Wash twice in washing buffer (C1) (300µL/well) | | |
| Incubation of the conjugate | 100µL | Ready to use Anti-IgG conjugate | 30 min Room temperature |
| Washing 2 | Wash 1 time in Dilution buffer (B1) (100µL/well) | | |
| Reading | <p>Add 100µL/well of Dilution buffer (B1) and proceed immediately with the reading in FIDIS™ Instrument</p> <p>The reading must be doing in the hour following adding the solution (B1). During this time, keep the plate at room temperature away from direct sunlight.</p> | | |

SYMBOLS USED

| | | | | | |
|---|---------------------------|---|-------------------------|---|-----------------------------------|
|  | Biological risk |  | Temperature limitation |  | Catalog Number |
|  | Read instructions for use |  | In Vitro Diagnostic Use |  | Lot Number |
|  | Number of tests |  | Use by |  | EC Declaration of Conformity |
|  | Microplate |  | Negative control |  | Microspheres |
|  | Specimen diluent |  | Positive control |  | IgG Conjugate |
|  | Calibrator |  | Wash Buffer |  | Microsphere reconstitution buffer |
|  | FIDIS test |  | Contains sodium azide |  | Reconstitute with |

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com