FIDIS™ Vasculitis

REF

MX 007

96

DÉFINITION

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BIBLIOGRAPHIE

HAGEN EC et al.

Development and standardization of solid-phase assays for the detection of anti-neutrophil cytoplasmic antibodies (ANCA) for clinical application: report of a large clinical evaluation study. Proceedings of the 6th International ANCA Workshop. Paris, 28 June - 1 July 1995.

Clin Exp Immunol, 101/1: 29, 1995.

JENNETTE JC et al.

Pathogenic potential of anti-neutrophil cytoplasmic auto-antibodies.

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RASMUSSEN N. et al.

An ELISA for the detection of antineutrophil cytoplasm antibodies (ANCA).

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Spectra Biologie, 93/3: 38-43, 1993.

SAXENA R et al.







Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome.

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96

INTENDED USE

The **FIDIS™ Vasculitis** kit is a quantitative homogeneous fluorescent-based microparticles immunoassay using flow cytometry readings. It is designed for the simultaneous detection of 3 autoantibodies specificities:

- **anti-neutrophil cytoplasm antibodies (ANCA)** directed against Myeloperoxidase (*MPO*) and Serine Proteinase 3 (*PR3*).
- **anti-glomerular basement membrane (GBM)** antibodies using a purified antigen based on type IV collagen.

SUMMARY AND EXPLANATION

• ANCA

ANCA antibodies form a family of autoantibodies of which the diagnostic and predictive value for the different necrotising vasculitis and 'overlap' systemic vasculitis, associated with glomerulonephritis or not, is now well established.

They are directed against antigens present in the primary azurophilic granules (alpha granules) of polynuclear neutrophils and in the lysosomes of monocytes.

ANCA antibodies are in the first place detected by means of indirect immunofluorescence using human ethanol-fixed neutrophil slides. The presence of ANCA is characterised in two distinct staining patterns: one cytoplasmic and the other perinuclear. These two types called c-ANCA and p-ANCA, are predominantly related respectively with antibodies to proteinase 3 (PR3) and myeloperoxidase (MPO).

The detection of ANCA with IFA only, may give problems in the differential diagnosis when one realizes that it is possible to have very similar patterns caused by other autoantibodies. This holds true as well in Wegener's granulomatosis and vasculitis, where the association of ANCA with other autoantibodies is described (rheumatoid factor, anti-nuclear antibodies, anti-smooth muscle antibodies, anti-thyroid, etc.) as in collagenosis, autoimmune hepatitis where, besides specific autoantibodies, also ANCA could develop.

The diagnostic sensitivity of ANCA is of particular importance with Wegener's granulomatosis (anti-PR3), Churg Strauss syndrome, microscopic periarteritis and idiopathic crescentic glomerulonephritis (anti-MPO).

These diseases represent a medical emergency and require a safe strategy for the identification of ANCA based on one characterisation in two ways:

- 1- By indirect immunofluorescence where two known patterns give a first orientation to the diagnosis.

- 2- A reliable identification of antibodies against the target antigens (MPO and PR3) with specific immunologic techniques using purified antigens.

• GBM

The measurement of anti-glomerular basement membrane (GBM) autoantibodies is intended as an aid in the diagnosis of Goodpasture's syndrome, based on clinical signs of rapidly progressive renal failure and lung haemorrhage. It concerns an autoimmune disease due to autoantibodies reacting with the components of the glomerular basement membrane (alpha chain 3 of type IV collagen).

This disease progresses rapidly and in a fatal manner in more than 75% of the cases without treatment. A quick diagnosis, based on the examination of anti-GBM autoantibodies as specific markers, is indispensable before the kidney and often the lung undergo important damage. Besides, the continuous presence of these autoantibodies implies a delay of renal transplantation.

The detection of anti-GBM antibodies allows to distinguish Goodpasture's syndrome from other nephropathies of immunologic origin (glomerulonephritis associated with anti-neutrophil cytoplasmic antibodies) who need a total different treatment. A negative result is a very important element for the clinician to exclude the diagnosis of Goodpasture's syndrome.

ASSAY PRINCIPLE

FIDIS™ Vasculitis kit is based on the use of distinct uniform size color-coded microspheres and a benchtop flow cytometer interfaced to digital signal processing hardware and software. A red diode laser beam in the flow cytometer classifies each set of microspheres on the basis of its unique fluorescence intensity (red to orange) thus identifying which analyte is being tested. At the same time, a green laser beam illuminates the external second molecule fluorescence to quantify the reaction related to the specific analyte.

Each antigen required for the assay is covalently coupled to an individual set of microspheres through its surface functional groups. The different antigen coupled microspheres are mixed together constituting the final microspheres reagent.

The **FIDIS™ Vasculitis** kit allows to detect three specific antibodies: anti-MPO, anti-PR3 and anti-GBM antibodies.

The test is performed in a 96 well microplate with a filtering membrane.

- In the first step, the sample is distributed in each well containing the microspheres mixture. If this sample contains one or more of the suspected antibodies, this(ese) antibody(ies) binds to the corresponding antigen(s) on the various sets of microspheres.
- After incubation, a wash step using a filtration process removes the unbound antibodies.

- A phycoerythrin labelled **anti-human IgG** conjugate is then added that binds to the previously bound antibodies.
- A final wash step allows to stop the reaction.
- The reaction is then directly measured by the flow cytometer, which differentiates each set of microspheres according to its fluorescence colour while simultaneously measuring the average fluorescence emitted by the conjugate.
- A calibration system allows the determination of the titer (AU/mL) of each sample by interpolation for each antigenic specificity.

The **FIDIS™ Vasculitis** kit could be used with **CARIS™** system (diluting and dispensing device).

SPECIMEN COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

REAGENTS SUPPLIED

96 wells microplate with filtering membrane and lid. MP	1 plate
Vial (A) of color-coded microsphere set of 3 sensitized by MPO (purified from human blood), PR3 (purified from human leucocytes) and GBM purified from bovine tissues) MICROSPHERES <u>Lyophilized</u> (to be reconstituted with the buffer named D)	Sq 6mL
Vial (B1) of sample dilution buffer (white vial) <u>Ready to use</u> DIL SPE	2 x 115mL
Vial of calibrator titrated in arbitrary unit for the specificities to be measured. <u>Ready to use</u> Each titer is printed on the vial label CAL	1 x 1,5mL
Vial of positive control concentrated. This control has a standard reactivity, which provides evidence of the proper reagents activity and proper assay performance. CONTROL + <u>To be diluted</u> Expected values are printed on the vial label.	1 x 250µL
Vial of negative control concentrate <u>To be diluted</u> CONTROL -	1 x 250µL
Vial of anti-human IgG coupled to phycoerythrin <u>Ready to use</u> CONJ IgG	1 x 12mL
Vial (C1) of washing buffer (black vial) <u>Ready to use</u> BUF WASH	1 x 100mL
Vial (D) of reconstitution buffer for the microsphere set <u>Ready to use</u> BUF MICROSPHERES	1 x 6mL

ADDITIONAL MATERIAL – NOT SUPPLIED

- Precision pipettes capable of accurately delivering 5µL to 1000µL
- Multichannel Pipettes or dispensers capable of accurately delivering 40µL to 300µL or 5µL to 2 mL
- Vortex mixer
- Laboratory timer to monitor incubation steps
- Ultrasonic bath
- Absorbent towel
- Serological pipettes
- Microplate sealing films

CAUTION

- Reagents in solution contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
- Avoid to use reagents if signs of contamination or other visible changes occur.
- The **FIDIS™ Vasculitis** kit has been developed according to CE Directives 67/548/EEC and 1999/45/EC regarding classification, packaging and labelling of dangerous preparations.
- The **FIDIS™ Vasculitis** kit has been optimized for use as described in this procedure. Do not substitute other manufacturer reagents. Dilution or alteration of these reagents may also alter the performance of the test. Follow thoroughly the test procedure to insure optimal performance.
- *Calibrators and controls are from human origin. The human sera used in the preparation of these products were tested and found non-reactive for antibodies to HIV-1, HIV-2, anti-HCV and Hepatitis B antigen. Because no test method can offer complete assurance that infectious agents are absent, handle as if capable of transmitting infectious diseases.*

STABILITY AND STORAGE

- Store reagents in their original packaging at +2°C to +8°C.
- Do not freeze reagents.
- Do not use kits beyond the expiration date.
- After use, store all components immediately back at +2°C to +8°C.

TEST SET-UP

FIDIS™ Washer

☞ Check the tubing connecting the pump to the stand and the manometer setting (wheel completely closed).

Using FIDIS™ Analyzer and MLX-BOOSTER Software

☞ See the User's Manual provided with the **FIDIS™ Instrument** for detailed instructions on running the equipment and for calculation. For additional information and/or trouble shooting problems, please contact bmd subsidiaries or distributors.

1. Start the run as described in the User's Manual.
NOTE: The FIDIS™ takes 30 minutes to warm up after being turned on. A new warm up is necessary after 4 hours of system inactivity.

2. Calibration and controls are described in the User's Manual. Calibration and controls should be routinely performed once per month and each time a new lot of Sheath fluid reagent is used in order to insure optimal instrument performance.

Calibration should also be performed when the temperature is out of the determined range shown on the run batch screen.

3. Programming a batch or a multibatch of test protocol on the FIDIS™ as described in the User's Manual.

4. Load the microplate into the plate holder of the FIDIS™ as described in the User's Manual.

5. Analyze the results according the User's Manual.

6. When finished for the day, perform the sanitizing and soak operations, prior to turning the Analyzer off according to the shutting down procedure described in the User's Manual.

Reagents Preparation

☞ Remove the individual components from storage, allow them to warm up to **room temperature** (+18°C to +25°C), and mix them well.

1. Microsphere preparation

- Reconstitute the microsphere solution by adding the vial **D** of reconstitution buffer to the vial **A** of lyophilized microspheres. Wait 5 minutes and vortex.
- Stable 2 months at +2°C to +8°C, after reconstitution.

2. Preparation of samples and controls.

- Dilute the samples and controls at **1:201** in dilution buffer (**B1**).
Ex.: 10µL sample in 2000µL dilution buffer (**B1**).
- **Vortex vigorously.**

3. Assay configuration

Use the work-sheet included in the kit to identify the location of the samples.

a. When setting up the test, systematically take into account the following well requirements: ⇒ See examples below

- ✗ 1 "reagent-blank" well
- ✗ 1 well for negative control
- ✗ 1 well for positive control
- ✗ 2 "calibrator" wells

b. Calculation of the correct number of wells necessary and their location.

In the following examples, 2 different configurations are described, according to the number of wells needed, and the availability of unused wells. The sample dispensing must be systematically carried out in a column. Do not leave empty wells.

The same microplate can be used for more than one serie of tests if unused wells are protected by microplate sealing film. Any single well can only be used once.

To avoid any error in using a well more than once, identify the used wells with a marker (see example 2 below).

c. Programming of the test protocol on the FIDIS™

☞ Refer to the FIDIS™ user's manual: "Creating a batch or a multi-batch".

Example 1:

- The microplate is totally blank (no well has yet been used): Use A1 as the 1st well.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	S4	S12	S20								
B	Ctrl-	S5	S13	S21								
C	Ctrl+	S6	S14	S22								
D	Cal	S7	S15	S23								
E	Cal	S8	S16	S24								
F	S1	S9	S17	S25								
G	S2	S10	S18									
H	S3	S11	S19									

Example 2:

- 16 wells have already been used: Use A3 as the 1st well.

	1	2	3	4	5	6	7	8	9	10	11	12
A			Blank	S4	S12	S20	S28	S36	S44	S52		
B			Ctrl-	S5	S13	S21	S29	S37	S45	S53		
C			Ctrl+	S6	S14	S22	S30	S38	S46			
D			Cal	S7	S15	S23	S31	S39	S47			
E			Cal	S8	S16	S24	S32	S40	S48			
F			S1	S9	S17	S25	S33	S41	S49			
G			S2	S10	S18	S26	S34	S42	S50			
H			S3	S11	S19	S27	S35	S43	S51			

ASSAY PROCEDURE

1. Microspheres dispensing

- Before to start the assay and if it is necessary, protect the unused wells with a microplate sealing film.
- **Vortex the microspheres reagent vigorously for 20s** and dispense 50µL in each well to be used.

PLEASE NOTE:

- Ⓢ If CARIS™ is used, the microspheres should be sonicated for 5 minutes.
- Ⓢ If a dispensing/diluting device is used (like CARIS™), the procedure steps 1 and 2 should be reversed: ⇒ dispense diluted sera in first and then, dispense microspheres

2. Addition and Incubation of the samples

- Dispense 100µL of sample dilution buffer (**B1**) for reagent blank.
- Dispense 100µL of prediluted controls, of ready to use calibrator and of prediluted samples.
- Cover the plate and incubate 30 minutes at room temperature, away from direct sunlight and without shaking.

3. Wash step 1

Wash the plate 2 times with washing buffer (C1) using filtration unit:

- a. Remove the microplate lid (verify if the vacuum break valve is closed) and place it on the washer.
- b. Start the pump. Stop it when the whole incubation liquid has been removed in totality and open the vacuum break valve. Close the vacuum break valve.
- c. Dispense 300µL of washing buffer (C1) in all the used wells.
- d. Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- e. Dispense 300µL of washing buffer (C1) in all the used wells.
- f. Start the pump. Stop it after 5 additional seconds when all the buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- g. Remove the microplate from the washer and remove the residual buffer by tapping vigorously the plate 10 times on an absorbent towel.
- h. Place again the microplate on the washer and start the pump for 5 seconds. Stop it and quickly open the vacuum break valve.
- i. Remove the microplate from the washer and remove the residual buffer by blotting vigorously the plate on an absorbent towel.
- j. Place the plate on a totally dry surface before starting the conjugate incubation.

4. Incubation of the conjugate

- Dispense 100µL of ready to use conjugate in each well used.
- Cover the plate and incubate 30 minutes at room temperature away from direct sunlight and without shaking. The incubation time starts after the conjugate has been added to all wells. If this timing is not followed, the results might be erroneous.

Note : the conjugate is photosensitive.
⇒ After using, close the vial carefully.

5. Wash step 2

Wash the plate using the filtration unit on 1 cycle with dilution buffer (B1).

*Use two different reagent reservoirs for the 2 washing steps.
Either clean and dry correctly the reservoir between both wash steps (avoid to mix both B1 and C1 buffers).*

- a. Remove the microplate lid (verify that the vacuum break valve is closed) and place it on the washer.
- b. Start the pump. Stop it when the whole incubation liquid has been removed and open the vacuum break valve. Close the vacuum break valve.

- c. Dispense 100µL of dilution buffer (B1) in all the used wells.
- d. Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- e. Remove the microplate from the washer and remove the residual buffer by blotting the plate on an absorbent towel.
- f. Dispense 100µL of dilution buffer (B1) in all the used wells. Place the plate on a totally dry surface at room temperature.
The reading must be doing in the hour following adding the solution (B1). During this time, keep the plate at room temperature away from direct sunlight.

6. Reading of the test

☞ Follow the **MLX-BOOSTER™** and **FIDIS™** user's manual: Processing batches.

VALIDATION CRITERIA OF THE RESULTS

Calibrator, negative and positive controls have to be run with every batch of samples to ensure that all reagents and procedures performed properly.

In order to validate the results, all the criteria listed below must be met. Otherwise, the test is invalid and must be repeated.

- a. The positive controls should show a value within the limits printed on the corresponding vial labels.
- b. The negative controls should be less than 20 AU/mL.

The microsphere set contains Internal standard beads allowing to verify the presence of serum in wells, as well as the good respect of the test procedure. Wells presenting a not corresponding fluorescent signal will be no valid with MLX BOOSTER™ software.

INTERPRETATION OF RESULTS

Arbitrary units (AU/mL)	< 20 AU/mL	20 – 25 AU/mL	> 25 AU/mL
Interpretation	Negative	Borderline ^(*)	Positive

() Borderline results should be repeated on a second sample and the interpretation of the results should be treated in the frame of additional testing and taking into account the clinical status of the patient.*

Each laboratory should establish and maintain its own references (normal) range values, based on the patient population and other local factors.

CALCUL OF RESULTS

☞ The results are automatically calculated by MLX-BOOSTER™ software and can be printed for each assay.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

The serums from Immunodeficient patients will give a no valid result.




QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities. Immunotrol I multiparametric control (bmd, Cat. n° HM036) contains human auto-antibodies directed against researched specificities. To be assayed in the same manner as the unknown samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

CORRELATION

The performances of **FIDIS™ Vasculitis**, lyophilized form were assessed on a population of 283 samples and compared to the performances of **FIDIS™ Vasculitis**, liquid form. A linear regression analysis of the two products showed that both are equivalent.

REPRODUCIBILITY / PRECISION (LYOPHILIZED FORM)

Antigenic specificity	Within-run (10 tests in the same run)		Between-run (5 tests in 5 different runs)	
	Mean value (AU/mL)	CV (%)	Mean value (AU/mL)	CV (%)
MPO	43	3.2	48	10.2
MPO	83	2.9	150	7.9
MPO	146	4.2	288	6.4
PR3	20	5.1	24	8.9
PR3	51	1.5	100	13.4
PR3	161	3.2	220	11.7
GBM	50	11.2	90	8.0
GBM	99	2.4	127	5.6
GBM	384	2.5	371	5.3

INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH OTHER METHODS

A comparative study was conducted to evaluate the performance characteristics of **FIDIS™ Vasculitis** test, relative to three other commercially available ELISA kits.

The study was performed on 182 samples from patients:

⇒ 117 positive samples for one parameter ANCA or GBM.

⇒ 65 negative samples and samples with potential biological interferences.

The evaluation was made on samples selected from a routine laboratory which were asked for markers tests of ANCA and GBM disorders by the physician. No information about age, gender and clinical status is available.

FIDIS™ Vasculitis and predicate devices use serum. All borderline results with **FIDIS™ Vasculitis** test or comparative ELISA tests were considered as negative.

MPO	ELISA	
	Positive	Negative
FIDIS™ Vasculitis Positive	31	3
FIDIS™ Vasculitis Negative	2	146

PR3	ELISA	
	Positive	Negative
FIDIS™ Vasculitis Positive	40	3
FIDIS™ Vasculitis Negative	7	132

GBM	ELISA	
	Positive	Negative
FIDIS™ Vasculitis Positive	21	0
FIDIS™ Vasculitis Negative	0	161

Performances for the selected population

Antigen	Positive percent agreement	Negative percent agreement	Overall agreement
MPO	93.9 % (31/33)	98.0 % (146/149)	97.3 % (177/182)
PR3	85.1 % (40/47)	97.8 % (132/135)	94.5 % (172/182)
GBM	100 % (21/21)	100 % (161/161)	100 % (182/182)

Discrepant results analysis

Antigen	FIDIS™+ / ELISA 1 – results	FIDIS™- / ELISA 1 + Results	Confirmed results with ELISA 2
MPO	3	/	2/3
	/	2	0/2
PR3	3	/	1/3
	/	7	5/7

⇒ a lot of these results are in limit of positivity, which can lead to no significant discrepancies.

Expected values

They were estimated from the 2 populations previously used for the comparative study:

⇒ 37 samples from blood donors

⇒ 28 samples selected for their potential biological interferences

The negative thresholds (20AU/mL) correspond to the 100th of percentile for MPO/GBM and 98.5th of percentile for PR3 for the populations studied.






BIBLIOGRAPHY

- HAGEN EC et al.
Development and standardization of solid-phase assays for the detection of anti-neutrophil cytoplasmic antibodies (ANCA) for clinical application: report of a large clinical evaluation study. Proceedings of the 6th International ANCA Workshop. Paris, 28 June - 1 July 1995.
Clin Exp Immunol, 101/1: 29, 1995.
- JENNETTE JC et al.
Pathogenic potential of anti-neutrophil cytoplasmic auto-antibodies.
lab Invest 70: 135-137, 1994.
- LESAVRE P. et al.
Les autoanticorps anti-cytoplasme des neutrophiles (ANCA). Signification clinique et rôle pathogène.
Polynucléaire Neutrophile et Biologie Clinique. Hôpital Rothschild PARIS. 7avril 1995.
- NOEL LH et al.
Les anticorps anticytoplasme des neutrophiles (ANCA): associations cliniques et aspects immunologiques.
Grüfeld JP, éd. Actualités Néphrologiques de l'hôpital Necker. Paris: Flammarion Médecine-Sciences, 223-253, 1992.
- RASMUSSEN N. et al.
An ELISA for the detection of antineutrophil cytoplasm antibodies (ANCA).
J Immunol Meth 127: 139-145, 1990.
- SANCHEZ-LALLOYER N.
Les anticorps anti-cytoplasme des polynucléaires neutrophiles.
Spectra Biologie, 93/3: 38-43, 1993.
- SAXENA R et al.
Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome.
J.INTERN. MED. 1995, 237, 143-152

SUMMARY OF THE TEST PROCEDURE

	VOLUME TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
Dispense 50µL of microspheres reagent in each well to be used			
Incubation of the samples	100µL	Dilution buffer (B1): reagent blank	30 min. Room temperature
	100µL	Diluted negative control <i>in dilution buffer (B1)</i>	
	100µL	Diluted positive control <i>in dilution buffer (B1)</i>	
	100µL	Ready-to-use Calibrator (in duplicate)	
	100µL	Diluted samples <i>in dilution buffer (B1)</i>	
Washing 1	Wash twice in washing buffer (C1) (300µL/well)		
Incubation of the conjugate	100µL	Ready to use Anti-IgG conjugate	30 min Room temperature
Washing 2	Wash 1 time in Dilution buffer (B1) (100µL/well)		
Reading	Add 100µL/well of Dilution buffer (B1) and proceed immediately with the reading in FIDIS™ Instrument <i>The reading must be doing in the hour following adding the solution (B1). During this time, keep the plate at room temperature away from direct sunlight</i>		

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Microplate		Negative control		Microspheres
	Specimen diluent		Positive control		IgG Conjugate
	Calibrator		Wash Buffer		Microsphere reconstitution buffer
	FIDIS test		Contains sodium azide		Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com