FIDIS™ APS

REF MX 012 IgG 

REF MX 013 IgM 

DEFINITION

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BAS DE LAAT H., et al.

Pathogenic anti-B2GP1 antibodies recognize domain I of B2GP1 only after a conformational change. *Blood* 2006 ; 107:1916-1924.

BEVER E.M.

Cofacteurs involved in the antiphospholipid syndrome. *Editorial. Lupus* 1992, 1 : 51-53

CHAMLEY L.W., et al

Conformationally altered B2GP1 is the antigen for anti-cardiolipin antibodies. *Clin Exp Immunol* 1999 ; 115:571-576

DE GROOT P.G., et al.

Anti-prothrombin antibodies and their relation with thrombosis and lupus anticoagulant. *Lupus* 7 : suppl2, S32, 1998

FONT J.

Anti-cardiolipin antibodies in patients with autoimmune diseases: isotype distribution and clinical associations. *Clin. Rheum.* 1989, 8 : 475-785

GALLI M., et al.

Antiprothrombin antibodies : detection and clinical significance in the antiphospholipid syndrom. *Blood* 1999 ; 93:2149.

HESS D.C.

Anti-cardiolipin antibodies: a study of frequency in TIA and stroke. *Neurology* 1991, 23 (suppl) : 129-132

ICHIKAWA K., et al.

A chimeric antibody with the human gamma1 region as a putative standard for assays to detect IgG B2GP1 dependent anti-cardiolipin and anti-B2GP1 antibodies. *Arthritis Rheum* 42 (1999) 2461-2470.

INFANTE RIVARD C. et al.

Lupus anticoagulant, anti-cardiolipin antibodies and fetal loss. A case control study. *N. Engl. J. Med.* 1991, 325 : 1063-1066

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




REF **MX 013** **IgM** 

INTENDED USE

The **FIDIS™ APS** ([bmd](#)) kit is a semi-quantitative homogeneous fluorescent-based microparticles immunoassay using flow cytometry readings. It is designed for the simultaneous detection of 3 autoantibodies specificities: cardiolipin, β_2 -glycoprotein I (β_2 -GPI) and prothrombin also called Factor II (FII).

✦ **FIDIS™ APS IgG** ([bmd](#)) kit is designed for the detection of human IgG isotype antibodies

✦ **FIDIS™ APS IgM** ([bmd](#)) kit is designed for the detection of human IgM isotype antibodies

The **FIDIS™ APS** kits can be used with a dispensing/diluting device, the **CARIS™ system**.

SUMMARY AND EXPLANATION

Anti-cardiolipin autoantibodies are present in clinical contexts comparable to those in which lupus anticoagulants (anti-prothrombinase) are found. Moreover, they are frequently associated. In such patients, it is not rare to observe a dissociated syphilitic serology with positive VDRL and negative TPHA. The most frequent associated pathologies are Systemic Lupus Erythematosus, drug-induced lupus, thromboembolic disease, thrombosis in general, recurrent abortion, autoimmune peripheral cytopenia, primary biliary cirrhosis, certain acute infections and HIV and HTLV retrovirus infections. With the exception of infections, in which anti-cardiolipin antibodies seem only transiently present, the antibodies would seem to have a predictive value with respect to thrombosis. They exacerbate the life threat in myocardial infarction and enhance the risk of graft rejection following bypass procedures. The term, "anti-cardiolipin antibody syndrome" has been suggested to designate the morbid signs and symptoms that are associated with the presence of the antibodies.

The detection of some anionic aPLs found in autoimmune diseases, by ELISA technique calls for the presence of a serum protein, β_2 -glycoprotein I (β_2 -GPI) or apolipoprotein H. Conversely, infections which are not generally accompanied by thrombotic elements, do not require its presence. β_2 -GPI is a 50 kDa plasma protein, occurring at a concentration of 200 μ g/ml, in the form of a polypeptide chain of 326 aminoacids.

Its great affinity for the phospholipids is the basis of its in vitro anticoagulant effect on coagulation and platelet aggregation.

Today most authors consider that the β_2 -GPI-dependent aPLs are probably directed against the protein which has been structurally modified by linkage with the PLs.

Antibodies targeting prothrombin alone have recently been described being highly associated with fetal loss in APS patients.

Prothrombin antibodies are the first marker for this serious complication as all other known antiphospholipid antibodies failed to correlate with fetal loss.

Laboratory criteria of anti-phospholipid syndrome (APS) classification were revised in 2006 and recommendations of the Sydney international consensus are :

1. Determination of both anti-cardiolipin and anti- β_2 -GPI antibodies (IgG and/or IgM) in serum or plasma.
2. Presence of anti-cardiolipin antibodies (IgG and/or IgM) in medium or high titer (**i.e. >40GPL or MPL, or > the 99th percentile**), on two or more occasions, at least 12 weeks apart.
3. Presence of anti- β_2 -GPI antibodies (IgG and/or IgM) in titer > the 99th percentile, on two or more occasions, at least 12 weeks apart.

ASSAY PRINCIPLE

FIDIS™ APS ([bmd](#)) kits are based on the use of distinct uniform size color-coded microspheres and a benchtop flow cytometer interfaced to digital signal processing hardware and software. A red diode laser beam in the flow cytometer classifies each set of microspheres on the basis of its unique fluorescence intensity (red to orange) thus identifying which analyte is being tested. At the same time, a green laser beam illuminates the external second molecule fluorescence to quantify the reaction related to the specific analyte.

Each antigen required for the assay is covalently coupled to an individual set of microspheres through its surface functional groups.

The different antigen coupled microspheres are mixed together, constituting the final microsphere reagent.

FIDIS™ APS IgG ([bmd](#)) kits allow the detection of 3 IgG specificities of: anti-cardiolipin, anti- β_2 -GPI and anti-prothrombin autoantibodies.

FIDIS™ APS IgM ([bmd](#)) kits allow the detection of 3 IgM specificities of: anti-cardiolipin, anti- β_2 -GPI and anti-prothrombin autoantibodies.

The test is performed in a 96 well microplate with a filtering membrane at the bottom of the wells.

- In the first step, the sample is distributed in each well containing the microspheres mixture. If this sample contains one or more of the suspected antibodies, this(ese) antibody(ies) bind to the corresponding antigen(s) on the various sets of microspheres.

- After incubation, a wash step using a filtration process removes the unbound antibodies.

- A phycoerythrin labeled **anti-human IgG or IgM conjugate** is then added that binds to the previously bound antibodies.

- A final wash step allows to stop the reaction.

- The reaction is then directly measured by the flow cytometer, which differentiates each set of microspheres according to its fluorescence color while simultaneously measuring the average fluorescence emitted by the conjugate.

- A calibration system allows of each sample by interpolation the determination of the titer (AU/mL) for β_2 -GPI and prothrombin specificities and in standardized unit (GPL/mL or MPL/mL) for cardiolipin specificity.

ADDITIONAL MATERIAL – NOT SUPPLIED

- Precision pipettes capable of accurately delivering 5 μ L to 1000 μ L
- Multichannel Pipettes or dispensers capable of accurately delivering 40 μ L to 300 μ L or 5 μ L to 2 mL
- Serological pipettes
- Laboratory timer to monitor incubation steps
- Ultrasonic bath - Absorbent towel
- Microplate sealing films - Vortex mixer

SPECIMEN COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

STABILITY AND STORAGE

- Store reagents in their original packaging at +2°C to +8°C.
- Do not freeze reagents.
- Do not use kits beyond the expiration date.
- After use, store all components immediately back at +2°C to +8°C.

REAGENTS SUPPLIED

96 well microplate with filtering membrane and lid. MP	1 plate
Vial (A) of color-coded microsphere set sensitized by synthetic cardiolipin, human purified β_2 -GPI and human purified prothrombin antigen(s). <u>Lyophilized</u> (to be reconstituted with the buffer named D) MICROSPHERES	sq 6 mL
Vial (B3) of sample dilution buffer (white vial). <u>Ready to use</u> DIL SPE	2 x 115 mL
Vial of calibrator* <u>Ready to use</u> CAL <i>Each titer is printed on the vial label</i>	1 x 1,5 mL
Vial of positive control* concentrated. This control has a standard reactivity, which provides evidence of the proper reagents activity and proper assay performance. <u>To be diluted</u> CONTROL + <i>Expected values are printed on the vial label.</i>	1 x 250 μ L
Vial of negative control concentrated. <u>To be diluted.</u> CONTROL -	1 x 250 μ L
Vial of anti-human IgG or IgM coupled to phycoerythrin <u>Ready to use</u> CONJ IgG CONJ IgM	1 x 12 mL
Vial (C3) of washing buffer (black vial) <u>Ready to use</u> BUF WASH	1 x 100 mL
Vial (D) of reconstitution buffer for the microsphere set <u>Ready to use</u> BUF MICROSPHERES	1 x 6mL

Calibrator titers are expressed in arbitrary units per ml (AU/ml) excepted for cardiolipin.*

Calibrator titers for cardiolipin are expressed in GPL/mL or MPL/mL referred to SAPPORO standard HCAL (IgG) and EY2C9 (IgM).*

CAUTION

- Reagents in solution contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
- Avoid to use reagents if signs of contamination or other visible changes occur.
- The **FIDIS™ APS** ([bmd](#)) kits have been developed according to CE Directives 67/548/EEC and 1999/45/EC regarding classification, packaging and labelling of dangerous preparations.
- The **FIDIS™ APS** ([bmd](#)) kits have been optimized for use as described in this procedure. Do not substitute other manufacturer reagents. Dilution or alteration of these reagents may also alter the performance of the test. Follow thoroughly the test procedure to insure optimal performance.
- *Calibrators and controls are from human origin. The human sera used in the preparation of these products were tested and found non-reactive for*

antibodies to HIV-1, HIV-2, anti-HCV and Hepatitis B antigen. Because no test method can offer complete assurance that infectious agents are absent, handle as if capable of transmitting infectious diseases.

FIDIS™ APS (bmd) cannot be used in multibatch with the other FIDIS™ kits.

TESTS SET-UP

FIDIS™ Washer

☞ Check the tubing connecting the pump to the stand and the manometer setting (wheel completely closed).

Using FIDIS™ Analyzer and MLX-BOOSTER Software

☞ See the User's Manual provided with the **FIDIS™ Instrument** for detailed instructions on running the equipment and for calculation. For additional information and/or trouble shooting problems, please contact bmd subsidiaries or distributors.

1. Start the run as described in the User's Manual.
NOTE: The FIDIS™ takes 30 minutes to warm up after being turned on. A new warm up is necessary after 4 hours of system inactivity.
2. Calibration and controls are described in the User's Manual. Calibration and controls should be routinely performed once per month and each time a new lot of Sheath fluid reagent is used in order to insure optimal instrument performance.
Calibration should also be performed when the temperature is out of the determined range shown on the run batch screen.
3. Programming a batch or a multibatch of test protocol on the FIDIS™.
4. Load the microplate into the plate holder of the FIDIS™ as described in the User's Manual.
5. Analyze the results according the User's Manual.
6. When finished for the day, perform the sanitizing and soak operations, prior to turning the Analyzer off according to the shutting down procedure described in the User's Manual.

Reagents Preparation

☞ Remove the individual components from storage, **allow them to warm up to room temperature (+18°C to +25°C)**, and mix them well.

1. Microsphere preparation

- Reconstitute the microsphere solution by adding the vial **D** of reconstitution buffer to the vial **A** of lyophilized microspheres. Wait 5 minutes and vortex.
- Stable 2 months at +2°C to +8°C, after reconstitution.

2. Preparation of samples and controls.

- Dilute the samples and controls at **1:101** in dilution buffer (**B3**).
Ex.: 10µL sample in 1000µL dilution buffer (**B3**).
- **Vortex vigorously.**

3. Assay configuration

Use the work-sheet included in the kit to identify the location of the samples.

a. When setting up the test, systematically take into account the following well requirements: ⇨ See examples below

- ✦ 1 "reagent-blank" well
- ✦ 1 well for negative control
- ✦ 1 well for positive control
- ✦ 2 "calibrator" wells

b. Calculation of the correct number of wells necessary and their location.

In the following examples, 2 different configurations are described, according to the number of wells needed, and the availability of unused wells. The sample dispensing must be systematically carried out in a column. Do not leave empty wells.

The same microplate can be used for more than one series of tests if unused wells are protected by microplate sealing film. Any single well can only be used once.

To avoid any error in using a well more than once, identify the used wells with a marker (see example 2 below).


c. Programming of the test protocol on the FIDIS™

☞ Refer to the FIDIS™ user's manual:

"Creating a batch or a multi-batch".


Example 1:

- The microplate is totally blank (no well has yet been used): Use A1 as the 1st well.

Example 2:

- 16 wells have already been used: Use A3 as the 1st well.

ASSAY PROCEDURE

1. *Microspheres dispensing*

- Before to start the assay and if it is necessary, protect the unused wells with a microplate sealing film.
- **Vortex the microspheres reagent vigorously for 20s** and dispense 50µL in each well to be used.

PLEASE NOTE:

- If **CARIS™** is used, the microspheres should be sonicated for 5 minutes.
- If a dispensing/diluting device is used (like **CARIS™**), the procedure steps 1 and 2 should be reversed:
⇒ dispense diluted sera in first and then, dispense microspheres

2. *Addition and Incubation of the samples*

- Dispense 100µL of sample dilution buffer (**B3**) for reagent blank.
- Dispense 100µL of prediluted controls, of ready to use calibrator and of prediluted samples.
- Cover the plate and incubate 30 minutes at room temperature, away from direct sunlight and without shaking.

3. *Wash step 1*

Wash the plate 2 times with washing buffer (**C3**) using filtration unit:

- Remove the microplate lid (verify if the vacuum break valve is closed) and place it on the washer.
- Start the pump. Stop it when the whole incubation liquid has been removed in totality and open the vacuum break valve.
Close the vacuum break valve.
- Dispense 300µL of washing buffer (**C3**) in all the used wells.
- Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve.
Close the vacuum break valve.
- Dispense 300µL of washing buffer (**C3**) in all the used wells.
- Start the pump. Stop it after 5 additional seconds when all the buffer has been removed and open the vacuum break valve.
Close the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by tapping vigorously the plate 10 times on an absorbent towel.
- Place again the microplate on the washer and start the pump for 5 seconds. Stop it and quickly open the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by blotting vigorously the plate on an absorbent towel.
- Place the plate on a totally dry surface before starting the conjugate incubation.

4. *Incubation of the conjugate*

- Dispense 100µL of ready to use conjugate in each well used.

- Cover the plate and incubate 30 minutes at room temperature away from direct sunlight and without shaking. The incubation time starts after the conjugate has been added to all wells. If this timing is not followed, the results might be erroneous.

Note : the conjugate is photosensitive.
⇒ **After using, close the vial carefully.**

5. *Wash step 2*

Wash the plate using the filtration unit on 1 cycle with dilution buffer (**B3**).

Use two different reagent reservoirs for the 2 washing steps.

*Either clean and dry correctly the reservoir between both wash steps (avoid to mix both **B** and **C** buffers).*

- Remove the microplate lid (verify that the vacuum break valve is closed) and place it on the washer.
- Start the pump. Stop it when the whole incubation liquid has been removed and open the vacuum break valve. Close the vacuum break valve.
- Dispense 100µL of dilution buffer (**B3**) in all the used wells.
- Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by blotting the plate on an absorbent towel.
- Dispense 100µL of dilution buffer (**B3**) in all the used wells. Place the plate on a totally dry surface at room temperature.

The reading must be doing in the hour following adding the solution (B3). During this time, keep the plate at room temperature away from direct sunlight.

6. *Reading of the test*

☞ Follow the **MLX-BOOSTER™** and **FIDIS™** user's manual: Processing batches.

VALIDATION CRITERIA OF THE RESULTS

Calibrator, negative and positive controls have to be run with every batch of samples to ensure that all reagents and procedures performed properly.

In order to validate the results, all the criteria listed below must be met. Otherwise, the test is invalid and must be repeated.

- The positive controls should show a value within the limits printed on the corresponding vial labels.
- The negative controls should be less than 10 AU/mL or 10 GPL/mL or 10MPL/mL.

The microsphere set contains Internal standard beads allowing to verify the presence of serum in wells, as well as the good respect of the test procedure. Wells presenting a not corresponding fluorescent signal will be no valid with **MLX BOOSTER™** software.

INTERPRETATION OF RESULTS

Arbitrary Units (UA/ml) Except cardiolipin	< 10	10 - 15	> 15
Result	Negative	Equivocal*	Positive

Standardized units (GPL/ml ou MPL/ml) for cardiolipin	< 10	10 - 15	> 15
Result	Negative	Equivocal*	Positive

(*) Equivocal results should be controlled on a second sample and the interpretation of the results should be done in the frame of additional testing and taking into account the clinical status of the patient.

Proposed Clinical classification

As per the Sydney international consensus meeting guidelines (S MIYAKIS, 2006); the positive threshold of the **FIDIS™ APS** (bmd), was established by the 99th percentile of the values obtained for a normal population (400 samples from « blood donors » and « biological interferences » populations).

Based on the evaluation of the **FIDIS™ APS** (bmd), bmd proposes the following classification; where a medium to high level has the most relevant clinical association.

≤15 GPL/mL or MPL/mL (99 th percentile)	16 to 39 GPL/mL or MPL/mL	40 to 80 GPL/mL or MPL/mL	> 80 GPL/mL or MPL/mL
Normal Level	Low Level	Medium Level	High Level

Each laboratory should establish and maintain its own references normal and clinical range values, based on the patient population and other local factors.

CALCULATION OF RESULTS

☞ The results are automatically calculated by MLX-BOOSTER™ software and can be printed for each assay.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

The serums from Immunodeficient patients will give a no valid result.

QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities. Immunotrol I and II multiparametric controls (bmd, Cat. n° HM036 and HM037) contain human auto-antibodies directed against researched specificities.

To be assayed in the same manner as the unknown samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

SENSITIVITY, SPECIFICITY AND CONCORDANCE COMPARED TO OTHERS METHODS

FIDIS™ APS (bmd) has been compared to 2 marketed ELISA kits:

The study was performed on 130 samples

Results

FIDIS APS IgG				
Antigenic specificity	Number samples	Relative sensitivity (%)	Relative specificity (%)	overall agreement (%)
Cardiolipin	85	95.5	88.9	94.1
β ₂ GPI	85	90.0	92.7	91.8
Prothrombin	82	100.0	94.8	95.1

FIDIS APS IgM				
Antigenic specificity	Number samples	Relative sensitivity (%)	Relative specificity (%)	overall agreement (%)
Cardiolipin	130	94.4	95.5	95.4
β ₂ GPI	128	90.0	96.6	96.1
Prothrombin	67	100.0	96.9	97.0

Method for determining threshold values

Calculated on the 400 samples from « blood donors » and « biological interferences » populations, the positive threshold (15UA/mL or 15GPL/mL or 15MPL/mL) corresponds to the 99th percentile of the values distribution for the 3 specificities.

The negative threshold (10UA/mL or 10GPL/mL or 10MPL/mL) corresponds to the 96th percentile for the 3 specificities.

Between these threshold, results are considered equivocal.

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





MIYAKIS S. et al.

International consensus statement on a update of the classification criteria for definite anti-phospholipid syndrome (APS). *J. Thromb. Haemost* 2006, 4:295-306.

SUMMARY OF THE TEST PROCEDURE

	VOLUME TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
Dispense 50µL of microspheres reagent in each well to be used			
Incubation of the samples	100µL	Dilution buffer (B3) : Reagent blank	30 min. Room temperature
	100µL	Diluted negative control 1:101	
	100µL	Diluted positive control 1:101	
	100µL	Calibrator (in duplicate)	
	100µL	Diluted samples 1:101	
Washing 1	Wash twice in washing buffer (C3) (300µL/well)		
Incubation of the conjugate	100µL	Ready to use Anti-IgG or anti-IgM conjugate	30 min Room temperature
Washing 2	Wash 1 time in Dilution buffer (B3) (100µL/well)		
Reading	<p>Add 100µL/well of Dilution buffer (B3) and proceed immediately with the reading in FIDIS™ Instrument</p> <p><i>The reading must be doing in the hour following adding the solution (B3). During this time, keep the plate at room temperature away from direct sunlight.</i></p>		

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Microplate		Negative control		Microspheres
	Specimen diluent		Positive control		IgG Conjugate
	Calibrator		Wash Buffer		Microsphere reconstitution buffer
	FIDIS test		Contains sodium azide		Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com