FIDIS™ Rheuma - RF

REF

MX 103

96

DÉFINITION

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REF **MX 103**

INTENDED USE

The **FIDIS™ Rheuma - RF** kit is a quantitative homogeneous fluorescent-based microparticles immunoassay using flow cytometry readings. It is designed for the simultaneous detection of IgM class rheumatoid factors directed against Fc determinants from human or animal G Immunoglobulins.

SUMMARY AND EXPLANATION

Rheumatoid factors (RF) form a population of autoantibodies reacting with the various antigenic sites located on the Fc fragment of the IgG class immunoglobulins. These antigen determinants may be specific to human IgG or common to the IgG of various animal species. Rheumatoid factors may belong to any immunoglobulin class, but the IgM class is dominant.

Rheumatoid factors are essentially associated to inflammatory rheumatic pathologies. The frequency of their presence is increased in the serum of patients suffering from rheumatoid polyarthritis. RF specificity is extended due to the fact that there is an associated positivity between human IgG and rabbit IgG.

Routinely, the diagnostic methods used are: latex agglutination test to identify IgM class RF directed against human antigen determinants and the Waaler-Rose haemagglutination test to identify the IgM class RF directed against rabbit antigen determinants.

However, EIA methods offer increased sensitivity and specificity as well as the advantage of a simultaneous determination of both specificities.

Apart from rheumatic pathologies, RF might appear in specific bacterial and viral infectious diseases or in inflammatory pathologies.

The prevalence is, then, ca. 2 %, the titer is low and the reactivity directed against human IgG solely. In a normal population, the RF prevalence is of 2 to 5%.

ASSAY PRINCIPLE

FIDIS™ Rheuma - RF kits are based on the use of distinct uniform size color-coded microspheres and a benchtop flow cytometer interfaced to digital signal processing hardware and software. A red diode laser beam in the flow cytometer classifies each set of microspheres on the basis of its unique fluorescence intensity (red to orange) thus identifying which analyte is being tested. At the same time, a green laser beam illuminates the external second molecule fluorescence to quantify the reaction related to the specific analyte.

Each antigen required for the assay is covalently coupled to an individual set of microspheres through its surface functional groups. The different antigen coupled microspheres are mixed together, constituting the final microsphere reagent.

FIDIS™ Rheuma – RF kits allow the detection of IgM class rheumatoid factors directed against **Fc determinants from human or animal G Immunoglobulins**.

The test is performed in a 96 well microplate with a filtering membrane at the bottom of the wells.

- In the first step, the sample is distributed in each well containing the microspheres mixture. If this sample contains one or more of the suspected antibodies, this(ese) antibody(ies) binds to the corresponding antigen(s) on the various sets of microspheres.
- After incubation, a wash step using a filtration process removes the unbound antibodies.
- A phycoerythrin labeled **anti-human IgM** conjugate is then added that binds to the previously bound antibodies.
- A final wash step allows to stop the reaction.
- The reaction is then directly measured by the flow cytometer, which differentiates each set of microspheres according to its fluorescence colour while simultaneously measuring the average fluorescence emitted by the conjugate.
- A calibration system allows the determination of the titer (IU/mL) of each sample by interpolation for each antigenic specificity.

The **FIDIS™ Rheuma - RF** kit could be used with **CARIS™** system (diluting and dispensing device).

SPECIMEN COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

REAGENTS SUPPLIED

96 wells microplate with filtering membrane and lid. MP	1 plate
Vial (A) of color-coded microsphere set sensitized by the Fc fragment of Human or Animal Immunoglobulins antigen(s). <u>Lyophilized</u> (to be reconstituted with the buffer named D) MICROSPHERES	Sq 6 mL
Vial (B2) of sample dilution buffer (white vial). <u>Ready to use</u> DIL SPE	2 x 115mL
Vial of calibrator* titered in international unit for the specificities to be measured. <u>Ready to use</u> CAL <i>Each titer is printed on the vial label</i>	1 x 1,5 mL
Vial of positive control concentrated. This control has a standard reactivity, which provides evidence of the proper reagents activity and proper assay performance. <u>To be diluted</u> CONTROL + <i>Expected values are printed on the vial label.</i>	1 x 250 µL
Vial of negative control concentrate <u>To be diluted</u> CONTROL -	1 x 250 µL
Vial of anti-human IgM coupled to phycoerythrin <u>Ready to use</u> CONJ IgM	1 x 12 mL
Vial (C2) of washing buffer (black vial) <u>Ready to use</u> BUF WASH	1 x 100 mL
Vial (D) of reconstitution buffer for the microsphere set <u>Ready to use</u> BUF MICROSPHERES	1 x 6mL

Calibrator* titers are printed on the vial label and are expressed in international units per ml (IU/ml) and standardized against the International Reference: 64/2 (1st British Standard for Rheumatoid Arthritis)

ADDITIONAL MATERIAL – NOT SUPPLIED

- Precision pipettes capable of accurately delivering 5µL to 1000µL
- Multichannel Pipettes or dispensers capable of accurately delivering 40µL to 300µL or 5µL to 2 mL
- Vortex mixer
- Laboratory timer to monitor incubation steps
- Ultrasonic bath
- Absorbent towel
- Serological pipettes
- Microplate sealing films

STABILITY AND STORAGE

- Store reagents in their original packaging at +2°C to +8°C.
- Do not freeze reagents.
- Do not use kits beyond the expiration date.
- After use, store all components immediately back at +2°C to +8°C.

CAUTION

- Reagents in solution contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
- Avoid to use reagents if signs of contamination or other visible changes occur.
- The **FIDIS™ Rheuma - RF** kit has been developed according to CE Directives 67/548/EEC and 1999/45/EC regarding classification, packaging and labelling of dangerous preparations.
- The **FIDIS™ Rheuma - RF** kit has been optimized for use as described in this procedure. Do not substitute other manufacturer reagents. Dilution or alteration of these reagents may also alter the performance of the test. Follow thoroughly the test procedure to insure optimal performance.
- *Calibrators and controls are from human origin. The human sera used in the preparation of these products were tested and found non-reactive for antibodies to HIV-1, HIV-2, anti-HCV and Hepatitis B antigen. Because no test method can offer complete assurance that infectious agents are absent, handle as if capable of transmitting infectious diseases.*

FIDIS™ Rheuma - RF cannot be used in multi-batch with the other FIDIS™ kits.

TEST SET-UP

FIDIS™ Washer

- ☞ Check the tubing connecting the pump to the stand and the manometer setting (wheel completely closed).

Using FIDIS™ Analyzer and MLX-BOOSTER Software

- ☞ See the User's Manual provided with the **FIDIS™ Instrument** for detailed instructions on running the equipment and for calculation. For additional information and/or trouble shooting problems, please contact bmd subsidiaries or distributors.

1. Start the run as described in the User's Manual.
NOTE: The FIDIS™ takes 30 minutes to warm up after being turned on. A new warm up is necessary after 4 hours of system inactivity.
2. Calibration and controls are described in the User's Manual. Calibration and controls should be routinely performed once per month and each time a new lot of Sheath fluid reagent is used in order to insure optimal instrument performance. Calibration should also be performed when the temperature is out of the determined range shown on the run batch screen.

3. Programming a batch or a multibatch of test protocol on the FIDIS™ as described in the User's Manual.
4. Load the microplate into the plate holder of the FIDIS™ as described in the User's Manual.
5. Analyze the results according the User's Manual.
6. When finished for the day, perform the sanitizing and soak operations, prior to turning the Analyzer off according to the shutting down procedure described in the User's Manual.

Reagents Preparation

☞ Remove the individual components from storage, allow them to warm up to **room temperature** (+18°C to +25°C), and mix them well.

1. Preparation of samples and controls.

- Dilute the samples and controls at **1:201** in dilution buffer (**B2**).
Ex.: 10µL sample in 2000µL dilution buffer (B2).
- **Vortex vigorously.**

2. Assay configuration

Use the work-sheet included in the kit to identify the location of the samples.

- a. *When setting up the test, systematically take into account the following well requirements:* ⇒ See examples below

- ✦ 1 "reagent-blank" well
- ✦ 1 well for negative control
- ✦ 1 well for positive control
- ✦ 2 "calibrator" wells

- b. *Calculation of the correct number of wells necessary and their location.*

In the following examples, 2 different configurations are described, according to the number of wells needed, and the availability of unused wells. The sample dispensing must be systematically carried out in a column. Do not leave empty wells.

The same microplate can be used for more than one serie of tests if unused wells are protected by microplate sealing film. Any single well can only be used once.

To avoid any error in using a well more than once, identify the used wells with a marker (see example 2 below).

- c. *Programming of the test protocol on the FIDIS™*

☞ Refer to the FIDIS™ user's manual: "Creating a batch".

Example 1:

- *The microplate is totally blank (no well has yet been used): Use A1 as the 1st well.*

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	S4	S12	S20								
B	Ctrl-	S5	S13	S21								
C	Ctrl+	S6	S14	S22								
D	Cal	S7	S15	S23								
E	Cal	S8	S16	S24								
F	S1	S9	S17	S25								
G	S2	S10	S18									
H	S3	S11	S19									

Example 2:

- *16 wells have already been used: Use A3 as the 1st well.*

	1	2	3	4	5	6	7	8	9	10	11	12
A			Blank	S4	S12	S20	S28	S36	S44	S52		
B			Ctrl-	S5	S13	S21	S29	S37	S45	S53		
C			Ctrl+	S6	S14	S22	S30	S38	S46			
D			Cal	S7	S15	S23	S31	S39	S47			
E			Cal	S8	S16	S24	S32	S40	S48			
F			S1	S9	S17	S25	S33	S41	S49			
G			S2	S10	S18	S26	S34	S42	S50			
H			S3	S11	S19	S27	S35	S43	S51			

ASSAY PROCEDURE

1. Microspheres dispensing

- Before to start the assay and if it is necessary, protect the unused wells with a microplate sealing film.
- **Vortex the microspheres reagent vigorously for 20s** and dispense 50µL in each well to be used.

PLEASE NOTE:

- ⓘ If CARIS™ is used, the microspheres should be sonicated for 5 minutes.
- ⓘ If a dispensing/diluting device is used (like CARIS™), the procedure steps 1 and 2 should be reversed: ⇒ dispense diluted sera in first and then, dispense microspheres.

2. Addition and Incubation of the samples

- Dispense 100µL of sample dilution buffer (**B2**) for reagent blank.
- Dispense 100µL of prediluted controls, of ready to use calibrator and of prediluted samples.
- Cover the plate and incubate 30 minutes at room temperature, away from direct sunlight and without shaking.

3. Wash step 1

Wash the plate 2 times with washing buffer (**C2**) using filtration unit:

- a. Remove the microplate lid (verify if the vacuum break valve is closed) and place it on the washer.
- b. Start the pump. Stop it when the whole incubation liquid has been removed in totality and open the vacuum break valve. Close the vacuum break valve.
- c. Dispense 300µL of washing buffer (**C2**) in all the used wells.
- d. Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- e. Dispense 300µL of washing buffer (**C2**) in all the used wells.
- f. Start the pump. Stop it after 5 additional seconds when all the buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- g. Remove the microplate from the washer and remove the residual buffer by tapping vigorously the plate 10 times on an absorbent towel.

- h. Place again the microplate on the washer and start the pump for 5 seconds. Stop it and quickly open the vacuum break valve.
- i. Remove the microplate from the washer and remove the residual buffer by blotting vigorously the plate on an absorbent towel.
- j. Place the plate on a totally dry surface before starting the conjugate incubation.

4. Incubation of the conjugate

- Dispense 100µL of ready to use conjugate in each well used.
- Cover the plate and incubate 30 minutes at room temperature away from direct sunlight and without shaking. The incubation time starts after the conjugate has been added to all wells. If this timing is not followed, the results might be erroneous.

Note : the conjugate is photosensitive.
 ⇒ After using, close the vial carefully.

5. Wash step 2

Wash the plate using the filtration unit on 1 cycle with dilution buffer (B2).

Use two different reagent reservoirs for the 2 washing steps.
 Either clean and dry correctly the reservoir between both wash steps (avoid to mix both B2 and C2 buffers).

- a. Remove the microplate lid (verify that the vacuum break valve is closed) and place it on the washer.
- b. Start the pump. Stop it when the whole incubation liquid has been removed and open the vacuum break valve.
Close the vacuum break valve.
- c. Dispense 100µL of dilution buffer (B2) in all the used wells.
- d. Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve.
Close the vacuum break valve.
- e. Remove the microplate from the washer and remove the residual buffer by blotting the plate on an absorbent towel.
- f. Dispense 100µL of dilution buffer (B2) in all the used wells. Place the plate on a totally dry surface at room temperature.

The reading must be doing in the hour following adding the solution (B2). During this time, keep the plate at room temperature away from direct sunlight.

6. Reading of the test

☞ Follow the **MLX-BOOSTER™** and **FIDIS™** user's manual: Processing batches.

VALIDATION CRITERIA OF THE RESULTS

Calibrator, negative and positive controls have to be run with every batch of samples to ensure that all reagents and procedures performed properly.

In order to validate the results, all the criteria listed below must be met. Otherwise, the test is invalid and must be repeated.

- a. The positive controls should show a value within the limits printed on the corresponding vial labels.
- b. The negative controls should be less than 25 IU/mL.

The microsphere set contains Internal standard beads allowing to verify the presence of serum in wells, as well as the good respect of the test procedure. Wells presenting a not corresponding fluorescent signal will be no valid with MLX BOOSTER™ software.

INTERPRETATION OF RESULTS

International Units (IU/ml) for the 2 specificities	< 25 IU/ml	25 – 30 IU/ml	> 30 IU/ml
Interpretation	Negative	Borderline ^(*)	Positive

(*) Borderline results should be repeated on a second sample and the interpretation of the results should be treated in the frame of additional testing and taking into account the clinical status of the patient.

Each laboratory should establish and maintain its own references (normal) range values, based on the patient population and other local factors.

CALCUL OF RESULTS

☞ The results are automatically calculated by MLX-BOOSTER™ software and can be printed for each assay.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

The serums from Immunodeficient patients will give a no valid result.



QUALITY CONTROL

It is recommended to use internally or externally sourced control material for the different specificities. HM050) contains human and animal antibodies directed against human specific IgM isotype Rheumatoid Factor (RF). This material is to be assayed in the same manner as the unknown samples.

CHARACTERISTICS AND PERFORMANCES OF THE TEST

CORRELATION

The performances of **FIDIS™ Rheuma-RF**, lyophilized form were assessed on a population of 265 samples and compared to the performances of **FIDIS™ Rheuma-RF**, liquid form. A linear regression analysis of the two products showed that both are equivalent.

REPRODUCIBILITY / PRECISION (LYOPHILIZED FORM)

Antigenic specificity	Within-run (10 tests in the same run)		Between-run (5 tests in 5 different runs)	
	Mean value (IU/mL)	CV (%)	Mean value (IU/mL)	CV (%)
FRh	65	0.8	92	7.1
FRh	103	2.8	201	9.9
FRh	407	2.0	492	3.7
FRa	82	3.6	127	15.0
FRa	143	3.1	154	11.0
FRa	443	1.4	504	5.0

INDICATION OF THE DIAGNOSTIC VALUE IN COMPARISON WITH OTHER METHODS

FIDIS™ Rheuma - RF kit has been compared to:

⇒ FR-LISA kit (bmd) based on ELISA method for the detection of rheumatoid factors from human and animal origin.

Confirmation methods for discrepant results:

⇒ ELISA (CLB) for human specificity and,
⇒ POLYARTITRE (Fumouze) based on haemagglutination method for animal specificity.

The study was performed on:

⇒ 81 samples related to rheumatoid diseases
⇒ 64 samples from blood donors
⇒ 78 samples selected for their potential biological interferences (cryoglobulinemia, hypergammaglobulinemia, IgG and IgM monoclonal immunoglobulins, complement, immune complexes).

Results

Human specificity		ELISA	
		Positive	Negative
FIDIS™	Positive	95	1
	Negative	0	127

Animal specificity		ELISA	
		Positive	Negative
FIDIS™	Positive	53	1
	Negative	1	168

Performances for the selected population

Antigen	Relative sensitivity (%)	Relative specificity (%)	Relative accuracy (%)
Human specificity	100	99.2	99.6
Animal specificity	98.2	99.4	99.1

Discrepant results analysis (3/223)

Antigenic specificity	FIDIS™ result	ELISA bmd	ELISA CLB or Hémagl	Associated specificities
FRh	positive (limit)	borderline	negative	Blood donor
FRa	borderline	positive	weakly positive	human
FRa	positive (limit)	negative	negative	human

The 3 results concern samples in limit of positivity, which can lead to no significant discrepancies.

Method for determining threshold values

Calculated on the 142 samples from « blood donors » and « biological interferences » populations and according standard 64/2.

The positive threshold (30IU/ml) corresponds to the 99.5th percentile of the values distribution for the 2 specificities.

The negative threshold (25IU/ml) corresponds to the 96.5th percentile for the two specificities.

Between these thresholds, results are considered borderline.

BIBLIOGRAPHY

WILLIAMS DG.
Rheumatoid arthritis : autoimmunity in rheumatoid arthritis.
In : Klippel JH, Dieppe PA, eds. Rheumatology, 1st edition.
St louis 1994 ;91 :9-11.

SHMERLING RH, DELBANCO TL.
The rheumatoid factor : an analysis of clinical utility.
Am J Med 1991 ;91 :528-534.

BORRETZEN M et al.
Rheumatoid Factors.
In : Autoantibodies. JB Peter and Y. Schoenfeld Eds.
Elsevier Sciences B.V. 1996 ;706-715.

ARNETT FC et al.
The American Rheumatism Association 1987 revised criteria for classification of rheumatoid arthritis.
Arthritis Rheum 1988 ;31 :315-324.






MOORE TT et al.
Rheumatoid Factors.
Clin Biochem 1993 ;26 :75-84.

CORPER AL et al.
Structure of human IgM rheumatoid factor Fab bound to its autoantigen IgG Fc reveals a novel topology of antibody-antigen interaction.
Nat Struct Biol. 1997 May;4(5):374-81.

SUMMARY OF THE TEST PROCEDURE

	VOLUME TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
Dispense 50µL of microspheres reagent in each well to be used			
Incubation of the samples	100µL	Dilution buffer (B2): reagent blank	30 min. Room temperature
	100µL	Diluted negative control <i>in dilution buffer (B2)</i>	
	100µL	Diluted positive control <i>in dilution buffer (B2)</i>	
	100µL	Ready-to-use Calibrator (in duplicate)	
	100µL	Diluted samples <i>in dilution buffer (B2)</i>	
Washing 1	Wash twice in washing buffer (C2) (300µL/well)		
Incubation of the conjugate	100µL	Ready to use Anti-IgM conjugate	30 min Room temperature
Washing 2	Wash 1 time in Dilution buffer (B2) (100µL/well)		
Reading	Add 100µL/well of Dilution buffer (B2) and proceed immediately with the reading in FIDIS™ Instrument The reading must be doing in the hour following adding the solution (B2). During this time, keep the plate at room temperature away from direct sunlight		

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Microplate		Negative control		Microspheres
	Specimen diluent		Positive control		IgM Conjugate
	Calibrator		Wash Buffer		Microsphere reconstitution buffer
	FIDIS test		Contains sodium azide		Reconstitute with

BioMédical Diagnostics SA

Office
Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com