Anti-R-TSH

REF RLTRE 96.3

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RÉFÉRENCES

B. Rees Smith et al

A new assay for thyrotropin receptor autoantibodies.

Thyroid 2004 14: 830-835

K. Kamijo et al

Clinical Evaluation of 3rd Generation assay for Thyrotropin Receptor Antibodies : The M22-biotin-based ELISA initiated by Smith

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Internet : www.bmd-net.com

Anti R-TSH

Quantitative determination of thyrotropin receptor autoantibodies in human serum

REF **RLTRE 96.3**

INTENDED USE

The **Anti R-TSH kit** () is intended for use by professional persons only for the quantitative determination of thyrotropin receptor autoantibodies in human serum. Hyperthyroidism in Graves' disease is due to the presence of autoantibodies to the TSH receptor and measurement of these autoantibodies can be useful in disease diagnosis and management.

ASSAY PRINCIPLE

In the ELISA, TSH receptor autoantibodies in patient sera, calibrators and controls are allowed to interact with TSH receptors coated onto ELISA plate wells.

After a 2 hour incubation, the samples are discarded leaving TRAb bound to the immobilised receptor.

A human monoclonal autoantibody to the TSH receptor labelled with biotin (M22-biotin) is added in a 2nd incubation step, where it interacts with immobilised TSH receptors, which have not been blocked by bound TRAb.

The amount of M22-biotin bound to the plate is then determined in a third incubation step by addition of streptavidin-peroxidase, which binds specifically to biotin.

Excess unbound streptavidin-peroxidase is then washed away and addition of 3,3',5,5'-Tetramethylbenzidine (TMB) resulting in the formation of a blue colour.

This reaction is stopped by addition of stop solution causing the well contents to turn from blue to yellow.

The absorbance of the yellow reaction mixture at 450nm is then read using an ELISA plate reader. A lower absorbance indicates the presence of TRAb in a test sample (as TRAb inhibit the binding of M22-biotin to TSH receptor coated plate wells).

MATERIALS REQUIRED AND NOT SUPPLIED

- ◆ Pipettes capable of dispensing 50µL, 75µL, 100µL and appropriate volumes for diluting SAPOD (F).
- ◆ Means of diluting concentrated wash (I).
- ◆ Pure water.
- ◆ ELISA plate reader suitable for 96 well formats and capable of measuring at 450nm.
- ◆ ELISA Plate shaker, capable of 500 shakes/min (not an orbital shaker).

STORAGE AND PREPARATION OF SERUM SAMPLES

- Sera to be analysed should be assayed soon after separation or stored, preferably in aliquots, at or below -20°C.
- 150µL is sufficient for one assay (duplicate 75µL determinations are recommended).
- Repeated freeze-thawing or increases in storage temperature must be avoided.
- Incorrect storage of serum samples can lead to loss of TRAb activity.
- Do not use lipaemic or haemolysed serum or samples containing particulates.
- Do not use plasma in the assay. When required, thaw test sera at room temperature and mix gently to ensure homogeneity.
- Centrifuge the serum prior to assay (preferably for 5 minutes at 10-15 000g in a microfuge) to remove any particulate matter. Please do not omit this centrifugation step for sera that are cloudy or contain particulates.

PREPARATION OF REAGENTS SUPPLIED

A TSH Receptor Coated Wells

12 breakapart strips of 8 wells (96 in total) in a frame and sealed in foil bag.

Allow to stand at room temperature for at least 20 minutes before opening.

MP

Ensure stripwells are fitted firmly into frame provided. After opening return any unused wells to the original foil packet and seal with adhesive tape. Place foil bag in the self-seal plastic bag with desiccant provided, store between +2°C and +8°C for up to 12 weeks.

B Start Buffer

10 mL – Coloured yellow

Ready to use

BUF **START**

C1-5 Calibrators

0.4, 1, 2.5, 10 and 30 u/L (units are NIBSC 90/672)

5 x 1.0 mL

Ready to use

CAL **C?**

D1 Negative control

1.0 mL

Ready to use

CONTROL **-**

- D2 Positive control**
(See label for concentration range)
1.0 mL
Ready to use
- E M22-Biotin**
15 mL – Coloured red
Ready to use
- F Streptavidin-Peroxidase (SAPOD)**
1 x 0.75 mL
Concentrated
Dilute 1 in 20 with diluent for SAPOD (G). For example, 0.5mL (F) + 9.5 mL (G).
Store at +2°C/+8°C after dilution for up to kit expiry date.
- G Diluent for SAPOD**
15 mL
Ready to use
- H Peroxidase Substrate (TMB)**
15 mL
Ready to use
- I Concentrated wash solution**
100 mL
Concentrated
Dilute to 1 litre with pure water before use. Store at +2°C/+8°C up to kit expiry.
- J Stop solution**
10 mL
Ready to use

CONTROL	+
---------	---

BIOTIN	M22
--------	-----

SA POD

DIL	SA POD
-----	--------

SUBS	TMB
------	-----

BUF	WASH	10x
-----	------	-----

SOLN	STOP
------	------

Store unopened kits and all kit components (A-J) at +2°C/+8°C.

ASSAY PROCEDURE

Allow all reagents and test samples to stand at room temperature (between +20°C and +25°C) for at least 30 minutes. A repeating pipette is recommended for steps 1, 5, 8, 10 & 11.

- Pipette **75µL** of **start buffer** (B) into each well to be used, leaving one well empty for blank (step 12).
- Pipette **75µL** of **calibrators** (C1-5), **controls** (D1 and D2) and **test sera** into respective wells (start with the 30u/L calibrator and descend down the plate to the negative control and then test sera) (except blank).
- Cover the frame and shake the wells for 2 hours at room temperature (+20°C/+25°C) on an ELISA plate shaker (500 shakes per min.).
- Aspirate well contents by use of a plate washing machine or discard by briskly inverting the frame of stripwells over a suitable receptacle. Wash the wells once with **350µL of diluted Wash Solution** (I), and aspirating the wash by use of a plate washing machine or discard the wash by briskly inverting the frame of stripwells over a suitable receptacle. Tap the inverted wells gently on a clean dry absorbent surface to remove excess wash solution (only necessary if washing plate by hand).
- Pipette **100µL of M22-Biotin** (E) into each well (except blank). Avoid splashing the material out of the wells during addition.
- Cover the plate, and incubate at room temperature for 25 minutes without shaking.
- Repeat wash step 4.

- Pipette **100µL of diluted Streptavidin-peroxidase** (F) into each well (except blank) and incubate at room temperature for 20 minutes without shaking.
- After incubation, discard the well contents by briskly inverting the frame of stripwells over a suitable receptacle. Wash the wells twice with **350µL of diluted Wash Solution** (I) followed by once with pure water (to remove any foam) and tap the inverted wells gently on a clean dry absorbent surface to remove excess wash solution. (if a plate washing machine is used, the plate can be washed 3 times with diluted wash solution (I) only).
- Pipette **100µL of TMB Substrate** (H) into each well (including blank) and incubate in the dark at room temperature for 30 minutes without shaking.
- Pipette **50µL Stop solution** (J) into each well (including blank) and shake the plate for approximately 5 seconds on a plate shaker. Ensure substrate incubation times are the same for each well.
- Read the absorbance of each well at 450nm using an ELISA plate reader, blanked against the well containing **100µL of TMB substrate** (H) and **50µL Stop solution** (J) only.

RESULTS ANALYSIS

A calibration curve can be established by plotting calibrator concentration on the x-axis (log scale) against the absorbance of the calibrators on the y-axis (linear scale).

The TRAb concentrations in patient sera can then be read off the calibration curve.


Other data reduction systems can be used. Results can also be expressed as inhibition (%I) of M22 binding calculated using the formula;

$$100 \times \left\{ 1 - \frac{\text{test sample absorbance at 450 nm}}{\text{Negative Control (D1) absorbance at 450 nm}} \right\}$$

TYPICAL RESULTS

(Example only, not for use in calculation of actual results)

Sample	Absorbance at 450nm (minus blank)	%I	u/L
Control D1	2.290	0	0
C1	1.987	13	0.4
C2	1.618	29	1
C3	1.108	52	2.5
C4	0.261	89	10
C5	0.079	97	30
Control D2	1.488	35	1.3

ASSAY CUT OFF

Cut off	u/L
Negative	< 0.4 u/L
Positive	≥ 0.4 u/L

CLINICAL EVALUATION

Clinical Specificity

139 samples from healthy blood donors were assayed in the 3rd generation **Anti R-TSH kit** ([bmd](#)). All 139 were found to be negative for TSH receptor autoantibodies.

Clinical Sensitivity

108 samples from patients with Graves's disease (treated and untreated patients) were assayed and 103 (95%) were identified as being positive for TSH receptor autoantibodies.

Lower Detection Limit

The kit negative control was assayed 50 times and the mean and standard deviation calculated. The lower detection limit at 2 standard deviations was 0.08 u/L.

Inter Assay Precision

Sample	u/L (n = 25)	CV (%)
1	5.5	8.7
2	1.6	8.8

Intra Assay Precision

Sample	u/L (n = 21)	CV (%)
1	1.3	5.5
2	5.1	4.2

Clinical Accuracy

Analysis of sera from patients with autoimmune diseases other than Graves' disease indicated no interference from autoantibodies to thyroglobulin; thyroid peroxidase; glutamic acid decarboxylase; 21-hydroxylase; acetylcholine receptor; dsDNA or rheumatoid factor.

Interference

No interference was observed when samples were spiked with the following materials; intralipid up to 30mg/mL; haemoglobin up to 5 mg/mL; bilirubin up to 0.2 mg/mL; human LH up to 10u/mL; hCG up to 160u/mL; human FSH up to 70u/mL and human TSH up to 3 u/L.

The data quoted in these instructions should be used for guidance only. It is recommended that each laboratory include its own panel of control samples in the assay. Each laboratory should establish its own normal and pathological reference ranges for TRAb levels.

SAFETY CONSIDERATION

- This kit is intended for *in vitro* use by professional persons only.
- Follow the instructions carefully.
- Observe expiry dates stated on the labels and the specified stability for reconstituted reagents. Refer to Materials Safety Data Sheet for more detailed safety information.
- Material of human origin used in the preparation of the kit have been tested and found non reactive for HIV1 and 2 and HCV antibodies and HBsAg but should, none the less, be handled as potentially infectious.
- Wash hands thoroughly if contamination has occurred and before leaving the laboratory. Sterilise all potentially contaminated waste, including test specimens before disposal.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy but these materials should be handled as potentially infectious.
- Some components contain small quantities of sodium azide as preservative.
- With all kit components, avoid ingestion, inhalation, injection and contact with skin, eyes and clothing.
- Avoid formation of heavy metal azides in the drainage system by flushing any kit component away with copious amounts of water.

SYMBOLS USED

	EC Declaration of conformity
	ELISA Test
	<i>In Vitro</i> Diagnostic Device
	Catalogue number
	Lot Number
	Expiry Date
	Number of test
	Consult Instructions
	Temperature limitation
	Biological risk
	Contains sodium azide
	Reconstitute

REFERENCES

B. Rees Smith et al

A new assay for thyrotropin receptor autoantibodies.
Thyroid 2004 14: 830-835

K. Kamijo et al

Clinical Evaluation of 3rd Generation assay for Thyrotropin Receptor Antibodies : The M22-biotin-based ELISA initiated by Smith
Endocrine Journal 2005 52: 525 – 529

ASSAY SCHEME

1. Allow all reagents and samples to reach room temperature (+20°C/+25°C) before use
2. Pipette: **75µL** Start buffer into each well (except blank)
3. Pipette: **75µL** Calibrators (starting with the highest concentration and descending to lowest), kit Controls, Patient Sera (except blank)
4. Incubate 2 hours at room temperature on an ELISA **plate shaker at 500 shakes/min**
5. Aspirate/Decant: Plate
6. Wash: Plate once **with 350 µL of diluted wash solution**
For automatic washing: Aspirate well contents by use of a plate washing machine
For manual washing: Invert and tap dry on absorbent material for manual washing.
7. Pipette: **100µL** M22-biotin into each well (except blank)
8. Incubate: 25 minutes at room temperature **without shaking**
9. Aspirate/Decant: Plate
10. Wash: Plate once **with 350 µL of diluted wash solution**
For automatic washing: Aspirate well contents by use of a plate washing machine
For manual washing: Invert and tap dry on absorbent material for manual washing.
11. Pipette: **100µL** SAPOD (diluted 1:20) into each well (except blank)
12. Incubate: 20 minutes at room temperature **without shaking**
13. Aspirate/Decant: Plate
14. For automatic washing: Wash Plate three times on automatic washer **with 350 µL of diluted wash solution.**
For manual washing: Wash twice **with 350 µL of diluted wash solution**, rinse once with pure water and dry on absorbent material
15. Pipette **100µL** TMB into each well (including blank)
16. Incubate 30 minutes at room temperature **in the dark without shaking**
17. Pipette **50µL** Stop Solution into each well (including blank) and shake for 5 seconds
18. Read absorbance at 450 nm

Do not perform the assay at temperatures above +25°C

BioMédical Diagnostics SA

Office

Actipole 25
4 bld de Beaubourg
77435 Marne la Vallée Cx2
France

Tel : 33 1 64 62 10 12
Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
Internet : www.bmd-net.com