Test QuikPac II *OneStep* Strep B

REF

SYN 10730

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Test QuikPac II *OneStep* Strep B

SYN 10730

Direct qualitative detection of group B Streptococcus (GBS) antigen from vaginal and cervical swab specimens

INTENDED USE

QuikPac II *OneStep* Strep B Test is a rapid diagnostic immunoassay for the direct qualitative detection of group B streptococcus (GBS) antigen from vaginal and cervical swab specimens, to aid in the diagnosis of streptococcus B infection. The assay may also be used as a confirmatory test for colonies of Group B streptococci grown on agar plates. This test is an *in vitro* diagnostic test for healthcare professional use only.

SUMMARY AND EXPLANATION OF THE TEST

Group B streptococci (GBS) are among the most frequent causes of life-threatening infections in neonates. Between 5% and 30% of all pregnant women are colonized with GBS¹. Several recent studies have shown that the intrapartum treatment of GBS-colonized women significantly reduces the incidence of GBS-caused sepsis.²⁻⁴ Therefore, screening for GBS is important. Standard culture methods require 24 to 48 hours, and the results may not be available soon enough for efficient treatment. Thus methods utilizing more rapid screening techniques are required.

QuikPac II *OneStep* Strep B Test is a simple, rapid and sensitive immunochromatographic assay for screening of GBS antigen from patient vaginal or cervical swab specimens. The test procedure takes less than 20 minutes and does not require special instrumentation.

PRINCIPLE OF THE TEST

The **QuikPac II *OneStep* Strep B Test** utilizes the chemical extraction of a GBS-specific antigen from the bacteria followed by the immunochromatographic assay for qualitative detection of GBS.

In the test procedure, a unique set of polyclonal antibodies are employed. One antibody is immobilized on a porous membrane while the other antibody is conjugated to dye particles as signalling component. A swab specimen from a patient is treated with Extraction Reagent to extract the antigen. The liquid extract is applied and allowed to migrate through the absorbent area and along the membrane. If GBS antigen is present, the labeled antibody-dye conjugate binds to it, forming antibody-antigen complex. As the mixture flows along the membrane, the complex is captured by the antibody immobilized in the test (marked "T") zone of the membrane, producing a visible rose-pink color band. Unbound dye-conjugated reagent is captured by the antibody immobilized in the control (marked "C") zone of the membrane.

A rose-pink line in the test zone indicates the presence of the GBS antigen. A rose-pink line in the control zone indicates that the test is working properly. When only a control line appears with no test line, the GBS antigen has not been detected and the test result is interpreted as negative.

The control line gives an added measure of quality control by demonstrating antibody recognition; assuring that the procedure was performed correctly and that the reagents are still chemically active. A desiccant is enclosed with the Testing Device to stabilize the incorporated biomolecular reagents.

KIT COMPONENTS

Testing Device TEST DEVICES	Contains dye-conjugated and membrane immobilized polyclonal antibodies
Sample Dropper DROPPERS	Disposable plastic dropper Sealed in the pouch with each Testing Device
Extraction Reagent (2x8 ml). SOLN EXTRACT	Detergent solution containing 0.1% sodium azide, in a dropper vial.
Positive Control (1 ml) CONTROL +	Suspension of non-viable group B streptococci, in a dropper vial. Contains 0.1% sodium azide.
Reaction/Extraction Cup.	Plastic w/Pour Spout, Natural
(Reaction) Cup Holder	8 Hole, White Plastic
Swab Dacron-Tipped	Plastic Shaft, Sterile-Wrapped
Product Insert	

MATERIALS REQUIRED BUT NOT PROVIDED

Clock (watch) or timer.

STORAGE AND STABILITY

Positive Control should be refrigerated between +2°C and +8°C. Other test components may be stored at room temperature (below +28°C). **DO NOT FREEZE**. Refer to expiration dates for stability.

WARNINGS AND PRECAUTIONS

1. Wear gloves while handling specimens.
2. Dispose of gloves and specimens using good microbiological practices.
3. Do not touch the swab tip at any time
4. Wash hands after performing the test.
5. Use only the sterile Dacron® swabs provided. Swabs from any other source may give erratic results.
6. Do not allow a sample swab to come in contact with any reagent vial dropper tip. Reagent or bacterial contamination may result.
7. Do not use reagents after their expiration dates.

SPECIMEN COLLECTION AND PREPARATION

1. The vaginal and cervical test specimens should be collected by a standard swab collection method. Grossly bloody samples should be avoided.
2. Patient swab may be stored or transported in a closed tube between +2°C and +8°C for up to 7 days before testing, either without media or with Modified Stuart's media (1 ml or less). Do not freeze the swabs.
3. To use QuikPac™ II OneStep Strep B Test as culture confirmation assay, collect two or three 24-hour colonies from an agar plate with a test swab. The test procedure is the same for all swab specimens.

TEST PROCEDURE

Note: Read all test instructions before running patient samples or controls.

Procedure Notes

1. Bring all specimens and controls to room temperature (between +15°C and +28°C) before testing.
2. Do not open the protective foil pouch until ready to perform the test.
3. Do not use commercial controls other than those included with the kit. They may contain additives which will interfere with test performance.

A. Extraction of GBS Antigen from Sample Swab

1. Label a Test Cup for each patient identification and place it in the Cup Holder.
2. Add **15 drops** of Extraction Reagent to the Test Cup. Place the specimen swab in the Test Cup and twirl briefly to mix ingredients. Leave the swab in the cup at room temperature for at least five minutes, but no longer than 30 minutes.
3. Twirl the swab vigorously for a few seconds, then expunge as much liquid as possible from the swab by pressing and rotating the fiber portion against the wall of the tube. Discard the swab.


The swab extract can be tested any time within the next 60 minutes by following the procedures outlined in the "Immunoassay of the Extract" section.

B. Control Testing Procedure

1. Add **8 drops** of Extraction Reagent and **one drop** of Positive Control to the Test Cup.
2. Swirl the cup gently to mix contents and leave for five minutes at room temperature.
3. Complete the test according to the "Immunoassay of the Extract" section below.


C. Immunoassay of the Extract

1. Remove a Testing Device and a Sample Dropper from the foil pouch. Discard the foil pouch and desiccant. Place the device on a level surface.
2. Using the Sample Dropper, add exactly **10 drops of the liquid extract** from the Extraction Cup to the well marked "S" of the Testing Device.
3. Read the test results at 10 minutes.

IMPORTANT: To avoid incorrect readings, do not interpret the test results after more than 10 minutes.

INTERPRETATION OF RESULTS

1. **Positive.** Two rose-pink color bands appear: one in the Control Zone ("C") and one in the Test Zone ("T"). GBS-specific antigen detected and the sample should be considered positive for GBS.
2. **Negative.** One rose-pink color band appears in the Control Zone ("C"), with no apparent color band in the Test Zone ("T"). GBS antigen is not present in the test sample at the level of sensitivity of the test.
3. **Invalid.** No color bands appear, or a band is visible in the Test Zone, but not in the Control Zone. The test is invalid. The specimen should be retested using new QuikPac II OneStep Strep B Test components.

Note: There is no meaning attributed to line color intensity and width.

QUALITY CONTROL

Internal Controls

QuikPac™ II OneStep Strep B Test contains built-in quality control features. Distinctive pink-rose band in control zone indicates that the test was performed correctly, and that all test components were in good functional condition.

External Controls

Good laboratory practice recommends the use of control reagents. Positive Control containing GBS antigen is provided in the test kit to verify test performance.

LIMITATIONS OF THE TEST

1. The test is limited to detection of GBS antigen in patient swab specimens or from colonies grown on agar plates.
2. For *in vitro* diagnostic use only.
3. Although the test is very accurate, a low incidence of false results may occur.
4. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
5. The sensitivity (limit of detection) value of the test is 5×10^5 CFU (colony forming units) per ml on the 50 μ l test sample. Positive results are obtained when the amount of GBS antigen in the test corresponds to this value or above. Conversely, negative results are obtained when the amount of the antigen is below the sensitivity of the test.

PERFORMANCE CHARACTERISTICS

Relative Sensitivity and Specificity

The performance of QuikPac II OneStep Strep B Test was compared to that of a conventional culture method in a study of 101 clinical specimens. Vaginal and cervical swab specimens were obtained during routine examinations at physician's offices and in outpatient clinical laboratories. The specimens were evaluated by assaying with QuikPac II OneStep Strep B Test and by a conventional culture method.

The results of the comparative study are shown below:

		QuikPac II OneStep Strep B	
		Positive	Negative
Culture	Pos. (50)	46	4
Method	Neg. (51)	5	46

Based on the comparison with a conventional culture method for detection of GBS from swab specimens, relative sensitivity of **QuikPac II OneStep Strep B** Test is 92.0% (46/50), and relative specificity is 90.2% (46/51).

Analytical Sensitivity

QuikPac II OneStep Strep B Test allows for the detection of GBS of serotypes Ia, Ib, Ic, II and III. A study of analytical sensitivity of the test was carried out by culturing serial dilutions of a logarithmic phase GBS broth culture, and by assaying the serial dilutions with three different lots of **QuikPac II OneStep Strep B** Test reagents. Comparing the colony counts of the cultures and the assay results yielded the sensitivity (limit of detection) value of 5×10^5 CFU (colony forming units) per ml, or 2.5×10^4 CFU/test, based on the 50 µl test sample.

Cross-Reactivity

Cross-reactivity studies were performed with a variety of non-GBS microorganisms. The microorganisms listed below tested negative at 1×10^8 CFU/test.

- *Candida albicans* (ATCC 14053)
- *Corynebacterium diphtheriae* (ATCC 9015)
- *Escherichia coli* (ATCC 25922)
- *Haemophilus influenzae* (ATCC 35056)
- *Klebsiella pneumoniae* (ATCC 13883)
- *Neisseria gonorrhoeae* (ATCC 9826)
- *Pseudomonas aeruginosa* (ATCC 27853)
- *Staphylococcus aureus* (ATCC 29213, 25923)
- *Staphylococcus epidermidis* (ATCC 12228)

- *Streptococcus, group A* (ATCC 19615)
- *Streptococcus, group C* (ATCC 12388)
- *Streptococcus, group D* (ATCC 12389)
- *Streptococcus, group F* (ATCC 12393)
- *Streptococcus, group G* (ATCC 12394)
- *Streptococcus pneumoniae* (ATCC 9163, 6303, 10015)

Reproducibility

In a study of reproducibility of the test, serial dilutions of a GBS broth culture in a logarithmic growth phase were used to prepare masked and coded specimen swabs. The dilutions were also assayed by colony counting. **QuikPac II OneStep Strep B** Test assays of the swabs were performed at three separate test sites, using three different lots of the test reagents. Specimens obtained from dilutions below the detection limit (5×10^5 CFU/ml) tested negative in all assays.

Specimens at or above the detection limit ($\geq 5 \times 10^5$ CFU/ml) tested positive in all assays. The above results demonstrate excellent reproducibility of the **QuikPac II OneStep Strep B** Test.

SYMBOLS USED

	Catalog number		Read instructions for use
	Batch code		EC Declaration of Conformity
	Number of tests		Use by
	<i>In vitro</i> diagnostic device		Storage temperature limitation
	Rapid test		Do not reuse
		Contains sodium azide	

BIBLIOGRAPHY

1. Finch, R.G., French, G.L., and Phillips, I. Group B streptococci in the female genital tract. *Br. Med. J.* 1: 1245-1247, 1976.
2. You, M.D., Mason, E.O., Leeds, L.J., Thompson, P.K., Clark, D.J., and Gardner, S.E. Ampicillin prevents intrapartum transmission of group B streptococcus. *JAMA* 241: 1245-1247, 1979.
3. Boyer, K.M., and Gotoff, S.P. Prevention of early-onset neonatal group B streptococcal disease with selective intrapartum chemoprophylaxis. *N. Engl. J. Med.* 314: 1665-1669, 1986.
4. Lim, D.V., Morales, W.J., Walsh, W.J., and Kazanis, D. Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemoprophylaxis. *J. Clin. Microbiol.* 23: 489-492, 1986.

QUICK OPERATING PROCEDURE

	Specimen	Positive control	Incubation time
Extraction	Add to the test cup: 15 drops of Extraction Reagent Immerse the specimen swab into the test cup Twirl briefly to mix and leave the swab in the cup <i>(no longer than 30 minutes)</i>	Add to the test cup: 8 drops of Extraction Reagent + 1 drop of Positive control Swirl the cup gently to mix	Leave for 5 minutes at room temperature
	Twirl the swab vigorously and press the swab against the wall of the tube to extract as much liquid as possible <i>(the swab extract can be tested within 60 minutes)</i>		
Immunoassay of the extract	Add 10 drops of the liquid extract from the extraction cup to the well marked "S" of the testing device		Leave for 10 minutes at room temperature
Reading	Interpret the test immediately		

BioMédical Diagnostics SA

Office
 Actipole 25
 4 bld de Beaubourg
 77435 Marne la Vallée Cx2
 France

Tel : 33 1 64 62 10 12
 Fax : 33 1 64 62 09 66

E-mail : support@bmd-net.com
 Internet : www.bmd-net.com